


## Thin-film Gold Electrodes

FMS-009



Electrochemical Technique	Suitable
Chronoamperometry	✓
Cyclic Voltammetry	✓
Square Wave Voltammetry	✓
Differential Pulse Voltammetry	✓
Open Circuit Potentiometry	✓
Electrochemical Impedance Spectroscopy	✓
Others	Enquire

## Description

The FMS-009 electrode features a sputtered gold working electrode and counter electrode on a polyester substrate. It is configured as a three-electrode system, incorporating a screen-printed Ag/AgCl reference electrode. These electrodes are designed for single use and are intended exclusively for research and development purposes.

Table 1: Technical specification for FMS -009 electrodes

Technical Specifications	
Working Electrode Material	Sputtered Gold
Counter Electrode Material	Sputtered Gold
Reference Electrode Material	Ag/AgCl
Substrate Material	PET
Conductive Track Material	Sputtered Gold
Working Area	12.6 mm <sup>2</sup>
Recommended Sample Volume	50 – 100 µL

## Sample Performance Data

The performance data presented in this document represents typical expected results for FMS-009 electrodes within the same batch. For detailed information regarding potential variations between batches, please contact the manufacturer.

### Method of Analysis


Analysis was performed by cyclic voltammetry using the settings outlined in table 2 with a mediator solution consisting of 5mM potassium ferricyanide/ferrocyanide in 10 mM PBS pH 7.4.

*Table 2: Analysis settings*

### Results

The voltammograms below show 3 successive scans on the same electrode, the first scan has been removed.

Setting	Value
E begin	0.1 V
E vertex1	0.6 V
E vertex2	-0.5 V
E step	0.01 V
Scan rate	50 mV/s
Number of scans	4



*Figure 1: 3 successive CV scans of mediator solution on a FMS-009 electrode.*

## Sample Performance Data

### Typical Electrode Performance

Table 3: Expected potassium ferricyanide / ferrocyanide oxidation & reduction peak values for FMS-009 electrodes.

Parameter	Peak Height (μA)	Peak potential (V)
Oxidation Peak	139.4	0.21
Reduction Peak	-129.4	0.10

### Repeatability

Table 4: Oxidation & reduction peak variability data from 20 electrodes within the same batch

Parameter	Intra – electrode Variability (%CoV)	Inter – electrode Variability (%CoV)
Oxidation Peak	0.9	2.2
Reduction Peak	0.6	2.4

For more information contact: [info@FlexMedical-Solutions.com](mailto:info@FlexMedical-Solutions.com)