




Version 5.11

User Manual



24/10/2024

If you have problems

First make sure to check the “Troubleshooting” section in this document and the Knowledge Base on our website: www.palmsens.com/knowledgebase/

This page contains support information on installation, software updates, and training.

Please make sure your software and firmware are up-to-date.



In case of persistent problems:

Use the contact form: www.palmsens.com/contact/

Give us a call: [+31 30 2459211](tel:+31302459211)

Or send an email: info@palmsens.com

Try to describe the problem as detailed as possible. Sending us the relevant method files, data files and screenshots can be helpful.

Please include your instrument model and serial numbers, as well as any applicable software and firmware version you are using.

Disclaimers

PalmSens BV cannot guarantee that its instruments will work with all computer systems, operating systems, and third-party software applications hardware/software. The information in this manual has been carefully checked and is believed to be accurate as of the time of compiling. However, PalmSens BV assumes no responsibility for errors that might appear.

See Appendix A for CE declarations of conformity.

Copyrights

PSTrace5 User Manual. Copyright © 2024, PalmSens BV, all rights reserved.

PSTrace, MultiTrace, PalmSens, EmStat, MultiEmStat and Sensit are trademarks of PalmSens BV.

No part of this document may be copied or reproduced in any form without the prior written consent of PalmSens BV.



Contents

1	Getting Started	1
1.1	Software requirements	1
1.2	Making a reliable USB connection	1
1.3	Which model do I have?	1
1.4	Installation	4
1.5	PSTrace basic principles	4
1.6	Configuring PSTrace	6
1.7	Connecting using a USB cable	12
1.8	Virtual EmStat Pico	14
1.9	Connecting using Bluetooth	14
2	First measurements	17
2.1	Which dummy cell do I have	17
2.2	First measurement on the TestSensor	19
2.3	First measurement on the Dummy Cell	20
2.4	First EIS measurement on the Dummy Cell	21
2.5	Readings	22
3	Measurements	25
3.1	Connecting a cell to the potentiostat	25
3.2	Setting up a measurement	29
3.3	Standard measurement sequence	32
3.4	Measuring versus OCP	33
3.5	Resolution and optimal range selection	34
3.6	Running a measurement	35
3.7	Peaks and levels	39
3.8	Multiplexer	39
3.9	Recording an auxiliary input	41
3.10	Recording the cell potential	42
3.11	Recording the WE potential	43
3.12	Recording the WE current	43
3.13	BiPotentiostat	44
3.14	Blank subtraction	45
3.15	Ohmic (iR) compensation	46
3.16	Potential and current limits	48
3.17	Override bandwidth filter	48
3.18	Digital triggers	49
3.19	Measuring the noise level of the instrument	49
3.20	Reducing noise	51
4	Electrochemical techniques	53
4.1	Linear Sweep Voltammetry (LSV)	54
4.2	Cyclic Voltammetry (CV)	57
4.3	Fast Cyclic Voltammetry (FCV)	61
4.4	AC Voltammetry (ACV)	63
4.5	Differential Pulse Voltammetry (DPV)	65

4.6	Square Wave Voltammetry (SWV).....	69
4.7	Normal Pulse Voltammetry (NPV)	72
4.8	Chronoamperometry (CA)	75
4.9	Zero Resistance Amperometry (ZRA)	78
4.10	MultiStep Amperometry (MA)	80
4.11	Fast Amperometry (FAM)	82
4.12	Pulsed Amperometric Detection (PAD)	84
4.13	Multiple Pulse Amperometry (MPAD).....	86
4.14	Open Circuit Potentiometry (OCP).....	88
4.15	Chronopotentiometry (CP)	89
4.16	Linear Sweep Potentiometry (LSP).....	91
4.17	MultiStep Potentiometry (MP).....	93
4.18	Stripping Chronopotentiometry (SCP)	96
4.19	Chronocoulometry (CC)	98
4.20	Mixed Mode (MM)	101
4.21	MethodSCRIPT Sandbox	107
5	Limitations and extra options for EmStat Pico	109
5.2	Extra options for EmStat Pico and Sensit Series	111
6	Electrochemical Impedance Spectroscopy	113
6.1	Introduction	113
6.2	Setting up an impedance measurement.....	117
6.3	Parameters for EIS.....	119
6.4	Parameters for Galvanostatic EIS	123
6.5	Parameters for Fast EIS and Fast GEIS.....	125
6.6	Running an EIS measurement	126
6.7	Optimizing EIS Parameters for Improved Results and Troubleshooting	128
6.8	Fitting measured EIS data in PSTrace	129
6.9	Export for analysis and circuit fitting	129
6.10	Limitations for EIS on EmStat Pico	132
7	Equivalent Circuit Fitting	135
7.1	Overview.....	135
7.2	Opening the circuit editor	136
7.3	Switching between Edit mode and Fit mode	136
7.4	Building a circuit in the graphical editor	137
7.5	Building a circuit with the CDC editor	140
7.6	Loading and saving circuits.....	145
7.7	Fitting or simulating a circuit	146
7.8	Batch-fitting a circuit on multiple curves.....	153
7.9	Overview of circuit components	153
7.10	Fitting Example	165
7.11	Fitting algorithm	172
8	Plot, curves and data	175
8.1	Handling curves.....	175
8.2	Session data	176
8.3	Switching between plots.....	178
8.4	Plot toolbar	179
8.5	Curve toolbar.....	184

8.6	Separating overlapping peaks	192
8.7	Curve operations	197
8.8	Curve calculations.....	203
8.9	Saving data.....	205
8.10	Exporting curves.....	207
9	Manual Control	209
9.1	Manual Control tab.....	209
10	Analytical mode	213
10.1	Supported techniques	213
10.2	Setting up parameters	214
10.3	Running a measurement.....	217
10.4	Result analysis	218
10.5	Example data files	220
11	Corrosion mode	221
11.1	Supported techniques	221
11.2	Running a corrosion measurement	222
11.3	Result analysis	223
11.4	Impedance Spectroscopy	229
11.5	Exporting results.....	232
11.6	Example data files	233
12	Files	235
12.1	File types.....	235
12.2	Loading and saving data	236
12.3	Exporting data to other file formats	240
13	Scripting	243
13.1	Features	243
13.2	Measure command.....	246
13.3	FastMode command.....	247
13.4	Cell command.....	247
13.5	SetCurrent command	248
13.6	ReadCurrent command	249
13.7	SetPotential command	250
13.8	ReadPotential command	250
13.9	Wait command.....	251
13.10	Repeat command	251
13.11	FindPeaks command.....	253
13.12	SetChannel command	253
13.13	NextChannel command.....	254
13.14	PrevChannel command	254
13.15	Stirrer command.....	254
13.16	SetAnalog command	254
13.17	ReadAuxiliary command	254
13.18	SetDigitalIO command	255
13.19	ReadDigitalIO command	255
13.20	Override parameter command.....	255

13.21	Running a script from the command line	257
14	MethodSCRIPT	259
14.1	About MethodSCRIPT	259
14.2	Generating MethodSCRIPT in PSTrace	259
14.3	Running a MethodSCRIPT	263
15	MUX8 and MUX16 multiplexers	265
15.1	MUX8 multiplexer	265
15.2	MUX16 multiplexer	270
16	MUX8-R2 multiplexer	275
16.1	MUX8-R2 functional diagram	276
16.2	Specifications	276
16.3	Configurations	276
16.4	Connections	277
16.5	MUX8-R2 Pin-outs	278
16.6	Stacking multiple MUX8-R2	280
17	PalmSens1 and PalmSens2	281
17.1	Description	281
17.2	Operating	281
17.3	Using a USB to serial converter	283
17.4	Auxiliary port pin-out	283
17.5	Cell connector pin-out	284
17.6	Specifications	284
17.7	BiPot specifications	285
17.8	Offset calibration	286
17.9	Battery replacement	286
18	PalmSens3	289
18.1	Description	289
18.2	Operating	289
18.3	Auxiliary port pin-out	291
18.4	Cell connector pin-out	292
18.5	Specifications	292
18.6	BiPot specifications for PalmSens3	293
18.7	Replacement of PalmSens1 or PalmSens2 battery	294
18.8	EIS Calibration for PalmSens3	297
18.9	Offset calibration	297
18.10	Performing a hard reset	297
18.11	Battery replacement	299
19	PalmSens4	303
19.1	Description	303
19.2	Operating	303
19.3	Display and power indicator	304
19.4	iR Compensation	305
19.5	Auxiliary port pin-out	305
19.6	Cell connector pin-out	307

19.7	Self Diagnostics.....	307
19.8	Performing a hard reset.....	309
19.9	Battery replacement.....	312
19.10	Specifications.....	312
20	EmStat instruments series	315
20.1	EmStat-based instruments overview	315
21	EmStat1 and EmStat2	319
21.1	Description	319
21.2	Operating	319
21.3	Status LED indicators	319
21.4	Cell connector pin-out	320
21.5	Specifications.....	320
21.6	Offset calibration	321
22	EmStat3 and EmStat3+	323
22.1	Description	323
22.2	Operating	323
22.3	Status LED indicators	324
22.4	Cell connector pin-out	324
22.5	Specifications.....	324
22.6	Offset calibration	325
23	EmStat3 Blue and EmStat3+ Blue	327
23.1	Description	327
23.2	Operating	327
23.3	Status LED indicators	328
23.4	Cell connector pin-out	328
23.5	Auxiliary port pin-out	329
23.6	Specifications.....	330
24	EmStat4S	331
24.1	Description	331
24.2	Operating	331
24.3	Status LED indicator.....	332
24.4	Cell connector pin-out	332
24.5	Specifications.....	333
25	EmStat4X	335
25.1	Description	335
25.2	Operating	335
25.3	USB power and charging.....	336
25.4	Display and status LED indicator	336
25.5	iR Compensation	337
25.6	Auxiliary port pin-out	338
25.7	Cell connector pin-out	339
25.8	SPE connector.....	339
25.9	Specifications.....	340

26	Sensit series	341
26.1	Description	341
26.2	Operating the Sensit Smart	342
26.3	Operating the Sensit BT	342
26.4	Sensit BT trigger button.....	343
26.5	Sensit Smart and Sensit BT specifications	344
26.6	Calibration	345
27	EmStat Pico based instruments	347
27.1	EmStat Pico specifications	347
27.2	Calibration	347
28	Troubleshooting	349
28.1	Running self-diagnostics	349
28.2	Verifying your potentiostat	349
28.3	Noise.....	350
29	Updating firmware	353
29.1	Firmware update.....	353
30	Maintenance and compliance	355
30.1	Lithium-Polymer batteries	355
30.2	Temperature compliance	355
30.3	Humidity compliance	355
30.4	Temperature drift	355
30.5	Atmospheric pressure	355
30.6	Cleaning.....	356
30.7	Periodic calibration and preventive maintenance	356
30.8	Service and repair	356
30.9	RoHS Compliance.....	356
Appendix A.	CE Declarations of conformity	357
Appendix B.	EU Waste Electrical and Electronic Equipment (WEEE) Directive	364
31	Index.....	365

1 Getting Started

1.1 Software requirements

System requirements for PSTrace:

- Windows 10 or 11 (64-bit)
- 2 GB RAM
- Minimum display resolution: 1366 × 768 pixels

Instrument requirements:

- a single channel or multiplexed instrument from PalmSens BV;
- or a multi-channel instrument from PalmSens BV (for controlling a single channel).

Required for test measurements:

- a sensor cable with crocodile clips and the PalmSens Dummy Cell.

1.2 Making a reliable USB connection

For a reliable and steady USB connection between your PC and the instrument, it is crucial to consistently employ the USB cable provided with the instrument by PalmSens BV. It is strongly advised against the use of USB cable extenders. Similarly, if feasible, it is recommended to refrain from utilizing a USB hub or docking station.



Please make sure to always use the USB cable that was supplied with the instrument by PalmSens BV.

1.3 Which model do I have?

The instruments listed below can be used with PSTrace.



Multi-channel instruments like MultiEmStat, MultiPalmSens and EmStat3-4WE are used with the MultiTrace software.



EmStat1
(sold until 2010)



EmStat2
(sold until 2013)



EmStat3 or EmStat3+
(sold until 2022)



EmStat4S LR or EmStat4S HR



EmStat3 or EmStat3+ Blue
(sold until 2023)



EmStat4X LR or EmStat4X HR



EmStatGo / EmStat4R



EmStatMUX8 or EmStatMUX16
(sold until 2017)



EmStatMUX8-R2



*PalmSens1 or 2
(sold until 2013)*



*PalmSens3
(sold until 2017)*



PalmSens4



Sensit Smart, based on EmStat Pico



Sensit BT, based on EmStat Pico



Sensit Wearable, based on EmStat Pico



EmStat Pico Development Board



EmStat4M Development Board

1.4 Installation

Install PSTrace by running the setup.exe from the media that came with the instrument or, when downloaded, from the location where the downloaded ZIP file was extracted.

After installation, a folder 'PSData' is created in your 'My Documents' folder. This folder contains some example method and data files. A PSTrace shortcut is placed in the Start menu and on the desktop.

1.5 PSTrace basic principles

The main window of PSTrace shows a Method Editor on the left-hand side which contains all measurement parameters and related information and a Plot area on the right-hand side which contains everything post-measurement related, including measurement results and analytical tools.




Figure 1 PSTrace main window

1.5.1 Methods

PSTrace uses 'Methods' as a starting point for a measurement. A Method contains all measurement parameters like which *Technique* (Linear Sweep Voltammetry, Square Wave Voltammetry, Electrochemical Impedance Spectroscopy, etc.) is used and information about post-measurement actions like data smoothing and peak searching. All these parameters can be edited in the Method Editor which is found at the left-hand side of the PSTrace window.

The parameters can be saved to and loaded from a '.psmethod' file using the menu: 'Method'. These files do not include any data and only contain the measurement parameters settings as shown in the Method Editor.



Switching between techniques or changes parameters in the Method Editor has no effect on anything that has already been measured and is displayed in the Plot.

1.5.2 Measurement data

As soon as a measurement is started a new 'Measurement' node appears in the Legend containing the default Curve for the corresponding technique. The Measurement contains the following information:

- Method parameters (as was defined in the Method Editor)
- Raw measurement data
- One or more curves

By clicking the Measurement in the Legend, the following window is shown:




Figure 2 Pop-up window shown when clicking a Measurement in the Legend

The window allows to change the Measurement title and notes directly. The notes shown here are a copy from the original Method used for this measurement. Also, the window allows to add new curves to the Plot based on the available Measurement data.

The "View method parameters" button shows the original Method used for this Measurement. This Method is a copy from the original Method defined in the Method Editor that was used with the selected Measurement.

1.5.3 Curves

A measurement can contain one or more Curves. By clicking a Curve in the Legend a window is shown with information about the Curve. It allows you to directly change the title of the Curve as shown in the Legend or change its appearance or view the Data used for this curve.




Figure 3 Pop-up window shown when clicking a Curve in the Legend

1.5.4 Saving your data

All available Measurement data and curves as well as the Method as shown in the Method editor can be saved to a single 'Session' file (.pssession). Use the menu 'Data' to save and load Sessions. Any titles changed or customised Curve appearances like colour and symbols used are saved as well.

1.6 Configuring PSTrace

Before you start using PSTrace it is important that all settings are determined. The Settings window is always shown at the first start of PSTrace.

The Settings window can also be found in the menu 'Tools' → 'General settings...'.




Figure 4 General settings in the Tools menu




Figure 5 Settings window

1.6.1 Program modes

PSTrace can run in three different modes.

The mode can also be set using the button in the main window at the left corner.




Figure 6 Mode-select buttons

1.6.1.1 Scientific mode

PSTrace operates primarily in Scientific mode, which serves as its default setting. This mode encompasses all the electrochemical techniques supported by the connected instrument and general tools designed for curve analysis. Users can engage in both automated and manual peak searches, facilitating the determination of peak potential, peak height, peak area, and peak width. Furthermore, users have the capability to conduct curve subtraction and apply linear regression to designated curve segments.

1.6.1.2 Analytical mode

The Analytical mode extends PSTrace for Voltammetric Analysis providing the possibility to do quantitative analysis by means of:

- Standard addition
- Using a calibration curve

See also section: [Analytical mode](#) on page 213.

1.6.1.3 Corrosion mode

The Corrosion mode of PSTrace translates the supported techniques to naming conventions generally used with Corrosion Analysis. Secondly a tab with analytical tools is added to the user interface for:

- Linear polarization, from which the polarization resistance is obtained;
- Tafel plots, from which the corrosion rate is obtained;
- and Impedance data analysis by means of equivalent circuit fitting for determining values like polarization resistance and corrosion rate.

See also section: [Corrosion mode](#) on page 221.

1.6.2 Hardware configuration

The tab 'Hardware' allows for setting the connected hardware configuration, so that the PSTrace interface shows the correct options.

If the connected device does not have the option or functionality you select, the option will appear, but no measurement / data will be possible with the respective tab. Thus, you can let it checked if you are using alternately more than one device with different hardware options.

1.6.2.1 Mains frequency

The Mains Frequency is used by PalmSens or EmStat to eliminate noise induced by electrical appliances.

1.6.2.2 BiPotentiostat module (BiPot)

If PalmSens is equipped with a BiPot module, a second WE (Working Electrode) can be used. Enabling the BiPot will add a new tab to the Method Editor for settings the BiPot parameters and it is supported in the Manual Control tab.

See also section: [BiPotentiostat](#) on page 44.

1.6.2.3 Stirrer

PalmSens or EmStat can be used with a stirrer which is activated during the conditioning and deposition phase of a measurement. It can also be controlled manually in the Manual Control tab. The stirrer is normally controlled using a Switchbox from PalmSens BV. In case a stirrer from BVT is used, its speed can be specified in the Manual Control tab.

See also section: [Manual Control tab](#) on page 209.

1.6.2.4 Multiplexer

Check the 'Multiplexer present' checkbox if you are using a MUX8, MUX16 or MUX8-R2 multiplexer. All multiplexers can either be in the form as an accessory with the potentiostat or have an integrated EmStat potentiostat.

Use the button 'Help select' to see which model you are using.

1.6.2.5 Comm options

Check the 'Bluetooth present' checkbox if you are using an instrument that is Bluetooth enabled. Enabling this checkbox adds a small button with a Bluetooth icon in the Connection box, found in the upper left corner of PSTrace.

See also section: [Connecting using Bluetooth](#) on page 14.

1.6.2.6 Auxiliary input

Some instruments have an auxiliary input port that can be used to measure an external signal simultaneously with a measurement. The standard accessories provide by PalmSens BV that make use of the auxiliary port are listed when clicking the 'Change aux. input' button.




Figure 7 The auxiliary input selection window

A custom auxiliary input can also be added, for example the output of a spectrometer. If the analog output of any external device is linear, it can be translated to any given unit by setting the slope, offset and unit information. Alternatively, the 'Calibrate' button for calibration can be used to determine the offset and slope.




Figure 8 Calibrate the selected auxiliary input

Calibration can be done by setting two points to determine both offset and slope of the linear relation, or by just adjusting the offset.

A temperature sensor for example can be calibrated with a high precision thermometer used in a low temperature and high temperature medium to enter two different values, e.g. room temperature and in boiling water at around 100 degree Celsius. By using the 'set point 1' and 'set point 2' buttons both the offset and slope will automatically be determined based on the given values.

1.6.3 Plot and data

The 'Plot and data' tab shows options on how to display and export measurement data.




Figure 9 Settings window – Plot and data tab

1.6.3.1 Preferred number format

In case the local decimal separator is a comma instead of a point, the checkbox 'Use local default decimal separator' can be used to override this.

1.6.3.2 Peak search options

- 'Automatic peak search': If checked, automatic peak search is done on measured data or when data is loaded from a file.
- 'Use alternative peak search algorithm': If checked, a different peak search algorithm is used which performs better on curves on a steep slope with no clear valleys.




Figure 10 Peaks found with the Alternative Peak Search option enabled.

1.6.3.3 Enable options for selecting reference electrodes

The checkbox 'Enable options for E corrections' adds options to specify the reference electrode used to the user interface of PSTrace. These new options can be found in the Method Editor and in the Plot toolbar:




Figure 11 Set reference electrode used in the Method Editor.




Figure 12 Extra toolbar added in plot with options to plot curves with potential values corrected for the reference electrode used.

1.6.3.4 Data visualization convention




Figure 13 Plot options for axes behavior

There are two common conventions for Voltammetric plots both of which are supported in PSTrace:

- The IUPAC (International Union of Pure and Applied Chemistry) makes sure that anodic currents are plotted in positive direction and cathod currents in negative direction.
- If the The Polarographic convention (also known as the North American, Polarographic, or Classic convention) is used, anodic currents will be plotted in negative direction and cathodic currents in positive direction.

1.6.3.5 Export data to ZView



Please take note that PSTrace offers an advanced built-in EIS spectrum analysis tool. The use of third-party tools such as ZView or EIS Analyser is entirely optional.

Measured EIS data can be exported with one click to ZView. The location for ZView is detected automatically. If for some reason PSTrace cannot detect the location of ZView automatically, browse to the location of the .exe file manually. Please note that some older versions of ZView may not be compatible. However, you can still export the data through the 'Data' menu.

See also section: [ZView](#) on page 130.

1.6.3.6 Export data to EIS Spectrum Analyser

Measured EIS data can be exported with one click to the free program EIS Analyser. To use the one-click-export button next to the plot, the location of this program needs to be set manually.

See also section: [EIS Spectrum Analyser](#) on page 129.

1.6.3.7 Export data to Origin

Measured data can be exported with one click to Origin. To use the one-click-export button next to the plot, the location of this program needs to be set manually. Please note that some older versions may not be compatible. However, you can still export the data through the 'Data' menu.

See also section: [Exporting Curves](#) on page 207.

1.7 Connecting using a USB cable

Make sure your instrument is powered on and connected to the PC using the USB cable that was supplied. Search for your instrument name on the displayed list and select it. However, note that certain first-generation PalmSens instruments may not be listed by their names but rather as a COM port. In such cases, choose the corresponding COM port associated with your instrument.





Figure 14 Connecting with a device

If a USB cable was plugged in or switched on after PSTrace was started, use the 'Refresh' button  to reinitialize the list.




PSTrace displays all available COM ports, inclusive of those not assigned to a PalmSens potentiostat. Attempting to establish a connection with unassigned COM ports or those associated with other devices may result in a connection error.

1.7.1 Using a USB to serial converter



This section only applies to PalmSens1 and PalmSens2 or an OEM product.

In case a USB – RS232 adapter is used; the driver of the adapter has to be installed according to the manual of the adapter. The COM number used can be found in Windows' Device Manager. This adapter requires a null-modem adapter (female-female) between the cable and PalmSens.

To open the Device Manager quickly, press the Windows key  and the Pause/Break key on the keyboard of your computer simultaneously. Click the Hardware tab, and then click Device Manager. The installed COM port and its number are listed under Ports.

Specify this number in the COM textbox in PSTrace if necessary.




Figure 15 Make sure the right COM port is used in PSTrace

1.8 Virtual EmStat Pico

This option is accessible in the instrument list, allowing for the development of MethodSCRIPT™ scripts without the need for a physical instrument. When connected, the software can generate MethodSCRIPT scripts for your methods, which can be viewed by clicking on the ‘Show MethodSCRIPT’ button below the Method Editor.

See also section: [MethodSCRIPT Sandbox](#) on page 107.

It's important to note that the virtual instrument lacks support for additional features available in our physical instruments, such as running measurements and controlling the cell. In instances where these actions do not result in an error message, the results are likely to be trivial.

1.9 Connecting using Bluetooth

To make a connection using Bluetooth, please follow these steps:

1. In case of a PalmSens1 or 2; make sure your PalmSens1 or 2 is Bluetooth ready. This can be checked by a marking with ‘BT’ at the bottom of the instrument.
2. Connect the provided Bluetooth dongle to the AUX port of the instrument or use an instrument with built-in Bluetooth capabilities, like the PalmSens4, Sensit BT or EmStat4X.
3. Make sure your PC is equipped with Bluetooth capabilities.
4. Make sure the ‘Bluetooth present’ is checked in the Settings window of PSTrace (see menu ‘Tools’ → ‘General settings...’)




Figure 16 “Bluetooth present” check box

5. Click the small button with the Bluetooth icon in the Connection box at the upper left corner of the main window.




Figure 17 Search Bluetooth devices-button in the Connection box



If the device is not listed, wait briefly, and click the 'Refresh' button again. It can take some time before Windows has initialized all available Bluetooth devices.

6. Select the Bluetooth instrument you want to connect to. The name consists of the letters PS, followed by the last four characters of the MAC address: PS-####. This address can usually be found on the instruments display or printed on the bottom of the housing.
7. Click 'Connect'. If this is the first time the PC is connecting to the Bluetooth extension, Windows will ask for pairing.

In Windows 7 the message appears "A Bluetooth device is trying to connect."

In Windows 8, 10 or 11, the message appears "Add a device".




Figure 18 Message appearing in Windows 8, 10 or 11

Click on this message and click "Allow". When prompted for a code, enter 1234.



Use pairing code: 1234

8. If the pairing procedure took too long, the connection might have timed out. In this case repeat the last few steps starting with step 5.
9. If the connection fails, try to remove the respective device PS-#### from the Bluetooth list, and start over. Pairing Bluetooth devices on Windows may present challenges, and it's not uncommon to encounter certain issues during this process. It's important to note that these challenges primarily stem from the Windows operating system itself and are not indicative of any problems with your device.

2 First measurements

To get acquainted with the instrument, a TestSensor (until late 2017) or Dummy Cell (since late 2017) is supplied to perform reproducible measurements.

2.1 Which dummy cell do I have

Before proceeding, please make sure to verify which dummy cell you have. See the next two images below for reference.




Figure 19 *PalmSens TestSensor A (shipped until late 2017)*




Figure 20 *PalmSens Dummy Cell*

2.1.1 PalmSens TestSensor A

The old PalmSens 'TestSensor A' has three stripped wires to connect to the electrodes using the provided croc clips. This test sensor simulates a RedOx circuit. The circuit contains two diodes in series in opposite direction with a cap in parallel and a series resistor. It simulates a non-diffusion limited RedOx system with a formal potential of 0 V. This circuit delivers *i*-E-curves, for example with Linear Sweep Voltammetry, that show an S-shape. That means the curve follows the Nernst equation. This circuit is prone to noise and thus can be used to investigate the noise in your measuring environment.

2.1.2 PalmSens Dummy Cell

The PalmSens Dummy cell offers three different circuits to check your system. To use the dummy cell, connect the croc clips to the corresponding connection pads on the dummy cell. The counter electrode (CE, black) and the reference electrode (RE, blue) are always connected to the same pads labelled CE and RE. The working electrode (WE, red) should be combined with the sense electrode (S, yellow/white) if present. For connecting the WE there are three different options.




Figure 21 Schematics of the PalmSens Dummy Cell circuits

- The WE A pad connects the working electrode to the same RedOx simulating circuit as found on the old TestSensor A. This circuit contains two diodes in series in opposite direction with a cap in parallel and a series resistor. It simulates a non-diffusion limited RedOx system with a formal potential of 0 V. This circuit delivers i-E-curves, for example with Linear Sweep Voltammetry, that show an S-shape. That means the curve follows the Nernst equation. This circuit is prone to noise and thus can be used to investigate the noise in your measuring environment.
- The WE B pad connects the working electrode to a 0.1% accurate 10 k Ω resistor. This resistor delivers a straight line as i-E-curves, for example as a result of a Linear Sweep Voltammetry measurement, following Ohm's law.
- The WE C pad connects the working electrode to a simplified Randles circuit. This circuit will show in a Nyquist plot of an Electrochemical Impedance Spectroscopy (EIS) experiment a semi-circle, which is characteristic for an RC system (resistor and capacitor in parallel).

2.2 First measurement on the TestSensor

This section describes how to do a first Linear Sweep measurement with any PalmSens BV instrument (EmStat or PalmSens-series) using the PalmSens TestSensor.

Please follow these steps if you have a PalmSens TestSensor as shown in the picture below.




Figure 22 PalmSens TestSensor A

1. Make sure the instrument is on and connected to the PC
2. Open PSTrace and connect to the instrument
3. Load the following method file (menu: 'Method' → 'Load'):
'User Name\Documents\PSData\PSDiffPulse.psmethod'
4. Connect the TestSensor using the croc clips as shown in the picture above.
5. Start measurement by clicking the start button in the toolbar.




Figure 23 The Start measurement button in the toolbar

6. See the picture below for the expected results.




Figure 24 PSTrace showing the expected results.



The TestSensor gives a current response in the order of nano amperes. This means that proper shielding is of significant importance. In case the curve shows a different shape than the one shown in the picture above, make sure to use a Faraday cage and connect the green GND lead to the Faraday cage.

See also section: [Reducing noise](#) on page 51.

2.3 First measurement on the Dummy Cell

This section describes how to do a first Linear Sweep measurement with any PalmSens BV instrument (EmStat or PalmSens-series) using the PalmSens Dummy Cell. Please follow these steps if you have a PalmSens Dummy Cell.

See also section: [Which dummy cell do I have](#) on page 17.

1. Make sure the instrument is on and connected to the PC
2. Open PSTrace and connect to the instrument
3. Load the following method file (menu: 'Method' → 'Load'):
'User Name\Documents\PSData\PSDummyCell_LSV_2V.psmethod'
4. Connect the Dummy Cell using the croc clips: CE, RE to the corresponding pads and WE to pad WE_B (10k resistor).
5. Start measurement by clicking start in the toolbar:




Figure 25 The Start measurement button in the toolbar

6. See the picture below for the expected results.




Figure 26 PSTrace in Scientific mode

2.4 First EIS measurement on the Dummy Cell

This section describes how to do a first impedance measurement with any PalmSens-series instrument using the PalmSens Dummy Cell.

Please follow these steps if you have a PalmSens Dummy Cell.

See also section: [Which dummy cell do I have](#) on page 17.

1. Make sure the instrument is on and connected to the PC
2. Open PSTrace and connect to the instrument
3. Load the data file (menu: 'Data' → 'Load data file...'): '[USER]\Documents\PSData\EIS examples\EIS on PS Dummy Cell WE_C.pssession'
4. Connect the Dummy Cell using the croc clips: CE, RE to the corresponding pads and WE to pad WE_C (Randles Circuit).
5. Select 'Overlay' in the drop-down next to the green start button:




Figure 27 Selecting 'Overlay' in the drop-down list

6. Start measurement by clicking start in the toolbar.
7. See the picture below for the expected results.




Figure 28 PSTrace in Scientific mode running EIS measurement

2.4.1 Fitting your data

Follow these additional steps to fit an equivalent circuit to your measurement data.

8. When the EIS measurement is finished, click the Circuit Fitting button at the top of the plot:




Figure 29 The 'Circuit Fitting' button for opening the 'Equivalent Circuit Analysis' window.

9. In the menu click 'Circuit → Load' and select the Simplified Randles circuit:




Figure 30 The 'Load circuit' window

10. Click the 'Load' button.
11. Now click the 'Fit Mode' button.
12. Click the 'Fit' button.
13. The found results should be similar to the results shown in the table below.

Table 1 Expected results when fitting the EIS data on a simplified Randles circuit

Element	Fitted Value	Min Value	Max Value	Unit	Error%
R 1	557.5	1.00E-6	1.00E+12	Ω	0.065
R 2	9966	1.00E-6	1.00E+12	Ω	0.035
C 1	0.033	1.00E-6	1.00E+3	μF	0.064
Chi-Squared:	4.60E-6	Iterations:	26		

See also section: [Electrochemical Impedance Spectroscopy](#) on page 113.

2.5 Readings

The status bar always shows the measured potential and the current as well as the noise if the cell is on.




Figure 31 Status bar

The measured potential, current (BiPot current if present), and noise readings are shown right at the bottom of each channel's Manual Control tab when the cell is set to ON.




Figure 32 Readings window

2.5.1 Underloads, overloads and oversteps

The box showing the status is used to note whether the measured values contain a current overload, current underload, voltage overload or timing overstep.

Measured current values are wrong when a current overload warning is shown in red. This occurs when the current is out of the range of the selected current range. A warning is given in orange for values that are exceeding 80% of the selected potential or current range. During a measurement with multiple ranges selected, the next current range, if available, will automatically be selected if this 80% threshold is exceeded.

If currents are below 4% of the selected current range, a current underload warning is given, since a lower current range can be applied. An underload will yield measurements with a low resolution and accuracy. Select lower current ranges if available to increase the available current resolution and accuracy.

See also section: [Resolution and optimal current range selection](#) on page 34.

The noise bar shows the noise level at the current range in use. In case the bar shows orange or red, it is advised to look for ways to limit the noise level.

See also section: [Reducing noise](#) on page 51.

Voltage overload means that the impedance between the counter and the reference electrode is too high. This can be found when:

- the counter or the reference electrode is not properly connected,
- the conductivity of the solution is too low, which is overcome by adding an electrolyte,
- an air bubble isolates the reference electrode from the solution.

An overstep occurs when the anticipated interval duration cannot be achieved, resulting in a time between two data points that surpasses the interval time. This can happen if there is a glitch in the communication, for example due to a bad Bluetooth connection. Or if too many instructions are added in the measurement loop of a MethodSCRIPT, when using the MethodSCRIPT Sandbox. Another reason can be setting an interval time that is lower than the interval time supported by the hardware.

3 Measurements

PSTrace allows to run all techniques supported by the connected instrument. This section explains how to set up a measurement. The options and limitations for each technique are described in the section: [Electrochemical techniques](#) on page 53.

For information about impedimetric measurements, see section: [Impedance Spectroscopy \(EIS\)](#) on page 113.

More theoretical background information about electrochemical techniques can be found in:

- Christopher M.A. Brett and Ana Maria Oliveira Brett, Electroanalysis (Oxford Chemistry Printers, 64) Oxford Science Publications, ISBN-13: 978-0198548164
- Joseph Wang, Analytical Electrochemistry 3rd ed, John Wiley & Sons, ISBN-13 978-0471678793

3.1 Connecting a cell to the potentiostat

PalmSens potentiostats are delivered with a cell cable terminated with 2-mm banana plugs as the lead connectors for the cell / electrodes / sample (excluding some models that come with SPE connectors). Each lead is accompanied by a crocodile (alligator) clip of the same color.



We advise keeping an unused lead connected with its crocodile clip. This practice prevents accidental electrical contacts with the exposed banana plug, minimizing the risk of unintended electrical contacts.

See the following table for PalmSens' standard cable lead colors.

Table 2 PalmSens standard cable lead colors

Connector color	Function	Label
Red	Working Electrode	W
Black	Counter / Auxiliary Electrode	C
Blue	Reference Electrode	R
Green	Ground	≡
Black with yellow or Black with white*	Sense lead	S
Yellow*	Working Electrode 2 / BiPot (if BiPot module is installed)	W2

*Note: these leads are not present in all models, refer to your potentiostat model specification for more details.

Most of PalmSens' instruments offer a 2- or 3-electrode cell setup. Some models (i.e. EmStat3+ and EmStat4 HR series) also are compatible with a 4-electrode setup.

Each electrode cable is individually shielded, and the entire cable features an external shielding mesh. The shielding extends up to the beginning of the banana plug. When you are using a Faraday cage, it is advisable to keep the entire banana plugs inside the Faraday cage housing to ensure the cage shields the banana plug.

See also section: [Reducing noise](#) on page 51.

3.1.1 2-electrode cell connection

For a 2-electrode setup, combine the counter and reference electrode leads. The counter electrode plug of PalmSens cables is stackable, allowing for convenient combination. Subsequently, connect the working electrode lead to one electrode and the combined counter and reference electrode leads to the other.

In models with a Sense lead, combine it with the working electrode, utilizing its stackable feature similar to the CE lead.




Figure 33 2-electrode cell connection

During a measurement in a 2-electrode setup, the potential will be applied between the two electrodes and the current will flow between the same two electrodes. This means the applied potential includes the potential drop across the working electrode's interface, the potential drop across the solution, and the potential drop across the counter electrode's interface. Any

changes at the counter electrode will influence the result. If your counter electrode is expected to be stable during the measurement due to low currents or a short duration of the measurement, you can use this setup almost equivalent to a 3-electrode setup. Otherwise, this setup is only used for applications where the full cell potential is an interesting parameter, for example, fuel cells, batteries, capacitors, dummy cells or electronic components.



Please avoid connecting a self-powered sample, such as batteries, fuel cells or supercapacitors, to an unpowered potentiostat, as this action can potentially damage the equipment. It is advisable to disconnect the sample before turning off the potentiostat to mitigate any risk of damage to the device.

3.1.2 3-electrode cell connection

The 3-electrode setup is the most common configuration for electrochemical measurements. In this arrangement, connect the corresponding leads to the respective electrodes. If your potentiostat model includes a Sense lead, combine it with the working electrode. The Sense plug of PalmSens cables is stackable, allowing for convenient combination.




Figure 34 3-electrodes cell connection

Similar to the 2-electrode setup, the working electrode (and sense if present) carries the current and potential. In this setup, the counter electrode also carries the current and polarizes versus the working electrode (WE). However, the potential is measured between the working (or sense if present) and reference electrodes, denoted as WE(+S) vs. RE. This potential is the one set and read in the plot results.

The potential WE vs. CE is referred to as the 'compliance potential,' and it is recorded in special cases.

See also section: [Recording the cell potential](#) on page 42.

In potential-controlled experiments, the potentiostat maintains a loop to ensure that the potential read (WE vs. RE) remains as programmed in the setup, consistently driving this potential towards the WE vs. CE polarization.

This has the advantage that the potential drop across the counter electrodes interface is compensated. Another advantage is that the reference electrode does not carry any current

and will thus keep a stable potential during the measurement. This means the potential of the working electrode is more accurate and reproducible.

3.1.3 4-electrode cell connection



This setup is only possible just with instruments equipped with a Sense lead.

In this configuration, you simply need to connect the corresponding leads to the electrodes, and no plug/lead combinations are required.




Figure 35 4-electrodes cell connections

The potential is measured between the sense lead and the reference electrode, while polarization drives the current flow between the working electrode and the counter electrode. This arrangement ensures that the potential drops across the counter or working electrodes' interfaces, and the associated reactions, do not contribute to the measurement. Processes specifically across the solution and obstacles, such as a membrane, between the sense and reference electrode, become the primary contributors to the measurement.

This setup is rarely used, with some examples including research solid-state electrolytes, membranes, or liquid-liquid interfaces.

3.1.4 Kelvin Connection

In addition to traditional 4-electrode cell connections, the "Kelvin Connection" serves as a specialized version of a 4-electrode setup, while functionally operating as a 2-electrode configuration. Notably, the sample is still a 2-electrode cell, typically with a very low impedance (high conductivity).

In the Kelvin Connection, like the 2-electrode connection, the CE lead is connected on the same side as the RE lead, while the WE lead is connected on the opposite side of the sample, alongside the Sense lead. The key distinction lies in its application when potential drops across cables and connectors must be meticulously considered. The primary objective of the Kelvin connection is to position the RE and Sense as close as possible to the sample. This proximity ensures that the impedance introduced by cables and connectors is effectively compensated through the functioning of the potentiostat.



Please avoid connecting a self-powered sample, such as batteries, fuel cells or supercapacitors, to an unpowered potentiostat, as this action can potentially damage the equipment. It is advisable to disconnect the sample before turning off the potentiostat to mitigate any risk of damage to the device.

3.2 Setting up a measurement

The 'Measurement' tab contains all the method parameters.

By default, the software displays the main parameters maximized, and secondary parameters minimized. Groups of parameters can be minimized and maximized by clicking on the section headers. The secondary parameters can be toggled using the '...' buttons.

The screenshot shows two expandable sections in the software interface. The first section, 'Pretreatment Settings', is expanded and shows a sub-section 'Linear Sweep Voltammetry Settings'. Below this, several parameters are listed with input fields and units:

t equilibration	2	s
E begin	-1.0	V
E end	1.0	V
E step	0.01	V
Scan rate	1.0	V/s

At the bottom of the 'Linear Sweep Voltammetry Settings' section, there is a button with three dots '...'. Below the main settings, there is another button with three dots '...'.

Figure 36 Section headers and the '...' button

With each change of parameters, the validation of the method is checked. Errors or incompatibilities are shown instantly at the bottom of the measurement tab.

The screenshot shows the 'Method Editor' window. It contains several input fields for parameters:

- E step: 0.01 V
- Scan rate: 1.0 V/s
- Number of scans: 1

Below these fields, there are several checkboxes and buttons:

- Post measurement (button)
- Cell on after measurement (checkbox)
- Auxiliary measurement (button)
- Aux. input (checkbox)
- Record BiPot (WE2) (checkbox)

At the bottom, there is a section for 'Expected duration: 00:00:03s' and '401 datapoints'. Below this, there is a red error message:

1 Error (for PalmSens4):
 [1]: Start potential must be between E vertex1 and E vertex2.

Figure 37 Method Editor showing an error.

Thus, the Method tab will prevent you from inserting a parameter that is outside the range of the connected instrument. If there is no instrument connected, PSTrace will allow you to exceed certain parameters, but warnings will appear as long as the device is connected.



If no instrument is connected, parameter validation is done on default values. The actual validity and limits can only be validated when the intended instrument is powered and connected.

If a parameter can be reached but there is an observation, a yellow warning will appear. This warning is advisory and will not prevent you from running the measurement.

3.2.1 Technique

The techniques list shows the currently selected technique and allows the user to change it. When an instrument is connected, the techniques list will be updated to show only the techniques supported by the connected instrument. Changing program mode will also filter out techniques that are not supported by the selected mode. If no instrument is connected and the default Scientific Mode is active, the list will show all the techniques supported by PSTrace.

See also section: [Program modes](#) on page 7.

The parameters for each technique are saved when switching to another technique. This allows the user to switch between techniques without the need to change parameters each time.

3.2.2 Sensor and sample

The Notes textbox can be used to describe information for example about the sensor used and sample that is measured.




Figure 38 The Notes box in the Method Editor.

3.2.3 Current range



Figure 39 Selection of allowed current ranges to be used during the measurement.

The current range selection buttons determine which range(s) are used for the measurement. The gray boxes refer to unselected ranges, while the blue boxes refer to selected ones. The small arrow hovering above one of the selected ranges determines at which current range the

measurement starts. The starting range can be changed by means of clicking above one of the other selected current ranges.

See also section: [Resolution and optimal current range selection](#) on page 34.

3.2.4 Method settings

Entered values for each method setting should be within the absolute limits any instrument by PalmSens BV can handle. If values are not supported by the instrument connected this is shown at the bottom of the method editor. See the next chapter for a description and allowed values of each parameter.

3.2.5 Validation

In case a value entered is invalid, an error message at the bottom is shown. Errors are shown in red and warnings are shown in yellow. In case of a warning, this can be ignored by the user and the measurement can be started. In case of an error, the values need to be changed first to a valid value to permit the start of the measurement.




Figure 40 Error shown in the Method Editor

A warning is shown in yellow. Warnings can be ignored but, on most occasions, this will result in a bad measurement due to wrong settings.




Figure 41 Warning shown in the Method Editor

In case a value for the step potential is given which cannot exactly be met by the instrument because of its hardware-determined resolution, the actual value is shown at the bottom as 'real E step'.




Figure 42 Real E step shown in the Method Editor

3.3 Standard measurement sequence

Voltammetric measurements have the following sequence.




Figure 43 Normal measurement sequence

The sequence of a voltammetric measurement (except for stripping chronopotentiometry) is:

1. In case "Measure vs OCP" option is enabled: determine the OCP.
2. Apply 'E condition', the conditioning potential, if 't condition' is not zero.
3. Apply 'E deposition', the deposition potential, if 't deposition' is not zero.
4. Apply 'E begin' and wait 't equilibration' seconds.
5. Start measurement at 'E start' and continue until 'E end' with steps 'E step', with the specified 'Scan rate'. In CV the scan is continued by reversing the scan direction. The current range is set automatically however with the constraints as specified.
6. Depending on the potentiostat model, the measurement is plotted during the measurement or when the measurement is finished.
7. If the cell must remain switched on, 'E standby' is applied; otherwise, the cell is switched off.

The sequence of an amperometric measurement:

1. In case "Measure vs OCP" option is enabled: determine the OCP.
2. Apply 'E condition' the conditioning potential, if 't condition' is not zero.
3. Apply 'E deposition' the deposition potential, if 't deposition' is not zero.
4. Apply E or 'E equilibration' and wait 't equilibration' seconds.
5. Start the measurement.

3.3.1 Pause and abort

During the measurement, the curve is shown on the screen in real-time if the scan rate allows for it. In that case also the 'Pause' button can be used to halt the scan until the same button is used again. This button is not available at higher scan rates.

It is possible to abort a measurement, by pressing the abort button above the plot.



Please note that when the 'Stop' or 'Pause' button is pressed, the potentiostat will complete the current measurement point. This can happen almost instantly in most cases, but it might take some time in situations like low-frequency impedance or very large time intervals.

3.3.2 Cell on after measurement

A standby potential can be set in the 'Post measurement' section in the Method Editor. The duration can be set using the 'For specified period' checkbox. This potential is applied after the measurement has concluded for the specified time duration or indefinitely if this parameter is not selected.

Post measurement

☒ Cell on after measurement

E standby V

☒ For specified period

t standby s

Figure 44 Using cell on after measurement.

3.4 Measuring versus OCP

As a standard practice, the potentials set in PSTrace refer to the values of the Working Electrode (WE) relative to the Reference Electrode (RE). This potential is also known as the "Polarization Potential". However, in many measurements, the potential can be specified with respect to the Open Circuit Potential (OCP) instead. This approach is commonly employed in corrosion and some energy-related applications.

In case one or more potentials are specified with respect to the OCP, this potential must be determined before the actual measurement is done. This OCP measurement requires a variable time, which is determined by the drift of the open circuit potential and the maximum time to measure the OCP value. The OCP value is set as soon as the drift is lower than the specified value for the 'Stability criterion' or when the 't Max. OCP' has elapsed. Thus, setting 'Stability criterion' to zero renders it ineffective.

☐ Stop when I < uA

☐ Stop when I > uA

☒ Measure vs OCP

☒ E begin versus OCP

☒ E end versus OCP

t Max. OCP s

Stability criterion mV/s

Figure 45 OCP parameters

3.5 Resolution and optimal range selection

The Method Editor shows one or two rows of buttons for the current and potential ranges. The potential ranges are only visible when a galvanostatic technique is selected.

With the range-buttons, the applicable current range(s) and potential range(s) during the measurement can be selected. If more than one button is selected, the instrument will select the most optimal range automatically (auto-ranging).

A measurement starts at the range with the small arrow above it.




Figure 46 Selection of applicable current and potential ranges shown for Potentiometry

The starting range can be changed by means of clicking above one of the other selected current ranges changes.

Double-clicking a range button will select only this range.

3.5.1 Automatic range selection / auto-ranging

During equilibration, the most appropriate range will be selected. If during a measurement the current or potential approaches a range limit, the instrument will switch automatically to a more suitable current or potential range (lower or higher). When a measurement exceeds the upper threshold (80% of the range), a higher range is selected within the range of selected ranges.

Likewise, when a measured current or potential falls repeatedly below the lower threshold (4% of the selected range), a lower range is selected if enabled by the user.

Because auto-ranging takes some time (in the order of milliseconds) this feature is only available at lower scan rates. For techniques at very high scan rates automatic ranging is not available. This is always shown as an error in the Method Editor if this is the case.

Table 3 The maximum measured currents for each instrument.

Instrument	Max. measured current
EmStat1, 2 and 3(+)	±2 * selected range
PalmSens1 and 2	±2 * selected range
PalmSens3	±3 * selected range
PalmSens4	±6 * selected range
EmStat4 LR and HR	±3 * selected range
EmStat Pico	±0.6 * selected range

An overload warning is given for measured values $\geq 80\%$ of the selected range.

For all instruments except for the EmStat Pico an underload warning is given for measured values $\leq 4\%$ of the selected range. For the EmStat Pico the underload warning is given for measured values $\leq 2\%$ of the selected range.

For instruments that support multiple potential ranges, the table below shows the available range for each instruments.

Table 4 The maximum measured potentials for each instrument.

Instrument	Max. measured potential
PalmSens4	± 10 * selected range
EmStat4 LR	± 3 * selected range
EmStat4 HR	± 6 * selected range

3.5.2 Spikes and jumps

Please note that auto-ranging can cause spikes in the measurement curve, especially during fast measurements and/or measurements with high Faradic currents. Spike or jumps may occur due to a short interruption or change in measurement interval when switching between the ranges.



If auto-ranging is not needed it is always better to enable a single current range to prevent spikes or jumps in your curves.

If you are uncertain about the current you will obtain, and the experiment requires manual range selection, start with the highest current range. You can then lower it in subsequent attempts if you still receive underload readings, even at the highest resulted currents.

3.6 Running a measurement

Measurements can be started and stopped with the buttons in the measurement tool strip or the 'Measure' menu.

The steps conditioning, deposition, and equilibration can be skipped each using the skip button.




Figure 47 Skip button

Next to the start button, you can choose to make a 'New' graph for the measurement. If 'Overlay' is selected, the measured curve will be added with another color to the existing curves in the plot. Or measure a 'Blank' curve as a background scan to subtract it later.




Figure 48 Options for the next measurement

If a blank curve is present an additional button is shown, see below. By enabling this button, the blank will be subtracted automatically after the measurement is finished.

See also section [Blank subtraction](#) on page 45.




Figure 49 Use Blank button

The progress and status of a measurement is shown in the status bar at the bottom.




Figure 50 Status bar

3.6.1 Autosaving measured data

If autosave is enabled, every single measurement will be saved automatically in the corresponding format after the measurement has finished or is aborted.

See also section: [Files](#) on page 235.



Figure 51 Enabling autosave

The autosave output can be configured in the Autosave settings window.




Figure 52 The Autosave Settings window

3.6.1.1 Output settings

The output file gets the prefix followed by a sequential number. This number is based on the number of existing files with the given prefix in the storage folder specified and therefore increments with each new file added to the folder.

3.6.1.2 Output data formats

The output file can be stored in different formats. The PSTrace format will store the file in a format that can be loaded in PSTrace again.

See also section: [File types](#) on page 235.

A CSV file is a comma-separated file and can easily be imported into third party programs like Excel and Origin.

As a third option, it is possible to store the data in a native Origin file format. The latter will also include a graph (except for EIS data).

3.6.1.3 EIS data export settings

The button 'EIS data export settings' allows which columns are exported to the non-PSTrace file formats.

3.6.1.4 Multiple curves output

Some measurements produce multiple curves, like a Cyclic Voltammetry measurement with multiple scans or a measurement done with a multiplexer. Enabling the checkbox will save the output of each separate curve to a single file instead of to a single file containing the data of all curves.

3.6.2 Recovery file




Figure 53 Recovery file settings

The 'Autosave settings' window also shows options for keeping a recovery file during measurements. This option can be useful for long-term measurements. In case of a PC crash or power outage, the most recent measurement data will still be available. The recovery file is overwritten for each new measurement.

3.6.3 On-device storage

The PalmSens4, EmStat4-series and Sensit BT are equipped with additional hardware for on-board data storage. In case the option 'Save on internal storage' is checked, the measurement will be stored on-board the instrument in a folder with the name of the day's date.



Figure 54 Save on internal storage in the Method Editor of PSTrace

The on-board data can be browsed and downloaded using the menu: 'Data' → 'Load data from internal storage...'.



For the PalmSens4 the following techniques are excluded for on-board storage:
EIS, MultiStep and MixedMode.

3.7 Peaks and levels





Figure 55 Peak/Level parameters tab (showing options for a technique as a function of time)

Table 5 Parameters in the Peak/Level parameters tab

Minimum Peak Width	<p>determines the minimum peak width. Peaks with a lower width may not be detected. The peak width is the Full Width at Half Maximum (FWHM).</p> 
Minimum Peak Height	<p>determines the minimum peak height. Peaks lower than this value are neglected.</p>
Smooth Window	<p>determines to which degree the measured curve is smoothed. The specified number is used to smoothen the data using the Savitzky-Golay method. Please note that too many points will influence the shape of the curve</p>
Oxidation or Reduction	<p>determines whether peaks or current levels are positive (Oxidation) or negative (Reduction).</p>
Current Level or Peaks	<p>determines whether peaks or current levels are measured</p>

Extra options for methods vs time:

3.8 Multiplexer

In the Multiplexer tab, it can be specified which multiplexer channels are measured and whether they should be measured consecutive (one by one) or scanned alternately (all channels simultaneously).

When the 'Consecutive' mode is chosen, any channel can be selected. In case 'Alternate' is checked, the manual selection of channels is limited to successive channels. So, it is possible to measure 2 to 8 or 2 to 16 successive channels, depending on the type of multiplexer in use.




Figure 56 Multiplexer parameters

3.8.1 Consecutive measurements

When doing Consecutive multiplexer measurements, the active multiplexer channel will switch to the next and be started after the entire measurement as defined in the Measurement tab has finished. This mode is available for every technique.

3.8.2 Alternate measurements

During the execution of alternating multiplexer measurements, all chosen channels will rapidly switch within each measurement interval. The complete switching process, including both switching and stabilization, for the multiplexer takes approximately ± 31 milliseconds. As a result, the lower limit of the interval time is defined as (number of channels selected) * 0.031 s. For instance, when using all 8 channels, the minimum interval time is 0.25 seconds.

The following techniques support the Alternate multiplexer mode:

- Chronoamperometry;
- Chronopotentiometry;
- Open Circuit Potentiometry;
- and (Galvanostatic) Electrochemical Impedance Spectroscopy.

3.8.3 MUX8 and MUX16 hardware settings

The multiplexer can be used in different configurations. For the older models, dip-switches are available on the underside of the multiplexer to change the configuration.

See also section: [MUX8 and MUX16 multiplexers](#) on page 265.

3.8.4 MUX8-R2 settings

The MUX8-R2 hardware configuration settings can be defined in software and will become visible in the 'Multiplexer' tab if the MUX8-R2 is either detected upon connection or selected in the 'General Settings' window.




Figure 57 Multiplexer MUX8-R2 hardware configuration settings

When switching techniques, the initial MUX8-R2 settings will be taken from the default settings as defined in the settings window of PSTrace. However, when saving a method to either a .psmethod file or .psession file the settings will be saved as present in the 'Multiplexer' tab.

See also section: [MUX8-R2 multiplexer](#) on page 275.

3.9 Recording an auxiliary input

Supported by all instruments with an auxiliary input port.

To record an auxiliary value, like an analog input voltage, temperature or cell potential, this needs to be enabled in the Method Editor in the section 'Record additional data':



Figure 58 Enabling recording on the auxiliary input in the Method Editor



Whether or not recording additional data is supported depends on a combination of the connected instrument and the selected technique.

If the checkbox 'Record Aux. input' is checked, the voltage of the corresponding analog input pin on the physical auxiliary port is recorded. The auxiliary input value can also be translated to another value, like temperature, RPM (For an RRDE) or pH, depending on what is connected. This can be set in the 'General Settings' window.

See also section: [Hardware configuration](#) on page 8.

Please refer to the instrument-specific sections in this manual for more information about the positions of the digital and analog I/O pins on your instrument's auxiliary port.

3.10 Recording the cell potential



Supported by a selection of techniques with PalmSens4, EmStat4 and EmStat Pico.

Recording the cell potential can be enabled in the Method Editor in the section 'Record additional data':




Figure 59 Enabling recording the cell potential in the Method Editor

The 'cell potential' is the potential between the working and the counter electrode and thus the potential applied to the whole cell. In other words, this is the working electrode's potential versus the counter electrode's potential. This potential is also called "Compliance Voltage", and as such, this specification refers to the maximum potential that can be applied between WE and CE.




Figure 60 Cell potential

Sometimes it is interesting to calculate the potential of the counter electrode versus the reference electrode. This can be easily calculated with the applied potential and the cell potential. As is visible from picture above the cell potential is calculated by

$$\text{Cell Potential} = E(\text{WE}) - E(\text{CE})$$

By rearranging the potential of the counter electrode $E(CE)$ can be calculated by

$$E(CE) = E(WE) - \text{Cell Potential}$$



For the EmStat Pico, as well as the Sensit series instruments, there is a limitation when recording cell potential. In certain situations, the recorded cell potential values may reach a certain limit below the actual value.

3.11 Recording the WE potential



Supported by a selection of techniques with (Multi)PalmSens4 and (Multi)EmStat4.

Recording the WE potential can be enabled in the Method Editor in the section 'Record additional data':

▲ Record additional data

☐ Record Aux. input

☐ Record cell potential

☒ Record WE potential

☐ Record BiPot (WE2)

Figure 61 Enabling recording the WE potential in the Method Editor

The 'WE potential' is the applied potential between the working and the reference electrode. This option can be enabled to measure the applied potential during supported Voltammetric and Amperometric techniques on supported devices and can be used to determine if the instrument was able to apply the requested potential.

3.12 Recording the WE current



Supported by a selection of techniques with PalmSens4 and EmStat4.

Recording the WE current can be enabled in the Method Editor in the section 'Record additional data':

▲ Record additional data

☐ Record Aux. input

☒ Record WE current

Figure 62 Enabling recording the WE current in the Method Editor

The 'WE current' is the applied current between the working and the reference electrode. This option can be enabled to measure the applied current during supported Potentiometric techniques on supported devices. The current will be measured in the applied current range. In the case of an OCP measurement, the current will be measured in the current range as specified in the 'Select current range(s) selection' bar.

3.13 BiPotentiostat

Supported by instruments with bipotentiostat functionality for using a second working electrode (WE2).

If a BiPot module is present in a PalmSens instrument and the 'BiPot present' setting is checked in the 'Settings' window, the second working electrode or WE2 can be used. Please note that the WE2 will be measured/polarized versus the same RE and CE used for the main WE1.

The Settings window can be found in the menu: 'Tools' → 'General Settings...'




Figure 63 Select BiPot present in Settings window

In case one of the channels has a BiPot module present, this will be detected automatically.

An additional tab 'BiPot' will become visible in the method editor. This tab contains the settings for the use of the second working electrode (WE2).




Figure 64 The additional BiPot tab in the Method Editor.

The E WE2 is the potential set on the WE2.

The BiPot can be used with the following methods:

- Linear Sweep Voltammetry
- Cyclic Voltammetry
- Amperometric Detection
- Multistep Amperometry
- Multistep Potentiometry

3.13.1 Plot $i(\text{WE2})$ vs $E(\text{WE1})$

The plot window will show the current obtained for WE1 as well as for WE2 or Auxiliary.

In case the checkbox 'Plot $i(\text{WE2})$ vs $E(\text{WE1})$ ' is checked the $i(\text{WE2})$ curve will always be shown as an overlay in the plot, also when WE2 was maintained at a fixed potential.

3.14 Blank subtraction

Blank subtraction can be used by first running a measuring as 'Blank' using the dropdown list next to the Start button. Or an existing curve can be designated as a Blank curve by clicking on this curve in the plot legend.




Figure 65 Select 'Blank' to measure a Blank curve for Blank subtraction before clicking the Start button.




Figure 66 By changing the curve function in the curve tool window, which is accessed by clicking on a curve in the plot legend.

When a Blank has become available, the 'Use Blank' button appears next to the measurement controls in the bar on the top of the screen.

When enabling the Use Blank button, the available Blank curve is subtracted immediately after the measurement from the measured curve. The Blank Curve is always shown in gray in the plot.




Figure 67 Button for Blank subtraction

3.15 Ohmic (iR) compensation

This setting is found under the button for advanced settings.



iR Compensation is optionally available for the PalmSens4 and EmStat4X.

The resistance between the reference electrode and the double layer of the specimen can cause a significant potential drop, decreasing the applied potential where it is required. This potential drop (E_{ohm}) is called "Ohmic drop" and it is directly proportional to the electrolyte resistance (R) and the actual current (i), thus $E_{ohm} = iR$.

The compensation module provides a "positive feedback" to compensate for the iR drop between the Reference electrode and the outside of the double layer of the electrochemical cell.

The resistance to compensate for can be entered directly in the Method Editor in PSTrace:




Figure 68 Resistance to compensate for as specified in the Method Editor

PalmSens recommends that instruments with iR-compensation module also have EIS capability, as it is the most common and straightforward technique for determining the resistance to be compensated. If the resistance of the electrolyte is accurately determined, then set the compensation to a maximum of 95% to prevent system instability.

Please note that if auto ranging is not allowed for the compensation used in combination with the selected current ranges, this is shown in the Method Editor:




Figure 69 Error shown in Method Editor

Make sure a single current range is selected in these circumstances. Moreover, each current range has a maximum resistance that can be compensated for.

See the following sections for more instrument-specific limitations:


- [iR compensation with PalmSens4](#) on page 305
- [iR compensation with EmStat4X](#) on page 337

3.15.1 Supported Techniques

The following techniques are supported for use with iR compensation:

- Linear Sweep Voltammetry
- Cyclic Voltammetry
- Square Wave Voltammetry
- Differential Pulse Voltammetry
- Normal Pulse Voltammetry
- Chronoamperometry
- Multistep Amperometry

3.16 Potential and current limits

The settings for current and potential limits are found under the  button for advanced settings.

The limits entered here apply to the entire measurement excluding the pretreatment stages. This parameter is also called “cut-off”.

If the referred measurement is inside a Script sequence, thus the Script will continue to the next step.

This resource is not available in some instruments.

3.17 Override bandwidth filter

This setting is found under the  button for advanced settings.

This function is only supported and visible when a MethodSCRIPT-enabled device like EmStat Pico or (Multi)EmStat4 is connected.




Figure 70 Bandwidth override setting in the Method Editor

Normally the filters of the potentiostat are set optimally for the required bandwidth. In some exceptional cases it is useful to override the automatically determined bandwidth setting.

When the ‘Override bandwidth’ checkbox is set, the setting will be overridden with the value given in Hertz.

3.18 Digital triggers

This setting is found under the  button for advanced settings.




Figure 71 Hardware trigger settings

In case “Trigger at ...” is set, the selected digital line(s) on the AUX port of the instrument will be set high when triggered. They will remain high until the end of the measurement.

See the instrument-specific sections in this manual for more information about the position of the digital pins on your instrument’s auxiliary port.

In case “Trigger at delay after start” is used, the delay will be rounded to the applicable interval time between each measured data point.

3.19 Measuring the noise level of the instrument

To determine the environmental noise levels, make sure the (green) test sensor is connected or the (blue) PS Dummy Cell is connected on pads RE, CE and WE_A.

1. Load method file ‘PSNoiseTest.psmethod’.
2. Start measurement.
3. Repeat the measurement but with a scan rate of 0.05 V/s. Note that the noise level is lower because the current sampling time is longer now. This decreases the measured noise level.




Figure 72 The curve shown in the figure above is a typical curve measured using the unmodified PSNoiseTest method and a PalmSens2. The sensor and connections are housed inside a grounded Faraday Cage.

In case the noise level is too high this might be due to a noisy power supply of the USB ports and therefore the use of a USB Hub with its own ac-adaptor between the PC or laptop and EmStat is advised.



Figure 73 The green (noisy) and red curves are both measured with the same conditions as before and with the same EmStat. The red curve is obtained when EmStat is connected to a USB hub instead of directly to a USB port of the PC.

The instrument can be controlled manually using the 'Manual Control' tab.

To evaluate noise levels in the Manual Control tab click 'Cell On'. Change the applied potential from $E = 0.000$ V to $E = 0.500$ V by manipulating the horizontal scrollbar or entering the value in the textbox. Select the current range of 10 nA (or 100 nA if 10 nA is not available). The measured values of the potential and current as well as the noise are shown.

If the noise bar shows orange or red, the noise level is high. You are advised to place the test sensor inside a Faraday cage with the green ground lead connected to the cage.

3.20 Reducing noise

Electrochemical sensors and cells are susceptible to noise. If noise is interfering with your measurement, the solving strategies are rather numerous, but the sources for noise are also numerous. Here we describe the most successful and common methods for noise reduction.

Your power grid is usually using an alternating current. This undulating current influences the measured currents. PSTrace and PSTouch have a filter setting for this mains frequency. Check in the Tools menu under General Settings if the mains frequency is set correctly.

Noise can also be introduced through the USB port of your computer. In that case using an external hub with its own power supply could reduce this noise.

Our environment is filled with electrical fields. Some of them are created by devices around us as side effects or in case of wireless communication on purpose. Although it is a bad idea to measure directly next to an electric arc furnace, it is usually not possible to have a workspace free of electrical fields, especially not during point-of-care measurements. A Faraday cage is usually sufficient to create a field-free environment. A metal box or cage out of metal mesh is a good Faraday cage. Even a shield out of aluminum foil can help. Place your electrochemical cell inside the Faraday cage and connect the cage to the ground lead (green) of the potentiostat. The cable delivered with your EmStat or PalmSens has an inbuilt shield and should protect your signal outside the Faraday cage. This is one of the most effective methods to reduce noise.



Be cautious to avoid accidental contact between electrodes or plugs and the Faraday cage. A short circuit between the Faraday cage (and thus the ground lead) and an electrode can lead to current leakage, resulting in incorrect responses and potential damage to the potentiostat in some cases.

Cables should not be excessively long, as they can inadvertently function as antennas that captures electromagnetic noise. However, the cables provided with your EmStat or PalmSens instruments include built-in shielding. When you use the original cable, concerns about noise induced by the cable are minimal. Nevertheless, avoid using extensions, but if they are necessary, they should be used within the Faraday cage.

Ground your measurement equipment. The best way to connect your equipment is star-shaped, that is all parts are connected with the ground at the same point. In an electrochemical lab that point is usually one small space of the Faraday cage. This way earth loops that induce noise are avoided.

Check if the contacts are corroded. If so, remove the stains, for example with sandpaper.

The Readings window (shown during the time the cell is on) displays the potential and the current and shows the noise level. The background color shows whether the noise might deteriorate the measurement. In case the noise level is higher than 0.1 times the selected current range, the bar will turn orange. In case the noise level exceeds 0.5 times the selected current range, the bar turns red. It is strongly advised to prevent measurements under such conditions.

The EmStats and Sensits are directly powered by the USB port. Some PC's or laptops have a noisy USB power supply. In this case, it is advised to use a USB hub with its own ac-adaptor between the PC or laptop and EmStat.

If you find that the suggestions to reduce noise haven't produced the desired results, it's possible that your data collection parameters need adjustment. If your results still appear noisy or jagged, consider increasing the 't interval,' 'E step,' or 'i step' parameters to achieve smoother data. Increasing these parameters will result in a more gradual change in time, potential, or current values, respectively, and can help in obtaining more refined and less fluctuating results. Experimenting with these settings may enhance the quality and precision of your data.

4 Electrochemical techniques

PSTrace supports the following techniques in the (default) Scientific Mode. Some instrument models only support a selection of these techniques. At the beginning of each technique-specific section is a list with the instrument models supporting the technique. All the above-listed techniques are described in this section, except EIS/GEIS.

See also section: [Electrochemical Impedance Spectroscopy](#) on page 113.

Table 6 All available techniques in PSTrace

Voltammetric techniques	
▪ Linear Sweep Voltammetry	LSV
▪ Cyclic Voltammetry	CV
▪ Fast Cyclic Voltammetry	FCV
▪ AC Voltammetry	ACV
Pulsed techniques	
▪ Differential Pulse Voltammetry	DPV
▪ Square Wave Voltammetry	SWV
▪ Normal Pulse Voltammetry	NPV
Amperometric techniques	
▪ Chronoamperometry	CA
▪ Zero Resistance Amperometry	ZRA
▪ MultiStep Amperometry	MA
▪ Fast Amperometry	FAM
▪ Pulsed Amperometric Detection	PAD
▪ Multiple-Pulse Amperometric Detection	MPAD
Galvanostatic techniques	
▪ Linear Sweep Potentiometry	LSP
▪ Chronopotentiometry	CP
▪ MultiStep Potentiometry	MP
▪ Open Circuit Potentiometry	OCP
▪ Stripping Chronopotentiometry	SCP (or PSA)

Coulometric techniques	
▪ Chronocoulometry	CC
Other	
▪ Mixed Mode	MM
▪ Electrochemical Impedance Spectroscopy	EIS
▪ Fast Electrochemical Impedance Spectroscopy	FEIS
▪ Galvanostatic Electrochemical Impedance Spectroscopy	GEIS
▪ Fast Galvanostatic Electrochemical Impedance Spectroscopy	FGEIS

4.1 Linear Sweep Voltammetry (LSV)

Supported instruments:

- PalmSens series
- EmStat series
- Sensit series

4.1.1 Description

In Linear Sweep Voltammetry a potential scan is performed from the begin potential, 'E begin', to the end potential 'E end'. The scan is not exactly linear, but small potential steps (E step) are made. The current is measured (sampled) during the last 25% interval period of each step. The number of points in a curve showing the current versus potential is

$$(E \text{ end} - E \text{ begin}) / E \text{ step} + 1.$$

The scan rate is specified in V/s, which determines the time between two steps and thus the sampling time. The interval time is equal to E step / scan rate. So, when E step is 0.005 V and the scan rate 0.1 V/s the interval time is 0.05 s.

4.1.2 Measuring

Increasing the 'E step' is suitable for rapid data collection and may improve the signal-to-noise ratio. On the other hand, decreasing the step potential is advantageous when you require enhanced resolution, sensitivity, and the ability to study complex electrochemical systems in detail. It's essential to strike a balance between scan rate, resolution, and sensitivity based on the characteristics of your electrochemical system and research goals. As a start point, it is common to use 'E step' between 0.003 and 0.01V for Scan rates in the range of 0.01 to 0.1 V/s.

In some applications, it is important that the current does not get too high. This might ruin the working electrode. If the potential at which this will occur is not known, it is possible to specify a maximum current value at which the scan stops. In this case, the end potential specified by the user may not be reached.

See also section: [Potential and current limits](#) on page 48.




Figure 74 Potential versus time during Linear Sweep Voltammetry

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

By default, the result shown is a graph of Current versus Potential.




Figure 75 Typical LSV plot for a non-diffusion limited redox system and the WE_A of the PalmSens Dummy Cell.

4.1.3 Technique specific parameters


E begin	The potential where the measurement begins at.																		
E end	The potential where the measurement stops at.																		
	<p>The applicable potential range of each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr> <td>PalmSens4</td><td>-10 V to +10 V</td></tr> <tr> <td>EmStat1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>EmStat3</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat3+</td><td>-4 V to +4 V</td></tr> <tr> <td>EmStat4 LR</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat4 HR</td><td>-6 V to +6 V</td></tr> <tr> <td>EmStat Pico</td><td>High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V</td></tr> </table> <p>See also section Limitations and extra options for EmStat Pico on page 109.</p>	PalmSens1 and 2	-2 V to +2 V	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V	EmStat1 and 2	-2 V to +2 V	EmStat3	-3 V to +3 V	EmStat3+	-4 V to +4 V	EmStat4 LR	-3 V to +3 V	EmStat4 HR	-6 V to +6 V	EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V
PalmSens1 and 2	-2 V to +2 V																		
PalmSens3	-5 V to +5 V																		
PalmSens4	-10 V to +10 V																		
EmStat1 and 2	-2 V to +2 V																		
EmStat3	-3 V to +3 V																		
EmStat3+	-4 V to +4 V																		
EmStat4 LR	-3 V to +3 V																		
EmStat4 HR	-6 V to +6 V																		
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V																		
E step	<p>The potential step size.</p> <p>The applicable step range for each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>1 mV to 250 mV</td></tr> <tr> <td>PalmSens3</td><td>0.15 mV to 250 mV</td></tr> <tr> <td>PalmSens4</td><td>0.075 mV to 250 mV</td></tr> <tr> <td>EmStat1 and 2</td><td>0.1 mV to 250 mV</td></tr> <tr> <td>EmStat3(+)</td><td>0.125 mV to 250 mV</td></tr> <tr> <td>EmStat4 LR</td><td>0.1 mV to 250 mV</td></tr> <tr> <td>EmStat4 HR</td><td>0.183 mV to 250 mV</td></tr> <tr> <td>EmStat Pico</td><td>High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV</td></tr> </table>	PalmSens1 and 2	1 mV to 250 mV	PalmSens3	0.15 mV to 250 mV	PalmSens4	0.075 mV to 250 mV	EmStat1 and 2	0.1 mV to 250 mV	EmStat3(+)	0.125 mV to 250 mV	EmStat4 LR	0.1 mV to 250 mV	EmStat4 HR	0.183 mV to 250 mV	EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV		
PalmSens1 and 2	1 mV to 250 mV																		
PalmSens3	0.15 mV to 250 mV																		
PalmSens4	0.075 mV to 250 mV																		
EmStat1 and 2	0.1 mV to 250 mV																		
EmStat3(+)	0.125 mV to 250 mV																		
EmStat4 LR	0.1 mV to 250 mV																		
EmStat4 HR	0.183 mV to 250 mV																		
EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV																		

Scan rate

The applied scan rate, also called sweep rate or sweep speed. The applicable range depends on the value of E step since the data acquisition rate is limited by the connected instrument.

The applicable scan rates for each instrument are:

PalmSens1 and 2	1 mV/s (1 mV step) to 25 V/s (5 mV step)
PalmSens3	0.02 mV/s (0.15 mV step) to 500 V/s (5 mV step)
PalmSens4	0.02 mV/s (0.075 mV step) to 500 V/s (10 mV step)
EmStat1 and 2	0.02 mV/s (0.1 mV step) to 5 V/s (5 mV step)
EmStat3(+)	0.025 mV/s (0.125 mV step) to 5 V/s (5 mV step)
EmStat4 LR	0.01 mV/s (0.1 mV step) to 500 V/s (200 mV step)
EmStat4 HR	0.01 mV/s (0.183 mV step) to 500 V/s (200 mV step)
EmStat Pico	0.01 mV/s to 10 V/s (10 mV step)

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

4.2 Cyclic Voltammetry (CV)

Supported instruments:

- PalmSens series
- EmStat series
- Sensit series

4.2.1 Description

In Cyclic Voltammetry, recurrent potential scans are performed between the potentials 'E vertex1' and 'E vertex2,' going back the number of times determined by the 'Number of scans.' The scan starts at the begin potential (E begin) which can be at one of these vertex potentials or anywhere in between. The experiment will always terminate at the same potential set as the 'E begin'.

The change of the potential per time is the scan rate. The scan rate is specified in V/s, which determines the time between two steps and thus the sampling time. The interval time is E step / scan rate. So, when E step is 0.005 V and the scan rate 0.1 V/s the interval time is 0.05 s.

The scan is not exactly linear, but an approximation by means of a staircase with small potential steps (E step). The current is measured (sampled) during the last 25% interval period of each step. The number of points per scan of the current versus potential curve is $2 * (E \text{ end} - E \text{ begin}) / E \text{ step}$.

The following graph shows the signal applied and the periods of current sampling.




Figure 76 Potential versus time during Cyclic Voltammetry

4.2.2 Measuring

As in LSV it is sometimes important that the current does not get too high. This might ruin the working electrode. If the potential at which this will occur is not known, it is possible to specify a maximum current value at which the scan direction will reverse. If the current set is reached, it will reverse the scan direction before the next vertex.

The user can also interactively determine where the scan direction is reversed, using the button 'Reverse'. In these cases, the specified vertex potential by the user may not be reached.

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

As in LSV, by default, the result shown is a graph of Current versus Potential. Successive cycles will be displayed in different colors.





Figure 77 Typical CV plot for a free diffusing redox species

4.2.3 Technique specific parameters

E start	The potential where the scan starts and stops at.
E vertex1	The first potential where the scan direction reverses.
E vertex2	The second potential where the scan direction reverses.
The applicable potential range of each instrument is:	
PalmSens1 and 2	-2 V to +2 V
PalmSens3	-5 V to +5 V
PalmSens4	-10 V to +10 V
EmStat1 and 2	-2 V to +2 V
EmStat3	-3 V to +3 V
EmStat3+	-4 V to +4 V
EmStat4 LR	-3 V to +3 V
EmStat4 HR	-6 V to +6 V
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V

See also section [Limitations and extra options for EmStat Pico](#) on page 109.

E step	The potential step size.	
	The applicable step range for each instrument is:	
	PalmSens1 and 2	1 mV to 250 mV
	PalmSens3	0.15 mV to 250 mV
	PalmSens4	0.075 mV to 250 mV
	EmStat1 and 2	0.1 mV to 250 mV
	EmStat3(+)	0.125 mV to 250 mV
	EmStat4 LR	0.1 mV to 250 mV
	EmStat4 HR	0.183 mV to 250 mV
	EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV
Scan rate	The applied scan rate, also called sweep rate or sweep speed. The applicable range depends on the value of E step.	
	The applicable scan rates for each instrument are:	
	PalmSens1 and 2	1 mV/s (1 mV step) to 25 V/s (5 mV step)
	PalmSens3	0.02 mV/s (0.15 mV step) to 500 V/s (5 mV step)
	PalmSens4	0.02 mV/s (0.075 mV step) to 500 V/s (10 mV step)
	EmStat1 and 2	0.02 mV/s (0.1 mV step) to 5 V/s (5 mV step)
	EmStat3(+)	0.025 mV/s (0.125 mV step) to 5 V/s (5 mV step)
	EmStat4 LR	0.01 mV/s (0.1 mV step) to 500 V/s (200 mV step)
	EmStat4 HR	0.01 mV/s (0.183 mV step) to 500 V/s (200 mV step)
	EmStat Pico	0.01 mV/s to 10 V/s (10 mV step)
Number of scans	The number of repetitions for this scan.	

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

4.2.4 Fast CV

A CV becomes a Fast CV if the scan rate in combination with E step results in a rate of over 2500 points / second (scan rate / E step > 2500).

See next section for more information.

4.3 Fast Cyclic Voltammetry (FCV)

Supported instruments:

- PalmSens series
- EmStat4 LR and HR

4.3.1 Description

(See also previous section: Cyclic Voltammetry)

In Cyclic Voltammetry a cyclic potential scan is performed between two vertex potentials E vertex1 and E vertex2. The scan can start (E begin) at one of these vertex potentials or anywhere in between.

A CV becomes a Fast CV if the scan rate in combination with E step results in a rate of over 2500 points / second (scan rate / E step > 2500).

Additional options for Fast CV compared to a normal CV include averaging multiple scans and using equilibration scans. This technique can be particularly useful for very noisy measurements. Of course, it would be preferable to eliminate noise at the source first.

Auto-ranging cannot be used at these high speeds, so only one current range can be selected.





Figure 78 Fast Cyclic Voltammetry with 'n scans'=2 and 'n averaged scans'=3

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

4.3.2 Technique specific parameters

E start	The potential where the scan starts and stops at.																		
E vertex1	The first potential where the scan direction reverses.																		
E vertex2	The second potential where the scan direction reverses.																		
	<p>The applicable potential range of each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr> <td>PalmSens4</td><td>-10 V to +10 V</td></tr> <tr> <td>EmStat1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>EmStat3</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat3+</td><td>-4 V to +4 V</td></tr> <tr> <td>EmStat4 LR</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat4 HR</td><td>-6 V to +6 V</td></tr> <tr> <td>EmStat Pico</td><td>High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V</td></tr> </table> <p>See also section Limitations and extra options for EmStat Pico on page 109.</p>	PalmSens1 and 2	-2 V to +2 V	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V	EmStat1 and 2	-2 V to +2 V	EmStat3	-3 V to +3 V	EmStat3+	-4 V to +4 V	EmStat4 LR	-3 V to +3 V	EmStat4 HR	-6 V to +6 V	EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V
PalmSens1 and 2	-2 V to +2 V																		
PalmSens3	-5 V to +5 V																		
PalmSens4	-10 V to +10 V																		
EmStat1 and 2	-2 V to +2 V																		
EmStat3	-3 V to +3 V																		
EmStat3+	-4 V to +4 V																		
EmStat4 LR	-3 V to +3 V																		
EmStat4 HR	-6 V to +6 V																		
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V																		
E step	<p>The potential step size.</p> <p>The applicable step range for each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>1 mV to 250 mV</td></tr> <tr> <td>PalmSens3</td><td>0.15 mV to 250 mV</td></tr> <tr> <td>PalmSens4</td><td>0.075 mV to 250 mV</td></tr> <tr> <td>EmStat1 and 2</td><td>0.1 mV to 250 mV</td></tr> <tr> <td>EmStat3(+)</td><td>0.125 mV to 250 mV</td></tr> <tr> <td>EmStat4 LR</td><td>0.1 mV to 250 mV</td></tr> <tr> <td>EmStat4 HR</td><td>0.183 mV to 250 mV</td></tr> <tr> <td>EmStat Pico</td><td>High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV</td></tr> </table>	PalmSens1 and 2	1 mV to 250 mV	PalmSens3	0.15 mV to 250 mV	PalmSens4	0.075 mV to 250 mV	EmStat1 and 2	0.1 mV to 250 mV	EmStat3(+)	0.125 mV to 250 mV	EmStat4 LR	0.1 mV to 250 mV	EmStat4 HR	0.183 mV to 250 mV	EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV		
PalmSens1 and 2	1 mV to 250 mV																		
PalmSens3	0.15 mV to 250 mV																		
PalmSens4	0.075 mV to 250 mV																		
EmStat1 and 2	0.1 mV to 250 mV																		
EmStat3(+)	0.125 mV to 250 mV																		
EmStat4 LR	0.1 mV to 250 mV																		
EmStat4 HR	0.183 mV to 250 mV																		
EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV																		
Scan rate	<p>The applied scan rate, also called sweep rate or sweep speed. The applicable range depends on the value of E step.</p> <p>The maximum scan rates for each instrument are:</p> <table> <tr> <td>PalmSens1 and 2</td><td>25 V/s (5 mV step)</td></tr> <tr> <td>PalmSens3</td><td>500 V/s (5 mV step)</td></tr> <tr> <td>PalmSens4</td><td>500 V/s (10 mV step)</td></tr> </table>	PalmSens1 and 2	25 V/s (5 mV step)	PalmSens3	500 V/s (5 mV step)	PalmSens4	500 V/s (10 mV step)												
PalmSens1 and 2	25 V/s (5 mV step)																		
PalmSens3	500 V/s (5 mV step)																		
PalmSens4	500 V/s (10 mV step)																		
Number of scans	The number of repetitions for this scan.																		

n averaged scans	The number of scans repetitions for averaging. In case 'Number of scans' is set to a value of more than 1, each scan in the plot is the result of an average of multiple scans, where the number of scans averaged is specified with this value.
n equil. scans	The number of equilibration scans. During these scans, no data is recorded.

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

As in CV, by default, the result shown is a graph of Current versus Potential. Successive cycles will be displayed in different colors.

4.4 AC Voltammetry (ACV)

Supported instruments:

- PalmSens series
- EmStat4 LR and HR

4.4.1 Description

In AC Voltammetry a potential scan is made with a superimposed sine wave which has a relatively small amplitude (normally 5 – 10 mV) and a frequency of 10 – 2000 Hz.

The AC signal superimposed on the DC-potential results in an AC response ($i_{ac\ rms}$). The resulting AC response is plotted against the potential.





Figure 79 Potential versus time during AC Voltammetry

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

4.4.2 Technique specific parameters

E begin	The potential where the measurement starts at.																
E end	The potential where the measurement stops at.																
	<p>The applicable potential range of each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr> <td>PalmSens4</td><td>-10 V to +10 V</td></tr> <tr> <td>EmStat4 LR</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat4 HR</td><td>-6 V to +6 V</td></tr> </table>	PalmSens1 and 2	-2 V to +2 V	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V	EmStat4 LR	-3 V to +3 V	EmStat4 HR	-6 V to +6 V						
PalmSens1 and 2	-2 V to +2 V																
PalmSens3	-5 V to +5 V																
PalmSens4	-10 V to +10 V																
EmStat4 LR	-3 V to +3 V																
EmStat4 HR	-6 V to +6 V																
E step	<p>The step potential size.</p> <p>The applicable step range for each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>1 mV to 250 mV</td></tr> <tr> <td>PalmSens3</td><td>0.15 mV to 250 mV</td></tr> <tr> <td>PalmSens4</td><td>0.075 mV to 250 mV</td></tr> <tr> <td>EmStat1 and 2</td><td>0.1 mV to 250 mV</td></tr> <tr> <td>EmStat3(+)</td><td>0.125 mV to 250 mV</td></tr> <tr> <td>EmStat4 LR</td><td>0.1 mV to 250 mV</td></tr> <tr> <td>EmStat4 HR</td><td>0.183 mV to 250 mV</td></tr> <tr> <td>EmStat Pico</td><td>High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV</td></tr> </table>	PalmSens1 and 2	1 mV to 250 mV	PalmSens3	0.15 mV to 250 mV	PalmSens4	0.075 mV to 250 mV	EmStat1 and 2	0.1 mV to 250 mV	EmStat3(+)	0.125 mV to 250 mV	EmStat4 LR	0.1 mV to 250 mV	EmStat4 HR	0.183 mV to 250 mV	EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV
PalmSens1 and 2	1 mV to 250 mV																
PalmSens3	0.15 mV to 250 mV																
PalmSens4	0.075 mV to 250 mV																
EmStat1 and 2	0.1 mV to 250 mV																
EmStat3(+)	0.125 mV to 250 mV																
EmStat4 LR	0.1 mV to 250 mV																
EmStat4 HR	0.183 mV to 250 mV																
EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV																
E ac	<p>The RMS amplitude of the applied sine wave.</p> <p><i>Applicable amplitude range: 1 mV to 250 mV</i></p>																
Frequency	<p>The frequency of the applied ac signal.</p> <p>The applicable frequency range for each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>0.12 Hz to 250 Hz</td></tr> <tr> <td>PalmSens3</td><td>1 Hz to 2000 Hz</td></tr> <tr> <td>PalmSens4</td><td>1 Hz to 2000 Hz</td></tr> <tr> <td>EmStat4 LR and HR</td><td>10 mHz to 50 kHz</td></tr> </table>	PalmSens1 and 2	0.12 Hz to 250 Hz	PalmSens3	1 Hz to 2000 Hz	PalmSens4	1 Hz to 2000 Hz	EmStat4 LR and HR	10 mHz to 50 kHz								
PalmSens1 and 2	0.12 Hz to 250 Hz																
PalmSens3	1 Hz to 2000 Hz																
PalmSens4	1 Hz to 2000 Hz																
EmStat4 LR and HR	10 mHz to 50 kHz																

Scan rate	The applied scan rate, also called sweep rate or sweep speed. The applicable range depends on the value of E step.	
	The applicable scan rates for each instrument are:	
	PalmSens1 and 2	1 mV/s (1 mV step) to 25 mV/s (5 mV step)
	PalmSens3	0.2 mV/s (1 mV step) to 50 mV/s (5 mV step)
	PalmSens4	0.2 mV/s (1 mV step) to 1 V/s (250 mV step)
	EmStat4 LR and HR	0.2 mV/s (1 mV step) to 1 V/s (250 mV step)
Measure DC current	When checked, the direct current (DC) will be measured separately and added to the plot as an additional curve.	

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

4.5 Differential Pulse Voltammetry (DPV)

Supported instruments:

- PalmSens series
- EmStat series
- Sensit series

4.5.1 Description

In Differential Pulse Voltammetry a potential scan is made using pulses with a constant amplitude of E pulse superimposed on the dc-potential. The amplitude is mostly in the range of 5 – 50 mV.




Figure 80 Potential versus time during Differential Pulse Voltammetry

The interval time between the pulses is equal to 'E step' / 'scan rate'.

The current is sampled twice in each step: one time just before applying the pulse and one time at the end of the pulse. The difference between these two current samples is plotted versus the potential.

The obtained current is proportional to the derivative of the curve obtained using Linear Sweep or Normal Pulse Voltammetry. A DPV thus has a peak-shaped curve. The peak height is (normally) proportional to the concentration in the solution. Be aware that the peak is not the redox potential!

$$E_{peak} = E_{1/2} - \frac{E_{pulse}}{2}$$

4.5.2 Measuring

As in Normal Pulse Voltammetry (NPV) the diffusion layer thickness increases with pulse time, thus the current will be lower when the pulse time is increased. However, a short pulse time will result in an increased capacitive current and therefore give a higher (non-linear) baseline.

In trace analysis, it is important to apply pulses with optimal pulse times. In general, the optimal value must be found by varying the pulse time.

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.




Figure 81 Typical DPV plot

4.5.3 Technique specific parameters

E begin	The potential where the measurement starts at.
E end	The potential where the measurement stops at.
	The applicable potential range of each instrument is:
PalmSens1 and 2	-2 V to +2 V
PalmSens3	-5 V to +5 V
PalmSens4	-10 V to +10 V
EmStat1 and 2	-2 V to +2 V
EmStat3	-3 V to +3 V
EmStat3+	-4 V to +4 V
EmStat4 LR	-3 V to +3 V
EmStat4 HR	-6 V to +6 V
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V

See also section [Limitations and extra options for EmStat Pico](#) on page 109.


E step	The potential step size.	
	The applicable step range for each instrument is:	
	PalmSens1 and 2	1 mV to 250 mV
	PalmSens3	0.15 mV to 250 mV
	PalmSens4	0.075 mV to 250 mV
	EmStat1 and 2	0.1 mV to 250 mV
	EmStat3(+)	0.125 mV to 250 mV
	EmStat4 LR	0.1 mV to 250 mV
	EmStat4 HR	0.183 mV to 250 mV
	EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV
Scan rate	The applied scan rate, also called as sweep rate or sweep speed. The maximum scan rate depends on the value of E step and t pulse. The scan rate must be $< (E \text{ step} / 2 / t \text{ pulse})$.	
	<i>For PalmSens1 and 2:</i> <i>In case the scan rate is so low that the time between two measured points is longer than approx. 0.05 s, the measured data points are displayed during the measurement. In other cases, the measurement is completed before the points are shown.</i>	
	The applicable scan rates for each instrument are:	
	PalmSens1 and 2	0.2 mV/s (0.1 mV step) to 50 mV/s (5 mV step)
	PalmSens3	0.02 mV/s (0.15 mV step) to 25 V/s (0.25 V step)
	PalmSens4	0.02 mV/s (0.07 mV step) to 25 V/s (0.25 V step)
	EmStat1 and 2	0.02 mV/s (0.1 mV step) to 5 V/s (5 mV step)
	EmStat3	0.025 mV/s (0.125 mV step) to 5 V/s (5 mV step)
	EmStat3+	0.025 mV/s (0.125 mV step) to 5 V/s (5 mV step)
	EmStat4 LR	0.1 mV/s (0.1 mV step) to 1 V/s (5 mV step)
	EmStat4 HR	0.1 mV/s (0.183 mV step) to 1 V/s (5 mV step)
	EmStat Pico	0.01 mV/s to 10 V/s (10 mV step)
E pulse	The pulse potential height.	
	The applicable pulse potential range for all instruments:	
	E pulse	1 mV to 250 mV

t pulse

The duration of the potential pulse. This duration needs to be set shorter than $0.5 \times \text{interval time}$ where the interval time is equal to $E \text{ step} / \text{Scan rate}$.

The applicable pulse time range for each instrument is:

PalmSens1 and 2	10 ms to 1 s
PalmSens3	0.2 ms to 1 s
PalmSens4	0.4 ms to 1 s
EmStat1, 2 and 3(+)	5 ms to 1 s
EmStat4 LR and HR	0.4 ms to 0.3 s
EmStat Pico	2 ms to 0.3 s

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

4.6 Square Wave Voltammetry (SWV)

Supported instruments:

- PalmSens series
- EmStat series
- Sensit series

4.6.1 Description

Square wave Voltammetry is in fact a special version of DPV.

DPV is SWV when t_{pulse} is equal to $\text{interval}/2$ (see DPV). The interval time is the inverse of the frequency (freq): $t_{\text{interval}} = 1 / \text{freq}$. As in DPV, the pulse amplitude is also normally in the range of 5 - 25 or 50 mV.




Figure 82 Potential versus time during Square Wave Voltammetry

4.6.2 Measuring

In trace analysis, it is important to apply a frequency with optimal values. As in DPV the optimal value must be found by varying the frequency.

SWV is sometimes used to measure the reaction rate of the electrode reaction. In this case, a plot of the currents observed in the positive potential pulses is plotted next to the currents measured in the negative potential pulses (forward and reverse currents).

The shape of these curves shows how reversible or how fast the electrode reaction is.

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.





Figure 83 Typical SWV plot

4.6.3 Technique specific parameters

E begin	Potential where scan starts.
E end	Potential where measurement stops.
	The applicable potential range of each instrument is:
PalmSens1 and 2	-2 V to +2 V
PalmSens3	-5 V to +5 V
PalmSens4	-10 V to +10 V
EmStat1 and 2	-2 V to +2 V
EmStat3	-3 V to +3 V
EmStat3+	-4 V to +4 V
EmStat4 LR	-3 V to +3 V
EmStat4 HR	-6 V to +6 V
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V

See also section [Limitations and extra options for EmStat Pico](#) on page 109.

E step	Step potential																
	The applicable step range for each instrument is:																
	<table><tr><td>PalmSens1 and 2</td><td>1 mV to 250 mV</td></tr><tr><td>PalmSens3</td><td>0.15 mV to 250 mV</td></tr><tr><td>PalmSens4</td><td>0.075 mV to 250 mV</td></tr><tr><td>EmStat1 and 2</td><td>0.1 mV to 250 mV</td></tr><tr><td>EmStat3(+)</td><td>0.125 mV to 250 mV</td></tr><tr><td>EmStat4 LR</td><td>0.1 mV to 250 mV</td></tr><tr><td>EmStat4 HR</td><td>0.183 mV to 250 mV</td></tr><tr><td>EmStat Pico</td><td>High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV</td></tr></table>	PalmSens1 and 2	1 mV to 250 mV	PalmSens3	0.15 mV to 250 mV	PalmSens4	0.075 mV to 250 mV	EmStat1 and 2	0.1 mV to 250 mV	EmStat3(+)	0.125 mV to 250 mV	EmStat4 LR	0.1 mV to 250 mV	EmStat4 HR	0.183 mV to 250 mV	EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV
	PalmSens1 and 2	1 mV to 250 mV															
	PalmSens3	0.15 mV to 250 mV															
	PalmSens4	0.075 mV to 250 mV															
	EmStat1 and 2	0.1 mV to 250 mV															
	EmStat3(+)	0.125 mV to 250 mV															
	EmStat4 LR	0.1 mV to 250 mV															
	EmStat4 HR	0.183 mV to 250 mV															
EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV																
Amplitude	<p>The amplitude of the square wave pulse. Values are half peak-to-peak.</p> <p><i>The amplitude range for all instruments: 0.1 mV to 250 mV</i></p>																
Frequency	<p>The frequency of the square wave.</p> <p>The applicable frequency range of each instrument is:</p>																
	<table><tr><td>PalmSens1 and 2</td><td>1 Hz to 400 Hz</td></tr><tr><td>PalmSens3</td><td>1 Hz to 1250 Hz</td></tr><tr><td>PalmSens4</td><td>1 Hz to 1250 Hz</td></tr><tr><td>EmStat1, 2 and 3(+)</td><td>1 Hz to 500 Hz</td></tr><tr><td>EmStat4 LR and HR</td><td>1 Hz to 2500 Hz</td></tr><tr><td>EmStat Pico</td><td>1 Hz to 500 Hz</td></tr></table>	PalmSens1 and 2	1 Hz to 400 Hz	PalmSens3	1 Hz to 1250 Hz	PalmSens4	1 Hz to 1250 Hz	EmStat1, 2 and 3(+)	1 Hz to 500 Hz	EmStat4 LR and HR	1 Hz to 2500 Hz	EmStat Pico	1 Hz to 500 Hz				
	PalmSens1 and 2	1 Hz to 400 Hz															
	PalmSens3	1 Hz to 1250 Hz															
	PalmSens4	1 Hz to 1250 Hz															
	EmStat1, 2 and 3(+)	1 Hz to 500 Hz															
	EmStat4 LR and HR	1 Hz to 2500 Hz															
EmStat Pico	1 Hz to 500 Hz																
Measure i forward/reverse	<p>If this option is enabled the plot will show two separate curves for the measured forward and reverse current.</p> <p><i>Not supported by PalmSens1 and PalmSens2.</i></p>																

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

4.7 Normal Pulse Voltammetry (NPV)

Supported instruments:

- PalmSens series
- EmStat series
- Sensit series

4.7.1 Description

In Normal Pulse Voltammetry (NPV) a potential scan is conducted in pulses by consistently increasing the pulse amplitude. With Normal Pulse Voltammetry the influence of diffusion

limitation on your i-E curve (Cottrel behavior) is removed. NPV is normally more sensitive than LSV, since the diffusion layer thickness will be smaller, resulting in a higher faradaic current.




Figure 84 Potential versus time during Normal Pulse Voltammetry

At the first potential step, the pulse is equal to 'E step', at next twice the value 'E step', until the end where the pulse is 'E begin' + n * 'E step' is equal to 'E end', where

$$n = (E_{\text{end}} - E_{\text{begin}}) / E_{\text{step}}$$

The pulse time t_{pulse} is specified by the user but must not exceed half the interval time.

So, the $t_{\text{pulse}} \leq E_{\text{step}} / (\text{scan rate} * 2)$.

4.7.2 Measuring

Since the diffusion layer thickness increases with time, the sampled current will be lower when the pulse time is increased. However, a short pulse time will result in an increased capacitive current and therefore give a higher (non-linear) baseline.

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.




Figure 85 Typical NPV plot

4.7.3 Technique specific parameters

E begin	The potential where the measurement starts at.
E end	The potential where the measurement stops at.
	The applicable potential range of each instrument is:
PalmSens1 and 2	-2 V to +2 V
PalmSens3	-5 V to +5 V
PalmSens4	-10 V to +10 V
EmStat1 and 2	-2 V to +2 V
EmStat3	-3 V to +3 V
EmStat3+	-4 V to +4 V
EmStat4 LR	-3 V to +3 V
EmStat4 HR	-6 V to +6 V
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V
See also section Limitations and extra options for EmStat Pico on page 109.	

E step	<p>The potential step size.</p> <p>The applicable step range for each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>1 mV to 250 mV</td></tr> <tr> <td>PalmSens3</td><td>0.15 mV to 250 mV</td></tr> <tr> <td>PalmSens4</td><td>0.075 mV to 250 mV</td></tr> <tr> <td>EmStat1 and 2</td><td>0.1 mV to 250 mV</td></tr> <tr> <td>EmStat3(+)</td><td>0.125 mV to 250 mV</td></tr> <tr> <td>EmStat4 LR</td><td>0.1 mV to 250 mV</td></tr> <tr> <td>EmStat4 HR</td><td>0.183 mV to 250 mV</td></tr> <tr> <td>EmStat Pico</td><td>High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV</td></tr> </table>	PalmSens1 and 2	1 mV to 250 mV	PalmSens3	0.15 mV to 250 mV	PalmSens4	0.075 mV to 250 mV	EmStat1 and 2	0.1 mV to 250 mV	EmStat3(+)	0.125 mV to 250 mV	EmStat4 LR	0.1 mV to 250 mV	EmStat4 HR	0.183 mV to 250 mV	EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV
PalmSens1 and 2	1 mV to 250 mV																
PalmSens3	0.15 mV to 250 mV																
PalmSens4	0.075 mV to 250 mV																
EmStat1 and 2	0.1 mV to 250 mV																
EmStat3(+)	0.125 mV to 250 mV																
EmStat4 LR	0.1 mV to 250 mV																
EmStat4 HR	0.183 mV to 250 mV																
EmStat Pico	High Speed mode: 0.395 mV to 250 mV Low Speed mode: 0.537 mV to 250 mV																
t pulse	<p>The duration of the potential pulse. This duration needs to be set shorter than $0.5 \cdot \text{interval time}$ where the interval time is equal to $E \text{ step} / \text{Scan rate}$.</p> <p>The applicable pulse time range for each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>10 ms to 1 s</td></tr> <tr> <td>PalmSens3</td><td>0.2 ms to 1 s</td></tr> <tr> <td>PalmSens4</td><td>0.4 ms to 1 s</td></tr> <tr> <td>EmStat1, 2 and 3(+)</td><td>5 ms to 1 s</td></tr> <tr> <td>EmStat4 LR and HR</td><td>0.4 ms to 0.3 s</td></tr> <tr> <td>EmStat Pico</td><td>2 ms to 0.3 s</td></tr> </table>	PalmSens1 and 2	10 ms to 1 s	PalmSens3	0.2 ms to 1 s	PalmSens4	0.4 ms to 1 s	EmStat1, 2 and 3(+)	5 ms to 1 s	EmStat4 LR and HR	0.4 ms to 0.3 s	EmStat Pico	2 ms to 0.3 s				
PalmSens1 and 2	10 ms to 1 s																
PalmSens3	0.2 ms to 1 s																
PalmSens4	0.4 ms to 1 s																
EmStat1, 2 and 3(+)	5 ms to 1 s																
EmStat4 LR and HR	0.4 ms to 0.3 s																
EmStat Pico	2 ms to 0.3 s																
Scan rate	<p>The applied scan rate. The maximum scan rate depends on the value of E step and t pulse.</p> <p>The scan rate must be $< (E \text{ step} / 2 / t \text{ pulse})$.</p> <p><i>For PalmSens2:</i> <i>In case the scan rate is so low that the time between two measured points is longer than approx. 0.05 s, the measured data points are displayed during the measurement. In other cases, the measurement is completed before the points are shown.</i></p>																

4.8 Chronoamperometry (CA)

Also known as Amperometric Detection, or Potentiostatic Step.

Supported instruments:

- PalmSens series
- EmStat series
- Sensit series

4.8.1 Description

The simplest, but widely used measurement technique is Chronoamperometry (or Amperometric Detection). Many sensors, like those for glucose or oxygen, require this technique. The instrument applies a constant DC-potential (E_{dc}) and the current is measured with constant interval times. The result is displayed in real time.




Figure 86 Potential versus time during Chronoamperometry

4.8.2 Measuring

Chronoamperometry is widely utilized in many applications, including highly sensitive electrochemical sensors (e.g., glucose biosensors), corrosion studies for assessing materials' resistance, as well as basic electrochemistry applications like studying electrode kinetics and reaction mechanisms, and electrodeposition processes for precise coating control in industries such as electronics, automotive, and aerospace.

The technique is also applied when electrochemical detection is used with a flow cell or flow injection cell (FIA). While batch measurements result in current levels, which depend on the concentration, a FIA setup shows peaks in the current.

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.




Figure 87 Typical CA curve obtained using capillary electrophoresis


4.8.3 Technique specific parameters

E dc	<p>The potential applied during the measurement. The applicable potential range of each instrument is:</p> <table border="1" data-bbox="544 1025 1353 1458"> <tbody> <tr> <td>PalmSens1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr> <td>PalmSens4</td><td>-10 V to +10 V</td></tr> <tr> <td>EmStat1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>EmStat3</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat3+</td><td>-4 V to +4 V</td></tr> <tr> <td>EmStat4 LR</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat4 HR</td><td>-6 V to +6 V</td></tr> <tr> <td>EmStat Pico</td><td>High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V</td></tr> </tbody> </table> <p>See also section Limitations and extra options for EmStat Pico on page 109.</p>	PalmSens1 and 2	-2 V to +2 V	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V	EmStat1 and 2	-2 V to +2 V	EmStat3	-3 V to +3 V	EmStat3+	-4 V to +4 V	EmStat4 LR	-3 V to +3 V	EmStat4 HR	-6 V to +6 V	EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V
PalmSens1 and 2	-2 V to +2 V																		
PalmSens3	-5 V to +5 V																		
PalmSens4	-10 V to +10 V																		
EmStat1 and 2	-2 V to +2 V																		
EmStat3	-3 V to +3 V																		
EmStat3+	-4 V to +4 V																		
EmStat4 LR	-3 V to +3 V																		
EmStat4 HR	-6 V to +6 V																		
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V																		
t interval	<p>The time between two current samples.</p> <p>The applicable time interval range of each instrument is:</p> <table border="1" data-bbox="544 1675 1353 1942"> <tbody> <tr> <td>PalmSens1 and 2</td><td>1 ms to 300 s</td></tr> <tr> <td>PalmSens3</td><td>0.2 ms to 300 s</td></tr> <tr> <td>PalmSens4</td><td>0.4 ms to 300 s</td></tr> <tr> <td>EmStat1, 2 and 3(+)</td><td>1 ms to 300 s</td></tr> <tr> <td>EmStat4 LR and HR</td><td>0.4 ms to 300 s</td></tr> <tr> <td>EmStat Pico</td><td>1 ms to 300 s</td></tr> </tbody> </table>	PalmSens1 and 2	1 ms to 300 s	PalmSens3	0.2 ms to 300 s	PalmSens4	0.4 ms to 300 s	EmStat1, 2 and 3(+)	1 ms to 300 s	EmStat4 LR and HR	0.4 ms to 300 s	EmStat Pico	1 ms to 300 s						
PalmSens1 and 2	1 ms to 300 s																		
PalmSens3	0.2 ms to 300 s																		
PalmSens4	0.4 ms to 300 s																		
EmStat1, 2 and 3(+)	1 ms to 300 s																		
EmStat4 LR and HR	0.4 ms to 300 s																		
EmStat Pico	1 ms to 300 s																		

t run

The total run time of the measurement.

The minimum and maximum duration of a measurement:
5 * interval time to 1,000,000 seconds (ca. 278 hours)

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

If you need to run faster steps use Fast Amperometry (FAM) instead. See details in section [Fast Amperometry \(FAM\)](#) on page 82.

4.9 Zero Resistance Amperometry (ZRA)

Supported instruments:

- PalmSens series
- EmStat series
- Sensit series



Use the technique Chronoamperometry in our software to perform ZRA measurements.

4.9.1 Description

All our potentiostats can be used as a Zero Resistance Ammeter (ZRA). A ZRA measures the current flowing through it without adding any resistance or polarization. This means the current is measured without the ZRA influencing the current.

4.9.2 Setup

If a potentiostat is used as a ZRA, in most setups it needs to be floating (Galvanically Isolated). The potentiostat's ground cable, is used for the ZRA setup and if the potentiostat is grounded additionally elsewhere the current could flow through that connection instead of the required path.

There are two options for a setup with a floating potentiostat.

1. Using Bluetooth and battery for communication and power. This means the potentiostat has no connection to ground.
2. Using a galvanic isolation dongle, which is inserted between the USB port of the computer and the potentiostat, makes the potentiostat float. This way you have a USB connection for communication and power, but no ground connection.

4.9.3 Connections



The reference electrode (blue) and counter electrode (black) connections of the potentiostat need to be short-circuited and left disconnected to anything else.

The current flows through the working electrode (red) and the ground (green). The Working Electrode is connected to the – input of the current follower and the + input is connected to Ground. The voltage difference between the Working Electrode and ground is zero and the current will be measured. See the following picture for a schematic representation.




Figure 88 Electrode connections for a ZRA measurement



Make sure that the current will not exceed the specification of your instrument model to avoid any damage. There will be no automatic interruption for current overload in this mode!

4.9.4 Recording ZRA currents

If you just want to see the current without recording, it is sufficient to go to the Manual Control tab and switch on the Cell.

If you want to record the values, you can start a Chronoamperometry measurement. The applied potential doesn't matter, because the reference and counter electrode of the ZRA is not part of the measured circuit.

4.10 MultiStep Amperometry (MA)

Supported instruments:

- PalmSens series
- EmStat series
- Sensit series

4.10.1 Description

Multistep amperometry simply allows the user to specify the number of potential steps they want to apply and how long each step should last. Each step works exactly as a Chronoamperometry step. The current is continuously sampled with the specified interval time. A whole cycle of steps can be repeated several times.




Figure 89 Potential versus time during Multistep Amperometry

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

4.10.2 Cottrell plot

A plot of i vs. $t^{-1/2}$ is often referred to as the “Cottrell plot”.

For more information see section: [Cottrell plot](#) on page 182.

4.10.3 Technique specific parameters

t interval	<p>The time between two current samples.</p> <p>The applicable time interval range of each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>1 ms to 300 s</td></tr> <tr> <td>PalmSens3</td><td>0.2 ms to 300 s</td></tr> <tr> <td>PalmSens4</td><td>0.4 ms to 300 s</td></tr> <tr> <td>EmStat1, 2 and 3(+)</td><td>1 ms to 300 s</td></tr> <tr> <td>EmStat4 LR and HR</td><td>0.4 ms to 300 s</td></tr> <tr> <td>EmStat Pico</td><td>1 ms to 300 s</td></tr> </table>	PalmSens1 and 2	1 ms to 300 s	PalmSens3	0.2 ms to 300 s	PalmSens4	0.4 ms to 300 s	EmStat1, 2 and 3(+)	1 ms to 300 s	EmStat4 LR and HR	0.4 ms to 300 s	EmStat Pico	1 ms to 300 s
PalmSens1 and 2	1 ms to 300 s												
PalmSens3	0.2 ms to 300 s												
PalmSens4	0.4 ms to 300 s												
EmStat1, 2 and 3(+)	1 ms to 300 s												
EmStat4 LR and HR	0.4 ms to 300 s												
EmStat Pico	1 ms to 300 s												
Cycles	The number of cycles (repetitions).												
Levels	<p>The number of potentials to apply within a cycle. Switching between levels adds an overhead time.</p> <p>The average overhead time for each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>~80 ms *</td></tr> <tr> <td>PalmSens3</td><td>~80 ms *</td></tr> <tr> <td>PalmSens4</td><td>~80 ms *</td></tr> <tr> <td>EmStat1, 2 and 3(+)</td><td>~80 ms *</td></tr> <tr> <td>EmStat4 LR and HR</td><td>~1 ms</td></tr> <tr> <td>EmStat Pico</td><td>~1 ms</td></tr> </table> <div> <p>* Other running processes in Windows and periodically saving to the recovery file (in case this setting is enabled) may add significantly to the overhead time.</p> </div>	PalmSens1 and 2	~80 ms *	PalmSens3	~80 ms *	PalmSens4	~80 ms *	EmStat1, 2 and 3(+)	~80 ms *	EmStat4 LR and HR	~1 ms	EmStat Pico	~1 ms
PalmSens1 and 2	~80 ms *												
PalmSens3	~80 ms *												
PalmSens4	~80 ms *												
EmStat1, 2 and 3(+)	~80 ms *												
EmStat4 LR and HR	~1 ms												
EmStat Pico	~1 ms												

E level [n]	The potential level at which the current is recorded.	
	The applicable potential range of each instrument is:	
	PalmSens1 and 2	-2 V to +2 V
	PalmSens3	-5 V to +5 V
	PalmSens4	-10 V to +10 V
	EmStat1 and 2	-2 V to +2 V
	EmStat3	-3 V to +3 V
	EmStat3+	-4 V to +4 V
	EmStat4 LR	-3 V to +3 V
	EmStat4 HR	-6 V to +6 V
	EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V
See also section Limitations and extra options for EmStat Pico on page 109.		
t [n]	The duration of the applied potential.	
	The minimum and maximum duration of a measurement: 5 * interval time to 1,000,000 seconds (ca. 278 hours)	

Additional options can be enabled using the  button:

Use limits for each level	Adds additional fields to each E level block for proceeding to the next level if a specified current is reached.
Select for which levels to record data	Adds additional 'Record' checkboxes to each E level block. If the 'Record' checkbox is not checked, the specified potential will be applied for the given time, without recording the current.
Use triggers	Adds additional digital trigger settings to each E level block. See for more information section Digital triggers.

If you need to perform actions which require more complexity, you can utilize Scripting.

See also section: [Scripting](#) on page 243.

4.11 Fast Amperometry (FAM)

Supported instruments:

- PalmSens series
- EmStat4 LR and HR

4.11.1 Description

Fast amperometry is a form of Amperometric Detection (Chronoamperometry) but with very high sampling rates or respectively very short interval times.




Figure 90 Potential versus time during Fast Amperometry



Please note that if you combine two or more FAM steps in a row using MethodSCRIPT, there will be an overhead time of approximately ~1 ms between them.

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

4.11.2 Technique specific parameters

E equilibration	Equilibration potential at which the measurement starts.						
E dc	Potential of the pulse. Note that this value is not relative to E equilibration, given above. The current is continuously sampled during this stage.						
	The applicable potential range for PalmSens series is: <table border="1"> <tr> <td>PalmSens1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr> <td>PalmSens4</td><td>-10 V to +10 V</td></tr> </table>	PalmSens1 and 2	-2 V to +2 V	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V
PalmSens1 and 2	-2 V to +2 V						
PalmSens3	-5 V to +5 V						
PalmSens4	-10 V to +10 V						
t run	Total run time of the measurement. Applicable run time: 1 ms to 30 s						

t interval	The time between two current samples.	
	Minimum interval times:	
	PalmSens1 and 2	0.25 ms
	PalmSens3	0.01 ms
	PalmSens4	0.02 ms
	EmStat4 LR and HR	0.001 ms

4.12 Pulsed Amperometric Detection (PAD)

Also known as Pulsed Amperometry

Supported instruments:

- PalmSens series
- EmStat series
- Sensit series

4.12.1 Description

With Pulsed Amperometric Detection a series of pulses (pulse profile) is periodically repeated. Pulsed Amperometric Detection can be used when higher sensitivity is required. Using pulses instead of constant potential might result in higher faradaic currents. PAD is also used when the electrode surface has to be regenerated continuously, for instance, to remove adsorbents from the electrode surface.





Figure 91 Potential versus time during Pulsed Amperometric Detection

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

4.12.2 Technique specific parameters

E dc	The dc or base potential.	
	The applicable potential range of each instrument is:	
	PalmSens1 and 2	-2 V to +2 V
	PalmSens3	-5 V to +5 V
	PalmSens4	-10 V to +10 V
	EmStat1 and 2	-2 V to +2 V
	EmStat3	-3 V to +3 V
	EmStat3+	-4 V to +4 V
	EmStat4 LR	-3 V to +3 V
	EmStat4 HR	-6 V to +6 V
	EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V
See also section Limitations and extra options for EmStat Pico on page 109.		
E pulse	Potential in pulse. Note that this value is not relative to E dc, given above.	
	The applicable pulse potential range of each instrument is:	
	PalmSens1 and 2	-2 V to +2 V
	PalmSens3	-5 V to +5 V
	PalmSens4	-10 V to +10 V
	EmStat2	-2 V to +2 V
	EmStat3	-3 V to +3 V
	EmStat3+	-4 V to +4 V
t pulse	The pulse time.	
	The applicable pulse time range for each instrument is:	
	PalmSens1 and 2	10 ms to 1 s
	PalmSens3	0.2 ms to 1 s
	PalmSens4	0.4 ms to 1 s
	EmStat1, 2 and 3(+)	5 ms to 1 s
	EmStat4 LR and HR	0.4 ms to 0.3 s
	EmStat Pico	2 ms to 0.3 s
Mode	DC:	i(dc) measurement is performed at potential E dc
	pulse:	i(pulse) measurement is performed at potential E pulse
	differential:	i(dif) measurement is i(pulse) - i(dc)

t run	Total run time of the measurement. The minimum and maximum duration of a measurement: 5 * interval time to 1,000,000 seconds (ca. 278 hours)
t interval	The time between two current samples.

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

4.13 Multiple Pulse Amperometry (MPAD)

Supported instruments:

- PalmSens series
- EmStat1, 2 and 3(+)

4.13.1 Description

The Multiple Pulse Amperometry (MPAD) technique involves applying a series of voltage pulses to an electrode immersed in a sample solution, and the resulting current of one of the pulses is measured.

This technique can be used when higher sensitivity is required. Using pulses instead of constant potential might result in higher faradaic currents. MPAD is used when the electrode surface must be regenerated continuously, for instance, to remove adsorbents from the electrode surface. Some examples include detection and quantification of specific analytes, particularly in liquid chromatography and flow-injection analysis.

Advantages of MPAD include its ability to reduce background noise and interference, making it especially useful for detecting trace-level analytes in complex sample matrices. Additionally, the choice of potential levels and pulse durations can be optimized for specific analytes, enhancing the method's sensitivity and selectivity. MPAD is commonly employed in techniques like high-performance liquid chromatography (HPLC) and capillary electrophoresis (CE) for the detection of a wide range of analytes, including pharmaceuticals, environmental contaminants, and biomolecules.




Figure 92 Potential versus time during Multiple Pulse Amperometry

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

4.13.2 Technique specific parameters

E1(measure)	The first potential level in which the current is recorded																		
E2	The second applied potential level																		
E3	The third applied potential level																		
	The applicable potential range of each instrument is:																		
	<table border="1"> <tr> <td>PalmSens1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr> <td>PalmSens4</td><td>-10 V to +10 V</td></tr> <tr> <td>EmStat1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>EmStat3</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat3+</td><td>-4 V to +4 V</td></tr> <tr> <td>EmStat4 LR</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat4 HR</td><td>-6 V to +6 V</td></tr> <tr> <td>EmStat Pico</td><td>High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V</td></tr> </table>	PalmSens1 and 2	-2 V to +2 V	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V	EmStat1 and 2	-2 V to +2 V	EmStat3	-3 V to +3 V	EmStat3+	-4 V to +4 V	EmStat4 LR	-3 V to +3 V	EmStat4 HR	-6 V to +6 V	EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V
PalmSens1 and 2	-2 V to +2 V																		
PalmSens3	-5 V to +5 V																		
PalmSens4	-10 V to +10 V																		
EmStat1 and 2	-2 V to +2 V																		
EmStat3	-3 V to +3 V																		
EmStat3+	-4 V to +4 V																		
EmStat4 LR	-3 V to +3 V																		
EmStat4 HR	-6 V to +6 V																		
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V																		
	See also section Limitations and extra options for EmStat Pico on page 109.																		
t1	The duration of the first applied potential																		

t2	The duration of the second applied potential
t3	The duration of the third applied potential
	Values t1, t2 and t3 can be 0.1 to 2 s
t run	Total run time of the measurement. The minimum and maximum duration of a measurement: 5 * interval time to 1,000,000 seconds (ca. 278 hours)

4.14 Open Circuit Potentiometry (OCP)

Supported instruments:

- PalmSens series
- EmStat series
- Sensit series

4.14.1 Description

For Open Circuit Potentiometry there is no polarization, and the so-called Open Circuit Potential (OCP) is measured and recorded with constant interval times. The result is a curve of Potential versus Time, displayed in real time. The OCP is also called Open Circuit Voltage (OCV).


In corrosion, it is referred to as the "Corrosion Potential" (E_{corr}), but in this context, it specifically denotes the potential of a metal or electrode when exposed to a corrosive environment.

If you just want to see the OCP value without recording a curve, the information is available in the status bar.

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

4.14.2 Technique specific parameters

t run	Total run time of the measurement. The minimum and maximum duration of a measurement: 5 * interval time to 1,000,000 seconds (ca. 278 hours)												
t interval	The time between two potential samples. The applicable time interval range of each instrument is: <table border="1"> <tr> <td>PalmSens1 and 2</td><td>1 ms to 300 s</td></tr> <tr> <td>PalmSens3</td><td>0.2 ms to 300 s</td></tr> <tr> <td>PalmSens4</td><td>0.4 ms to 300 s</td></tr> <tr> <td>EmStat1, 2 and 3(+)</td><td>1 ms to 300 s</td></tr> <tr> <td>EmStat4 LR and HR</td><td>0.4 ms to 300 s</td></tr> <tr> <td>EmStat Pico</td><td>1 ms to 300 s</td></tr> </table>	PalmSens1 and 2	1 ms to 300 s	PalmSens3	0.2 ms to 300 s	PalmSens4	0.4 ms to 300 s	EmStat1, 2 and 3(+)	1 ms to 300 s	EmStat4 LR and HR	0.4 ms to 300 s	EmStat Pico	1 ms to 300 s
PalmSens1 and 2	1 ms to 300 s												
PalmSens3	0.2 ms to 300 s												
PalmSens4	0.4 ms to 300 s												
EmStat1, 2 and 3(+)	1 ms to 300 s												
EmStat4 LR and HR	0.4 ms to 300 s												
EmStat Pico	1 ms to 300 s												

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

4.15 Chronopotentiometry (CP)

Supported instruments:

- PalmSens series
- EmStat4 LR and HR

4.15.1 Description

Chronopotentiometry (CP), also called Galvanostatic Step, is an electrochemical technique that requires a galvanostat instead of a potentiostat. It is supported by many PalmSens instruments. In this method, a constant current is applied, and the resulting potential (voltage) is continuously recorded over time in a definite time interval. This technique is particularly useful for studying electrochemical reactions, kinetics, and processes under non-steady-state conditions, offering valuable insights into how the electrode potential evolves in response to the applied current. It serves as a fundamental tool in various electrochemical investigations and analyses.

This technique is a valuable tool in situations where a Reference Electrode (RE) is not available. When combined with a known electrode surface area, it offers a reliable and reproducible alternative for quantitative and comparative electrochemical studies. Its ability to maintain a constant current density, even with changes in the system's resistance, allows for highly reproducible experiments, ensuring precise control over current conditions and minimizing the impact of resistance variations. This feature makes it particularly valuable for quantitative analyses, such as determining reaction rate constants and charge transfer coefficients.

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.




Figure 93 Current versus time during Chronopotentiometry

4.15.2 Technique specific parameters

i applied	<p>The current to apply. The unit of the value is the selected current range, so if 10 μA is selected and 1.5 is given as a value, the applied current will be 15 μA.</p> <p>The applicable applied currents for each instrument are:</p> <table border="1"> <thead> <tr> <th></th><th>Applied</th><th>Available Ranges (CR)</th></tr> </thead> <tbody> <tr> <td>PalmSens2</td><td>-2 to 2 * CR</td><td>1 μA to 10 mA</td></tr> <tr> <td>PalmSens3</td><td>-3 to 3 * CR</td><td>1 μA to 10 mA</td></tr> <tr> <td>PalmSens4</td><td>-6 to 6 * CR</td><td>1 nA to 10 mA</td></tr> <tr> <td>EmStat4 LR</td><td>-3 to 3 * CR</td><td>10 nA, 1 μA, 100 μA, 10 mA</td></tr> <tr> <td>EmStat4 HR</td><td>-3 to 3 * CR</td><td>1 μA, 100 μA, 10 mA and 100 mA</td></tr> </tbody> </table>			Applied	Available Ranges (CR)	PalmSens2	-2 to 2 * CR	1 μ A to 10 mA	PalmSens3	-3 to 3 * CR	1 μ A to 10 mA	PalmSens4	-6 to 6 * CR	1 nA to 10 mA	EmStat4 LR	-3 to 3 * CR	10 nA, 1 μ A, 100 μ A, 10 mA	EmStat4 HR	-3 to 3 * CR	1 μ A, 100 μ A, 10 mA and 100 mA
	Applied	Available Ranges (CR)																		
PalmSens2	-2 to 2 * CR	1 μ A to 10 mA																		
PalmSens3	-3 to 3 * CR	1 μ A to 10 mA																		
PalmSens4	-6 to 6 * CR	1 nA to 10 mA																		
EmStat4 LR	-3 to 3 * CR	10 nA, 1 μ A, 100 μ A, 10 mA																		
EmStat4 HR	-3 to 3 * CR	1 μ A, 100 μ A, 10 mA and 100 mA																		
t run	<p>Total run time of the measurement.</p> <p>RuntimeRange</p>																			

t interval	The time between two potential samples.	
	The applicable time interval range of each instrument is:	
	PalmSens1 and 2	1 ms to 300 s
	PalmSens3	0.2 ms to 300 s
	PalmSens4	0.4 ms to 300 s
	EmStat1, 2 and 3(+)	1 ms to 300 s
	EmStat4 LR and HR	0.4 ms to 300 s
	EmStat Pico	1 ms to 300 s

4.16 Linear Sweep Potentiometry (LSP)

Supported instruments:

- PalmSens4
- EmStat4 LR and HR

4.16.1 Description

In Linear Sweep Potentiometry is very similar to Linear Sweep Voltammetry, but in this case, the current is controlled, instead of the potential. It is also known as Galvanodynamic Step. A current scan is performed from the begin current, 'i begin', to the end current, 'i end'. The scan is not really linear, but small current steps ('i step') are made. The potential is measured (sampled) during the last 25% interval period of each step. So, the number of points in the curve of the current versus potential is $(i_{\text{end}} - i_{\text{begin}}) / E_{\text{step}} + 1$.

The scan rate is specified in current per time (i.e. $\mu\text{A/s}$), which determines the time between two steps and thus the sampling time. The interval time is 'i step' / 'Scan rate'. I.e.: when 'i step' is 0.005 mA and the scan rate 0.1 mA/s the interval time is 0.05 s.

4.16.2 Measuring

As other Galvanostatic techniques, The LSP is particularly useful in systems where a Reference Electrode is not possible, but the Electrochemical Surface Area is well known. As a controlled current is applied to the electrode, this ensures a stable and consistent rate of reaction, which is advantageous when dealing with complex or multi-step electrochemical reactions. Galvanodynamic experiments also allow for the determination of rate constants and reaction mechanisms directly from the current-time data, facilitating quantitative analysis. Moreover, they are less susceptible to artifacts caused by double-layer charging effects and potential drift, often encountered in potentiodynamic experiments.

In some applications, it is important that the potential does not get too high. This might ruin the working electrode. If the current at which this will occur is not known, it is possible to specify a maximum potential value at which the scan stops. In this case, the end current specified by the user may not be reached.

See for more information section [Potential and current limits](#) on page 48.




Figure 94 Current versus time during Linear Sweep Potentiometry

During the measurement, the curve is shown on the screen in real-time. It is possible to abort the measurement, by pressing the abort button above the plot.


During a measurement, the use of the 'Pause' button will halt the scan until the same button is used again. This button is not available at higher scan rates.

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

4.16.3 Technique specific parameters

Applied current range	<p>The applied current range. If a range of 10 μA is selected and 1.5 is given as an applied current value, the applied current will be 15 μA.</p> <p>The applicable applied currents for each instrument are:</p> <table border="1"> <thead> <tr> <th></th><th>Applied</th><th>Available Ranges (CR)</th></tr> </thead> <tbody> <tr> <td>PalmSens2</td><td>-2 to 2 * CR</td><td>1 μA to 10 mA</td></tr> <tr> <td>PalmSens3</td><td>-3 to 3 * CR</td><td>1 μA to 10 mA</td></tr> <tr> <td>PalmSens4</td><td>-6 to 6 * CR</td><td>1 nA to 10 mA</td></tr> <tr> <td>EmStat4 LR</td><td>-3 to 3 * CR</td><td>10 nA, 1 μA, 100 μA, 10 mA</td></tr> <tr> <td>EmStat4 HR</td><td>-3 to 3 * CR</td><td>1 μA, 100 μA, 10 mA and 100 mA</td></tr> </tbody> </table>			Applied	Available Ranges (CR)	PalmSens2	-2 to 2 * CR	1 μA to 10 mA	PalmSens3	-3 to 3 * CR	1 μA to 10 mA	PalmSens4	-6 to 6 * CR	1 nA to 10 mA	EmStat4 LR	-3 to 3 * CR	10 nA, 1 μA , 100 μA , 10 mA	EmStat4 HR	-3 to 3 * CR	1 μA , 100 μA , 10 mA and 100 mA
	Applied	Available Ranges (CR)																		
PalmSens2	-2 to 2 * CR	1 μA to 10 mA																		
PalmSens3	-3 to 3 * CR	1 μA to 10 mA																		
PalmSens4	-6 to 6 * CR	1 nA to 10 mA																		
EmStat4 LR	-3 to 3 * CR	10 nA, 1 μA , 100 μA , 10 mA																		
EmStat4 HR	-3 to 3 * CR	1 μA , 100 μA , 10 mA and 100 mA																		
i begin	Current in the selected range where scan starts at.																			
i end	Current in the selected range where scan stops at.																			

i step	The current step size.	
	The applicable minimum 'i step' for each instrument is:	
	PalmSens1 and 2	$0.001 * CR$
	PalmSens3	$0.00015 * CR$
	PalmSens4	$0.000075 * CR$
	EmStat4 LR	$0.0001 * CR$
	EmStat4 HR	$0.000183 * CR$
Scan rate	The applied scan rate. The applicable range depends on the value of i step since the data acquisition rate is limited by the connected instrument.	

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

4.17 MultiStep Potentiometry (MP)

Supported instruments:

- PalmSens series
- EmStat4 LR and HR

4.17.1 Description

MultiStep Potentiometry simply allows the user to specify the number of current steps they want to apply and how long each step should last. The potential response is continuously sampled with the specified interval.

A whole cycle of steps can be repeated several times.

(Note: if only one current step is used, this technique is identical to Potentiometry which provides a larger range of measurement rates.)





Figure 95 Current versus time during Multistep Potentiometry

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.


4.17.2 Technique specific parameters

t interval	<p>The time between two potential samples.</p> <p>The applicable time interval range of each instrument is:</p> <table border="1"> <tbody> <tr> <td>PalmSens1 and 2</td><td>1 ms to 300 s</td></tr> <tr> <td>PalmSens3</td><td>0.2 ms to 300 s</td></tr> <tr> <td>PalmSens4</td><td>0.4 ms to 300 s</td></tr> <tr> <td>EmStat1, 2 and 3(+)</td><td>1 ms to 300 s</td></tr> <tr> <td>EmStat4 LR and HR</td><td>0.4 ms to 300 s</td></tr> <tr> <td>EmStat Pico</td><td>1 ms to 300 s</td></tr> </tbody> </table>	PalmSens1 and 2	1 ms to 300 s	PalmSens3	0.2 ms to 300 s	PalmSens4	0.4 ms to 300 s	EmStat1, 2 and 3(+)	1 ms to 300 s	EmStat4 LR and HR	0.4 ms to 300 s	EmStat Pico	1 ms to 300 s
PalmSens1 and 2	1 ms to 300 s												
PalmSens3	0.2 ms to 300 s												
PalmSens4	0.4 ms to 300 s												
EmStat1, 2 and 3(+)	1 ms to 300 s												
EmStat4 LR and HR	0.4 ms to 300 s												
EmStat Pico	1 ms to 300 s												
Cycles	The number of repetitions.												

Levels	The number of potentials to apply within a cycle.		
	The average overhead time for each instrument is:		
	PalmSens1 and 2	~80 ms *	
	PalmSens3	~80 ms *	
	PalmSens4	~80 ms *	
	EmStat1, 2 and 3(+)	~80 ms *	
	EmStat4 LR and HR	~1 ms	
	EmStat Pico	~1 ms	
<div>* Other running processes in Windows and periodically saving to the recovery file (in case this setting is enabled) may add significantly to the overhead time.</div>			
i level [n]	The current level at which the potential is recorded.		
	The applicable applied currents for each instrument are:		
		Applied	Available Ranges (CR)
	PalmSens2	-2 to 2 * CR	1 µA to 10 mA
	PalmSens3	-3 to 3 * CR	1 µA to 10 mA
	PalmSens4	-6 to 6 * CR	1 nA to 10 mA
	EmStat4 LR	-3 to 3 * CR	10 nA, 1 µA, 100 µA, 10 mA
	EmStat4 HR	-3 to 3 * CR	1 µA, 100 µA, 10 mA and 100 mA
t [n]	The duration of the applied current.		
	The minimum and maximum duration of a measurement: 5 * interval time to 1,000,000 seconds (ca. 278 hours)		

Additional options can be enabled using the  button:

Use limits for each level	Adds additional fields to each i level block for proceeding to the next level if a specified potential is reached.
Select for which levels to record data	Adds additional 'Record' checkboxes to each i level block. If the 'Record' checkbox is not checked, the specified current will be applied for the given time, without recording the potential.
Use triggers	Adds additional digital trigger settings to each i level block. See for more information section Digital triggers.

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

4.18 Stripping Chronopotentiometry (SCP)

Also known as Potentiometric Stripping Analysis (PSA)

Supported instruments:

- PalmSens series

4.18.1 Description

Chronopotentiometric Stripping or Stripping chronopotentiometry is a sensitive analytical technique.

The sequence of a stripping chronopotentiometry measurement:

1. Apply E condition, the conditioning potential, if t condition is not zero.
2. Apply E deposition, the deposition potential, if t deposition is not zero.
3. Apply E deposition and wait t equilibration seconds.
4. If Stripping Current = 0 then the cell is switched off, otherwise, the specified constant current is applied. The measurement with a rate of 40 kHz starts. The measurement stops when either the measured potential is below 'E end' or the run time is exceeded.




Figure 96 Stripping chronopotentiometry (SCP)

4.18.2 Measuring

The actual measurement can be done in two modes:

- Chemical stripping, using a chemical oxidant (or reductant)
- Using the instrument as a galvanostat and applying a constant stripping current.

In both cases, the potential versus time (E vs t) is measured. The resulting curve is re-calculated to the inverse derivative, so dt/dE (in s/V) vs E.

In case a component was deposited at the electrode surface, it requires current to oxidize or reduce. The higher the amount of absorbed component, the more electrical charge (time integral of the current) it requires, so the longer it takes to change the electrode potential.

The plot of dt/dE vs E therefore will show a peak the potential where the oxidation or reduction occurs.

The direction of the stripping current depends on the end potential (E end) minus deposition potential (E deposition). In case no deposition time was used, the sign is determined by the end potential (E end).

4.18.3 Technique specific parameters

E end	<p>Potential where measurement stops.</p> <p>The applicable potential range for PalmSens series is:</p> <table border="1" data-bbox="544 1037 1318 1176"> <tr> <td>PalmSens1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr> <td>PalmSens4</td><td>-10 V to +10 V</td></tr> </table>	PalmSens1 and 2	-2 V to +2 V	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V
PalmSens1 and 2	-2 V to +2 V						
PalmSens3	-5 V to +5 V						
PalmSens4	-10 V to +10 V						
Stripping Current	<p>If specified as 0, the method is called chemical stripping otherwise it is constant current stripping.</p> <p>The direction of the stripping current is determined by E end minus E deposition.</p> <p>The applicable stripping current and potential range of each instrument is:</p> <table border="1" data-bbox="544 1496 1353 1720"> <tr> <td>PalmSens1 and 2</td><td>±0.001 μA to ±2 mA. -2 V to +2 V.</td></tr> <tr> <td>PalmSens3</td><td>±0.001 μA to ±30 mA. -5 V to +5 V.</td></tr> <tr> <td>PalmSens4</td><td>±1 nA to ±30 mA. -10 V to +10 V</td></tr> </table>	PalmSens1 and 2	±0.001 μ A to ±2 mA. -2 V to +2 V.	PalmSens3	±0.001 μ A to ±30 mA. -5 V to +5 V.	PalmSens4	±1 nA to ±30 mA. -10 V to +10 V
PalmSens1 and 2	±0.001 μ A to ±2 mA. -2 V to +2 V.						
PalmSens3	±0.001 μ A to ±30 mA. -5 V to +5 V.						
PalmSens4	±1 nA to ±30 mA. -10 V to +10 V						
Measurement time	<p>The maximum measurement time. This value should always exceed the required measurement time. It only limits the time of the measurement. When the potential response is erroneously and E end is not found within this time, the measurement is aborted.</p>						

4.19 Chronocoulometry (CC)

Supported instruments:

- PalmSens3 and PalmSens4
- EmStat series
- Sensit series

4.19.1 Description

Chronoamperometry (CA) and Chronocoulometry (CC) have the same potential waveform but are used for different purposes. In CA, the current is monitored as a function of time, whereas in CC, the charge is monitored as a function of time. This means that in essence, the potentiostat performs the same experiment as in CA, keeping the potential constant and registering the current over time. However, in CC, instead of displaying the current, the software displays the charge *versus* time. The charge is determined by integrating the current.

While the measured current in the CA shows a linear correlation with the reaction rate, the charge measured in CC shows a linear correlation with the number of converted reactants. This is described in Faraday's law. This makes CC useful during deposition or electrochemical synthesis. For the characterization of batteries and capacitors, CC is also helpful.

In analytical chemistry, CC is employed to determine the adsorbed number of active species in a solution of freely diffusing active species. Multiplying the Cottrell equation by t yields that the charge Q of the freely diffusing species' reaction is proportional to $t^{3/2}$. The total charge also includes the contribution of the charge stored in the electrochemical double layer and the charge due to reactions of adsorbed species. These two effects are a lot faster than the reaction of free diffusing species. Plotting Q versus $t^{1/2}$ (Anson plot) delivers (ideally) a jump in charge followed by a linear increase. When the linear part is extrapolated, the intersection with the charge axis delivers the contribution of the double layer and adsorbed species. A previous blank measurement allows the determination of the double layer contribution and thus the calculation of the adsorbed species contribution.

To ensure that the charge step at the beginning is recorded, this measurement is conducted in two steps. The first step at a reaction-free potential followed by the step where the reaction is initiated, just like the classic Cottrell experiment.




Figure 97 Potential versus time during Chronocoulometry

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

4.19.2 Anson plot

For more information about the Anson plot, see section [Anson plot](#) on page 183.

4.19.3 Technique specific parameters


t interval

The time between two current samples.

The applicable time interval range of each instrument is:

PalmSens1 and 2	1 ms to 300 s
PalmSens3	0.2 ms to 300 s
PalmSens4	0.4 ms to 300 s
EmStat1, 2 and 3(+)	1 ms to 300 s
EmStat4 LR and HR	0.4 ms to 300 s
EmStat Pico	1 ms to 300 s

E step 1	<p>Potential applied during the first potential step.</p> <p>The applicable potential range of each instrument is:</p> <table border="1" data-bbox="544 331 1353 763"> <tr><td>PalmSens1 and 2</td><td>-2 V to +2 V</td></tr> <tr><td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr><td>PalmSens4</td><td>-10 V to +10 V</td></tr> <tr><td>EmStat1 and 2</td><td>-2 V to +2 V</td></tr> <tr><td>EmStat3</td><td>-3 V to +3 V</td></tr> <tr><td>EmStat3+</td><td>-4 V to +4 V</td></tr> <tr><td>EmStat4 LR</td><td>-3 V to +3 V</td></tr> <tr><td>EmStat4 HR</td><td>-6 V to +6 V</td></tr> <tr><td>EmStat Pico</td><td>High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V</td></tr> </table> <p>See also section Limitations and extra options for EmStat Pico on page 109.</p>	PalmSens1 and 2	-2 V to +2 V	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V	EmStat1 and 2	-2 V to +2 V	EmStat3	-3 V to +3 V	EmStat3+	-4 V to +4 V	EmStat4 LR	-3 V to +3 V	EmStat4 HR	-6 V to +6 V	EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V
PalmSens1 and 2	-2 V to +2 V																		
PalmSens3	-5 V to +5 V																		
PalmSens4	-10 V to +10 V																		
EmStat1 and 2	-2 V to +2 V																		
EmStat3	-3 V to +3 V																		
EmStat3+	-4 V to +4 V																		
EmStat4 LR	-3 V to +3 V																		
EmStat4 HR	-6 V to +6 V																		
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V																		
t 1	<p>Duration of the first potential step.</p> <p>The minimum and maximum duration of a measurement: 5 * interval time to 1,000,000 seconds (ca. 278 hours)</p>																		
E step 2	<p>Potential applied during the second potential step.</p> <p>The applicable potential range of each instrument is:</p> <table border="1" data-bbox="544 1167 1353 1599"> <tr><td>PalmSens1 and 2</td><td>-2 V to +2 V</td></tr> <tr><td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr><td>PalmSens4</td><td>-10 V to +10 V</td></tr> <tr><td>EmStat1 and 2</td><td>-2 V to +2 V</td></tr> <tr><td>EmStat3</td><td>-3 V to +3 V</td></tr> <tr><td>EmStat3+</td><td>-4 V to +4 V</td></tr> <tr><td>EmStat4 LR</td><td>-3 V to +3 V</td></tr> <tr><td>EmStat4 HR</td><td>-6 V to +6 V</td></tr> <tr><td>EmStat Pico</td><td>High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V</td></tr> </table> <p>See also section Limitations and extra options for EmStat Pico on page 109.</p>	PalmSens1 and 2	-2 V to +2 V	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V	EmStat1 and 2	-2 V to +2 V	EmStat3	-3 V to +3 V	EmStat3+	-4 V to +4 V	EmStat4 LR	-3 V to +3 V	EmStat4 HR	-6 V to +6 V	EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V
PalmSens1 and 2	-2 V to +2 V																		
PalmSens3	-5 V to +5 V																		
PalmSens4	-10 V to +10 V																		
EmStat1 and 2	-2 V to +2 V																		
EmStat3	-3 V to +3 V																		
EmStat3+	-4 V to +4 V																		
EmStat4 LR	-3 V to +3 V																		
EmStat4 HR	-6 V to +6 V																		
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V																		
t 2	<p>Duration of the second potential step.</p> <p>The minimum and maximum duration of a measurement: 5 * interval time to 1,000,000 seconds (ca. 278 hours)</p>																		

Additional options can be enabled using the  button. See for more information the applicable sub-sections in section [Measurements](#) on page 25.

4.19.4 Overhead time between steps

The delay between the execution of the two potential steps varies, depending on the instrument model. The following table gives an indication of the average overhead time for each model.

The average overhead time for each instrument is:

PalmSens1 and 2	~80 ms *
PalmSens3	~80 ms *
PalmSens4	~80 ms *
EmStat1, 2 and 3(+)	~80 ms *
EmStat4 LR and HR	~1 ms
EmStat Pico	~1 ms



* Other running processes in Windows and periodically saving to the recovery file (in case this setting is enabled) may add significantly to the overhead time.

4.20 Mixed Mode (MM)

Supported instruments:

- PalmSens series
- EmStat1, 2 and 3(+) (limited)
- EmStat4 LR and HR

4.20.1 Description

Mixed Mode is a flexible technique that allows for switching between potentiostatic, galvanostatic, and open circuit measurements during a single run.

See section [Standard Measurement Sequence](#) on page 32 for information about pretreatment and post measurement settings.

The Mixed Mode uses different stages similar to the levels during Multistep Amperometry or Potentiometry, but each stage can be galvanostatic or potentiostatic independent of the previous stage.

The available stage types are 'ConstantE', 'ConstantI', 'SweepE', 'OpenCircuit' and 'Impedance'. 'SweepE' offers a potential linear sweep (ramp), so as a regular LSV step. During an Impedance stage the impedance is measured by applying a small AC potential superimposed with a DC potential. This corresponds to an EIS single frequency step (Scan type = Default, Frequency type = Fixed).

Each stage can use the previous stage's potential as a reference point, for example, a constant current is applied for a fixed period and afterward, the reached potential is kept constant for a fixed period.

Furthermore, each stage can end, because a fixed period has elapsed, or certain criteria are met. Available criteria include reaching a maximum current, minimum current, maximum potential, and minimum potential.

These modes are useful especially for energy conversion and storage research, i.e.: battery, solar cell, or supercapacitors research. A classic test for batteries involves charging and discharging them through multiple cycles. A constant current is applied, and the potential change is recorded. If a certain potential is reached, the next stage is triggered, which is usually applying the inverted constant current. This is repeated for multiple cycles. While this method could be performed with Multistep Potentiometry as soon as further steps are introduced the Mixed Mode is necessary. For example, to determine the capacity you would like to discharge a battery, but you need to take care that the terminal potential isn't crossed. First the battery is charged. You can apply a constant current until a set potential is reached, e.g.: the termination potential. After that, the termination potential is kept constant until a minimum current limit is reached or the time has elapsed. Then this process is repeated with a negative current to discharge the battery until a minimum potential is reached. Impedance steps can be inserted between to determine the parameters called State of Charge (SoC) and State of Health (SoH) of the battery.




Figure 98 Scheme of the potential and current profile during the above-described example.

Mixed Mode Settings

t interval: 60.0 s

Cycles: 3

Stage 1: Constant I

I applied: 1.0 * 100 mA

t run: 86400.0 s

☐ Proceed when E < 0.0 V

☒ Proceed when E > 1.7 V

Stage 2: Constant E

E dc: 1.7 V

t run: 3600.0 s

☐ E vs previous E

☒ Proceed when |I| < 50000.0 μ A

☐ Proceed when |I| > 0.0 μ A

Stage 3: Constant I

I applied: -1.0 * 100 mA

Figure 99 Part of the Method needed to perform the above-described method

Another example is to study how well a supercapacitor stores charge. First, the capacitor is charged with a fixed current followed by observing the OCP of the supercapacitor.

4.20.2 Technique specific parameters

t interval	The time between two current samples.												
	The applicable time interval range of each instrument is:												
	<table> <tr> <td>PalmSens1 and 2</td><td>1 ms to 300 s</td></tr> <tr> <td>PalmSens3</td><td>0.2 ms to 300 s</td></tr> <tr> <td>PalmSens4</td><td>0.4 ms to 300 s</td></tr> <tr> <td>EmStat1, 2 and 3(+)</td><td>1 ms to 300 s</td></tr> <tr> <td>EmStat4 LR and HR</td><td>0.4 ms to 300 s</td></tr> <tr> <td>EmStat Pico</td><td>1 ms to 300 s</td></tr> </table>	PalmSens1 and 2	1 ms to 300 s	PalmSens3	0.2 ms to 300 s	PalmSens4	0.4 ms to 300 s	EmStat1, 2 and 3(+)	1 ms to 300 s	EmStat4 LR and HR	0.4 ms to 300 s	EmStat Pico	1 ms to 300 s
PalmSens1 and 2	1 ms to 300 s												
PalmSens3	0.2 ms to 300 s												
PalmSens4	0.4 ms to 300 s												
EmStat1, 2 and 3(+)	1 ms to 300 s												
EmStat4 LR and HR	0.4 ms to 300 s												
EmStat Pico	1 ms to 300 s												
Cycles	The number of repetitions.												

4.20.3 Parameters for Constant E stage

E dc	<p>Potential during measurement.</p> <p>The applicable potential range of each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr> <td>PalmSens4</td><td>-10 V to +10 V</td></tr> <tr> <td>EmStat1 and 2</td><td>-2 V to +2 V</td></tr> <tr> <td>EmStat3</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat3+</td><td>-4 V to +4 V</td></tr> <tr> <td>EmStat4 LR</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat4 HR</td><td>-6 V to +6 V</td></tr> <tr> <td>EmStat Pico</td><td>High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V</td></tr> </table> <p>See also section Limitations and extra options for EmStat Pico on page 109.</p>	PalmSens1 and 2	-2 V to +2 V	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V	EmStat1 and 2	-2 V to +2 V	EmStat3	-3 V to +3 V	EmStat3+	-4 V to +4 V	EmStat4 LR	-3 V to +3 V	EmStat4 HR	-6 V to +6 V	EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V
PalmSens1 and 2	-2 V to +2 V																		
PalmSens3	-5 V to +5 V																		
PalmSens4	-10 V to +10 V																		
EmStat1 and 2	-2 V to +2 V																		
EmStat3	-3 V to +3 V																		
EmStat3+	-4 V to +4 V																		
EmStat4 LR	-3 V to +3 V																		
EmStat4 HR	-6 V to +6 V																		
EmStat Pico	High Speed mode: -1.7 V to +2 V Low Speed mode: -1.25 V to +2 V																		
t run	<p>The run time of the stage.</p> <p>The minimum and maximum duration of a measurement: 5 * interval time to 1,000,000 seconds (ca. 278 hours)</p>																		

4.20.4 Parameters for Constant I stage

i applied	<p>The current to apply. The unit of the value is the selected current range at the top. So if 10 uA is selected and 1.5 is given as a value, the applied current will be 15 uA.</p> <p>The applicable applied currents for each instrument are:</p> <table><tr><th></th><th>Applied</th><th>Available Ranges (CR)</th></tr><tr><td>PalmSens2</td><td>-2 to 2 * CR</td><td>1 µA to 10 mA</td></tr><tr><td>PalmSens3</td><td>-3 to 3 * CR</td><td>1 µA to 10 mA</td></tr><tr><td>PalmSens4</td><td>-6 to 6 * CR</td><td>1 nA to 10 mA</td></tr><tr><td>EmStat4 LR</td><td>-3 to 3 * CR</td><td>10 nA, 1 µA, 100 µA, 10 mA</td></tr><tr><td>EmStat4 HR</td><td>-3 to 3 * CR</td><td>1 µA, 100 µA, 10 mA and 100 mA</td></tr></table>		Applied	Available Ranges (CR)	PalmSens2	-2 to 2 * CR	1 µA to 10 mA	PalmSens3	-3 to 3 * CR	1 µA to 10 mA	PalmSens4	-6 to 6 * CR	1 nA to 10 mA	EmStat4 LR	-3 to 3 * CR	10 nA, 1 µA, 100 µA, 10 mA	EmStat4 HR	-3 to 3 * CR	1 µA, 100 µA, 10 mA and 100 mA
	Applied	Available Ranges (CR)																	
PalmSens2	-2 to 2 * CR	1 µA to 10 mA																	
PalmSens3	-3 to 3 * CR	1 µA to 10 mA																	
PalmSens4	-6 to 6 * CR	1 nA to 10 mA																	
EmStat4 LR	-3 to 3 * CR	10 nA, 1 µA, 100 µA, 10 mA																	
EmStat4 HR	-3 to 3 * CR	1 µA, 100 µA, 10 mA and 100 mA																	
t run	<p>The run time of the stage.</p> <p>The minimum and maximum duration of a measurement: 5 * interval time to 1,000,000 seconds (ca. 278 hours)</p>																		

4.20.5 Parameters for Open Circuit stage

t run	The run time of the stage. The minimum and maximum duration of a measurement: 5 * interval time to 1,000,000 seconds (ca. 278 hours)
--------------	--

4.20.6 Parameters for Sweep E stage

E begin	Potential where scan starts.																
E end	Potential where measurement stops.																
	PotentialRanges																
E step	Step potential StepRanges																
Scan rate	<p>The applied scan rate. The applicable range depends on the value of E step since the data acquisition rate is limited by the connected instrument.</p> <p>The applicable scan rates for each instrument are:</p> <table> <tr> <td>PalmSens1 and 2</td><td>1 mV/s (1 mV step) to 25 V/s (5 mV step)</td></tr> <tr> <td>PalmSens3</td><td>0.02 mV/s (0.15 mV step) to 500 V/s (5 mV step)</td></tr> <tr> <td>PalmSens4</td><td>0.02 mV/s (0.075 mV step) to 500 V/s (10 mV step)</td></tr> <tr> <td>EmStat1 and 2</td><td>0.02 mV/s (0.1 mV step) to 5 V/s (5 mV step)</td></tr> <tr> <td>EmStat3(+)</td><td>0.025 mV/s (0.125 mV step) to 5 V/s (5 mV step)</td></tr> <tr> <td>EmStat4 LR</td><td>0.01 mV/s (0.1 mV step) to 500 V/s (200 mV step)</td></tr> <tr> <td>EmStat4 HR</td><td>0.01 mV/s (0.183 mV step) to 500 V/s (200 mV step)</td></tr> <tr> <td>EmStat Pico</td><td>0.01 mV/s to 10 V/s (10 mV step)</td></tr> </table>	PalmSens1 and 2	1 mV/s (1 mV step) to 25 V/s (5 mV step)	PalmSens3	0.02 mV/s (0.15 mV step) to 500 V/s (5 mV step)	PalmSens4	0.02 mV/s (0.075 mV step) to 500 V/s (10 mV step)	EmStat1 and 2	0.02 mV/s (0.1 mV step) to 5 V/s (5 mV step)	EmStat3(+)	0.025 mV/s (0.125 mV step) to 5 V/s (5 mV step)	EmStat4 LR	0.01 mV/s (0.1 mV step) to 500 V/s (200 mV step)	EmStat4 HR	0.01 mV/s (0.183 mV step) to 500 V/s (200 mV step)	EmStat Pico	0.01 mV/s to 10 V/s (10 mV step)
PalmSens1 and 2	1 mV/s (1 mV step) to 25 V/s (5 mV step)																
PalmSens3	0.02 mV/s (0.15 mV step) to 500 V/s (5 mV step)																
PalmSens4	0.02 mV/s (0.075 mV step) to 500 V/s (10 mV step)																
EmStat1 and 2	0.02 mV/s (0.1 mV step) to 5 V/s (5 mV step)																
EmStat3(+)	0.025 mV/s (0.125 mV step) to 5 V/s (5 mV step)																
EmStat4 LR	0.01 mV/s (0.1 mV step) to 500 V/s (200 mV step)																
EmStat4 HR	0.01 mV/s (0.183 mV step) to 500 V/s (200 mV step)																
EmStat Pico	0.01 mV/s to 10 V/s (10 mV step)																

4.20.7 Parameters for EIS stage

Frequency	The applied frequency in Hz.
E dc	The dc potential applied during the EIS scan.
E ac	The amplitude of the E ac signal has a range of 0.0001 V to 0.25 V (RMS). In many applications, a value of 0.010 V (RMS) is used. The actual amplitude must be small enough to prevent a current response with considerable higher harmonics of the applied ac frequency.

t run	The total run time of a scan.
t. Min. sampling	Each measurement point of the impedance spectrum is performed during the period specified by minimum sampling time 't Min sampling'. This means that the number of measured sine waves is equal to t Min sampling * frequency. If this value is less than 1 sine wave, the sampling is extended to 1 / frequency. So, for a measurement at a frequency, at least one complete sine wave is measured. Reasonable values for the sampling are in the range of 0.1 to 1 s.
t. Max equilibration	The impedance measurement requires a stationary state. This means that before the actual measurement starts, the sine wave is applied during 't Max equilibration' only to reach the stationary state. The maximum number of equilibration sine waves is however 5. The minimum number of equilibration sines is set to 1, but for very low frequencies, this time is limited by 't Max equilibration'. The maximum time to wait for a stationary state is determined by the value of this parameter. A reasonable value might be 5 seconds. In this case, this parameter is only relevant when the lowest frequency is less than 1/ 5 s so 0.2 Hz.

4.20.8 Overhead time between steps

The delay between the execution of the two potential steps varies, depending on the instrument model. The following table gives an indication of the average overhead time for each model.

The average overhead time for each instrument is:

PalmSens1 and 2	~80 ms *
PalmSens3	~80 ms *
PalmSens4	~80 ms *
EmStat1, 2 and 3(+)	~80 ms *
EmStat4 LR and HR	~1 ms
EmStat Pico	~1 ms



* Other running processes in Windows and periodically saving to the recovery file (in case this setting is enabled) may add significantly to the overhead time.

4.21 MethodSCRIPT Sandbox

Supported instruments:

- EmStat4 LR and HR
- EmStat Pico
- Sensit series



MethodSCRIPT™ is a scripting language designed for researchers and developers.

4.21.1 Description

The MethodSCRIPT Sandbox allows you to write your own MethodSCRIPT and run them conveniently in PSTrace.

The MethodSCRIPT language allows for programming a human-readable script directly into the potentiostat. The simple script language makes it easy to combine different measurements and other tasks.

See for more information section [MethodSCRIPT](#) on page 259.

5 Limitations and extra options for EmStat Pico

The EmStat Pico module has some limitations for each technique which are explained in this section. This section also applies to instruments based on the EmStat Pico like the Sensit Smart and Sensit BT.



Figure 100 Sensit Smart and Sensit BT based on EmStat Pico module



Figure 101 EmStat Pico Development Board

5.1.1 EmStat Pico limitations for each mode

The EmStat Pico can be used in 3 different modes; Low Speed, High Speed or Max Range:

- Low Speed mode: for scan rates up to 1 V/s or a bandwidth of 100 Hz.
- High Speed mode: for high scan rates and frequencies.
- Max Range mode: a combination of the Low and High Speed modes for optimal dynamic dc-potential range

Each mode comes with distinct limitations regarding potential ranges. The following table illustrates the applicable ranges for each mode, along with the main performance differences between them.

Table 7 Limitations for each EmStat Pico mode

	Low Speed mode	High Speed mode	Max Range mode
Full dc-potential range	-1.2 to +2 V	-1.7 to +2 V	-1.7 to +2 V
Dynamic dc-potential range ¹	2.2 V	1.2 V	2.6 V
Compliance voltage	-2.0 to +2.3 V ²		
Maximum current	±3 mA		
Max. acquisition rate (datapoints/s)	100	1000	100
Supports FRA/EIS	NO	YES	NO
Applied dc-potential resolution	537 µV	395 µV	932 µV
Applied potential accuracy	< 0.2%	< 0.5%	< 0.5%
Available current ranges	100 nA, 2 µA, 4 µA, 8 µA, 16 µA, 32 µA, 63 µA, 125 µA, 250 µA, 500 µA, 1 mA, 5 mA	100 nA, 1 µA, 6 µA, 13 µA, 25 µA, 50 µA, 100 µA, 200 µA, 1 mA, 5 mA	100 nA, 1 µA, 6 µA, 13 µA, 25 µA, 50 µA, 100 µA, 200 µA, 1 mA, 5 mA
Current accuracy	< 0.5% of current ±0.1% of range	< 1% of current ±0.1% of range ³	
Measured current resolution	0.006% of range (5.5 pA on 100 nA range)		

¹ The dynamic range is the range that can be covered during a single scan within the full potential range. For example; a linear scan can start at -1.5 V and end at 1.1 V or vice versa, covering 2.6 V dynamic range.

² The compliance voltage is the maximum potential between Working and Counter electrode and depends on the selected mode.

³ Channel 2 has an uncompensated series resistor (typical 110 Ω) in series with the WE2 signal. This additional resistance must be taken into account.

5.2 Extra options for EmStat Pico and Sensit Series

If an EmStat Pico is connected additional functions become visible when expanding the advanced options for a technique in the Method Editor

☒ **Override potentiostat mode**
 Potentiostat mode: Low Speed
 Potentiostat channel: Channel 0
☒ **Override E range**
 Min E range override: 0.0
 Max E range override: 0.0

Figure 102 Advanced options for EmStat Pico

5.2.1 Override potentiostat mode

By default, the 'Override potentiostat mode' option is disabled, causing the potentiostat mode to be automatically selected based on the predetermined parameters. While the automated selection is generally efficient, there are instances where a manual override may be preferable. An example of this scenario is when dealing with a cell that demands a higher compliance voltage. In such cases, selecting the 'Max Range' option manually proves to be the most suitable, in combination with the 'Override E range'. See also section "Override E range" below.

5.2.2 Potentiostat channel

The EmStat Pico is equipped with two channels each with their own RE, CE and WE. Here you can specify on which channel the measurement should be performed.

Please be aware that in Sensit BT SNS version and Sensit Smart, the Reference Electrode (RE) and Counter Electrode (CE) of the second channel (Channel 1) are in an open circuit configuration. Therefore, if you simply select this channel, you will always obtain an open circuit response.

5.2.3 Override E range

To achieve the desired potential between the Working Electrode (WE) and Reference Electrode (RE), a potentiostat applies a potential between the WE and Counter Electrode (CE). This potential is generally higher than the requested potential between the WE and RE. The maximum potential that the potentiostat can apply on the CE vs. WE is called the compliance potential.

See also section [Sensit Smart and Sensit BT specifications](#) on page 344 for a more detailed overview of the limitations of the EmStat Pico concerning compliance potential.

On the EmStat Pico, the compliance potential range can be manually set by enabling the 'Override E range' option and configuring the values for 'Min E range override' and 'Max E range override'. This specified potential range governs the potentials applied during the measurement, directly influencing the available compliance potential. Typically, PSTrace automatically selects the potential range based on the measurement potentials, which are defined as WE vs. RE.

An instance where this feature proves useful is when dealing with a cell that requires a higher compliance voltage on the positive end. In automatic mode, there might be some compliance range allocated for the lower potential end that goes unused. By setting 'Min E range override' and 'Max E range override' towards more positive potentials, the potentiostat can provide additional potential compliance when required.

A way to verify that the compliance potential is exceeded is to check the box 'Record Cell potential' in the Method Editor.

An alternative way to verify if the compliance potential is exceeded is to measure the CE lead (black) versus the GND lead (green). If the potential reaches 2.8V vs GND, the lower (negative) limit of the potential range should be extended. If the potential reaches 0.2V vs GND, the upper (positive) limit of the potential range should be extended.

6 Electrochemical Impedance Spectroscopy

Electrochemical Impedance Spectroscopy (EIS) is an electrochemical technique to measure the impedance of a system in dependence of the AC potentials frequency.

In classic Electrochemical Impedance Spectroscopy (EIS) measurements, they are typically conducted in "potentiostatic mode." In this mode, a sine wave is superimposed on the specified DC potential, resulting in the equation $E = E_{dc} + E_{ac}$. This configuration produces a DC current with an additional AC component, expressed as $I = I_{dc} + I_{ac}$. From these two waves the potentiostat calculates total impedance (Z), phase angle (Φ), real impedance (Z') and complex / imaginary impedance (Z''). The spectrum is made by measuring these parameters for potential waves with different frequencies.

Indeed, the term "spectroscopy" implies a frequency sweep, which is the most common measurement. However, our instruments provide the flexibility to set the frequency and vary other parameters, such as DC potential and time.

Alternatively, some of our instruments offer the option of EIS in galvanostatic mode (GEIS), where the current is controlled while the potentials are measured.

An introduction can be found in for instance the textbook:

Electrochemical Methods: Fundamentals and Applications, written by Allen J. Bard and Larry R. Faulkner, ISBN-13: 978-0471043720

Supported instruments:

- PalmSens3
- (Multi)PalmSens4
- (Multi)EmStat4 series
- EmStat Pico and Sensit series

The EIS option is an optional feature, except for the Sensit Series. Please verify if it was included in your purchase order.

6.1 Introduction

PSTrace provides different modes of EIS measurements:

- a frequency scan at a fixed dc-potential (default EIS)
- frequency scans at each dc-potential in a potential scan
- frequency scans at specified time intervals (time scan)
- a single frequency applied at each dc potential in a potential scan (Mott-Schottky)
- a repeated single frequency at specified time intervals

For Galvanostatic EIS (GEIS) the modes are:

- a frequency scan at a fixed dc-current
- frequency scans at each current in a current scan
- frequency scans at specified time intervals (time scan)
- a single frequency applied at each current in a current scan
- a single frequency at specified time intervals

The measured data can be presented in different plot formats.

For fitting EIS data on an equivalent circuit, see section: [Equivalent Circuit Fitting](#) on page 135.

The measurements yield the impedance defined by $Z(\omega) = E_{ac}(\omega) / I_{ac}(\omega)$.

The impedance is a complex number $Z(\omega) = Z'(\omega) - jZ''(\omega)$, where Z' is the real part and Z'' the imaginary part of the impedance and $\omega = 2\pi f$, where f is the applied frequency.

The phase shift φ is defined as $\tan(\varphi) = -Z'' / Z'$. (Note: we do not use (ω) from here anymore)

The impedance data are often presented in the Nyquist plot, with Z'' vs Z' or in a Bode plot, with the decimal logarithm of the magnitude of Z and phase shift versus the logarithm of the frequency.

The magnitude Z is defined as

$$|Z|^2 = (Z')^2 + (-Z'')^2.$$




Figure 103 Example of a Nyquist plot




Figure 104 Example of a Bode plot

The inverse of the impedance is the admittance, which is also a complex number $Y = Y' + jY''$, where $Y = 1/Z$, $Y' = Z' / (Z'^2 + Z''^2)$ and $Y'' = Z'' / (Z'^2 + Z''^2)$.

Sometimes impedance measurements are represented in a plot of Y vs $\log(f)$ or Y'' vs Y' .




Figure 105 Example of an admittance plot

The example shown in the plots above represents a measurement of an ideal working electrode at which a faradaic reaction occurs. The equivalent circuit for such a simple system is shown in the figure below.




Figure 106 Equivalent circuit for a simplified Randles circuit

Where R_{Ω} is the ohmic resistance of the electrolytic solution, C_{dl} the capacitance of the double layer of the working electrode and R_{ct} is the resistance due to the electron transfer of the faradic reaction.

At high frequencies, the impedance is dominated by the value of R_{Ω} , so $Z = R_{\Omega}$, with a phase shift of 0° .

At low frequencies the impedance of the double layer is so high that it can be neglected, so the impedance is equal to $R_{\Omega} + R_{ct}$ and the phase shift is again 0° .

At intermediate frequencies the influence of the capacitance of the double layer (which has an impedance of $Z = -j / (\omega C_{dl})$) results in a phase shift and a measured impedance between the two limits of R_{Ω} and $(R_{\Omega} + R_{ct})$.

6.2 Setting up an impedance measurement

The Measurement tab contains all the method parameters. The Notes textbox can be used to describe the sensor used and sample that is measured.

With each change, the validation of the method is checked. Errors or incompatibilities are shown instantly at the bottom of the measurement tab.

6.2.1 Current ranges



Figure 107 Buttons for selecting the applicable current ranges

With the current range buttons at the top, the applicable current range(s) during the measurement can be selected. If more than one button is selected (blue), the instrument will select the most optimal current. The highest current to be enabled is determined by the lowest impedance value and equal to $1/Z$ (lowest value) as well as the highest obtained value for the dc-current. The lowest enabled current range is determined by the highest impedance to be measured and equal to $1/Z$ (highest) and again the lowest value to the dc-current.

So, if the lowest impedance is in the order of 100 Ω or lower, enable 10 mA and if the highest impedance is in the order of M Ω , enable current ranges down to 1 μ A or 100 nA. Nevertheless, it's important to be aware that for high frequencies, the autoranging routines consider the maximum bandwidth when determining how to establish optimal gain settings. This involves a combination of switching to the optimal range and applying a gain of 10 or 100 on the active current range. Consequently, if you configure the potentiostat to operate with a lower Current Range (CR), for instance, by deselecting the higher CR option, you may experience inaccuracies in responses. The chosen CR might lack the necessary bandwidth to adequately sample certain high frequencies.

In general, it is best to enable all current ranges for EIS. In a general way, for EIS experiments the instruments can select the proper current range automatically without any issues. Note that maybe some ranges can be unavailable according to other measurements parameters (mainly 'Max. frequency' and 'E ac').

For a typical Electrochemical Impedance Spectroscopy (EIS) measurement, which involves a frequency scan starting from a high frequency (i.e., >1kHz), data collection is rapid, typically taking only a few seconds for each frequency measurement. In this scenario, the potentiostat efficiently searches for a suitable range, starting from the highest one.



An EIS measurement starts at the highest selected current range.

It's worth noting that if you start with a low frequency (i.e., >1Hz), each measurement will take longer, and the potentiostat may require additional time on the first point to select the proper range. In such cases, situations where manual range selection can be more straightforward may arise. Another example where the manual range may be suitable is a time scan at low frequencies.

6.2.2 Auto current ranging in Galvanostatic EIS

Since the current follower is in the controlled loop, current range switching is not possible while measuring. Auto current ranging is therefore not supported in galvanostatic mode.

6.2.3 Pretreatment settings

Pretreatment is always done in potentiostatic mode, also for Galvanostatic EIS.

Before a measurement starts, the working electrode or sensor can be pretreated according to the specified values for E condition1 and 2 and related times.

Each time a new dc-potential is applied, the equilibration time is used before the impedance measurements start.

6.2.4 Frequency scan settings

A frequency scan starts at the 'Max. frequency' and ends at 'Min. frequency'. The frequency distribution can be logarithmic, linear or custom. If frequency distribution is logarithmic, a fixed number of frequencies per decade are applied. Typically, EIS experiments encompass a broad range of frequencies, spanning from kilohertz to millihertz. Consequently, the use of logarithmic frequency distribution is more suitable, and a common practice is to set 10 frequencies per decade (order of magnitude). For instance, in a scan ranging from 10 kHz to 100 Hz, the total number of frequencies will be 2 decades times 10 + 1, resulting in 21 frequencies. The actual applied number of frequencies can be entered or the number of frequencies per decade.

n frequencies	21	=	10	/ dec.	Edit
Max. frequency	10000.0			Hz	
Min. frequency	100.0			Hz	

Figure 108 Parameters for defining the frequency scan in the Method Editor




Figure 109 'Edit applied frequencies' window

The list of frequencies is shown by clicking the 'Edit' button in the method editor. It opens a dialog which allows the user to choose a mode (Logarithmic, linear, or custom) to generate the frequency list. Any edits made to the frequency list turns it into a custom list which is written back to the method. The modified list can also be saved in the application settings to be used later if needed.

6.3 Parameters for EIS

Scan type

Default: this is the setting for the regular EIS experiment. The dc-potential 'E dc' will remain fixed at a determinate value.

E dc scan: multiple measurements will be conducted with different dc-potential levels, starting from 'E begin' and ending at 'E end' in increments determined by the 'E step'. At each dc-potential level, a single fixed frequency is applied, or a frequency scan is made, according to the set 'Frequency type'.

Time scan: repeated measurements at a defined 'E dc' will be conducted and recorded versus time.

6.3.1.1 Parameters for Scan type = E dc scan

E begin	<p>The dc-potential at which the measurement starts.</p> <p>The applicable potential range for each instrument is:</p> <table> <tr> <td>PalmSens3</td><td>-5 V to +5 V</td></tr> <tr> <td>PalmSens4</td><td>-10 V to +10 V</td></tr> <tr> <td>EmStat4 LR</td><td>-3 V to +3 V</td></tr> <tr> <td>EmStat4 HR</td><td>-6 V to +6 V</td></tr> <tr> <td>EmStat Pico</td><td>-1.25 V to +2 V</td></tr> </table> <p>See also section: Limitations and extra options for EmStat Pico on page 109.</p>	PalmSens3	-5 V to +5 V	PalmSens4	-10 V to +10 V	EmStat4 LR	-3 V to +3 V	EmStat4 HR	-6 V to +6 V	EmStat Pico	-1.25 V to +2 V
PalmSens3	-5 V to +5 V										
PalmSens4	-10 V to +10 V										
EmStat4 LR	-3 V to +3 V										
EmStat4 HR	-6 V to +6 V										
EmStat Pico	-1.25 V to +2 V										
E step	<p>The potential step size. Sets the increment to be used between 'E begin' and 'E end'.</p> <p>The applicable step range for each instrument is:</p> <table> <tr> <td>PalmSens3</td><td>0.15 mV to 250 mV</td></tr> <tr> <td>PalmSens4</td><td>0.075 mV to 250 mV</td></tr> <tr> <td>EmStat4 LR</td><td>0.1 mV to 250 mV</td></tr> <tr> <td>EmStat4 HR</td><td>0.183 mV to 250 mV</td></tr> <tr> <td>EmStat Pico</td><td>0.537 mV to 250 mV</td></tr> </table>	PalmSens3	0.15 mV to 250 mV	PalmSens4	0.075 mV to 250 mV	EmStat4 LR	0.1 mV to 250 mV	EmStat4 HR	0.183 mV to 250 mV	EmStat Pico	0.537 mV to 250 mV
PalmSens3	0.15 mV to 250 mV										
PalmSens4	0.075 mV to 250 mV										
EmStat4 LR	0.1 mV to 250 mV										
EmStat4 HR	0.183 mV to 250 mV										
EmStat Pico	0.537 mV to 250 mV										
E end	<p>The potential at which the scan ends. The applicable range is same as 'E begin'</p>										

6.3.1.2 Parameters for Scan type = Time scan

t run	<p>The total run time of a scan.</p> <p>Run time for each instrument: <1 s to 1,000,000 s</p>
t interval	<p>The minimum interval time between each data point (Frequency type=Fixed) or between each frequency scan (Frequency type=Scan). We recommend a time higher than the required time to measure the data point or perform the frequency scan + overhead time. While it's possible to use a shorter time, doing so may lead to incorrect impedance calculations.</p>

6.3.1.3 Parameters for all scan types

E dc	<p>The dc-potential applied during the EIS scan. Also called "DC Bias" or "level". The most common setting for this parameter is 0 V vs. OCP. By default, the option 'Measure vs. OCP' comes checked.</p> <p>See also section: Measuring versus OCP on page 33.</p>
E ac	<p>The amplitude of the E ac signal has a range of 0.001 V to 0.25 V (RMS). In many applications, a value of 0.010 V (RMS) is used.</p>

Frequency type	<p>Fixed: a single frequency is applied for the given duration or at each potential step or time interval.</p> <p>Scan: a frequency scan is performed starting at the given 'Max. frequency' to the given 'Min. frequency'. This is the setting for the regular EIS experiment.</p>
Frequency	The applied frequency (in case of Frequency type=Fixed).
Show time domain data	Only visible for EmStat4-series instruments. When enabled a Scope window is shown during the measurements to view the time domain data as sinewaves or in a Lissajous plot.

6.3.1.4 Parameters for Frequency type = Scan

n frequencies	The number of frequencies to apply between the given Max. frequency and Min. frequency.
Max. frequency	The frequency to start the frequency scan on.
Min. frequency	The frequency to end the frequency scan on.
Measure vs OCP	If enabled, the Open Circuit Potential will be determined first as a reference point for applied dc potentials.

See also:

- “Frequency scan settings” in the section: [Setting up an impedance measurement](#) on page 117.
- Section: [Measuring versus OCP](#) on page 33.

6.3.1.5 Advanced parameters under [...] button

Pretreat each scan	The potentials set in the “Pretreatment Settings” will be applied before each frequency scan is started
Force max stability filter	Applies to PalmSens4 only. Overrides the frequency-dependent stability filter to its maximum value, can prevent a shift in impedance due to 100uA to 10uA CR switch.
t. Min. sampling	Each measurement point of the impedance spectrum is performed during the period specified by minimum sampling time ‘t Min sampling’. This means that the number of measured sine waves is equal to t Min sampling * frequency. If this value is less than 1 sine wave, the sampling is extended to 1 / frequency. So, for a measurement at a frequency, at least one complete sine wave is measured. Reasonable values for the sampling are in the range of 0.1 to 1 s.
t. Max equilibration	The impedance measurement requires a stationary state. This means that before the actual measurement starts, the sine wave is applied during ‘t Max equilibration’ only to reach the stationary state. The maximum number of equilibration sine waves is however 5. The minimum number of equilibration sines is set to 1, but for very low frequencies, this time is limited by ‘t Max equilibration’. The maximum time to wait for a stationary state is determined by the value of this parameter. A reasonable value might be 5 seconds. In this case, this parameter is only relevant when the lowest frequency is less than 1/ 5 s or 0.2 Hz. Increasing this parameter can be beneficial when measuring unstable or noisy samples at low frequencies, although it will extend the overall measurement time.

6.3.1.6 Parameters applicable for PalmSens3 only

Allow AC coupled measurements > 200 Hz	Enables the use of ac-coupled measurements when reading the current during measurements above 200 Hz. This option increases accuracy, but also increases sensitivity to noise and might introduce a ‘jump’ around 200 Hz.
Disable use of High Stability mode	The High Stability Mode is enabled by default for measurements at frequencies below 400 Hz. This mode filters out high-frequency noise and increases the stability of the measurement. If a small jump in phase shift is observed in a measurement around 400 Hz, it is advised to disable this mode.

Sensitivity mode

Determines the maximum gain used and therefore resolution in the measurements. If measured signals in the scope window seem too noisy, a lower sensitivity might be helpful.

6.4 Parameters for Galvanostatic EIS

In Galvanostatic-EIS (GEIS) the current through the cell is controlled rather than the potential across the cell as with standard EIS. In some cases, it's preferred to perform the EIS measurement in galvanostatic mode e.g. when measuring low impedance cells. To prevent large currents on low impedances low amplitudes can be used in potentiostatic mode, however, the cell impedance may change during the measurement giving unexpected change in currents. In Galvanostatic mode, the current amplitude is controlled and set by the user.

Scan type

Default: this is the setting for the regular GEIS experiment. The dc-current 'i dc' will remain fixed at a determinate value.

i dc scan: multiple measurements will be conducted with different dc-current levels, starting from 'i begin' and ending at 'i end' in increments determined by the 'i step'. At each dc-current level, a single fixed frequency is applied, or a frequency scan is made, according to the set 'Frequency type'.

Time scan: repeated measurements at a defined 'i dc' will be conducted and recorded versus time.

6.4.1.1 Parameters for Scan type = i dc scan

Applied current range	The range in which the specified current values (such as 'i begin', 'i end' and 'i step') will be applied.	
	Applicable applied currents for each instrument:	
	Applied	Available Ranges (CR)
	PalmSens3	-3 to 3 * CR 1 µA to 10 mA
	PalmSens4	-6 to 6 * CR 1 nA to 10 mA
	EmStat4 LR	-3 to 3 * CR 10 nA, 1 µA, 100 µA, 10 mA
	EmStat4 HR	-3 to 3 * CR 1 µA, 100 µA, 10 mA and 100 mA
i begin	Current at which the scan starts expressed in the Applied current range.	

i step	<p>The current step size. Sets the increment to be used between 'i begin' and 'i end'.</p> <p>The applicable minimum 'i step' for each instrument is:</p> <table> <tr> <td>PalmSens1 and 2</td><td>0.001 * CR</td></tr> <tr> <td>PalmSens3</td><td>0.00015 * CR</td></tr> <tr> <td>PalmSens4</td><td>0.000075 * CR</td></tr> <tr> <td>EmStat4 LR</td><td>0.0001 * CR</td></tr> <tr> <td>EmStat4 HR</td><td>0.000183 * CR</td></tr> </table>	PalmSens1 and 2	0.001 * CR	PalmSens3	0.00015 * CR	PalmSens4	0.000075 * CR	EmStat4 LR	0.0001 * CR	EmStat4 HR	0.000183 * CR
PalmSens1 and 2	0.001 * CR										
PalmSens3	0.00015 * CR										
PalmSens4	0.000075 * CR										
EmStat4 LR	0.0001 * CR										
EmStat4 HR	0.000183 * CR										
i end	Current at which the scan ends expressed in the Applied current range.										

6.4.1.2 Parameters for Scan type = Time scan

t run	<p>Total run time of the scan.</p> <p>Run time for each instrument: <1 s to 1,000,000 s</p>
t interval	The minimum interval time between each data point (Frequency type=Fixed) or between each frequency scan (Frequency type=Scan). The interval time cannot be lower than the required time to measure the data point or perform the frequency scan + overhead time.

6.4.1.3 Parameters for all scan types

i dc	The dc current applied during the EIS scan is expressed in the Applied current range.
i ac	The amplitude is in the RMS value as a factor of the used current range. The maximum factor is 0.4 per current range except for the 10mA current the max factor is 0.2 x 10 mA.
Frequency type	<p>Fixed: a single frequency is applied for the given duration or at each potential step or time interval.</p> <p>Scan: a frequency scan is performed starting at the given 'Max. frequency' to the given 'Min. frequency'. This is the setting for the regular GEIS experiment.</p>
Frequency	The applied frequency (in case of Frequency type=Fixed).
Show time domain data	Only visible for EmStat4-series instruments. When enabled a Scope window is shown during the measurements to view the time domain data as sinewaves or in a Lissajous plot.

6.4.1.4 Parameters for Frequency type = Scan

n frequencies	The number of frequencies to apply between the given Max. frequency and Min. frequency.
----------------------	---

Max. frequency	The frequency to start the frequency scan on. In galvanostatic mode, the bandwidth of the current follower affects the stability of the system. For this reason, the maximum frequency is limited to 100 kHz for the Palmsens4.
Min. frequency	The frequency to end the frequency scan on.

See also “Frequency scan settings” in the section: [Setting up an impedance measurement](#) on page 117.

6.4.1.5 Advanced parameters under [...] button

t. Min. sampling	Each measurement point of the impedance spectrum is performed during the period specified by minimum sampling time ‘t Min sampling’. This means that the number of measured sine waves is equal to t Min sampling * frequency. If this value is less than 1 sine wave, the sampling is extended to 1 / frequency. So, for a measurement at a frequency, at least one complete sine wave is measured. Reasonable values for the sampling are in the range of 0.1 to 1 s. Increasing this parameter can be beneficial when measuring unstable or noisy samples at high frequencies.
t. Max equilibration	The impedance measurement requires a stationary state. This means that before the actual measurement starts, the sine wave is applied during ‘t Max equilibration’ only to reach the stationary state. The maximum number of equilibration sine waves is however 5. The minimum number of equilibration sines is set to 1, but for very low frequencies, this time is limited by ‘t Max equilibration’. The maximum time to wait for a stationary state is determined by the value of this parameter. A reasonable value might be 5 seconds. In this case, this parameter is only relevant when the lowest frequency is less than 1/ 5 s or 0.2 Hz. Increasing this parameter can be beneficial when measuring unstable or noisy samples at low frequencies, although it will extend the overall measurement time.

6.5 Parameters for Fast EIS and Fast GEIS



This technique is only supported by instruments from the EmStat4-series.

The Fast EIS and Fast EIS techniques provide a form of EIS where the latency between each measured datapoint is reduced to a minimum of 1 ms, starting at a frequency of 10 kHz. This allows for observing rapid impedance changes on a cell at a fixed frequency.

See the previous section for more detailed information about the allowed ranges for each technique parameter.

6.6 Running an EIS measurement

The sequence of the measurement depends of course on the specified parameters.

1. In case Measure vs OCP is enabled: determine the OCP.
2. If 't condition 1' > 1 s, apply E condition1 and wait the specified time
3. If 't condition 2' > 1 s, apply E condition2 and wait the specified time
4. Apply the dc-potential at which the impedance has to be measured and wait 't equilibration' seconds.
5. Apply sine wave and wait for the stationary state, using the parameters as described above.
6. Sample the impedance value by measuring at least one complete sine wave, but during at least the specified 't Min. sampling'.
7. Perform the necessary calculations and present them in the 'EIS Plot' window.

The next measurement is now done by stepping back to step 1, step 3 or step 4:

8. Goto step1: if the checkbox 'Pretreat each frequency' or 'Pretreat each scan'
9. Goto step 3: if the next dc-potential has to be applied
10. Goto step 4: if another frequency has to be applied and the checkbox 'Pretreat each frequency' has NOT been checked.

In case a 'Time scan' is made, the program waits until the 't interval' has been completed.

6.6.1 The Scope window

During the measurement the measured time domain response is presented in the Scope window. For PalmSens4, this data is always shown. For other instruments, the recording of time domain data needs to be enabled in the Method Editor. The EmStat Pico and derived Sensit-series instruments are not capable of recording time domain data.



The EmStat Pico does not support the Scope window shown during measurements.




Figure 110 The Scope window

The Scope window is used to check that the response is not distorted by higher harmonics or noise (especially 50/60 Hz from mains) and whether the current is not measured in an overload situation.

By utilizing the selection buttons, you have the ability to present time domain data in two distinct forms: sinusoidal waveforms or Lissajous plots. In the latter representation, the plotted relationship between measured current and potential provides valuable insights. Importantly, data is stored for each frequency, allowing separate viewing after measurements are done.



Closing the Scope window clears the time domain data from memory, which is not saved with the measurement file.

6.6.2 Measurement results

The presented numerical results are:

- dc current
- ac current, both in units of the applied current range
- Impedance amplitude Z and phase shift in degrees
- Real or in-phase impedance Z' and imaginary or out-of-phase impedance Z''

The measured data can graphically be presented in different ways.




Figure 111 Impedance plot options

For a frequency scan:

- $-Z''$ vs Z' , also called the “Nyquist plot” or “Complex plot”
- $\log |Z|$ vs $\log f$
- $\log |Z|$ and $-\text{Phase (shift)}$ vs $\log f$, also called the “Bode plot”
- Y'' vs Y' , also called the “Admittance Nyquist plot”
- $\log Y$ vs $\log f$, also called the “Admittance Bode plot”
- Y' and Y'' vs $\log f$
- C' vs C'' , also called “Complex Capacitance plot”

For a time or potential scan:

- Z' and Z'' vs time (t) or potential (E)
- Z and $-\text{Phase (shift)}$ vs t or E
- Y' and Y'' vs t or E
- Y vs t or E
- C vs t or E
- C' and C'' vs t or E

Please note that if a plot of Y vs potential is made, the plot in fact is an ac-voltammogram.

The numerical values of impedance, phase shift, and admittance are shown in the ‘Data’ tab.

The series capacitance C_s is $1 / (j * f * 2 * \pi * Z'')$ where the imaginary component j can be ignored as there is no real part to contribute to the magnitude. The real and imaginary capacitance are calculated as follows $C' = Z'' / (f * 2 * \pi * |Z|^2)$ and $C'' = Z' / (f * 2 * \pi * |Z|^2)$.

RRCRC_0V

Copy to clipboard

Index	Frequency/Hz	Z/Ohm	ZRe/Ohm	ZIm/Ohm	-Phase/Degr	Cs/F	CR
1	10000	669.882	661.115	-108.025	9.28	1.473E-07	1 mA
2	7943.28	769.209	707.219	-302.53	23.16	6.623E-08	1 mA
3	6309.57	909.844	750.994	-513.639	34.37	4.911E-08	1 mA
4	5011.87	1101	800.139	-756.39	43.39	4.198E-08	1 mA
5	3981.07	1354	873.853	-1034	49.8	3.866E-08	1 mA
6	3162.28	1680	971.615	-1371	54.67	3.672E-08	1 mA
7	2511.89	2080	1104	-1763	57.93	3.595E-08	1 mA
8	1995.26	2585	1319	-2223	59.33	3.588E-08	1 mA
9	1584.89	3206	1641	-2754	59.22	3.646E-08	1 mA
10	1258.93	3962	2122	-3346	57.61	3.779E-08	1 mA

Figure 112 Data tab showing the raw measured data

6.7 Optimizing EIS Parameters for Improved Results and Troubleshooting

EIS measurements can be very sensitive to electromagnetic noise. To address and mitigate potential noise issues, please refer to section “[Reducing noise](#)” on page 51” for detailed guidance. Here are the primary parameters you can adjust to troubleshoot the most common issues encountered when conducting impedance measurements.

E dc: as the EIS data may be challenging to interpret, setting the DC potential (E_{dc}) to match a stable OCP is a common practice. This approach helps initiate the experiment as close to a ‘Steady State’ as possible, minimizing variations in other system properties throughout the

frequency sweep. Therefore, it is advisable to ensure a stable OCP curve before starting an EIS measurement, unless your experiment necessitates a different configuration.

E_{ac}, also called “AC amplitude”, plays a critical role in EIS experiments. To ensure accurate and reliable data, it's important to strike the right balance. Keeping the amplitude small is essential to prevent the appearance of significant higher harmonics of the applied AC frequency, which can complicate data interpretation and introduce redox reactions that are typically beyond the scope of common EIS experiments. To assess the suitability of your amplitude, if your instrument offers a 'Scope window', you can examine the shape of the current wave. If it deviates substantially from a sine wave, it may indicate that the system is operating outside its linearity region, suggesting the need for a lower amplitude.

Conversely, there are instances where increasing the amplitude becomes necessary. If your system exhibits high impedance results, often in excess of megaohms, and you still encounter noise or inconsistent curves, adjusting the amplitude upwards may be beneficial. In such cases, even with an increased potential amplitude, the net current remains relatively low, allowing you to maintain linearity while improving signal quality.

Max. frequency: when conducting EIS measurements, it's essential to be mindful of the maximum frequency setting. For the most common 3-electrode electrochemical cells, like those with metal electrodes in aqueous solutions, frequencies exceeding 10 kHz are often unnecessary. The reference electrode's response time is a limiting factor, and any variations observed at higher frequencies may be attributable to artifact signals. While high frequencies offer faster data collection, they may require time-consuming data filtering or exclusion. To avoid this, consider initiating your measurements with a lower frequency setting.

However, if your application demands higher frequencies, such as when working with solid electrolytes, be aware of potential inductance issues arising from cables and connections. These can affect data accuracy, so proper consideration of cable quality and connections becomes essential for precise measurements.

t. Min sampling: when confronted with inconsistent data in the high-frequency region, consider increasing the minimum sampling time. This adjustment can help stabilize data collection and improve the reliability of measurements, especially when dealing with rapid changes or fluctuations in the system's response at higher frequencies.

6.8 Fitting measured EIS data in PSTrace

PSTrace comes with an advanced Equivalent Circuit fitting tool.

See also section: [Equivalent Circuit Fitting](#) on page 135.

6.9 Export for analysis and circuit fitting

PSTrace comes with an advanced Equivalent Circuit fitting tool. However, PSTrace also allows exporting EIS data to Scribner's ZView and the free program EIS Spectrum Analyser. Measured data in PSTrace can be opened or added to ZView with a one-click action.

6.9.1 Free EIS Spectrum Analyser

Data can be saved in a file that is loaded in the third-party program (free of charge) EIS Spectrum Analyzer, written by Alexander S. Bondarenko and Genady A. Ragoisha and is available from:

<http://www.palmsens.com/eisspectrumanalyser>.

This program allows users to fit the measured data to specific equivalent circuits and for instance obtain the best values of the elements in the equivalent circuit, for instance, R_{Ω} , R_{ct} and C_{dl} .

The file can be saved using the button in the toolbar left from the plot:



Or from the EIS Plot menu:




Figure 113 EIS plot menu

For detailed instructions consider the manual of EIS Spectrum Analyser.

6.9.2 ZView

If ZView 3.3 or later is installed, PSTrace will detect this automatically. A one-click-export button will appear next to the plot. The measurement selected in the Legend will be exported to ZView.



Figure 114 ZView one-click-export button




Figure 115 Equivalent circuit model in ZView

Multiple measurements can be run from PSTrace and fitted on the active circuit model in ZView with a single click from PSTrace.

For detailed instructions consider the manual of ZView.

6.9.3 Origin

Origin from OriginLab is scientific graphing and data analysis software widely used at universities.

See for more information: <http://www.originlab.com/>

If Origin is installed on the PC where PSTrace is running on an extra button appears in the toolbar next to the plot:



Figure 116 Origin one-click-export button

When clicking this button, Origin will be opened showing the selected EIS data directly in a Work sheet.

If the button is not showing and Origin is installed, open the General settings window in the menu 'Tools' → 'General settings...' and check in the 'Plot and data' tab if the location for Origin is properly set.




Figure 117 Location for Origin specified in the Settings window

To change the exported columns for EIS data, click the button 'Set columns for exporting EIS data' on the left-hand side of the Settings window.

6.10 Limitations for EIS on EmStat Pico



This section also applies to instruments based on the EmStat Pico potentiostat module, including the Sensit BT and Sensit Smart.

Electrochemical Impedance Spectroscopy (EIS) measures the change of the complex impedance over a frequency range (spectrum). The common way to calculate the complex impedance is by dividing the complex voltage by the complex current. The ADUCM355 has only 1 ADC so it is not possible to measure both signals (voltage and current) at the same time. Thus the proposed solution to enable EIS measurement using only the current signals is a 3-stage measurement. The total EIS measurement for one single frequency point is split up into three measurements on three different impedances:

1. Zcell + Rload
2. Rload (load resistor in the WE signal path)
3. Rcall (calibration reference resistor)

Since all three measurements are performed under the same conditions and Rcall is a reference resistor of known value, the final complex impedance of Zcell can be calculated having only the complex currents of the three measurements.

6.10.1 Complex voltage correction

For the plain ADuCM355 the conditions for the 3 measurements are equal, for the Emstat Pico the measurement on the Zcell+Rload is performed using external RE-buffers introducing a complex transfer function for the applied AC voltage on the Zcell+Rload. This complex voltage transfer function is modeled by an electronic circuit simulation of the gain and phase behavior of the transfer function.

6.10.2 Time-domain sinewave

The Emstat Pico integrates a DFT calculation block enabling onboard complex current calculations. In contrast to other PalmSens devices having the EIS feature, the raw ADC signal is not available and therefore the time-domain signal cannot be shown in a host application (PSTrace, PStouch)

6.10.3 Measurement duration

The accuracy of the complex current depends on the number of the applied frequency cycles presented to the DFT process. For higher frequencies, the time to measure multiple cycles is relatively short in contrast to the lower frequencies. For example, measuring 8 cycles of 1 Hz takes 8 seconds resulting in a 24 seconds measurement duration for a complete 3-stage measurement.

6.10.4 Current ranges

Since the conditions must be the same, all 3 measurements must be performed using the same current range. Rload (~100 Ω) and the Rcall (1 k Ω or 100 k Ω) are fixed values, however the complex current measurements are frequency depended due to parasitic effects of the signals

path. The changes of Z_{cell} vs frequency can be so large that it cannot be covered by staying in the same current range. Auto current ranging dynamically changes the current range (in combination with the PGA) to cover the frequency range in the EIS measurement.




Figure 118 DFT signal path block diagram for EmStat Pico

More details can be found in the “ADuCM355 Hardware Reference Manual- UG-1262”, pages 137-143.

7 Equivalent Circuit Fitting

The Circuit Editor can be used to build, simulate, and fit equivalent circuits on your Electrochemical Impedance Spectroscopy data. The following chapters present an overview of the circuit editor's options and features.

For more detailed information on fitting equivalent circuits electrochemical impedance spectroscopy measurement the following literature is recommended:

- Electrochemical Impedance Spectroscopy. Chapter 23: An Integrated Approach to Impedance Spectroscopy. Mark E. Orazem & Bernard Tribollet, ISBN: 978-0-470-04140-6.
- Electrochemical Impedance Spectroscopy and its Applications. Chapter 14: Modeling of Experimental Data. Andrzej Lasia, ISBN: 978-1-4614-8932-0.

7.1 Overview

The circuit editor has two different modes:

- Edit mode: draw the circuit or type CDC circuit. You typically start with this mode, but afterwards you can return to it whenever you need to redraw the circuit.
- Fit mode: fit the EIS data on the drawn circuit or simulate a response of the circuit.

7.1.1 Circuit Edit Mode




Figure 119 Circuit Editor main window in Edit Mode.

7.1.2 Circuit Fit Mode




Figure 120 Circuit Editor main window in Fit Mode.

7.2 Opening the circuit editor

The circuit editor can be opened from the EIS Plot tab in PSTrace. To fit a circuit the measurement to fit on must be selected from the Session Manager first. When no measurement is selected it is only possible to simulate circuits.




Figure 121 Button for opening Circuit Fitting window

7.3 Switching between Edit mode and Fit mode

The circuit editor has two modes.

1. In Edit Mode circuits can be built, edited, loaded and saved.
2. In Fit Mode the circuit built in the Edit Mode can be simulated or fit on the electrochemical impedance spectroscopy data. (Note: Fitting is only possible when a measurement was selected in the Session Manager while opening the circuit editor.)

To switch from the edit mode to the switch mode either press the Fit Mode button or select Fit Mode from the view menu.




Figure 122 Switching between the two modes

To switch from the fit mode to the edit mode either press the Edit Mode button or select Edit Mode from the view menu.




Figure 123 Switching between the two modes

7.4 Building a circuit in the graphical editor

To build a circuit the circuit editor must be in edit mode. Circuits can be built graphically using the circuit components and the graphical circuit editor.

7.4.1 Graphical circuit editor

Circuits can be built graphically using the circuit components and the graphical circuit editor.

7.4.2 Selecting a circuit component

By default, the resistor is selected. To select another circuit component click on it. The eraser is used to remove components from the circuit.




Figure 124 Circuit component selection buttons

7.4.3 Adding a component to the circuit

Components can be added in series or in parallel with other components. The first component is added by hovering the mouse cursor over the open connector in the graphical circuit editor and clicking the left mouse button.




Figure 125 The open connector

When hovering the mouse cursor over a connector it will turn green and a grey preview of the component is displayed.




Figure 126 Preview of the component

By clicking the left mouse button the component is added to the circuit.




Figure 127 The inserted component.

To add components in serial the mouse cursor must hover over a connector. Clicking on the left mouse confirms the placement of the component.




Figure 128 Adding components in series

Components can be added in parallel to other components by hovering the mouse cursor over the component it must be placed parallel to. When hovering the mouse cursor over another component it will turn green to indicate that it is selected and a preview parallel component is shown in grey. Click the left mouse button to confirm the placement of the component in the circuit.




Figure 129 Preview of a parallel component




Figure 130 Component added in parallel

When creating larger circuits it is useful to know that components can be added in parallel within or over a parallel group in the circuit.

To place a component in parallel to a component within a parallel group directly hover the mouse cursor over the other component and click on the left mouse button.




Figure 131 Placing a component in parallel to a component within a parallel group

To place a component in parallel to the entire group hover the mouse cursor over the group. This will select the entire group (i.e. all its components will turn green). Clicking on the left mouse button will confirm the placement of the component.




Figure 132 Placing a component in parallel to the entire parallel group

7.4.4 Removing components from the circuit

Components can be removed from the circuit using the eraser tool. To remove a component, select the eraser tool.




Figure 133 Eraser tool

To remove a component or a group of parallel components hover the mouse cursor over it. The component(s) will turn green indicating that they are selected. Clicking the left mouse button will remove the selected components from the circuit.




Figure 134 Removing a single component




Figure 135 Removing a group of components

7.5 Building a circuit with the CDC editor

Every equivalent circuit for EIS data fitting can be represented by a CDC or Circuit Description Code. The CDC circuit editor is located at the bottom of the circuit editor. To build a circuit using the CDC circuit editor click on it and type the circuits CDC code. If the CDC code is valid the circuit is updated directly in the graphical circuit editor. The CDC code is described by Boukamp in:

Boukamp BA (1986) A package for impedance/admittance data analysis. *Solid State Ionics* 18 & 19, 136-40.












Figure 136 Field for entering CDC directly

7.5.1 Supported CDC characters

The following table shows which characters are supported by the CDC circuit editor.

Table 8 Supported character for writing the CDC







Key	Function
	Inserts a resistor.
	Inserts a capacitor.
	Inserts an inductor.
	Inserts a Warburg impedance.
	Inserts a constant phase element.
	Inserts a parenthesis defining the beginning of a parallel group.
	Inserts a parenthesis defining the end of a parallel group.
	Inserts a bracket defining the beginning of a series group.
	Inserts a bracket defining the end of a series group.





7.5.2 Supported CDC controls

The following table shows which keyboard control keys are supported by the CDC circuit editor.

The mouse cursor places the text cursor in the CDC editor by clicking the left mouse button. It can also be used to select text in the CDC editor by holding down the left mouse button and dragging the cursor over the CDC characters it should select.

Table 9 Supported character for writing the CDC

Key	Function
 	The arrow keys move the text cursor left and right respectively. In combination with Shift, they select CDC code.
 	Home and end move the cursor to the beginning or end of the CDC code respectively. In combination with Shift, they select CDC code.
	Backspace removes the CDC character to the left.
	Delete removes the CDC character to the right. Backspace cannot be used to remove a selection of CDC characters.

	Control + A selects all of the CDC code.
	Control + C copies the selected CDC code to the clipboard.
	Control + X cuts the selected CDC code and sends it to the clipboard.
	Control + V pastes the contents of the clipboard to the location of the text cursor.

7.5.3 Adding components

Components are added to the circuit by typing the keys of the supported characters.




Figure 137 Pressing R will add a resistor to your circuit

To add components in parallel parenthesis must be used. After opening the parenthesis the CDC editor will turn red indicating the CDC code is invalid at the moment. Ignore this and add the components that should be placed parallel to each other.




Figure 138 Entering CDC code for components parallel to each other

As soon as the parenthesis are closed the CDC code will become valid again and the graphical circuit editor will display the circuit as well.




Figure 139 Graphical representation of the (RC) circuit in the graphical circuit editor




Figure 140 Example of a CDC code with three components parallel to each other

Within a group of components parallel to each other brackets must be used to add a component in series to one of the components. After opening the brackets the CDC editor will turn red indicating the CDC code is invalid at the moment. Ignore this and add the components that should be placed series to each other.



Figure 141 Entering CDC code for components in series within a parallel group




Figure 142 Graphical representation of the ([RW]C) circuit in the graphical circuit editor

In some cases it is easier to copy and paste parts of the CDC code to quickly build a large circuit. This can be done using the supported CDC controls.



Figure 143 Copying and pasting in the CDC editor

7.5.4 Removing components

Components can be removed from the circuit with the CDC circuit editor by placing the text cursor next to the component that should be removed and using the supported CDC circuit editor controls. Pressing backspace or delete removes the CDC character to the left or right of the text cursor respectively. The changes to the circuit are directly visible in the graphical circuit editor.




Figure 144 Backspace removes the constant phase element to the left of the text cursor




Figure 145 Delete removes the capacitor to the right of the text cursor

Removing a parenthesis or bracket from a valid CDC code will invalidate it. The CDC circuit editor will indicate this by turning red. This happens typically when you want to remove a parallel or series group from a circuit.




Figure 146 Removing a parallel group from a circuit using backspace

In this example the parentheses around the second resistor are removed as they no longer have any meaning.

A selection of CDC characters can be removed by cutting or deleting it with control + x, backspace or delete. The selection can also be replaced by pasting data from the clipboard, control + v, or replacing it with any of the supported characters.



7.5.5 Invalid CDC code

In the case of an invalid CDC code the CDC circuit editor indicates this by turning red. The graphical circuit editor is not updated until the CDC code is valid again, this is indicated by the CDC circuit editor turning white again.

When adding or removing a parallel or series group from a circuit using the CDC circuit editor ignore that the code is invalid, as it will be valid again once the entire parallel or series group has been added or removed.



When using the graphical circuit editor in combination with the CDC circuit editor hovering the mouse over the graphical representation will remove any invalid CDC code in the CDC circuit editor.

R(C[R)W]

Figure 147 Invalid placement of closing parenthesis

Here the invalid CDC code can be resolved by removing the closing parenthesis between the resistor and Warburg impedance and placing it at the end of the CDC code.

R(C[RW])RC)

Figure 148 Unbalanced amount of opening and closing parentheses or brackets

Here the invalid CDC code can be resolved by adding an opening parenthesis or removing a closing parenthesis, as the amount of opening and closing parenthesis is unbalanced in this example.

R(RCfewfe)

Figure 149 Unsupported CDC characters have been pasted into the CDC circuit editor

This code is invalid because unsupported CDC characters were pasted into the CDC circuit editor from the clipboard. To resolve this, the unsupported CDC characters must be removed.

7.6 Loading and saving circuits

Circuits can be stored in and retrieved from the circuit database which is stored locally in a folder on the computer. The load and save options are in the Circuit menu.




Figure 150 “Circuit” menu

To load or save a circuit the circuit editor must be in edit mode.

7.6.1 Loading a circuit

To load a circuit, open the load circuit screen by selecting load in the Circuit menu. In this screen the available circuits stored on your computer are listed in the circuit browser on the left. A preview of the selected circuit is presented in the preview panel on the right.




Figure 151 “Load circuit” window

The selected circuit can be loaded by pressing on the Load button or double-clicking on the circuit.

7.6.2 Saving a circuit

Circuits can be saved using the save circuit screen. To open the save circuit screen select save in the Circuit menu. The option to save a circuit is only available if there is a circuit containing at least one component in the editor. The save circuit screen contains a circuit browser and a field to input the name of circuit that is being saved.




Figure 152 "Save circuit" window

To save the circuit press the save button. Optionally a name can be assigned to the circuit in the circuit name input field.

7.7 Fitting or simulating a circuit

After an equivalent circuit has been built or loaded into the circuit editor in "Edit Mode" it can be simulated or fit onto the data of the electrochemical impedance spectroscopy measurement.

7.7.1 Simulating a circuit

The impedance of circuits can be simulated using the circuit editor. When simulating the fit options are disabled. Check the title of the circuit editor or the view menu to see whether the circuit editor is set to simulation mode.




Figure 153 Circuit editor set to simulation mode and fit mode. In fit mode the name of the measurement is visible in the title bar and the simulation is not checked in the view menu.

Values of the components can be changed in the circuit table and the graphical circuit editor. Changing the value of a component in either the graphical circuit editor or the table will automatically update the value in the other. The effects of these changes will appear directly in the preview plot.




Figure 154 Changing the value of the Warburg Coefficient for a Warburg Impedance in the graphical circuit editor.

Element	Fitted Value	Unit	Error%
R 1	1000	Ω	✓
R 2	1000	Ω	✓
W 1	1000	σ	✓
C 1	1.000	μF	✓

Figure 155 Changing the value of the Warburg Coefficient for a Warburg Impedance in the circuit table.

For example, changing the value of the Warburg Coefficient to 100 will lower the impedance at the lower frequencies.




Figure 156 Simulation of the absolute impedance over frequency of a Randle's circuit using the default values.




Figure 157 Simulation of the absolute impedance over frequency of a Randle's circuit with the value of the Warburg Coefficient lowered to 100.

Using the plot selector, it is possible to view the effects of lowering the Warburg Coefficient to 100 on the negative phase angle over frequency and the Nyquist plot as well.




Figure 158 Plot selector

7.7.2 Setting the simulation's frequency range

The frequency range of the simulation can be set in the "Advanced Options" window which can be accessed by clicking on the 'Advanced Options' button.




Figure 159 The 'Advanced Options' button

The simulations frequency range (Hz) can be specified in the simulation settings, by default this is 0.01 to 100000 Hz.




Figure 160 Advanced Simulation Settings

7.7.3 Fitting a circuit

Before fitting the circuit, it is important to check whether the circuit editor is in simulation or fit mode. Check the title of the circuit editor or if the simulation option is checked in the view menu.




Figure 161 Circuit editor set to simulation mode and fit mode. In fit mode the name of the measurement is visible in the title bar and the simulation is not checked in the view menu.

To ensure a good fit on a measurement it is important to build the right equivalent circuit. It is often possible to achieve a good fit with different but similar circuits. Although these other circuits can provide a good fit, they are not necessarily an accurate representation of the cell that was measured. Selecting appropriate initial values for the components in the models is also recommended, as in some cases the fitting algorithm will get stuck in a local minimum and the quality of the fit is not optimal.

As the complexity of the proposed circuit increases, achieving optimal results with the algorithm becomes more challenging without reasonable initial values. In certain instances, specific elements or sets of elements may tend towards values close to zero or infinite after a run (i.e.: resistor with less than a few microohms). This indicates that the algorithm is almost ignoring this part of the circuit, and is stuck called "local minima". To address this, it becomes

necessary to supply more realistic initial values. Additionally, you can also manually edit the minimum and maximum values in the table to further refine the simulation parameters.

For more details, see the example of Maxwell fit: [Fitting Example](#) on page 165.

7.7.4 Fitting on a specific frequency range

The frequency over which you desire to fit the circuit can be specified either by clicking in the plot or by entering the values in the corresponding textboxes.

Note: Specifying the frequency range over which to fit the data is only possible when in fit mode. To specify the frequency for a simulation please refer to the help section on simulating a circuit. Specifying the fitting frequency range is only possible in the plots that also show the measured data (i.e., the blue dots).

Specifying fitting frequency with the mouse is done by moving your mouse over the sample from where you would like to start or end your fit, then double-clicking on that sample (blue dot). Then move your mouse to the data point where you would like to respectively end/start your fit and click on it once. The frequency range selected for fitting is indicated by the solid blue dots.




Figure 162 Specifying the frequency range for fitting. Left: First move your mouse to the data point from which you would like to start the fit and click on it. Middle: Click on it again, now it will be the only data point selected (in this case the sample at 11.51 Hz). Right: Then move your mouse to the data point where you would like fitting to end and click on it. The selected frequency range is now indicated by the solid blue dots and the values in the text boxes below the plot.

7.7.5 Excluding outliers

To exclude outliers from a fit, hold down the Ctrl-button on your keyboard and drag a box around the outlier with your mouse to deselect it. Selecting that point again by dragging a box while holding down Ctrl, will reselect it. While dragging an area using the Ctrl- and the left mouse button, it is possible to (de)select all points in the specified area.




Figure 163 Deselecting specific points by holding Ctrl-button and dragging the left mouse button.

7.7.6 Fixate a component's value

When necessary, it is possible to fixate the value of any given component in the model by checking its checkbox in the Fixed column.

Results						
Fixed	Element	Fitted Value	Min Value	Max Value	Unit	Error%
<input checked="" type="checkbox"/>	R 1	540.0	1.00E-6	1.00E+12	Ω	▼
<input type="checkbox"/>	R 2	9982	1.00E-6	1.00E+12	Ω	▼
<input type="checkbox"/>	C 1	0.033	1.00E-6	1.00E+3	μF	▼
Chi-Squared:		0.0001	Iterations:		6	


Figure 164 Fixating a value.

7.7.7 Undo/Redo changes to component values

The undo/redo buttons allow you to revert unwanted changes to the model's component values either due to an incorrect fit or a mistake when entering the values.

Results

Fixed	Element	Fitted Value	Min Value	Max Value	Unit	Error%
<input type="checkbox"/>	R 1	707.4	1.00E-6	1.00E+12	Ω	22.25
<input type="checkbox"/>	R 2	1.000E+12	1.00E-6	1.00E+12	Ω	3.004E+9
<input type="checkbox"/>	C 1	1000	1.00E-6	1.00E+3	μF	37.75
	Chi-Squared:	0.8852	Iterations:	51		

Algorithm: Levenberg-Marquardt

Results

Fixed	Element	Fitted Value	Min Value	Max Value	Unit	Error%
<input type="checkbox"/>	R 1	556.3	1.00E-6	1.00E+12	Ω	28.30
<input type="checkbox"/>	R 2	9970	1.00E-6	1.00E+12	Ω	3.013E+17
<input type="checkbox"/>	C 1	0.033	1.00E-6	1.00E+3	μF	1.142E+6
	Chi-Squared:	0.8852	Iterations:	51		


Algorithm: Levenberg-Marquardt

Figure 165 In the case of a poor fit the values of a previous fit can be recovered by clicking on the undo button.

In the next sections, brief instructions are given on selecting the components for the circuit. However, further reading on the topic is recommended:

- <http://www.consultrsr.net/resources/eis/index.htm>
- Electrochemical Impedance Spectroscopy. Chapter 23: An Integrated Approach to Impedance Spectroscopy. Mark E. Orazem & Bernard Tribollet, ISBN: 978-0-470-04140-6.
- Electrochemical Impedance Spectroscopy and its Applications. Chapter 14: Modelling of Experimental Data. Andrzej Lasia, ISBN: 978-1-4614-8932-0.

7.8 Batch-fitting a circuit on multiple curves

When multiple EIS curves are available, a single equivalent circuit can be fitted on multiple curves.




Figure 166 Panel for selecting multiple EIS curves.

When this panel is visible, the 'Fit' button applies to the checked curves only. The table showing the fit results and the results shown in the equivalent circuit apply to the curve selected in the same list.

After fitting, the checkboxes can be used to show or hide the results in the Plot with the EIS data curves and fitted curves.

The initial component values used for fitting can be copied to all curves in the list, using the button 'Copy values to selected items'.

7.9 Overview of circuit components

The following components can be used to build equivalent circuits in the circuit editor. The two chapters after this chapter provide two examples of selecting the components for an equivalent circuit, setting/adjusting their values to obtain a good fit, and interpreting the quality of the fit.

7.9.1 Resistor

$$Z_R = R$$

The impedance of a resistor is independent of frequency, and it only contributes to the real component of impedance. Hence it does not affect phase shift and it is represented by a single dot in the Nyquist plot. By default, resistors are 1000 Ω in the circuit editor.

When modeling an electrochemical impedance measurement on a cell a resistor can be used to model the solution resistance (resistance between the working electrode's surface and the tip of the reference electrode).




Figure 167 Left: The absolute impedance of a resistor plotted against frequency. Top right: Phase shift of a resistor over frequency. Bottom right: In the Nyquist plot a resistor is seen as a single dot (little red dot in the center of the plot).

7.9.2 Capacitor

$$Z_c = 1/j\omega C$$

A capacitor's effect on impedance pertains to its imaginary component and decreases with increasing frequency. On its own a capacitor causes a phase shift -90° independent of frequency (by default the negative phase is represented upwards). As it only affects the imaginary component of impedance, it is represented by a vertical line in the Nyquist plot. The default value of a capacitor is $1 \mu\text{F}$ in the circuit editor.




Figure 168 Left: The absolute impedance of a capacitor plotted against frequency. Top right: Phase shift of a capacitor over frequency. Bottom right: Nyquist plot of capacitor.

When a capacitor is placed in series with a resistor the circuit's phase shift becomes frequency dependent. At low frequencies, it approaches 90° and at high frequencies it approaches 0° . Thus, at high frequencies the capacitor's effect on the circuit's impedance and phase shift

becomes negligible and the circuit effectively behaves as a single resistor. The Nyquist plot below shows that the real component of the impedance remains 1 k Ω regardless of frequency.

A resistor and capacitor in series model the impedance of an ideally polarized liquid electrode / blocking electrode, i.e., an electrode that does not transfer charge with the surrounding solution. In this model, the resistor represents the solution resistance and the capacitor represents the product of the real surface of the electrode and its double layer capacitance. With increasing frequency there is less time for the electrochemical double layer to charge and its effects on impedance and phase shift diminish.

Another application of a resistor and capacitor in series is modeling an ideal coating. Again, the resistor represents the solution resistance. But in this model, the capacitor represents the coating capacity. A coating's capacity depends on the coating's thickness, surface area, and dielectric constant.

For more information on modeling corrosion, we recommend reading:

Electrochemical Impedance Spectroscopy and its Applications. Chapter 11: Coatings and Paints. Andrzej Lasia, ISBN: 978-1-4614-8932-0.




Figure 169 Left: The phase shift of a resistor and capacitor in series plotted against frequency. Top right: Absolute impedance over frequency. Bottom right: Nyquist plot.

7.9.3 Inductor

$$Z_L = j\omega L$$

An inductor is the reciprocal of a capacitor its effect also pertains to the imaginary component of impedance. However, in contrast to a capacitor, the impedance of an inductor increases with frequency and its phase shift is +90° frequency (by default the positive phase is represented downwards). It is also represented by a vertical line in the Nyquist plot, but in the opposite direction. In the circuit editor, the default value of an inductor is 100 μ H.

An example of a model where an inductor is used is a model of a faradaic reaction involving adsorption of one or more species. In certain cases, on the Nyquist plot inductive loops can be observed together with capacitive loops and it is possible to model inductance as the product

of the resistance of the charge transfer squared and the absorption pseudocapacitance. For further reading see:

Electrochemical Impedance Spectroscopy and its Applications. Chapter 5: Impedance of the faradaic reactions in the presence of adsorption. Andrzej Lasia, ISBN: 978-1-4614-8932-0.




Figure 170 Left: The absolute impedance of an inductor plotted against frequency. Top right: Phase shift of an inductor over frequency. Bottom right: Nyquist plot of Inductor.




Figure 171 Nyquist plot of a circuit that can model faradaic reactions involving adsorption and subsequent desorption of a species.

7.9.4 Warburg element

$$Z_W = \frac{\sigma}{\sqrt{\omega}} - j \frac{\sigma}{\sqrt{\omega}}$$

A Warburg element is a component used to model the transfer of charge between the electrode and a redox species in the solution and the depletion of the diffusion layer's inner layer. This Warburg element assumes there is a semi-infinite linear diffusion layer. Both the real and imaginary components of impedance increase equally with frequency. The equal decrease in real and imaginary impedance can be seen in the Nyquist plot and results in a constant phase

shift of 45° . The unit of the Warburg coefficient, σ , is Ohm to the power -0.5. The default value of the Warburg coefficient in the circuit editor is 1000.




Figure 172 Left: The absolute impedance of a pure Warburg element plotted over frequency. Top right: Phase shift of a pure Warburg element over frequency. Bottom right: Nyquist plot of a pure Warburg element.

The Randle's circuit is often used to model the transfer of charge between the electrode and a redox species solution. The Randle's circuit consists of two resistors, a capacitor, and a Warburg element. The first resistor is in series and represents the solution resistance. The electrode's double-layer capacitance is represented by a capacitor parallel to a resistor that represents the charge transfer resistance and the Warburg element. At high concentrations of the oxidant, the semicircle and the diffusion are separated and at low concentrations, they overlap.




Figure 173 Nyquist plot of a Randle's circuit with a large value for the Warburg coefficient of the Warburg element. This represents the transfer of charge between the electrode and redox species with a low concentration of the oxidant (overlap of semicircle and diffusion).




Figure 174 Nyquist plot of a Randle's circuit with a small value for the Warburg coefficient of the Warburg element. This represents the transfer of charge between the electrode and redox species with a high concentration of the oxidant (no overlap of semicircle and diffusion).



Note that some authors refer to the Randle Circuit with just a R(CR) set, so without the Warburg element. This document refers to this circuit as the "Simplified Randle's Circuit"

7.9.5 Warburg Short/Open

In contrast to the Warburg element, the Warburg short/open do not assume semi-infinite diffusion.

$$Z_O = \frac{\sigma B}{(B\sqrt{j\omega})^\phi} \tanh(B\sqrt{j\omega})^\phi$$

$$Z_T = \frac{\sigma B}{(B\sqrt{j\omega})^\phi} \coth(B\sqrt{j\omega})^\phi$$

Respectively the Warburg short (Z_O) and Warburg open (Z_T) model the transmissive and reflective boundaries of a diffusion layer. Some authors use the terms "O element" and "T element" to describe these elements. Nevertheless, it's important to note that these letters refer to the hyperbolic cotangent and tangent, respectively, but in the admittance formulas (thus, the opposite). "B" represents the thickness (l) in meters (m) and the diffusion coefficient (D) of the diffusion layer in square meters per second (m^2/s):

$$B = l/\sqrt{D}.$$

The experimental parameter (ϕ) has a maximum value of 1 when diffusion is uniform and is smaller when diffusion is nonuniform. The default values for the Warburg coefficient, B, and the experimental parameter are 1000 σ , 1 vs and 0.5 respectively.




Figure 175 Left: The absolute impedance of a Warburg open plotted against frequency. Top right: Phase shift of a Warburg open over frequency. Bottom right: Nyquist plot of Warburg open.




Figure 176 Left: The absolute impedance of a Warburg short plotted against frequency. Top right: Phase shift of a Warburg short over frequency. Bottom right: Nyquist plot of Warburg short.

The figure below shows various circuits illustrating the differences in the transfer of charge between the electrode and a redox species a solution, assuming finite diffusion transmissive/reflective boundary or semi-infinite diffusion.




Figure 177 Nyquist plot demonstrating the differences in diffusion at low frequencies for a Randle's circuit (yellow line and bottom circuit), an electrode with a reflective boundary (blue line and top circuit) and an electrode with a transmissive boundary (red line and middle circuit).

7.9.6 Gerischer element

$$Z_G = \frac{Z_0}{\sqrt{k + j\omega}}$$

The effects of a Gerischer element are similar to those of a Warburg open. However, the parameters of the model are different, Z_0 is the magnitude of the impedance at $\omega = 1$ rad/s and k is a reaction rate constant. The Gerischer element introduces a phase lag that is related to the kinetics of charge transfer processes. It is designed to model the effect of an electrochemical species reacting with something to form an inactive substance or absorbed species during the diffusion process.




Figure 178 Nyquist and Bode plots showing response of the Gerischer element

7.9.7 Constant Phase element

$$Z_{CPE} = 1/T(j\omega)^\Phi$$

The Constant Phase element (CPE) is referred to by some authors as “Element Q”. It deserves to describe a frequency dispersion present in most real electrochemical double-layers. The CPE equation has two parameters: T, that is related to its capacitance and Φ , that is the constant phase exponent which is related to the deviation from an ideal capacitor. At a Φ of 1 the behavior of a constant phase element resembles that of a regular capacitor, for a Φ of 0.5 that of a Warburg element, and for a Φ of 0.0 that of a resistor. The unit of T is Siemens per second to the power Φ . The value of the constant phase exponent should be between 0 and 1; it regulates the constant phase which is -90 times Φ . The phase shift of a constant phase element is independent of frequency, the slope of the line in the Nyquist plot corresponds to the

magnitude of the phase shift specified by Φ . The default values of the constant phase element in the circuit editor are $1\ \mu\text{T}$ and $1\ \Phi$.




Figure 179 Left: The absolute impedance of a constant phase element with a constant phase exponent of $1\ \Phi$ (i.e. a capacitor) plotted against frequency. Top right: Phase shift of a constant phase element over frequency. Bottom right: Nyquist plot of constant phase element.




Figure 180 Left: The absolute impedance of a constant phase element with a constant phase exponent of $0.6\ \Phi$ plotted against frequency. Top right: Phase shift of a constant phase element over frequency. Bottom right: Nyquist plot of constant phase element.

The constant phase element is used instead of a regular capacitor to model the double layer capacitance when an electrode displays a frequency-dependent dissipation/dispersion of energy (typical for nonuniform electrode surfaces). When studying a redox reaction without diffusion limitations the dissipation/dispersion of energy presents as a depressed semicircle in the Nyquist. See the figure below.




Figure 181 Nyquist plot of a redox reaction without diffusion limitations, with a solution resistance of 10 Ω , a charge transfer resistance of 1 $k\Omega$ and a constant phase element that represents the double layer capacitance of an ideally polarizable liquid electrode, i.e. a constant phase element with a phase exponent of 1 Φ (or a capacitor).




Figure 182 Nyquist plot of a redox reaction without diffusion limitations, with a solution resistance of 10 Ω , a charge transfer resistance of 1 $k\Omega$ and a constant phase element that represents the double layer capacitance of a solid electrode with dispersion/dissipation of energy. Here the semicircle is shifted downwards as a result of the dispersion/dissipation, modelled as a constant phase element with a phase exponent smaller than 1 Φ .

In the Bode plot, the most notable effect of the CPE compared to a regular capacitor is the reduction of the (negative) maximum phase shift.

7.9.8 Bisquert Short/Open

$$Z_N = \sqrt{\zeta\chi} \tanh\left(L \sqrt{\frac{\chi}{\zeta}}\right)$$

$$Z_M = \sqrt{\zeta\chi} \coth\left(L \sqrt{\frac{\chi}{\zeta}}\right)$$

where $\zeta \equiv \frac{1}{\frac{1}{Z_{R_k}} + \frac{1}{Z_{Q_m}}} = \frac{1}{\frac{1}{R_k} + T(j\omega)^\phi}$ and $\chi \equiv Z_{R_m} = R_m$

The Bisquert Short (N) and Bisquert Open (M) can be used to model porous electrodes. These electrodes have a large surface area but the pores themselves contain a limited volume affecting the rate with which ions enter and react in the pores. L is a fixed parameter that represents the depth of the pore. ζ relates to the impedance of the active porous part of the electrode, R_k models the reaction resistance and Q_m the diffusion / double-layer capacitance (a constant phase element is used to compensate for inhomogeneities of the pores). Finally, the χ relates to the impedance of the electrolyte within the pore. Both elements assume that reactions can only take place at the active porous part of the electrode. The Bisquert Short models an absorbing boundary condition where the contribution of the base electrode to the reaction can short the porous film. The Bisquert Open models a reflective boundary condition where the base electrode is insulating and cannot contribute to the reaction. For both elements the parameters are in the respective order; pore impedance (R_m), reaction resistance (R_k), diffusion of active porous area modeled by a constant phase element (T, ϕ) and pore depth (L).




Figure 183 Nyquist plot demonstrating the differences in impedance of a porous electrode with an absorbing boundary. The blue line ($R_m = 100 \Omega$, $R_k = 1e20 \Omega$, $CPE = (1 \mu T, 1 \phi)$, $L = 1$) represents a situation where the diffusion reaction Q_m is dominant, impedance of $R_k \gg$ than R_m and Q_m . The red line ($R_m = 100 \Omega$, $R_k = 0.1 \Omega$, $CPE = (1 \mu T, 1 \phi)$, $L = 1$) where there is a homogeneous reaction with the absorbing boundary condition and a diffusion reaction.




Figure 184 Nyquist plot demonstrating the differences in impedance of a porous electrode with a reflecting boundary with different relations between reaction rate and diffusion rate. The blue line ($R_m = 20 \Omega$, $R_k = 1e20 \Omega$, $CPE = (0.1 \mu T, 1 \varphi)$, $L = 1$) represents a situation where there is diffusion with a reflective boundary, impedance of $R_k \gg$ than R_m and Q_m . The red line ($R_m = 20 \Omega$, $R_k = 100 \Omega$, $CPE = (0.1 \mu T, 1 \varphi)$, $L = 1$) where the reaction at the active porous part is slower than the diffusion and the yellow line ($R_m = 20 \Omega$, $R_k = 10 \Omega$, $CPE = (0.1 \mu T, 1 \varphi)$, $L = 1$) where the reaction is faster than the diffusion.

For more information, please refer to:

Theory of the Impedance of Electron Diffusion and Recombination in a Thin Layer. Juan Bisquert. (2002) *The Journal of Physical Chemistry B*. 106 (2), 325-333

7.10 Fitting Example

This example demonstrates the fitting of an equivalent circuit on a measurement with two distinct time constants (semi-circles). The measurement used in this example was performed on a PalmSens EIS Dummy Cell described in the figure below.

This file is installed with PSTrace, by default in the folder:

"<username>\Documents\PSData\EIS Examples"

The figure below shows the dummy cell's circuit, the values of its components, and a simulated of a Nyquist plot of the dummy cell.

This example aims to demonstrate that three similar equivalent circuits can all fit the data accurately. Although the quality of the fit is similar the values of the components differ significantly emphasizing the importance of choosing the right equivalent circuit to interpret a measurement.




Figure 185 Nyquist plot of a simulation of the PalmSens EIS Dummy Cell

7.10.1 Designing equivalent circuits

There are two distinct time constants (semi-circles) visible in the Nyquist plot. Time constants can be modeled by placing a capacitor or constant phase element in parallel to a resistance (and optionally a Warburg element). For this example, three different circuits will be used, (1) a circuit with a Voigt structure (the same as the circuit inside the dummy cell), (2) a circuit with a ladder structure and (3) a circuit with a Maxwell structure. Voigt circuits are typically used for modeling the redox process on an electrode. Ladder circuits can be used when there are one or more adsorbed species, or to fit many types of corrosion-protective coatings. Maxwell circuits have been used to study dielectric phenomena.




Figure 186 Three different circuits that can model measurements with two-time constants. 1: A Voigt circuit (the actual model of the dummy cell). 2: A Ladder circuit. 3: A Maxwell circuit.

7.10.2 Fitting the circuits

After the circuit has been drawn the Fit mode can be opened. Often the default values of the components are sufficiently close to their actual values and the fitting algorithm will fit the circuit correctly. However, in some cases circuit fitting algorithms end up in so-called local

minima, this inherits to mathematical optimization algorithms. The fitting of the three equivalent circuits on the measurement is demonstrated below. In the fit of the Maxwell circuit, some tips are also given on how to avoid getting stuck in local minima.

7.10.3 Voigt

The following image shows the fitting tab with at the top a Nyquist plot of the measurement (blue dots) and a simulation of the Voigt circuit with the circuit editors default values (red line). Clicking on the fit button will fit the equivalent circuit on the measurement. Assuming the measurement was performed correctly this circuit will directly give an accurate fit with the default values of the components.




Figure 187 Fitting of the Voigt circuit. Left: the default values of the components result in the small red semi-circle. Right: After pressing the fit button the circuits semi-circles (time constants) no longer overlap and fit the measurement accurately.

7.10.4 Ladder

Just as the Voigt circuit the ladder circuit directly fits on the measured data without first getting stuck in a local minimum. The fitted values of the components are even quite close to those of the Voigt model. However, in the Maxwell circuit the resistor R3 is used to obtain the value of the pseudocapacitance and not for modeling the charge transfer resistance.




Figure 188 Fitting of the ladder circuit. Left: the default values of the components result in the small red semi-circle. Right: After pressing the fit button the circuit's semi-circles (time constants) no longer overlap and fit the measurement accurately.

7.10.5 Maxwell – local minima example

As seen below the Maxwell circuit will not fit directly on the measurement, as it gets stuck in a local minimum. Therefore, it is necessary to refine the initial values of the circuit's components to avoid ending up in this local minimum.




Figure 189 Fitting of the Maxwell circuit. Left: the default values of the components are far from the actual values as the overlapping red semi-circles (time-constants) are much smaller than the ones that were measured. Right: clicking on the fit button will result in the fitting algorithm getting stuck in a local minimum, the model of the semi-circle clearly does not overlap with the measured data also the values several components are unrealistic. Thus, values of the components must be adjusted to obtain an accurate fit.



A general tip for avoiding local minima is providing a rough initial estimation in which the red semicircle(s) (i.e., the modelled time constants) approaches that of the measurement, preferably set the initial values so that the circuit's semicircles are slightly larger than the measured semi circles.

Based on the Nyquist plot the real component (x-axis) of the impedance is close to $25\text{ k}\Omega$ at the lowest frequencies. In the Maxwell circuit resistor R2 is parallel to the second time constant (C2 and R3) and the current will primarily flow through this resistor at low frequencies, thus increasing R2 to $25\text{ k}\Omega$ should increase the total size of the overlapping semicircles. Before we will attempt to fit again, we will change the size of the capacitors to decrease the overlap of the resistors. During the last fit, the value of C2 got stuck at 48 kF which is unrealistically large, therefore, we will change it back to $1\text{ }\mu\text{F}$. The value of R3, $349\text{ k}\Omega$, is also too high, so we will set it somewhere between 0 and $25\text{ k}\Omega$ (this is where the imaginary component (y-axis) of the impedance approaches zero). Then, click on the fit button again and the fit should now be accurate.




Figure 190 Left: the value of resistor R2 is increased to match the real impedance of the circuit at low frequencies ($25\text{ k}\Omega$ as seen on the x-axis on the plot above), as a result, the semicircle is now almost as wide as the measured semi-circles. Right: The unrealistically high capacitance of C2 is lowered back to $1\text{ }\mu\text{F}$, this reduces the overlap of model's the semi-circles.




Figure 191 Left: The very high resistance of R3 is lowered to 10 kΩ, a value between 0 and 25 kΩ (i.e. the range of the real component of the impedance (x-axis)). This changes how the model's semi-circles overlap. Right: with the values of the components as specified in the left image the fitting algorithm now finds an accurate fit.

7.10.6 Inspecting the quality of the fit

The results of the fit are displayed in the table after fitting the circuit on the data. The chi-squared test is a general indication of the fit's quality. The lower the chi-squared value, the better the fit. The square root of the chi-squared statistic gives the average error of the fit. For these fits, this is 0.4% which is very good. The individual errors of the components are also given in the table.

Achieving errors and chi-square values as low as these is uncommon in a real electrochemical cell, where numerous other variables can significantly influence the impedance response. The complexities and variability inherent in practical electrochemical systems often contribute to higher levels of uncertainty and deviations from idealized model fits. Therefore, these exceptionally low error and chi-square values may be challenging to replicate in real-world electrochemical scenarios where multiple factors come into play.

Although the quality of the fit is good, the values of the components for the ladder and Maxwell model cannot be interpreted correctly as the specific processes of the cell that these components represent do not comply with the actual measurement, thus it is demonstrated that blindly assuming that your model is an accurate representation, just because it fits, might lead to circuits lacking actual physical meaning.

Table 10 Fitting results using Voigt model

Voigt			
Element	Fitted Value	Unit	Error%
R 1	201.5	Ω	0.250
C 1	1.028	μF	0.357
R 2	14.98	$\text{K}\Omega$	0.164
C 2	32.54	nF	0.097
R 3	9.997	$\text{K}\Omega$	0.124
Chi-Squared:	1.33E-5	Iterations:	17

Table 11 Fitting results using Ladder model

Ladder			
Element	Fitted Value	Unit	Error%
R 1	201.5	Ω	0.250
C 1	31.54	nF	0.095
R 2	10.63	$\text{K}\Omega$	0.113
C 2	1.042	μF	0.357
R 3	14.34	$\text{K}\Omega$	0.170
Chi-Squared:	1.33E-5	Iterations:	20

Table 12 Fitting results using Maxwell model

Maxwell			
Element	Fitted Value	Unit	Error%
R 1	201.5	Ω	0.250
C 1	31.54	nF	0.095
R 2	24.97	$\text{K}\Omega$	0.092
C 2	0.343	μF	0.343
R 3	18.52	$\text{K}\Omega$	0.200
Chi-Squared:	1.33E-5	Iterations:	10

Visual inspection of the fit's quality is also possible by changing the plot to either the error in the magnitude of the impedance over frequency, the error in the real component of impedance over frequency or the error in phase shift over frequency. These plots represent the difference in magnitude of impedance, real impedance, and phase shift between the fit of the equivalent circuit and the measured data. As seen in the figure below there are subtle discrepancies in impedance at the lower frequencies and phase shift at higher frequencies between the equivalent circuit and the measurement.




Figure 192 Error plots displaying the differences in impedance and phase shift between the equivalent circuit and the measurement.

7.10.7 Adjusting the advanced fit options

In some cases, you may wish to adjust the fitting options to obtain the type of fit that you want. In the Advanced Options window of the circuit editor, several advanced options can be set. These options include selecting a frequency range over which you would like to fit your equivalent circuit and adjustments to the default parameters of the fitting algorithm.




Figure 193 Advanced Options window for fitting EIS data

7.11 Fitting algorithm

There are two different fitting algorithms available, which can be selected in the “Advanced Options” window in the circuit editor: the Levenberg-Marquardt algorithm and the Nelder-Mead algorithm.

7.11.1 Levenberg-Marquardt

The Levenberg-Marquardt algorithm is available for finding the minimum complex non-linear sum of squares. The advantages of this algorithm are that it is robust and generally quick at finding a solution. A disadvantage of the Levenberg-Marquardt algorithm is that it was not designed to work with specified minimum and maximum values for components. As a result, the algorithm sometimes ends up in local minima or saddle points with unrealistic values for the components (i.e., extremely large/small or negative values). Local minima and saddle

points are inherent to mathematical optimization algorithms and in some cases, a good fit will require changing the default values of the components before performing the fit. Typically, when fitting equivalent circuits the Levenberg-Marquart will provide better results than the Nelder-Mead algorithm.

7.11.2 Nelder-Mead

The Nelder-Mead algorithm, also known as the simplex algorithm, is an iterative optimization method commonly used for nonlinear optimization problems, including curve fitting. It is particularly useful when the objective function (the function being minimized or maximized) is not differentiable or when gradients are hard to obtain. Nelder-Mead is best suited for fitting equivalent circuits with a few parameters, for equivalent circuits with more than five parameters we highly recommend the Levenberg-Marquardt algorithm.

Note that the 'Error%' column in the 'Results' table will remain empty when using the Nelder-Mead algorithm.

7.11.3 Maximum Iterations

Generally, it will not be necessary to increase the maximum number of iterations for fitting. As the Levenberg-Marquardt algorithm usually requires few iterations to achieve a fit. However, when the minimum step size parameters are lowered this could be necessary. When you do not achieve a good fit and see that the number of iterations was 500 in the table of the fit tab increasing the maximum number of iterations may help. Usually, it is better to carefully review the design of your circuit and enter more appropriate values for its components before fitting again.

Typically, the Nelder-Mead algorithm requires more iterations to converge than the Levenberg-Marquardt. The number of steps also depends on the convergence tolerance, the lower this is set the more steps it will take to converge.

7.11.4 Minimum Step Size

The minimum delta sum of squares and minimum parameter step size are the desired stopping conditions for the fitting algorithm.

The Levenberg-Marquardt algorithm attempts to reduce the sum of squares with each step, the minimum delta sum of squares specifies when further improvements become negligible. This can be lowered to achieve a more precise fit for the value of the components. For this to have an effect you will usually have to lower the minimum parameter step size as well. It is questionable whether lowering the minimum delta sum of squares will improve your fit as generally the accuracy of the fit is limited by the fact you are fitting a model which is a simplification of the actual ongoing processes in the cell and that will never fit the measurement perfectly.

As the Levenberg-Marquardt algorithm approaches the global minimum (i.e. the desired fit) or a local minimum or saddle point (i.e. an undesired fit) the differences in the values of the components become smaller each step. Lowering this value is recommended if the expected component in the circuit has a value in the range of 1-10 pico Farad/Henri/Ohm etc.

Minimum step size does not apply to the Nelder-Mead algorithm.

7.11.5 Lambda

The Levenberg-Marquardt algorithm uses lambda as a damping/scaling parameter. The lambda factor increases/decreases the change in the values of the components for the next step. In the event of a successful step (i.e. the sum of squares decreases) the lambda factor decreases and becomes increasingly smaller as the algorithm approaches a minimum. In the event of an unsuccessful step, it becomes larger allowing the algorithm to approach a minimum with fewer steps.

The starting lambda value can be lowered in case you want to further optimize a fit that was already good or increased to find a minimum in fewer steps. The same holds for the lambda scaling factor reducing it will increase the number of iterations required but may be useful if further optimization is required and increasing it could decrease the number of iterations required to find a minimum or it could make the algorithm unstable.

7.11.6 Convergence Tolerance

The Nelder-Mead will stop when the simplex shrinks as it approaches a global/local minimum or a saddle point. Lowering the convergence tolerance will let it shrink further and could help it move along a saddle point and converge to a local minimum. For equivalent circuits with many parameters, convergence is difficult to achieve.

8 Plot, curves and data

In this section, all available functions related to measured curves are explained. Impedance data is presented differently from regular curves. For information about working with impedance data, see section:

[Electrochemical Impedance Spectroscopy](#) on page 113.

8.1 Handling curves

There are two Plot tabs in PSTrace; impedance measurement curves (EIS curves) are shown in the 'EIS Plot' tab, all other curves are shown in the 'Plot' tab.




Figure 194 PSTrace window showing the Plot tab.

The legend next to the Plot tab shows all available measurements (session data) with corresponding curves.

A 'Session' can contain multiple measurements with curves.

See also section [Files](#) on page 235.

8.2 Session data

When a measurement is started the resulting data is added in the Session data list or Legend next to the plot window and the default curve is added.




Figure 195 The legend showing session data

The curve selected in the curves list determines which curve is active. The active curve is used when saving a single curve, showing data in the Data tab, marking peaks or slopes, etc.

The asterisk (*) next to the titles indicates the measurement or curve is not saved yet. They will disappear when the data is saved to a session file.

When a measurement is selected in the legend, a window appears next to the legend showing all data that was recorded during the measurement together with some other information about the measurement:

8.2.1 Automatic naming of curves

In the case of a new measurement, the curve is named with the technique acronym followed by the respective plotted variables "y vs x", i.e.: "LSV i vs E". In the case of a Cyclic Voltammetry measurement, the curve receives an addition in the name 'Scan n', where n is a sequential number. In the case of a multiplexer measurement or measurement using an auxiliary channel, it receives the name of the channel used.

8.2.2 Measurement Tool window

Selecting a Measurement item in the Legend opens the following tool window for this measurement.

Selected measurement

Title: **RM Au**

[Click here to add notes...](#)

Technique: Cyclic Voltammetry
 File: C:\...\RM Au electrode.psession
 Date: 22/02/2017
 Start time: 17:09
 Device used: Unknown, EmStat2 (PS4A16A000010)
 Datapoints: 239

Generate new curve

Select data for Y axis:

- time/s
- potential/V
- Scan 1/ μ A
- charge/ μ C**
- Scan 2 current/ μ A
- Scan 3 current/ μ A

vs.

Select data for X axis:

- time/s
- potential/V
- Scan 1/ μ A
- charge/ μ C
- Scan 2 current/ μ A
- Scan 3 current/ μ A

Add Q vs t curve

Preview plot: Charge/ μ C vs Time/s. The plot shows a cyclic voltammogram with a peak at approximately 6 seconds.

Buttons:

- Delete:** Delete the recorded Measurement including all related data and information (notes, Method parameters used).
- Add Q vs t curve:** Add curves generated from the data recorded during the measurement.
- View method parameters:** View the method parameters used for this measurement. Also useful for changing automatic peaks detection settings for the curves in this measurement.
- View raw data:** View the raw data recorded for this measurement in the 'Data' tab.

Annotations:

- The measurement title as shown in the legend
- Notes as specified in the method editor before starting the measurement. Can be in this window as well.
- Preview plot of the new curve generated from the selected data.

Figure 196 Window with measurement information shown next to the selected measurement in the legend.

8.2.3 Curve Tool window

Selecting a Curve item in the Legend opens the following Curve tool window.




Figure 197 Windows with tools for the selected curve shown next to the legend after selecting a curve

8.3 Switching between plots

In case the legend contains curves with different units on the X and/or Y axis, a plot selection bar appears on top of the legend that allows for switching between different plots.




Figure 198 The plot selection toolbar





When starting a measurement, the appropriate plot will automatically be selected, or added if necessary. The automatic selection of the active plot can be disabled in the 'Plot visual settings window'.

See also section: [Plot toolbar](#) on page 179.

8.4 Plot toolbar

The plot toolbar at the left-hand side of the plot contains the following buttons.

8.4.1 Scaling and selection




	Automatic scaling to show all data points.
	Use pointer to select a datapoint.
	Use the pointer to drag an area that defines the new scaling boundaries.
	Enter minimum and maximum values for both axes.

8.4.2 Clear helplines



Removes all existing helplines (for example LLS or integration lines).

8.4.3 Plot options

	Toggle smart scaling to automatically round up / down minimum and maximum values for both axes.
	Toggle symbol visibility on each datapoint in the curve(s).
	Opens the 'Plot visual settings' window.

8.4.4 Plot visual settings window




Figure 199 The 'Plot visual settings' window

The 'Plot visual settings' window allows for changing the default appearance of the plot.

The 'Default Symbol for datapoint' can be set here for all new measurements. This setting is saved upon closing PSTrace. The button 'Apply to all existing curves' makes sure that the new setting is also applied on any existing curves in the plot.

The 'Plot area' appearance can be customized by toggling the checkboxes for

- 'Background gradient':
Default set to checked. When unchecked, the light blue gradient in the plot is replaced with a plain white background.
- 'Grid':
Default set to checked. When unchecked, the grid lines in the plot are no longer showing.
- 'Show axes origin lines':
Default set to checked. When unchecked, the blue horizontal and vertical lines at the axis 0 point are no longer showing.
- 'Autoscale during measurement':
Default set to checked. When unchecked, the plot will no longer automatically rescale when new added datapoints during a measurement are out of the axis bounds. This can be helpful to reduce CPU usage when running (very) long duration measurements with many datapoints.
- 'Autoswitch active plot':
Default set to checked. In case curves with different units are present in the legend, a plot selection toolbar is shown above the legend. If the checkbox is set, the relevant plot will automatically be selected when a new measurement is started.
See also section: [Switching between plots](#) on page 178.

- 'Use polarographic conventions (default IUPAC)':
There are two common conventions for Voltammetric plots both of which are supported in PSTrace:
 - The IUPAC (International Union of Pure and Applied Chemistry) makes sure that anodic currents are plotted in a positive direction and cathodic currents in a negative direction.
 - If the Polarographic convention (also known as the North American, Polarographic, or Classic convention) is used, anodic currents will be plotted in a negative direction and cathodic currents in a positive direction.
- 'Draw peak search windows (in Analytical mode)':
If enabled the potential windows that are assigned to an analyte are made visual in the plot. This helps the user to visually check the chosen potential windows.
- 'Convert area to charge':
By enabling this checkbox, the pop-up shown when clicking a peak label in the plot shows the charge for the peak next to the label for area:










Figure 200 The peak label showing peak charge

8.4.5 Axes options

The data shown in the plot can be presented in different ways, depending on the data available. The table below shows which options for the axes are available.




Table 13 Buttons for different plot options

	Plots the derivative of the curve using the centered three-point method.
	Plots the second derivative of the curve.
	Plots the integral of the curve.
	Plots the logarithm (Log ₁₀) of the Y values.
	Shows the Cottrell plot.
	Shows the Anson plot.

8.4.6 Plot export options



The data that is presented in the plot at any given time can be exported to third-party software. The table below describes the buttons found in the plot area for exporting the data.

Table 14 Options for exporting the data that is shown in the plot

	Takes a snapshot of the plot area and copies the result to the clipboard for pasting in other software. The snapshot settings can be changed in the 'Plot visual settings' window.
	This button is only visible when exporting to Excel is enabled. See also section Exporting curves on page 207.
	This button is only visible when exporting to Origin is enabled. See also section Exporting curves on page 207.

The data in the plot can be presented with additional information using the buttons described in the table below.

Table 15 Print buttons in the Plot toolbar

	Shows a print preview and printing options for the plot
	Opens the general printer settings window, allowing you to select the default printer for printing in PSTrace and changing the default printer settings.

8.4.7 Cottrell plot



The option for displaying the Cottrell plot is only visible for current measured versus time (multistep amperometry).

The Cottrell Equation is well known in electrochemistry. It describes how the current I due to a reaction of a free diffusing species develops over time t :

$$I = zFAc^* \sqrt{\frac{D}{\pi t}}$$

Where D is the species' diffusion coefficient, z the number of electrons transferred per molecule, A the active area of the electrode, c^* the bulk concentration and F the Faraday constant.

If a potential step is performed from a reaction free potential to a potential where a reaction is no longer potential controlled, but diffusion controlled, we expect the current to look like in the figure below.




Figure 201 a) current response to a potential step (Cottrell experiment) b) linear form of a)

In the Cottrell plot the current I is plotted versus $t^{1/2}$ and according to the Cottrell equation a linear plot is expected as shown in the figure above (b). The slope m of this plot is

$$m = zFAc^* \sqrt{\frac{D}{\pi}}$$

Depending on which factors are known the Cottrell plot can be used to determine for example D , c^* or A . The challenge is to choose the right part of the Cottrell plot to perform the linear fit.

In the initial and final segments of the curve, non-linearity is observed. The non-linearity at the curve's end is attributed to the diffusion layer's inability to expand infinitely due to convection or other influencing factors. Two factors contribute to the non-linearity at the curve's beginning and end. Firstly, achieving an ideal potential step is unattainable as potentiostats have rise times in the order of nanoseconds or microseconds. Despite reaching the set potential swiftly, it remains slower than perfectly instantaneous. Secondly, the capacitive current or capacitive charging current arises from ions forming a capacitor at the electrochemical double layer. These combined factors introduce deviations from the ideal potential step, influencing the curve's shape at its initial and final points.

8.4.8 Anson Plot



The option for displaying the Anson plot is only visible for charge measured versus time (chronocoulometry).

While measured current in the MultiStep Amperometry (MSA) shows a linear correlation with the reaction rate, the charge measured in Chronocoulometry (CC) shows a linear correlation with the amount of converted reactants. This is described in Faraday's law. This characteristic renders Chronocoulometry (CC) a valuable tool for quantitative analysis, especially in applications such as electrodeposition or electrochemical synthesis. The linear relationship it exhibits enables precise measurement and analysis of the quantity of reactants transformed during the electrochemical process. Additionally, Chronocoulometry proves beneficial in the characterization of batteries, as it can be employed to accurately determine their capacity.

In analytical chemistry CC is used to determine the adsorbed amount of active species in a solution of free diffusing active species. Multiplying the Cottrell equation with t delivers that the charge Q_{diff} of the free diffusing species' reaction is proportional to $t^{3/2}$.

$$Q_{diff} = zFAc^* \sqrt{\frac{Dt}{\pi}}$$

Where D is the species' diffusion coefficient, z the number of electrons transferred per molecule, A the active area of the electrode, c^* the bulk concentration and F the Faraday constant.

The total charge Q during a potential step like the Cottrell experiment (see Cottrell plot) also includes the contribution of the charge stored in the electrochemical double layer Q_{dl} and the charge due to reactions of adsorbed species. Q_{ads} (see figure below)

$$Q = Q_{diff} + Q_{ads} + Q_{dl}$$

The cause of Q_{ads} and Q_{dl} are a lot faster than the reaction of free diffusing species. Plotting Q versus $t^{1/2}$, known as Anson plot, delivers (ideally) a jump in charge followed by a linear increase (see figure below).




Figure 202 a) current response to a potential step (Cottrell experiment) b) linear form of a)

When the linear part is extrapolated the intersection with the charge axis delivers $Q_{ads} + Q_{dl}$. A previous blank measurement allows the determination of Q_{dl} and thus the calculation Q_{ads} .

Analog to the Cottrell plot the Anson plot won't show linear behavior at the beginning and the end for real experiments for the same reasons.

To make sure that the charge step, in the beginning, is recorded CC is performed in 2 steps. The first step at a reaction-free potential is followed by the step where the reaction is initiated, just like the classic Cottrell experiment.

8.5 Curve toolbar

This section describes the tools found in the curve toolbar.



Figure 203 The Curve toolbar

8.5.1 Marking peaks manually

The buttons to manually mark peaks in the plot can be found in the curve toolbar under the button 'Manual peak'.

8.5.1.1 Fixed baseline



Draw a baseline fixed on the curve for peak detection.




Figure 204 Baseline on curve

8.5.1.2 Free baseline



Draw a baseline that is not attached to the curve.




Figure 205 Free baseline

8.5.1.3 Three points



Mark three points on the curve to determine the left and right sides of the baseline and the top of the peak.




Figure 206 Three points determine baseline and corresponding peak

8.5.1.4 Extrapolate baseline



A linear slope (using linear regression) is marked on the curve from which the baseline is extrapolated.




Figure 207 Extrapolate from slope

As soon as the slope for the baseline line has been drawn, the Curve calculations window is shown:




Figure 208 Calculations window showing the properties of the slope

The values shown in the window can be changed manually and applied again on the curve.

8.5.1.5 Non-linear baseline



A polynomial with a specified order is fitted through the selected points on the curve.

The following buttons are shown during selection:




Figure 209 Buttons to set the degree for the polynomial baseline

The selected 'Max. order' defines the maximum (n) order, or degree, applied for the polynomial:

$$a_0 + a_1x + a_2x^2 + a_3x^3 \dots + a_nx^n$$

The number in red represents the current order 'n' applied. When a polynomial is set, the button 'Apply' shows a preview of the peak using the polynomial as the baseline. Each time the polynomial is changed; the 'Apply' button should be clicked to mark a new peak using the changed polynomial as the baseline. The 'Close' button should be clicked when a satisfactory peak is found.




Figure 210 Marking a peak in a CV curve

The baselines can also be used for subtraction. This can be done in the Curve Operations Window. See section [Curve operations](#) on page 197.




Figure 211 Curve operations window

8.5.1.6 Add notes



Annotate a feature of a curve by marking a point and adding a note

The following buttons are shown during selection:

Notes: Add Note Finish

Figure 212 Input field and buttons to enter and set a note to mark feature on a curve

Mark a point to annotate by clicking in the plot with the left mouse button and enter the note you want to add in the input field.




Figure 213 Marking a point on a curve to annotate

Press the 'Add Note' button to add the note to the curve. Labels in the plot marked with 'i' symbol contain notes.




Figure 214 A plot label containing a note

8.5.1.7 Manual level



Add a level to an amperometric or potentiometric curve with an optional note and/or concentration

The following buttons are shown during selection:

Concentration: 0 mol/L Notes: Add Level Finish

Figure 215 Input fields and buttons to add a manual level

Mark a level by clicking in the plot with the left mouse button and enter the note and/or concentration you want to add in the input fields.




Figure 216 Marking a level manually

8.5.2 Peaks toolbar buttons

Detect peaks



Automatically detect peaks in all curves or a selected curve. A single curve can be selected using the small drop-down arrow next to the button. The settings in the [Peaks tab](#) on page 39 are used for peak finding. Results are shown in the Peaks Data window. Note: This button is replaced by the 'Find levels' button if the last measurement was versus time. This can be reversed in the Peaks tab.

8.5.2.1 Remove peaks



Remove all peaks or remove the peaks of a specific curve. Peaks in a single curve can be removed using the small drop-down arrow next to the button.

8.5.2.2 Find levels



Available if measured as a function of time and *Current Level* is checked in the Peaks tab. The settings in the Peaks tab are used for the level finding. Results are shown in the Peaks Data window.

8.5.2.3 Remove levels



Remove all levels or remove the levels of a specific curve. Peaks in a single curve can be removed using the small drop-down arrow next to the button.

8.5.2.4 Windows



Overlapping peaks separation window.




Figure 217 Window showing separation of overlapping peaks in a CV measurement.

The [overlapping peaks separation window](#) on page 192 allows for the separation of two or more overlapping peaks. This tool window is only accessible for LSV and CV measurements.



Toggle peaks or levels data window.





Figure 218 Windows showing peak (left) or level (right) data

The table shows the values for the peak potential, peak height, peak area, peak width at half-height, the value at the baseline as well as the maximum values for the ascending and descending slopes, and the sum of both values for each curve in a measurement. These results can be copied to excel or the clipboard for the selected or all measurements. The measurement's automatic peak detection parameters can be adjusted, and the automatic peak detection can be performed again. In the case of a curve plotted versus the time it is possible to switch between peak and level detection.

8.6 Separating overlapping peaks

Overlapping peaks in a cyclic or linear sweep voltammetry measurement can be separated using the peak separation window. The peak separation window can be opened from the toolbar above the plot . Peak separation works best for reversible peaks, the amount of overlap and number of overlapping peaks determine how well the peaks can be separated from each other (separating peaks is not guaranteed to work properly in all cases).






Figure 219 Overview of the peak separation window

First, select the measurement that contains the curve that you would like to separate in the legend of the Plot tab in the main window. Then the peaks or levels data window  should be opened and the  separate overlapping peaks button clicked. Select the anodic or cathodic curve that you would like to separate peaks in from the dropdown menu in the top left corner of the peak separation window.




Figure 220 The dropdown menu for selecting the curve

Before you can separate the overlapping peaks, they should be marked in the plot to the right. Click on mark peaks, this will change the view in the plot to the differintegral view. The peaks are still at approximately the same potential, however, in this view, they are more distinct.




Figure 221 Button for the first step; marking the peaks to separate




Figure 222 Left: Anodic curve with several overlapping peaks. Right: Differintegral view of the same curve.

Once you selected the peaks you want to separate, click on 'Apply' Tip: you can zoom in by scrolling the mouse and pan using the right mouse button. Selecting peaks that do not overlap will not work and selecting more than three peaks takes much longer to compute.




Figure 223 Marking the peaks in a curve

Once the peaks have been marked, they can be separated by clicking on the respective button. This will change the plot back to the regular view and display the results.




Figure 224 Results of the peak separation. Blue curve: original curve. Red and green curves: curves of the peaks that were separated. Yellow curve: baseline curve (both peaks removed). Purple curve: sum of the separated curves.

The results are also shown in the table on the left. This table can be used to hide/show and show the peaks corresponding curves, select curves to add to your measurement and to delete curves (by clicking on the red cross).

Results		
	Curve	
<input checked="" type="checkbox"/>	Curve: Anodic, Peak 1	✕
<input checked="" type="checkbox"/>	Curve: Anodic, Peak 2	✕
<input checked="" type="checkbox"/>	Curve: Anodic, Baseline	✕
<input checked="" type="checkbox"/>	Curve: Anodic, Reconstructed curve	✕

Add selected curves
Clear all

Figure 225 The list with the curves as a result of peak separation

The results can be added to your measurement in the main Session manager (the legend in the main window showing all available session data) using the add button below the table. The clear button removes all the results from the table and the plot. It is also possible to export the curves to excel, this can be useful when you need to export the fractional differintegral (the checkbox below the plot).

8.6.1 Advanced peak separating settings

In some cases, the separation of the peaks can be improved by adjusting the default settings in the advanced settings.




Figure 226 Advanced settings for peak separation

The order of the differintegral is specified in the fractional derivative order text box. By default, the 0.5 semi-derivative is used to separate the peaks from each other. The peaks are best separated when they do not or barely overlap in the differintegral view and when the beginning and end of the peaks are on the same level / the baseline is not slanted.




Figure 227 Results of the peak separation with the default settings. Notice that the baseline (purple curve) is slanted in the differintegral view (right), this often leads to sub optimal results (the purple curve deviates from the original blue curve).





Figure 228 Increasing the order of the fractional derivative to 0.61 improves the separation of the peaks significantly.

The show as differintegral check box toggles between the regular and differintegral view. Toggling back to the differintegral view is useful when optimizing the separation of the peaks.

The window size of the moving average filter is used to prevent noise from interfering with the peak separation. The separation works best with smooth curves, however, there is a trade-off between large and small window sizes. Large window sizes can result in an underestimation of the peak height and small window sizes can result in errors when separating the peaks.

Finally, depending on the type of reaction switching between the reversible and irreversible peak model can improve the peak separation. The reversible peak model works best with peaks that are symmetrical in the differintegral view and the irreversible peak model works best with asymmetrical peaks (peaks where the second half is steeper).

8.7 Curve operations

The curve operations window is available in the Curve toolbar and can be opened with the button .




Figure 229 Curve operations window


1. Select one or multiple curves to apply operations on.
2. Choose the mode of operation.

New curves will be generated and added to the corresponding measurement, based on the selected operation.

For more information about the Baseline buttons, see the following sections:

- [Linear baseline subtraction](#) on page 197
- [Non-linear baseline subtraction](#) on page 198
- [Moving average baseline](#) on page 200

8.7.1 Linear baseline subtraction

To do a linear baseline subtraction, click on the 'Subtract linear baseline' button in the Curve Operations window. This is opened using the  in the plot toolbar. Or in the menu 'Peaks' → 'Curve operations'.




Figure 230 The button for starting linear baseline subtraction in the “Perform curve operations” window.

Clicking on this button will open a new window with a plot of the selected curve. Set two markers on the curve. As soon as the second marker is set, a preview is shown. A third click on the plot will remove the preview, reset the markers and place the first one again. If the preview seems satisfactory, click ‘Accept’.





Figure 231 Plot showing a preview of a linear baseline subtraction.

8.7.2 Non-linear baseline subtraction

To do a non-linear baseline subtraction, click on the ‘Subtract non-linear baseline’ button in the Curve Operations window. This is opened using the  button in the plot toolbar.




Figure 232 The button for starting non-linear baseline subtraction in the “Perform curve operations” window.

The following buttons are shown during selection:




Figure 233 Buttons to set the degree for the polynomial baseline

The selected ‘Max. order’ defines the maximum (n) order, or degree, applied for the polynomial:

$$a_0 + a_1x + a_2x^2 + a_3x^3 \dots + a_nx^n$$

A polynomial with a specified order is fitted through the selected points on the curve. A polynomial requires a minimum of n+1 markers. So, a 3rd order polynomial requires three markers.

The number in red represents the current order n applied. When a polynomial is set, the button ‘Apply’ shows a preview of the subtracted curve using the polynomial as the baseline. Each time the polynomial is changed; the ‘Apply’ button should be clicked to update the preview. The ‘Accept’ button should be clicked when the preview is satisfactory.





Figure 234 Setting markers on a DPV curve.

8.7.3 Moving average baseline subtraction

To do a moving average baseline subtraction, click on the 'Moving average baseline' button in the Curve Operations window. This is opened using the  in the plot toolbar.



Moving average baseline subtraction is not applicable for Cyclic Voltammograms.




Figure 235 The 'Moving average baseline' button in the Curve Operations window

By clicking on the 'Moving average baseline' button, a preview window opens. In this window, the window size and the maximum number of sweeps can be adjusted to obtain the desired result. The grey line represents a preview of the baseline correction, it is updated instantly after changing the window size and the maximum number of sweeps.




Figure 236 Moving average baseline settings

When performing a moving average baseline correction, it is recommended to change the window size first.

Increasing the window size will increase the number of data points that are used to determine the new value of the data points with the moving average correction. The ideal window size depends on the number of data points in your measurement and the width of the peaks. A window size that is too large will result in unwanted rounding errors.

When the desired result cannot be obtained by increasing the window size the maximum number of sweeps (i.e., the iterations of the moving average correction applied to the curve) can be increased to achieve the desired effect.

To see a preview of the moving average baseline correction, click on 'Apply'. Clicking on 'Accept' will add the curve to the session manager.




Figure 237 Plot with result preview of the baseline correction

8.7.4 Discarding data points

Using this option, it is possible to create a new curve from a specific selection of a curve. Specifying the selection of the curve to keep with the mouse by moving your mouse over the sample from where you would like to start or end your selection, then clicking on that sample (blue dot) twice. Then move your mouse to the data point where you would like to respectively end/start your selection and click on it once. The selected range for the new curve is indicated by the solid markers.




Figure 238 Discarding datapoints in a curve.

Alternatively, one or more curve points can be deselected by holding down both the Ctrl-button and the left mouse button and dragging an area around the outlier(s) with your mouse. Repeating this step on the points again will reselect them.





Figure 239 Selecting a series of data points in a curve for discarding.

Click 'Accept' to add the new curve based on the selection to the current session.

8.8 Curve calculations

Use the button  to open the Curve calculations window.




Figure 240 Curve calculations window

The values obtained by using the pointer are shown in the 'Use manual values:' window so they can be fine-tuned and re-applied. These values are saved when PSTrace is closed.

Using the button 'Copy to clipboard' the results can be pasted (usually using Ctrl+V) as plain text into any other program.

The 'Match X range with selected curve' button selects the entire range of the curve.

8.8.1 Cyclic Voltammograms

In case the targeted curve is a Cyclic Voltammogram, a selection tool appears below the plot toolbar, allowing the use to select the Anodic or Cathodic part of the curve for doing the calculations on. Use the 'Finish' button in the toolbar when done setting the markers.




Figure 241 Toolbar for selecting the Anodic or Cathodic part of a CV curve for calculations.

8.8.2 Extrapolate



See for more information regarding Extrapolate baseline button section: [Marking peaks manually](#) on page 184.

The buttons for integration and marking a linear regression line can be found under the *Curve operations* button.




Figure 242 Calculations on the curve

8.8.3 Linear regression



Mark begin and end on the curve for linear regression using linear least squares (LLS).




Figure 243 Calculation of a slope using LLS.

As soon as the line has been drawn, the Curve calculations window with the results is shown:




Figure 244 Calculations window showing the properties of the slope.

The values shown in the window can be changed manually and applied again on the curve.

8.8.4 Integration



Draw a baseline for integration on a curve.




Figure 245 Result of drawing a baseline for integration.

The integration result is shown in the Curve calculations window:




Figure 246 Calculations window showing the integration results in two different units.

The values shown in the window can be changed manually and applied again for integration of part of the curve.

The “Match X range with selected curve” button selects the entire range of the curve. This can be useful when determining the total charge of a curve for example.

8.9 Saving data

All measurements and curves that are present in the Plot and EIS Plot tabs and the method present in the Method Editor can be saved to a single ‘.psession’ file.




Figure 247 A PSession (.psession) file contains basically all the data available at a given moment.

Data is saved via the menu: 'Data' → 'Save data...'




Figure 248 Data menu




Figure 249 Save Data dialog for saving to a PSession file.

In the tree at the right side of the Save Data dialog window data can be selected that will be included in the PSSession file.

8.9.1 CSV

Files can also be exported to a CSV (Comma Separated Values) format. This is a common file format supported by many applications like Excel, OpenOffice Calc and Origin.

8.9.2 OriginLab

Files can be exported to a native .OPJ file for use in Origin, from OriginLab Corporation.

See also section: [Files](#) on page 235.

8.9.3 Legacy curve files

Since the release of PSTrace v5 and MultiTrace v4, the following file formats can no longer be used for saving data:

.pss, .pst, .mux, .psd and .pseisdata.

The new “.pssession” file format has been introduced to consolidate and replace the need for these various formats. It's important to note that although the old file formats cannot be used for saving, they can still be loaded in PSTrace v5 and MultiTrace v4.

See also section [Files](#) on page 235.

8.10 Exporting curves

PSTrace allows curves to be exported to third-party software with a single click.




Figure 250 Export options shown in the Plot menu.

8.10.1 Excel

If Excel version 2007 or newer is installed an 'Open in Excel' menu item is added to the Plot menu and an extra button is added at the left-hand side of the Plot. Clicking the menu item or the button will open an Excel window and automatically export the data of all visible curves to the spreadsheet. A graph will be generated for all exported data. The template for the graph can be found in the installation folder of PSTrace.

8.10.2 Origin

If Origin (from OriginLab Corporation) is installed an 'Open in Origin' menu item is added to the Plot menu and an extra button is added at the left-hand side of the Plot. Clicking the menu item or the button will open an Origin window and automatically export the data of all visible curves to a separate book and graphs will be generated for each exported curve.

If the button is not showing and Origin is installed, open the General settings window in the menu ('Tools' → 'General settings...') and check in the 'Plot and data' tab if the location for Origin is properly set.




Figure 251 Location for Origin specified in the Settings window.

9 Manual Control

The Manual Control tab enables direct control of the instrument, allowing users to set conditions in real time without the necessity of configuring a method.

9.1 Manual Control tab

The 'Manual Control' tab shows the following options to control the connected potentiostat directly.




Figure 252 'Manual Control' tab

- Cell On is used to establish a connection for polarizing the Working Electrode (WE) to a specific potential or current. By default, the cell is in the "Off" state, signifying that the Counter Electrode (CE) lead is unplugged, and thus there is only a potential reading. If the BiPot module is present, WE2 can also be switched on or off.
- Control E or i is used to set the mode of PalmSens. Selecting E is to set the potential and to read the current response, i.e. the instrument is used as a potentiostat.

Selecting *i* will allow to control the current and to read the potential, i.e. the instrument is set as galvanostat. In short: Potentiostat (E) or Galvanostat (i).

- Potential or Current allows to set either one of these values, depending on the set mode with 'Control E or i'. The scrollbar and textbox are used to specify the potential or the current.
- Current Range or Potential Range specifies the sensitivity of the measurement.

Table 16 The maximum measured currents for each instrument.

Instrument	Max. measured current
EmStat1, 2 and 3(+)	±2 * selected range
PalmSens1 and 2	±2 * selected range
PalmSens3	±3 * selected range
PalmSens4	±6 * selected range
EmStat4 LR and HR	±3 * selected range
EmStat Pico	±0.6 * selected range

An overload warning is given for measured values $\geq 80\%$ of the selected range.

For all instruments except for the EmStat Pico an underload warning is given for measured values $\leq 4\%$ of the selected range. For the EmStat Pico the underload warning is given for measured values $\leq 2\%$ of the selected range.

For instruments that support multiple potential ranges, the table below shows the available range for each instrument.

Table 17 The maximum measured potentials for each instrument.

Instrument	Max. measured potential
PalmSens4	±10 * selected range
EmStat4 LR	±3 * selected range
EmStat4 HR	±6 * selected range

- Multiplexer Channel is enabled if the use of a multiplexer is enabled in the Settings window (Tools → General settings...). In this box the channel which is used for control and measurement can be selected.
- The WE2 Potential and WE2 Current Range frames are visible if BiPot is enabled in the Configuration window.
- Checking Stirrer On sets the digital ports d0 to high and d1 to low (and vice versa) of the auxiliary port of the instrument. See the instrument-specific sections in this manual for the auxiliary port pin-out.
- Analog In reads the voltage applied on the auxiliary analog input port of the connected instrument.
See the instrument-specific sections in this manual for the auxiliary port pin-out. The function of the Analog In is determined by the Auxiliary input settings found under the menu: 'Tools' → 'General settings...' → 'Auxiliary input'
- Analog Out set a specific potential between on the auxiliary analog output port of the connected instrument.
See the instrument-specific sections in this manual for the auxiliary port pin-out.
- In case the 'Speed control' for the stirrer is enabled, Analog Out controls the stirrer speed.
- Digital Output allows for setting the digital output lines on the auxiliary port manually.

- See the instrument-specific sections in this manual for the auxiliary port pin-out.
- Digital Input allows for reading the digital output lines on the auxiliary port manually.
See the instrument-specific sections in this manual for the auxiliary port pin-out.

10 Analytical mode



The analytical extension of program PSTrace for Voltammetric Analysis provides the possibility to do quantitative analysis by means of Standard addition or using a Calibration curve.

10.1 Supported techniques

The supported electrochemical techniques for voltammetric analysis are shown in the Technique list:

- Linear Sweep voltammetry
- Differential Pulse voltammetry
- Square Wave voltammetry
- Chronopotentiometric Stripping




Figure 253 List showing available techniques for the Analytical mode

The other available techniques for PalmSens or EmStat are not listed when PSTrace is set in Analytical mode. These techniques are identical to those found in the Scientific Mode, employing the same signal and adhering to the same specifications. The only distinction lies in

the outcome, as the data results are directed to calculation tools for concentration determination in the Analytical Mode.

10.2 Setting up parameters

When the program is set to 'Analytical mode' in the Settings window or the top left of the main window, new tabs called 'Analysis' and 'Analytical result' are shown. The 'Analysis' tab in the Method Editor is used to enter analysis parameters beforehand for a series of measurements. The 'Analytical result' tab is used to show the standard addition or calibration curve(s) and the calculated peak heights as well as the obtained concentration. In this tab, it is also possible to edit the analysis parameters of existing measurements. The parameters of the standard solution as well as the analyte(s) or component(s) to be determined can be defined in the 'Analysis settings' group in the top left of the tab. The program allows measurement of up to four analytes using up to four standard solutions.

PSTrace allows subtraction of a separately measured blank voltammogram. This blank voltammogram can be measured before the other curves or can be loaded as data file.

See for more information about using blank curves, section: [Blank subtraction](#) on page 45.

Analysis parameters can be entered *before* a measurement in the 'Analysis' tab of the Method Editor or edited for existing measurements in the 'Analysis settings' group of the 'Analytical result' tab.



There is no use in changing the parameters in the 'Analysis' tab of the Method Editor *after* the measurement have been carried out.




Figure 254 'Analysis' tab in Method Editor (left) and the 'Analysis settings' group in the 'Analytical result' tab (right).

The mode of the determination is either Standard addition or Calibration. The content of the table is either Added volumes or Concentration in the cell.

Analytes
Concentration unit: ppm ▼

☒ Added volumes in μl ☐ Concentration in cell

	Analyte	Analyte	Analyte	Analyte
► Id	Cu	Pb	Cd	
Sol nr.	1	2	3	4
Conc	10	10	10	10000
St. 1	50	100	100	10
St. 2	100	200	200	10
St. 3	100	100	10	10
St. 4	10	10	10	10

Figure 255 Volumes added for each analyte

To specify the concentration unit, choose from the provided list. It's important to clarify that the selection of the concentration unit does not impact the calculations derived from the curves. The software does not execute a unit conversion between the calibration and the result. Therefore, the chosen unit primarily serves as a format for display purposes and does not alter the numerical outcomes.

The column Id is used to specify which analyte(s) or component(s) has to be determined, e.g. Cd, Pb, Cu, antibiotics, pesticides. If an Id field is empty, the column is marked in grey and the specified values are not used.

The Id's of the peaks are also given in the voltammetric plots in the peak values, *after* the 'Recalculate' button was used.




Figure 256 Peaks found for each Id showing in the Plot tab

10.2.1 Specifying added volumes

If the Table content is Added volume in μl then the fields St. n of the table contains the added volumes of standard solutions in μl for up to 4 additions. In the conc line the concentration given in the specified concentration unit of each analyte in the standard solution must be set. In case more than 1 analyte is measured, a standard solution may contain 1 or more analytes if they have the same concentration. If an analyte is added from another solution or has another concentration than the previous, a new Solution number (Sol nr.) has to be specified.

So if for example only 1 analyte is determined, the Sol nr. can be 1. In case 4 analytes are determined and four different standard solutions are used, the Sol. nr values can be different so for instance 1, 2, 3 and 4. If a standard solution contains more than 1 analyte, the next row is automatically set to the first value.

10.2.2 Specifying concentrations in cell

If concentration in cell was chosen in the Analytes section, the absolute concentration of analyte in the cell according to the specified unit needs to be put into the table assuming that the original sample contained no analyte.

The cell volume is used to correct peak heights or peak areas for dilution effects due to addition of standard solutions. So the cell volume is the total volume in which the measurement is performed, which can for instance be the sample + separate base electrolyte. If the cell is filled with sample solution only, then the cell volume and sample volume are equal. This dilution factor due to the standard additions is shown in the 'Analytical result' window.

10.2.3 Peak settings

Three different peak features can be used for determination:

- Peak height
- Peak area
- And the sum of the peak's slopes

Peak Settings

Use value: ☒ height ☐ area ☐ sum of slopes

	E peak	Auto	E left	E right
► Cu	-0.3	<input checked="" type="checkbox"/>	9999	9999
Pb	-0.6	<input checked="" type="checkbox"/>	9999	9999
Cd	-0.8	<input checked="" type="checkbox"/>	9999	9999
		<input type="checkbox"/>		

Auto peak search window around E peak (V):

For minimum peak width and height: see tab 'Peaks'

Figure 257 Peak locations for each analyte with options for manual baselines

The column 'E peak' gives the expected peak potential. If a peak is found in the peak window E peak +/- Peak window, it will be assigned to the given Id.

If a manual baseline should be used for the specific analyte, E left and E right can be specified by entering values in the corresponding cells. As soon as the values are entered the peaks will be updated in the plot window. It is also possible to draw a baseline manually using the manual peak options in the Plot toolbar.

See also section: [Marking peaks manually](#) on page 184.

The 'Auto peak search window' determines the maximum width of the peak. The peak for each analyte will be searched in E peak +/- Peak window.

If a peak is not assigned, the peak potential is wrong, or the peak window is too narrow; change the peak parameters in the Peaks tab on the left and click the Mark peaks button to apply the new peak parameters.



Figure 258 Mark peaks button

The sample volume is needed in order to calculate the concentration of the analyte(s) in the sample. The program first calculates the concentration of the analyte(s) in the cell. These values are shown in the tables of Analytical result under cCell. The value of the concentration in the sample is simply found by multiplying cCell by the factor 'Cell volume / Sample volume'. Thus, the sample volume cannot exceed the cell volume.


Please note that in case the addition of the standard solutions have a negligible effect on the cell volume, the final cell volume is not relevant.

All results

	cSample/ ppm	cCell/ ppm	Corr R ²	Slope
▶ Cu	193.98	19.398	0.9893	0.0623
Pb	1146.2	114.62	0.9969	0.2140
Cd	1189.6	118.96	0.9993	0.1050

Figure 259 Table showing results of each analyte

10.3 Running a measurement

Before starting a new determination, click the New button  in the toolbar, or select 'New analysis' in the Analysis menu.




Figure 260 Analysis measurement toolbar

For a determination, a 'Sample' and at least one 'Standard' curve is required. The type of curve measured after pressing the green 'Start' button is shown as selected item in the list next to this button. If blank correction is required, the 'Blank' solution also has to be measured or the blank curve has to be loaded from file. After a blank curve has been measured, it can be saved as blank file by using 'Save blank' in the Analysis menu.

The 'Plot' tab shows the voltammograms and the 'Analytical results' tab gives the calibration or standard additions plots. Please note that the correlation coefficient for linear least squares lines is given. Also, the slope of the calibration or standard addition curves are shown. This value of the slope is an important value for the sensitivity of the electrode or sensor.




Figure 261 Analytical Result tab

10.3.1 Standard addition measurement

A determination by means of standard addition requires at least two curves: sample and one standard addition. As soon as both curves are available, the analytical results can be obtained in the window 'Analytical result'. It is advised to use at least two additions in order to check the linearity of the curve. The standard additions have to be done before a standard curve is measured.

10.3.2 Measurement by using a calibration curve

The calibration curve is measured or loaded before a sample is measured. If only one standard is measured, the calibration curve is calculated with $c=0$ and $i_{peak}=0$ as second point.

After the sample has been measured, the calibration plot will show lines from which the concentration is obtained.

10.3.3 Using Blank subtraction

See section [Blank subtraction](#) on page 45.

10.4 Result analysis

After measurements of the curves, parameters in the Analysis window can be changed if required. New calculations of the analytical results are shown after pressing the button

'Recalculate'. Secondly the peak labels are updated with their corresponding analyte *after* pressing the 'Recalculate' button.




Figure 262 The 'Recalculate' button

It is possible correct for a blank curve and also to undo blank correction afterwards by using the 'Use blank' box.

Results of measurements can be stored by using the 'Analysis' menu. Results can be printed with the 'Print' button in the 'Analysis Result' window or copied to another program (e.g. Excel) by using the button 'All to clipboard'.

If a point in the calibration or standard addition curve has to be neglected, the corresponding checkbox in the 'Results for each analyte' window can be used.

Cu					
	Use	cCell/ppm	E _{Peak} /V	i _{Peak} /μA	Total Added Volume μl
sample	<input checked="" type="checkbox"/>	19.398	-0.3300	1.3948	0
st. 1	<input type="checkbox"/>	OK	-0.3000	3.9484	0
st. 2	<input checked="" type="checkbox"/>	100.00	-0.2800	7.6224	0

Figure 263 Uncheck to remove this value from the linear plot

10.4.1 Standard Addition

The Standard Addition Method is performed in the sample's matrix, thus the influences of other species on the sensitivity are already taken into account. A known volume of the sample and a known volume of the supporting electrolyte are added to the cell. The measurement is performed. A fixed volume of standard solution with a known concentration is added. The measurement is repeated. These two steps are repeated until a sufficient amount of standard additions have been measured. From these values a linear regression is calculated. The zero point of this function is the analyte's concentration in the cell.

To create the points for the linear regression, the concentration of the artificially added analyte in the cell c_{Cell} must be known after each addition of standard. For the measurement of the sample no standard was added, so the value is 0. For the other measurements c_{Cell} is calculated as follows:

$$c_{Cell} = \frac{\text{Total Added Volume in } \mu\text{L} * \text{Conc}}{\text{Cell Volume} + \text{Total Added Volume in } \mu\text{L}}$$

Conc is the concentration of the standard solution.

After the I_{Peak} is plotted versus c_{Cell} the slope of the linear regression is calculated with the general equation:

$$Slope = \frac{\sum(x - \bar{x})(y - \bar{y})}{\sum(x - \bar{x})^2}$$

The IPeak values are y and \bar{y} is the average of all known y values. The cCell values are x and \bar{x} is the average of all known x. The intersection of the linear regression is:

$$Intersection = \bar{y} - Slope * \bar{x}$$

With these two parameters you can calculate the null of the curve. The concentration of cCellSample is the null * (-1). This means:

$$cCellSample = \frac{Intersection}{Slope}$$

For the concentration of the Analyte in the Sample cSample, the dilution by the supporting electrolyte needs to be considered:

$$cSample = \frac{cCellSample * Cell\ volume}{Sample\ volume}$$

10.5 Example data files

The program comes with two example files. These files are stored in the default PS Data folder, which is typically:

"My Documents\PSData\Analytical mode examples".

You can load these in PSTrace using the menu: 'Data → 'Load data file...'

The file "Standard Addition Example.psession" shows a determination of some heavy metal ions in a nickel galvanic bath. This determination is based on standard. The measurement was done by adding 1 ml sample to a base electrolyte with a final cell volume of 10 ml and using 3 different standard solutions.

The file "Calibration Example.psession" shows a determination by using a calibration curve. In this file the table 'Analysis' contains the concentration values of the analytes in the cell.

11 Corrosion mode



The Corrosion mode in PSTrace allows the user to run specific corrosion measurement techniques and use analytical tools for corrosion research on the measurement data.

11.1 Supported techniques

The available electrochemical measurement techniques in the Corrosion mode are:

- Potentiostatic (apply constant potential, equivalent to Chronoamperometry)
- Galvanostatic (apply constant current, equivalent to Chronopotentiometry)
- Linear Polarization (potential sweep, equivalent to Linear Sweep Voltammetry)
- Cyclic Polarization (bi-directional potential sweep, equivalent to Cyclic Voltammetry)
- Corrosion Potential (open circuit potential, equivalent to Open Circuit Potentiometry)
- Impedance Spectroscopy

These techniques are identical to those found in the Scientific Mode, employing the same signal and adhering to the same specifications. The only distinction lies in the display and outcome, as the data results are directed to calculation tools for corrosion determination.




Figure 264 List showing available techniques



The other standard available electrochemical techniques that do not apply to corrosion research are not listed when PSTrace is set in the Corrosion mode.

11.2 Running a corrosion measurement

The sequence of a corrosion measurement is:

1. If $t_{\text{Cond1}} > 0$ s then E_{Cond1} is applied for t_{Cond1} s.
2. If $t_{\text{Cond2}} > 0$ s then E_{Cond2} is applied for t_{Cond2} s.
3. The cell is switched off if the measurement of E_{oc} or OCP is required.
4. Now the value of E_{oc} is continuously measured until either 'Max. OCP time' is reached OR until the stability criterion is met OR the button 'Force continue' is pressed.
5. Now the cell is switched on at the value of E or E_{begin} for ' t_{eq} ' s, after which the actual measurement starts.




Figure 265 Waiting for the OCP to meet the stability criterion.

11.2.1 Measure versus OCP

Corrosion measurements can be done by specifying the potential scan with respect to the Open Circuit Potential OCP or with absolute values versus the reference electrode.

In case one or more potentials are specified with respect to the OCP, the open circuit potential must be determined before the actual measurement is done. This OCP measurement requires a variable time, which is determined by the drift of the open circuit potential and the maximum time to measure the OCP value. The OCP value is set as soon as the drift is lower than the specified value for the 'Stability criterion' or when the ' $t_{\text{Max. OCP}}$ ' has elapsed.

☒ **Measure vs OCP**
☒ **E begin versus OCP**
☒ **E end versus OCP**

t Max. OCP s
Stability criterion mV/s

Figure 266 OCP parameters

It is possible to condition the electrode or the corrosion sample before this OCP is obtained. For this purpose two potentials can be applied: 'E cond1', during a period of 't cond1' and 'E cond2' during a period of 't cond2'. In case the value of 't cond' is set to 0 s, the corresponding conditioning potential is not applied. So, if both values for 't cond' are 0, the electrode is not conditioned at all.

11.3 Result analysis

The 'Corrosion' tab in the plot window is used for analysis of linear polarization and electrochemical impedance spectroscopy measurements. The supported measurements are shown in the legend on the right and the material settings of the selected measurement are shown on the left.




Figure 267 The Corrosion analysis tab

11.3.1 Setting up material parameters

The material parameters required for analysis can be entered before the measurement in the Method Editor or after the measurement in the 'Corrosion' tab.

Measurement Material	
Surface area (cm ²):	0.001
Equiv. Weight (g/mol):	0.0
Density (g/cm ³):	0.1
β anodic (V/dec):	0.1
β cathodic (V/dec):	0.0
Data smoothing (Savitzky-Golay):	0: spike rejection

Figure 268 Tab with Material properties in the Method Editor

The values set in the 'Material' tab of the Method Editor are copied to the tab for Corrosion analysis when a measurement is started.

In the 'Corrosion' tab, the material parameter input fields are found beneath the plot.




Figure 269 Material parameters used for calculations in the Corrosion tab

The material parameters are required to determine the corrosion rate. Linear polarization measurements require surface area, equivalent weight and density and impedance measurements also require the anodic and cathodic Tafel slopes, β_{anodic} and $\beta_{cathodic}$.

- Surface area is the area of the sample in cm².
- Equivalent weight is the equivalent mass of one mole of the sample material in g/mol.

- Density is the density of the sample material in g/cm³.
- β_{anodic} and $\beta_{cathodic}$ are the Tafel slope parameters for the sample material, these can be determined from a linear polarization measurement or from literature.

11.3.2 Linear polarization

Linear polarization is typically used to study the corrosion response of metallic coatings. The following analysis techniques are supported for the estimation of the corrosion rate based on linear polarization measurements:

- Auto Butler-Volmer fit: Fitting the Butler-Volmer model over an automatically detected range.
- Manual Butler-Volmer fit: Fitting the Butler-Volmer model over a manually selected range.
- Manual Tafel slope fit: Fitting Tafel slopes in the linear regions of the anodic and cathodic slopes.

Note: To achieve an accurate estimation of the corrosion rate it is recommended to use a measurement with at least one linear Tafel slope that ranges over one decade in current density. Additionally, the distance between the Tafel slope and the corrosion potential should at least be 50 mV.

The results of these analysis techniques are presented in the Tafel fit results table.


Tafel fit results		
Butler-Volmer		
Manual	Auto	
E corrosion, V	-0.319	-0.319
I corrosion, μA	13.17	13.99
I corrosion per surface area, A/cm ²	2.635E-5	2.797E-5
Polarization Resistance, Ω	7926	7891
Anodic β Tafel constant, V/decade	0.347	0.358
Cathodic β Tafel constant, V/decade	0.784	0.879
Corrosion rate, mm/year	0.129	0.137
Tafel slopes		
Manual		
E corrosion, V	-0.312	
I corrosion, μA	5.874	
I corrosion per surface area, A/cm ²	1.175E-5	
Polarization Resistance, Ω	1.006E+4	
Anodic β Tafel constant, V/decade	0.212	
Cathodic β Tafel constant, V/decade	0.380	
Corrosion rate, mm/year	0.058	
Measured Eoc value:	0.000 V	
		 Copy to clipboard

Figure 270 Table with Tafel fit results

- The corrosion potential (E corrosion) is the potential at which the anodic and cathodic reaction rates are equal. The measured current approaches zero at the corrosion potential, because all the electrons released by the dissolving of the metal are consumed by reduction reactions. For the Butler-Volmer techniques the corrosion potential is determined as the potential where the log the measured current is the smallest. In the Tafel slope method the corrosion potential is the potential where the anodic and cathodic Tafel slopes intersect.

- The corrosion current ($I_{corrosion}$) is a measure of the rate of corrosion, measuring it directly is not possible. The corrosion current can be estimated the current at which the Tafel slopes intersect or by fitting the Butler-Volmer equation on a linear polarization measurement.
- $$I = I_{corrosion} \left(e^{\frac{2.303(E - E_{corrosion})}{\beta_{anodic}}} - e^{\frac{-2.303(E - E_{corrosion})}{\beta_{cathodic}}} \right)$$
- β_{anodic} and $\beta_{cathodic}$ are the Tafel slopes, these represent the change in Volts per decade of current in the Tafel plot.
- When plotting current over potential a linear slope is equals resistance (i.e. $R=U/I$). The slope close to the corrosion potential is approximately linear, this is referred to as the polarization resistance. The polarization resistance is inversely proportional to the corrosion current, assuming that the Tafel slopes are constants (Stern-Geary equation).
- $$I_{corrosion} = \frac{\beta_{anodic} \cdot \beta_{cathodic}}{R_{polarization}(\beta_{anodic} + \beta_{cathodic})}$$
- The corrosion rate in mm/year can be calculated according to the standard practice described in the ASTM Standard G 102. To calculate this an estimation of corrosion current is needed as well as the following material parameters: equivalent weight (EW) in g/mol, the density (d) in grams/cm³, and the sample area (A) in cm² of the study sample. Combined with a constant (K) defined by the ASTM (3272 mm/(amp*cm*year)) this information is used to determine the corrosion rate in mm/year.
- $$Corrosion\ Rate = \frac{I_{corrosion} \cdot K \cdot EW}{d \cdot A}$$

11.3.3 Selecting a curve for corrosion analysis

On the right-hand side of the screen is a legend which contains all compatible corrosion measurements. To perform a Butler-Volmer or Tafel slope fit select a Linear Polarization measurement from the legend. The checkbox in front of the measurements indicate whether they are also visible in the plot window.




Figure 271 The legend to select measurement data for performing a model fit

11.3.4 Performing an automatic Butler-Volmer fit

The auto Butler-Volmer fit is automatically applied after selecting a Linear Polarization curve from the legend and after changing the value of one of the material parameters. Alternatively, clicking on the Auto Butler-Volmer fit button in the toolbar above the plot also applies the fit.




Figure 272 Toolbar in the 'Corrosion' tab with buttons for data fitting

11.3.5 Performing a manual Butler-Volmer fit

By selecting the Manual Butler-Volmer fit from the toolbar above the plot the range for the Butler-Volmer can be set manually.

Usage: click two points on the curve to mark start and end point




Figure 273 Marking two points for a logarithmic fit

11.3.6 Performing a Tafel slope fit

To perform a Tafel slope fit click on the Manual Tafel slope fit button in the toolbar above the plot. The Tafel slopes (Evan's diagram) are drawn in the plot, they are either fit automatically or loaded from the previous fit. The ranges of the slopes can be adjusted in two ways. The range can be specified by clicking twice on a linear section of either the anodic or cathodic part of the corrosion measurement.

Usage: mark two points for the anodic slope line and/or two points for the cathodic slope line. After selecting two points a green line is drawn representing the Tafel slope.




Figure 274 Setting the anodic and cathodic slope lines for Tafel plot analysis.

Alternatively, the ranges of the slopes can be adjusted using the Tafel slope settings. Either by adjusting the ranges of the slopes using sliders or by entering the potential ranges in the text boxes.




Figure 275 Adjustments of the Tafel slope settings

To confirm the fit of the Tafel slopes click on Confirm in the top left corner of the plot.

11.3.7 Correcting for the potential of the reference electrode

In the corrosion mode it is possible to apply a reference electrode correction to the measured potential. Click on the downward arrow of the Show E corrected button in the toolbar above the plot and select Determine E offsets.




Figure 276 From this menu you can select a correction for the potential of the reference electrode.

11.4 Impedance Spectroscopy

Electrochemical Impedance Spectroscopy (EIS) can be used to study corrosion and the effects of a wide range of coatings. For example, anodized coatings (anodized aluminium), conversion coatings (Chromate conversion coating), or organic coatings (paint). The corrosion rate and the pitting/disbanding of coatings are studied by fitting equivalent circuit models on the EIS measurement.





Figure 277 Fitting an equivalent circuit on measured corrosion data

11.4.1 Selecting a curve for corrosion analysis

On the right-hand side of the screen is a legend which contains all compatible corrosion measurements. To perform an equivalent circuit fit select a EIS measurement from the legend. The checkbox in front of the other EIS measurements indicate whether they are also visible in the plot window.



Toggle visibility in the plot window

11.4.2 Selecting an equivalent circuit

To change the selected equivalent circuit either click on the 'Load corrosion model presets' button above the plot window to open the circuit library or click on the downward arrow of the button to quickly select a circuit.




Figure 278 Corrosion circuit library.




Figure 279 Corrosion circuits quick access.

Loading a circuit from the circuit library or the quick select menu will update the circuit model settings. A graphical representation of the selected equivalent circuit is displayed, below the image is a table with circuits corresponding parameters. When selecting one of the four corrosion rate circuits the corrosion parameters are added to the table: polarization resistance

(R polarization in Ohm (Ω)), corrosion current ($I_{\text{corrosion}}$ in Ampere (A)) and corrosion rate in mm/year.

Note: for the corrosion parameters to be calculated all material parameter values must be greater than 0.

11.4.3 Fitting the equivalent circuit

Before fitting a circuit, the frequency range on which to perform the fit can be selected. By default, a fit is performed over the entire measurement. To select a frequency enter the desired range in frequency range text boxes below the plot.

Fit from frequency (Hz): to

Figure 280 The frequency range to fit EIS data on.

Alternatively, it is possible to specify the frequency range in the plot by selecting the 'Select fit range' mouse pointer button from the toolbar left to the plot and clicking in the plot.




Figure 281 Selecting a frequency range for data fitting

Selecting the fit frequency range. Double-click on a point in the plot. Then select the point to where the fit should be performed.

To fit the model on the measurement, click on the 'Fit model' button. The result of the fit is stored in your measurement, but to store the fit result it is necessary to save your measurement. Simple circuits with few components will often fit directly, but more complex circuits will require fine tuning of the circuit parameters before a good fit is obtained. To fine tune the model, you can either use your prior knowledge of the cell or you can estimate certain parameter values from the Bode and Nyquist plots. For more information on circuit fitting please refer to the help section on [circuit fitting](#) on page 165.


11.5 Exporting results

There are three ways to export the results from your Tafel or circuit fit:


Butler-Volmer	Manual	Auto
E corrosion, V	-0.319	-0.319
I corrosion, μA	16.23	11.35
I corrosion per surface area, A/cm^2	3.246E-5	2.270E-5
Polarization Resistance, Ω	7902	8102
Anodic β Tafel constant, V/decade	0.403	0.311
Cathodic β Tafel constant, V/decade	1.103	0.665
Corrosion rate, mm/year	0.159	0.111

Tafel slopes	Manual	
E corrosion, V	-0.317	
I corrosion, μA	6.627	
I corrosion per surface area, A/cm^2	1.325E-5	
Polarization Resistance, Ω	9854	
Anodic β Tafel constant, V/decade	0.236	
Cathodic β Tafel constant, V/decade	0.416	
Corrosion rate, mm/year	0.065	

Measured Eoc value:

 Copy to clipboard

Advanced fit options

 Fit model

Element	Fitted Value	Min Value	Max Value	Unit	Error%
Rs	965.9	1.00E-6	1.00E+12	Ω	7.645
CPEdl	0.060	1.00E-6	1.00E+3	μT	4.940
CPEdl	0.708	0	1.00	ϕ	0.696
Rct	0.304	1.00E-12	1.00E+6	$\text{M}\Omega$	13.13
WD	4.983	1.00E-12	1.00E+6	$\text{M}\sigma$	1.893
WD	0.254	1.00E-12	1.00E+6	\sqrt{s}	3.751
WD	0.992	0	1.00	φ	4.573
Chi-Squared:	0.0002	Iterations:	27		
R polarization:	5.287E+6	I corrosion:	1.709E-8		
Corrosion rate:	4.149E-8				

Figure 282 Left: Copying Tafel fit results to clipboard. Right: Copying circuit fit results to clipboard.

1. Exporting a Tafel plot or a circuit fit curve to Excel (Bode or Nyquist) by pressing the export to Excel button on the left-hand side of the plot.
2. Printing/saving a report as a .rtf file by clicking on the print report button or selecting the save report button in the toolbar above the plot.




Figure 283 Left: Tafel fit report. Right: Circuit fit report.

11.6 Example data files

The program comes with an example file stored in the default PS Data folder:
"My Documents\PSData\Corrosion mode examples\Corrosion.psession".

You can load this file in PSTrace using the menu 'Data → 'Load data file...'

12 Files

PSTrace uses a number of different file formats. This chapter describes which file types are used by PSTrace.

12.1 File types

PSTrace creates a number of different files in ASCII format.

The files are default stored in the folder:

My Documents\PSData

The following files are currently supported by PSTrace and Multitrace

Table 18 Supported file types in PSTrace and Multitrace

Method	.psmethod
Raw measurement data, curves and methods	.psession
Script	.psscript
MultiTrace project file	.psproject

12.1.1 Obsolete file formats

The following list of file formats has become obsolete. They can still be loaded in latest versions of PSTrace, but no longer be saved in the particular format.

Table 19 Obsolete file types

Method (old style)	.pms (scans) or pmt (vs time)
Curve	.pss (scans) or .pst (vs time)
EISdata	.pseisdata
Multiple curves	.mux
Analysis curves	.pds
Info	.rtf

12.1.2 Method files

Methods can be saved separately or are automatically saved with a curve.

12.1.3 Curve files (obsolete)

If a curve (.pss, .pst) is saved or multiple curves (.mux, .pds) are saved, a method file with the same name and corresponding extension is always saved. This method file also contains the last plot view window used (zoom).

The curve file contains measured data of potential and current or time and current.

12.1.4 EIS data files (obsolete)

Electrochemical Impedance Spectroscopy data is saved in an .pseisdata file. These files are plain text and contain values for frequency, current range, Z, phase, Idc and Iac.

12.1.5 Session files

Session files can contain all the data available at any time in a PSTrace session.

All available measurements, corresponding curves and method, including the active method in the method editor can be stored to a single .psession file format.




Figure 284 A PSession (.psession) file contains basically all the data available at a given moment.

12.1.6 PSScript files

The .PSScript files are plain ASCII files containing a script that can be loaded and run in the Script window of PSTrace. (See section [Scripting](#) on page 243)

12.2 Loading and saving data

Methods and measurements can be loaded and stored from and to a hard-drive or other storage solution.

12.2.1 Loading methods

There are three ways to load a method (.psmethod) file:

1. Selecting the load method file button from the toolbar;



2. Selecting load method from the method menu;




Figure 285 The 'Method' menu

3. Selecting a recently loaded method from the recent methods list;




Figure 286 Recent methods in the 'Method' menu

4. Drag and drop a method file to the PSTrace window from the Windows File Explorer.

12.2.2 Saving methods

Methods can be saved as psmethod files by:

1. Selecting the save method file button from the toolbar;



Figure 287 Button for saving methods

2. Selecting save method from the method menu;




Figure 288 Menu item for saving methods

12.2.3 Loading measurements

Measurements and sessions can be loaded or added to your current session, i.e. loading a pssession/pss/pst/pseisdata/mux/pds file will replace all other data in your current session.

There are different ways to load measurements:

1. Selecting the load data button from the toolbar.



Figure 289 Button for loading methods.

2. Selecting load data file from the 'Data' menu.




Figure 290 'Load data file' menu item in the 'Data' menu.

3. Selecting a recently loaded data file from the recent data files list.




Figure 291 The recent data files list in the 'Data' menu

The 'Load data' dialog shows a 'Preview' tab and a 'Details' tab. The 'Preview' tab contains some additional settings that influence the assignment of colors for the loaded curves. The available settings are:

- 'Default': use the setting stored with each curve.
- 'Do not auto assign': ignores automatically assigning of colors and assigns last color stored with the curve.
- 'Auto assign all': ignores the colors stored with the curves and auto assign new colors for all curve.




Figure 292 The 'Load data' dialog with curve color settings in the 'Preview' tab.

Data files can be added to your current session by:

1. Selecting the add curve button from the toolbar.



Figure 293 Button for adding existing data to the current session

2. Selecting load data file from the data menu;




Figure 294 Menu item for adding existing data to the current session

3. Dragging and dropping data files to the PSTrace window from the Windows File Explorer.

12.2.4 Saving measurements

The measurements in your current session can be saved in a psession file. In the save dialog you can specify the measurements and curves that should be stored in the psession file. The save dialog can be opened by:

1. Pressing CTRL+s on the keyboard;
2. Selecting the save data button from the toolbar;



Figure 295 Button for saving data

3. Selecting save data from the data menu;




Figure 296 Menu item for saving data

On the right-hand side of the save dialog you can (un)check the data that you would like to save, only checked measurements/curves are saved.




Figure 297 Save data window with list to select the data to include

12.3 Exporting data to other file formats

Measured data can be exported to different file formats. The export options are available in the 'Curve' menu.




Figure 298 Curve menu showing export options

12.3.1 CSV

The CSV (Comma Separated Values) is a popular format in plain ASCII, supported by many applications like Excel, OpenOffice Calc and Origin. To change the values exported for EIS data, open the 'Settings' window (menu: 'Tools' → 'General settings...') and click the 'Plot and data' tab. Then click the button 'Set columns for exporting EIS data'.

12.3.2 Origin

Origin from OriginLab is scientific graphing and data analysis software widely used at universities.

See for more information: <http://www.originlab.com/>

To change the values exported for EIS data, open the 'Settings' window (menu: 'Tools' → 'General settings...') and click the 'Plot and data' tab. Then click the button 'Set columns for exporting EIS data'.

Origin does not need to be installed on the PC to export curves to the Origin format.

13 Scripting

PSTrace provides a straightforward and user-friendly interface for executing individual measurements. The Script menu extends this functionality, allowing users to perform more complex workflows and include additional commands. One of the versatile applications of the Script feature is the ability to execute a sequence of diverse methods, establish loops, introduce wait periods, and define specific parameters.

An example script file is available in the default data folder, typically located at 'My Documents\PSData.'

In the default data folder (normally 'My Documents\PSData') a script file is available as example.

13.1 Features




Figure 299 Scripting main window

13.1.1 File menu




Figure 300 The 'File' menu

The table below explains each menu item.

Table 20 The menu items in the 'File' menu.

Menu item	Action
New	Clears the script present in the editor.
Load script	Loads a saved script.
Save script	<p>Saves a script to a specific location. All the file references to method or curve files in the script will be made relative to the location of the script file. For example, if the script file is in location: C:\directory A\script.psscript</p> <p>and the script refers to files in location C:\directory A\directory B\</p> <p>the file references in the script file will be relative: directory B\file1.psmethod directory B\file2.psmethod and not to the absolute path: C:\directory A\directory B\file1.psmethod C:\directory A\directory B\file2.psmethod</p> <p>So, if the directory A including script and sub-directory is copied to another location or PC, the references will still be intact.</p>
Export all to single folder	<p>Copies all files that are referred to and the actual script file to the folder chosen to save the new script file into. So, if a script file is at location: C:\directory A\script.psscript</p> <p>and file references in the script are for example to: C:\directory B\file1.psmethod C:\directory C\file2.psmethod</p> <p>All files will be copied to a single folder and all file references will be changed to this folder. This allows you to save all relevant files to a single location for back-up purposes.</p>

Export all and pack to ZIP

This function does the same as ‘Export all to single folder’ as described above, but instead of targeting a folder all files are saved to a single ZIP file.

13.1.2 Composing a script

The script contains of a sequence of commands which is run from the top down. Add commands using the list of available commands shown at the left of the script.

The order of the commands in the script can easily be changed by dragging and dropping them. Click and hold on the text of a command and drag the command to a location where you want. The black arrow shows where the command will be dropped.

In the picture below the ‘Wait’ command is dragged above the ‘Measure’ command:




Figure 301 Inserting a ‘Wait’ command.

From PSTrace version 5.10, the Scripting validates the compatibility of the current command with the connected device. If the command is not compatible, it will be highlighted in red, and an accompanying error message, also displayed in red, will provide an explanation for the encountered incompatibility.

13.1.3 Changing parameters

The parameters for each command can be changed in the panel at the right side of the script. This panel shows the settings for the selected command. Click the ‘?’ button to get help for the selected command.

13.1.4 Using the output file

The output file can be specified by typing in the textbox or using the ‘Browse’ button.




Figure 302 Setting the script output file.

If one of the following commands is used in the script and the output checkbox is checked, the corresponding results are appended to the text file provided.

- ReadCurrent
- ReadPotential

- FindPeaks
- ReadAnalog
- ReadDigitalIO

The file can be read using any text editor like Notepad, or they can be loaded in Excel. The columns are separated using tabs.

13.2 Measure command

This command loads a Method file and runs it, like a standard run in the main PSTrace window. This command is designed to operate as a standalone instruction, eliminating the need to incorporate additional commands for a complete method run. Other commands such as 'Cell', 'Read current', 'Set potential', etc., are intended for users who need to create alternative workflows or customize the measurement sequence. However, if the objective is to run the measurement as a typical operation in the main PSTrace, there is no requirement to include additional commands.

Measurement ?

☒ Clear plot before start measurement

Load method from specified file:

C:\Users\Youname\PSData\PSDiffPulse ...

Create new method... Edit...

☒ Subtract Blank

☒ Load existing Blank curve:

C:\Users\Youname\PSData\PSDiffPulse. ...

Create new method... Edit...

☐ Measurement using method:

C:\Users\Niels\Documents ...

Create new method... Edit...

☐ Measurement Blank first

☒ Save result

Prefix: dpvData

Format: PSTrace format

Save to folder: C:\Users\Youname\Documents\PSDa ...

Figure 303 Panel with options for the 'Measure' command.

- Clear plot before start measurement: If checked; the plot and legend will be cleared before the measurement is started. In case a Blank is measured first or afterwards, both the Blank and the actual measurement curve will be displayed in the plot.
- Load method from specified file: The specified method will be loaded for this measurement. A new method can be created and saved instantly using the 'Create new method...' button. The 'Edit...' button allows to make changes in the method file and save (overwrite) these changes to the specified file.



If the same method file is assigned in multiple scripts, any alterations made will be applied since the method file is loaded at the initiation of each measurement. If you need to modify a single parameter from a method, please refer to 'Override' command.

- Subtract Blank: If subtract blank is checked a Blank curve will be subtracted from the measured curve. This curve can be an existing curve ('Load existing Blank curve') or can be measured ('Measure using method').
- Measure Blank first: When this option is checked, the method designed for the Blank will be loaded initially and employed for measurement. Subsequently, the actual measurement will take place. Upon completion, the system will automatically subtract the Blank curve from the measured data.
- Save result: When checked, the measured curve (or curves in case a Blank was used) will be saved using the prefix, followed by a number. This number automatically increments if multiple curves are saved with the same prefix. If curves are not saved from within the Script, all curves in the PSTrace plot can be saved manually to a single session file (.psession) using the 'Data' menu as well.
- Cell on after measurement: In case a method is used with the setting 'Cell on after measurement' together with 'cell off after n seconds', the latter setting is ignored in the script so the cell can be turned off after a specific period in the script.

To keep the cell on after a measurement for a period, you should use the 'Wait' command in combination with the 'Cell' command.

13.3 FastMode command

When a potentiostat is idle, it normally sends an idle status package with voltage, current and auxiliary readings every second. If an instrument receives a command while measuring the values for the idle status package it will not be able to process the received command immediately. Setting 'FastMode' on disables the idle packages being sent every second. This way the response time of the instrument is always optimal.

13.4 Cell command



Figure 304 Panel with options for the 'Cell' command.

Turns the cell on or off. By default, the potentiostat is initialized with the cell turned off. In this state, only the Working Electrode (WE) or Sense (S) and Reference Electrode (RE) leads are connected to the potentiostat, enabling the reading of its Open Circuit Potential (OCP). Upon turning the cell on, the Counter Electrode (CE) lead is connected, resulting in polarization of the cell.

It's important to note that this command is not necessary for preceding or succeeding 'Measurement' commands, as it automatically manages the cell state as needed, turning it on or off accordingly.

Example:




Figure 305 Example showing the use of the 'Cell' command.

13.5 SetCurrent command




Figure 306 Panel with options for the 'SetCurrent' command.

The 'SetCurrent' sets the connected device in galvanostatic mode and applies a constant current at the specified current range. It is important to note that this command is not supported by certain models without galvanostatic capabilities (i.e., EmStat1/2/3 series, EmStat Pico, and Sensit series).

13.6 ReadCurrent command




Figure 307 Panel with options for the 'ReadCurrent' command.

The 'ReadCurrent' command should be used when the cell is on. It reads the current and stores the value in the text file specified at the bottom of the script window.

13.6.1 Auto ranging

If the 'Use auto ranging' checkbox is used, the current range will be evaluated at each interval time and changed if necessary. If the correct current range has been found or the minimum or maximum range is set, the current will be read.

The first current range to evaluate is the highest current range.

If 'Use auto ranging' is not used the 'ReadCurrent' command will leave the current range unchanged. The current range can be changed in the script using the 'Current range' command. See next section.

13.6.2 Output

If there is no column in the text file yet, it will add the line with columns first and then the line with values. The same columns are used for potentials.

The values are all separated by a TAB and can be imported in Excel.

Example output:

```
Script output 12/07/2012 - 17:42:43:
Time                Potential in V                Current in uA
12/07/2012 17:42:45      -3.313E-5
12/07/2012 17:42:46      -7.501E-6
12/07/2012 17:42:47       5.250E-5
12/07/2012 17:42:48      -4.250E-5
```

13.7 SetPotential command




Figure 308 Panel with options for the 'SetPotential' command.

The 'SetPotential' command sets the connected device in potentiostatic mode and sets the potential to the given value. The cell state (on or off) is not changed.

13.8 ReadPotential command




Figure 309 Panel for the 'ReadPotential' command.

The 'ReadPotential' command outputs the potential to the output file. Depending on the potentiostat model, certain potential ranges or other options may become available. These options are equivalent to those encountered in the OCP (Open Circuit Potentiometry) technique.

13.8.1 Output

If there is no column in the text file yet, it will add the line with columns first and then the line with values. The same columns are used for currents.

The values are all separated by a TAB and can be imported in Excel.

Example output:

```
Script output 12/07/2012 - 17:45:59:
Time                Potential in V          Current in uA
12/07/2012 17:46:01    -1.501
12/07/2012 17:46:02    -1.401
12/07/2012 17:46:03    -1.299
12/07/2012 17:46:04    -1.200
12/07/2012 17:46:05    -1.098
```


13.9 Wait command

Figure 310 Panel with options for the 'Wait' command.

The 'Wait' command places the script in an idle state for a specified period or until a specified event occurs.

- **Wait for:** The script continues after the specified number of seconds.
- **Resume today at:** The script continues at the specified time of day (this can go on for multiple days)
- **Wait until user resumes:** The script shows the message in the script result box (see Features) and waits until the user clicks the active Pause button to resume.
- **Wait until digital line:** The script waits until the IO input D0 on the auxiliary port is set HIGH or LOW. This feature can be used to let an external device determine if the script can continue. When waiting for the digital line to be HIGH or LOW, the pin status is polled about every 20 ms.

In case a quicker response is required using MethodSCRIPT™ might be preferable.

Please refer to the instrument-specific sections in this manual for more information about the positions of the digital and analog I/O pins on your instrument's auxiliary port.

13.10 Repeat command

Figure 311 Panel with options for the 'Repeat' command.

The 'Repeat' command allows for the repetition of a set of commands. To repeat specific commands, attach them to the 'Repeat' instruction in the list. This can be done by clicking and dragging the desired commands into the 'Repeat' command or by adding a new command with the 'Repeat' option selected to include it in the repetition. The commands earmarked for repetition must be presented as a branch within the 'Repeat' command.

Example:




Figure 312 An example script involving a nested 'Repeat' command.



Since PSTrace 5.10 a Repeat branch can also contain nested Repeat branches, for creating a 'loop-in-a-loop'.

13.11 FindPeaks command

Figure 313 Panel with options for the 'FindPeaks' command.

The 'FindPeaks' command performs a peak search over the last measured curve available. This instruction is therefore logically preceded by the Measure command.

The peaks found are written to the output file specified.

- All peaks: All peaks are searched using the parameters as specified in the method editor. This command always uses the default algorithm for peak search and may give different results compared to using the auto-detect peaks button in the main window.
- Use window: The number of peaks must be specified and for each peak the left and right of the search window must be provided.
- Manual baseline: The provided left and right values for each peak are used as left and right of the peak baseline, therefore forcing the finding of a peak.

13.11.1 Output

If there are no columns set yet in the output file, a line with columns will be added first. All values are separated by a TAB and can be exported to Excel.

13.12 SetChannel command

Figure 314 Panel with options for the 'SetChannel' command.

The 'SetChannel' command sets the channel of a CH8, MUX or MUX8-R2 multiplexer.

13.13 NextChannel command

The 'NextChannel' command sets the connected CH8, MUX or MUX8-R2 multiplexer to the next channel in sequence. If this command is executed when the multiplexer is already at the last available channel, it will wrap around and set to the first channel.

13.14 PrevChannel command

The 'PrevChannel' command sets the connected CH8, MUX or MUX8-R2 multiplexer to the previous channel in sequence. If this command is executed when the multiplexer is already at the first available channel, it will wrap around and set to the last channel.

13.15 Stirrer command




Figure 315 Panel with options for the 'Stirrer' command.

A stirrer can be controlled by any instrument from PalmSens BV that has an auxiliary port. This command sets the stirrer on or off by using the D0 and D1 digital lines.

13.16 SetAnalog command



Figure 316 Panel with options for the 'SetAnalog' command.

The 'SetAnalog' command sets a potential on the Analog output of the auxiliary port of the connected instrument. The output range depends on the connected instrument.

Please refer to the instrument-specific sections in this manual for more information about the positions of the digital and analog I/O pins on your instrument's auxiliary port.

13.17 ReadAuxiliary command

The 'ReadAuxiliary' command reads the value for the auxiliary input port of the connected device and appends the result to the output file. The input range depends on the instrument connected and/or the auxiliary input type connected.

Please refer to the instrument-specific sections in this manual for more information about the positions of the digital and analog I/O pins on your instrument's auxiliary port.

See for setting the type of auxiliary input: [Configuring PSTrace](#) on page 6

Output example for a temperature sensor:

```
Script output 18/08/2023 - 10:00:32:
18/08/2014 10:00:32  Aux. In:      T = 29.276 °C
18/08/2014 10:01:32  Aux. In:      T = 30.275 °C
18/08/2014 10:02:32  Aux. In:      T = 30.388 °C
18/08/2014 10:03:32  Aux. In:      T = 30.497 °C
18/08/2014 10:04:32  Aux. In:      T = 30.613 °C
```

13.18 SetDigitalIO command



Figure 317 Panel with options for the 'SetDigitalIO' command.

Set the digital IO ports of the connected device. Selected and highlighted buttons stand for a 1 otherwise 0, so in the example above the digital lines set are: 0101

Please refer to the instrument-specific sections in this manual for more information about the positions of the digital and analog I/O pins on your instrument's auxiliary port.

13.19 ReadDigitalIO command

Reads the digital input states of the 3 or 4 digital lines of the connected device and appends the result to the output file.

Output example:

```
Script output 12/07/2023 - 18:15:56:
12/07/2012 18:15:56  Digital lines: 0 0 0 0
12/07/2012 18:15:57  Digital lines: 0 0 0 1
12/07/2012 18:15:58  Digital lines: 0 0 1 0
```

Please refer to the instrument-specific sections in this manual for more information about the positions of the digital and analog I/O pins on your instrument's auxiliary port.

13.20 Override parameter command

The 'Override parameter' command is always attached to a Measure command. The value provided can either be 'Fixed' or an 'Added value'.




Figure 318 Panel with options for the 'Override parameter' command.

The command can be found under the 'advanced' tab. This command is helpful, if you want to repeat the same methods in one experiment several times with different parameters each time.

- Parameter: The method parameter to override.
- Original value: Refers to the parameter value retrieved from the loaded method. This value is presented in a muted grey shade for clarity and differentiation from user-modified values.
- New value: Refers to the value intended to replace the original value when the option 'New Value is Fixed' is selected. Alternatively, if the option 'New Value is Added' is selected, it represents the quantity/amount to be added (summed) to the original value. In case 'New Value is Added' is selected, and the measure command is used within a Repeat block as shown in the example below (Pseudo Polarography), this means that the result of the previous iteration will be used for each new value. If the resulting value is not valid, because it exceeds the parameter limits, the script will continue using the last valid value.
- Skip first iteration: This option is designed for use when there is an Override command within a Repeat and is visible only when 'New Value is Added' is chosen. When checked, it ensures that the first measurement (iteration) applies the original value instead of adding (summing). If no Repeat is associated, selecting this option will render the Override without effect.

13.20.1 Using the override parameter for Pseudo Polarography

The override parameter can be used to change the deposition potential for each iteration in a repeat loop. In the example below the deposition potential (E_{dep}) is read as 1.000 V in the method provided for the Measure command. The Override parameter command is attached to the Measure command. The Override parameter command adds 0.1 V with each iteration to the deposition potential. This results in a deposition potential in the first iteration of -0.4 V, in the second of -0.3 V, etc.




Figure 319 Using the 'override parameter' for pseudo polarography.

13.21 Running a script from the command line

Running PSTrace from the command line provides the following options:

```
> pstrace.exe ["path to .psscript file"] [-q] [-d#] [-h]
```

If a script file is specified, PSTrace will do the following automatically at startup:

- Connect to the first known device (PalmSens or EmStat) found.
If -d# is used it will connect to the device on position # as shown in the devicelist . So -d1 will connect to the first device in the list, -d2 to the second, etc.
(If -h is used in combination if -d#, PSTrace and the script window will remain invisible to the user.)
- Open the script window.
- Load the specified script file.
- Start the script.
- If -q is provided then PSTrace will be closed when the script has finished.



Make sure to use quotes around the file path if it has spaces.

Example:

The following two command lines could be included in a batch file (.BAT) or run from an external application to automate a measurement sequence:

```
cd "C:\Program Files (x86)\PalmSens BV\PSTrace 5.9"
pstrace.exe "C:\Users\John\Documents\PSData\script.psscript" -q
```


14 MethodSCRIPT



14.1 About MethodSCRIPT

MethodSCRIPT™ is the language our latest generation of potentiostats speaks. MethodSCRIPT gives you full control over the instrument or potentiostat module.

The script language is in plain text, uses SI units where possible and has a shallow learning curve.

The script can be sent directly to our instruments through PSTrace, a terminal, or virtually any software written in any programming language. Once the instrument receives the script, it compiles and executes it on-board, with no pre-processing required beforehand.

MethodSCRIPT allows you to customize and combine measurement techniques or to create custom techniques. It also allows for programming an instrument with a script to function stand-alone, without any host (PC or smartphone) intervention. MethodSCRIPT supports actions for:

- storing data on the internal storage;
- going to a sleep mode;
- waking up on time, digital or physical triggers;
- controlling external peripherals;
- doing calculations;
- using conditional statements;
- doing data analysis (e.g. peak detection);
- and more.

14.2 Generating MethodSCRIPT in PSTrace

PSTrace can be used to generate MethodSCRIPT from the parameters set in the Method Editor. This can be done by clicking the “Show MethodSCRIPT” button, which is found underneath the Method Editor.




Figure 320 The Show MethodSCRIPT button

14.2.1 MethodSCRIPT button visibility

In case the MethodSCRIPT button is not showing, there can be several reasons:

1. There is no instrument connected
2. The connected instrument does not support MethodSCRIPT
3. A setting needs to be changed in the General Settings window.

To change the setting for showing the MethodSCRIPT button, open the Settings window via the menu 'Tools' → 'General Settings ...'. The first tab 'Hardware' shows a checkbox for "Show MethodSCRIPT Editor".




Figure 321 Setting in General Settings window for showing the MethodSCRIPT button

14.2.2 Adding more generated MethodSCRIPT

When the button is clicked it will open the MethodSCRIPT Generator window. This window allows to add some additional MethodSCRIPT for certain functionality.




Figure 322 The MethodSCRIPT window showing generated script

The following functionality can be added to the generated script:

- “Send data during pretreatment”
The pretreatment part (present if $t > 0$) consists of maximum three ‘meas_loop_ca’ commands, which is basically a chronoamperometry measurement. One for conditioning, one for deposition and one for equilibration. During each pretreatment stage, the measured current data can optionally be sent by the potentiostat. If the checkbox is set, PSTrace will generate curves for each pretreatment stage.
- “Add time data”
When enabled, each ‘meas_loop_’ command will receive an extra variable for time. In PSTrace an extra column in the data tab will appear with the time for each datapoint, starting at 0 s for the first datapoint.
- “Hibernate and repeat”
This checkbox is particularly useful when using the Sensit Wearable. It creates an infinite loop around the measurement where it first sets the device in hibernation at the start of each iteration. As soon as the device is woken up it runs the measurement and resumes hibernation for the next iteration. The device will wake up by any form of communication by default. Connecting the device to USB and trying to connect for instance will wake the device up from hibernation.
The following options are available for waking the device from hibernation.
 - “Wake on double tap”

If the device's hardware, like the Sensit Wearable, has 'tap detection' this checkbox will make it wake up from hibernation when double tapped (physically) by the user.

- “Wake on pin/button”
Will wake up the device if the device’s hardware has a push button or a dedicated IO pin for waking up. Supported devices include, Sensit BT, Sensit Wearable and the EmStat Pico Development Board.
- “Wake after ... s”
Will make the device wake up as soon as the given time has passed.

“Store script in permanent memory”

This will add a command at the end of the script to place the script in non-volatile memory making it remain there when the device is powered off. In combination with the “AUTO-RUN” function it will allow the device to work independently and run the script always when it is powered on.



Be careful with this option and make sure to test your script first to see if it can be interrupted, before placing it in the permanent memory and enabling the AUTO-RUN option.

A faulty script in permanent memory in combination with the AUTO-RUN function enabled can brick your device!

To ensure that a hibernating device can always be woken up, the first bit of the ‘hibernate’ command must be enabled. You can easily verify this by checking if the argument passed to the hibernate command is an odd number. The MethodSCRIPT generator in PSTrace automatically handles this by default.

14.3 Running a MethodSCRIPT

It is not possible to run scripts directly from the MethodSCRIPT generation window. The generated script first needs to be transferred to the MethodSCRIPT Editor so it can be ran as a “MethodSCRIPT Sandbox” measurement. To do this, simply click on the “MethodSCRIPT Editor” button.




Figure 323 The “Open in Editor” button which opens the generated script in the Editor, allowing it to be edited and run

Now that the script is in the Editor, it can be ran by clicking the ‘Start Measurement’ green triangular button at the top of the main window.

PSTrace parses and stores the variables sent in ‘pcks’. Curves are generated automatically for each meas_loop that defines a pck with two or more variables, scripts with multiple meas_loops will generate multiple curves. The first variable in the pck will be set as the X-axis and a curve is created for each subsequent variable in the pck. Please note that to plot data versus time you will need to add a variable with the time to the pck.

For altering the script please refer to the MethodSCRIPT documentation which can be found under the Help button and on our website.

See for more information and documentation and example scripts:

www.palmsens.com/methodscript

15 MUX8 and MUX16 multiplexers



Figure 324 MUX multiplexer

The MUX8 and MUX16 are used to switch electrodes to the EmStat or PalmSens potentiostat.

Please note that a potentiostat with a multiplexer is not the same as a multi-channel potentiostat.

The EmStat-modules up to 3 (and 3+) were offered in combination with the MUX8-R2 multiplexer integrated in a single housing as the model “EmStat MUX”.

For the PalmSens and EmStat Blue the multiplexers come optional in a separate housing. The multiplexers require an available auxiliary port (AUX) DE-15. These products were discontinued and replaced by the MUX8-R2.

15.1 MUX8 multiplexer

The MUX8 multiplexer is used with 2- or 3- electrode sensors or cells up to 8 channels.

15.1.1 Functional diagram

The figure below shows the functional diagram of the MUX multiplexer.




Figure 325 Functional diagram of the MUX multiplexer

15.1.2 Cell connector pin-out

The connector on the MUX8 circuit board and the DSUB connector have the pinout as shown in the figure below.




Figure 326 MUX8 multiplexer cell connector pin-out

15.1.3 Specifications

The following table shows the main specifications of the MUX8 multiplexer.

Table 21 MUX8 multiplexer specifications

Number of channels	2-8
--------------------	-----

Multiplexer modes	Sensorarrays with up to eight working, reference and counter electrodes
	Sensorarrays with up to eight working and eight combined reference/counter electrodes
	Sensorarrays with up to eight working electrodes all sharing the same reference and the same counter electrode
	Sensorarrays with eight working electrodes all sharing the same combined reference/counter electrode
On resistance	WE, CE and RE channels: 2 Ω typical
Leakage current	10 pA typical at 25 °C
Charge injection	20 pC typical
Connections:	Shielded flat cable, with stripped end leads or by means of the MUX8 Connection Terminal (in shielded housing)
	8x WE
	8x CE
	8x RE
	one CE used when all WE's share one counter electrode
	one RE used when all WE's share one reference electrode
	analog ground for shielding

15.1.4 Sensor configurations

The multiplexer can be used in different modes. Each mode is set by several jumpers on the board or by means of the corresponding switches found at the bottom of the instrument.




Figure 327 Jumpers on the MUX8 circuit board




Figure 328 Switches at the bottom of the MUX8

Possible sensor configurations are:

- 1 Sensorarrays with (up to) eight working, reference and counter electrodes
- 2 Sensorarrays with eight working and eight combined reference/counter electrodes
- 3 Sensorarrays with eight working electrodes sharing a reference and a counter electrode
- 4 Sensorarrays with eight working electrodes sharing a combined reference/counter electrode

In all configurations the sensors can be multiplexed, leaving the not-selected sensors or cells at open circuit. Sensor configurations 2, 3 and 4 have the possibility to leave not-selected sensors or cells at open circuit or to apply the same potential to all sensors or cells.

15.1.5 Jumpers / switches

J1: RE to CE

Set ON when the sensor has a combined reference and counter electrode. This jumper therefore connects RE to CE.

J2: CE Direct

Set ON if the sensor array has more than one working electrode, but only one shared counter electrode. The CE from the potentiostat is switched directly to pin 1 and pin 2 of the cell connector.

J3: RE Direct

Set ON if the sensor array has more than one working electrode, but only one shared reference reference. The RE from the potentiostat is switched directly to pin 13 and 14 of the cell connector.

J4: CE MUX

Set ON if CE must be multiplexed. This is the case when each of the sensors has its own counter electrode (using pin 3 to 10 of the cell connector).

J5: RE MUX

Set ON if RE has to be multiplexed. This is the case when each of the sensors has its own reference electrode (using pin 15 to 22 of the cell connector).

J6: Unselected WE to GND:

Set ON if all unselected working electrodes or sensors must remain polarized at the potential as set by the potentiostat. If J6 is OFF only the selected channel is polarized, and all the idle channels will be in open circuit.

In case CE and RE are multiplexed as in Conf 1, so when J4 and J5 are ON, this jumper/switch is not relevant since only the selected channel's WE, CE and RE are polarized.

J7: Unselected WE to standby potential

Always leave it OFF when using the PSTrace.

15.1.6 Recommended settings for the available configurations

- Conf. 1: Sensor array with up to eight working, eight reference and eight counter electrodes:

- o Jumpers to be ON: J4 and J5

The potential is only applied to the selected channel. All other(idle) channels will remain at open circuit. (See remark below *)

- Conf. 2a: Sensor array with up to eight working and eight combined reference/counter electrodes:

- o Jumpers to be ON: J1, J4 and J5

Note: all leads CE1-8 and RE1-8 are connected together and this combined lead is connected to all eight combined reference/ counter electrodes.

The potential is only applied to the selected channel. All channels NOT selected are at open circuit. (See remark below *)

When J6 is also ON, the potential is applied to all working electrodes continuously.

- Conf. 2b: Sensor array with up to eight working and eight combined reference/counter electrodes:

- o Jumpers ON: J1, J2 and J3

Note: the combined reference/counter electrodes are connected to the leads CE Direct and/or RE Direct

When J6 is also ON, the potential is applied to all working electrodes continuously.

- Conf. 3: Sensor array with up to eight working electrodes all sharing one reference and one counter electrode:

- o Jumpers ON: J2 and J3

Note: the reference and counter electrodes are connected to RE Direct and CE Direct respectively.

When J6 is also ON, the potential is applied to all working electrodes continuously.

- Conf. 4: Sensor array with up to eight working electrodes all sharing one combined reference/counter electrode:

- o Jumpers ON: J1, J2 and J3

Note: the reference/counter electrode is switched to RE Direct and/or CE Direct.

When J6 is also ON, the potential is applied to all working electrodes continuously.



*** IMPORTANT REMARK**


It is not possible to apply a potential simultaneously to more than one sensor or cell each with three electrodes. This requires a multi-channel potentiostat, one potentiostat for each channel. This is however possible with two electrode sensors or cells, so when combined counter and reference electrodes are applied.

15.2 MUX16 multiplexer

The MUX16 multiplexer is used with 2 electrode sensors or cells up to 16 channels. It can have a shared counter and reference electrode, or each working electrode with each its own combined reference/counter electrode.

15.2.1 Functional diagram

The figure below shows the functional diagram of the MUX multiplexer.



15.2.2 Cell connector pin-out

The connector on the MUX16 circuit board and the DSUB connector have the pinout as shown in the figure below.




Figure 329 MUX16 multiplexer cell connector pin-out

15.2.3 Specifications

The following table shows the main specifications of the MUX16 multiplexer.

Table 22 MUX16 multiplexer specifications

Number of channels	2-16
Multiplexer modes	Sensorarrays with up to 16 working electrodes all sharing the same counter and same reference electrode
	Sensorarrays with up to 16 eight working and the same combined reference/counter electrode
	Sensorarrays with up to 16 working electrodes, each with its own combined reference/counter electrode
On resistance	WE channels: 2 Ω typical
Leakage current	10 pA typical at 25 °C
Charge injection	20 pC typical
Connections	Shielded flat cable, with stripped end leads or by means of the MUX16 Connection Terminal (in shielded housing)
	16x WE
	1x CE
	1x RE
	RE and CE can be connected together and in this case connected to up to 16 combined CE/RE electrode
	analog ground for shielding

Dimensions of PCB

76 x 74 mm

15.2.4 Sensor configurations




Figure 330 Jumpers on the MUX16 circuit board




Figure 331 Switches at the bottom of the MUX16

Possible sensor configurations are:

- 1 Sensorarrays with up to 16 working and 16 combined reference/counter electrodes
- 2 Sensorarrays with up to 16 working electrodes all sharing a reference and a counter electrode
- 3 Sensorarrays with up to 16 working electrodes sharing a combined reference/counter electrode

So the sensor configuration with 16 WE's all having their own reference and separate counter electrode is NOT applicable with the MUX16 multiplexer.

In these configurations the sensors or electrodes can be multiplexed, leaving the not-selected sensors at open circuit when J3 is not ON.

J3 is ON when the potential has to be applied to all working electrodes simultaneously.

15.2.5 Jumpers / switches

J1: RE to CE

Connects the CE and RE leads when ON. Use J1 when only a combined CE and RE is applied. An alternative is to leave J1 OFF and always connect both the CE and RE leads to your RE/CE electrode.

J2: Active

Leave ON

J3: (unselected WE) to AGND

Keeps all WE's at applied potential when ON. If OFF only the selected WE has an applied potential leaving not-selected WE's at open circuit.

J4: (unselected WE) to standby potential

Leave OFF

15.2.6 Recommended settings for the available configurations

- Conf. 1: Sensor array with up to 16 working and 16 combined reference/counter electrodes:
 - o Jumpers to be ON: J1 and J2

When J3 is also ON, the potential is applied to all working electrodes continuously.
(Use pin 1, 3, 5 to 31 of the cell connector to connect all combined RE/CE electrodes.)
- Conf. 2: Sensor array with up to 16 working electrodes all sharing one reference and one counter electrode:
 - o Jumper to be ON: J2

When J3 is also ON, the potential is applied to all working electrodes continuously.
Use pin 1, 5, 9, 13, 17, 21, 25 or 29 of the cell connector to connect RE.
Use pin 3, 7, 11, 15, 19, 23, 27 or 31 of the cell connector to connect CE.
- Conf. 3: Sensor array with up to 16 working electrodes all sharing one combined reference/counter electrode:
 - o Jumpers to be ON: J1 and J2

When J3 is also ON, the potential is applied to all working electrodes continuously.
Use pin 1, 3, 5..... to 31 of the cell connector to connect RE/CE.



* IMPORTANT REMARK

It is not possible to apply a potential simultaneously to more than one sensor or cell each with three electrodes. This requires a multi-channel potentiostat, one potentiostat for each channel. This is however possible with two electrode sensors or cells, so when combined counter and reference electrodes are applied.

16 MUX8-R2 multiplexer

The MUX8-R2 multiplexer can be used to expand a PalmSens3, PalmSens4, EmStat4X, EmStat3 Blue or EmStat3+ Blue potentiostat.

The multiplexer allows to increase productivity by automatically switching between multiple electrochemical cells each with their own WE, RE and CE electrodes. Additionally, it can be employed to automatically switch the active Working Electrode (WE) within a cell containing multiple WEs, while utilizing a shared RE/CE.

The MUX8-R2 model has replaced both the MUX8 and MUX16 models. This new version incorporates a new feature of electronic switching, enabling sensor configurations (such as the combination or merging of channel leads) through software, eliminating the need for manual switchers or jumpers.



Figure 332 PalmSens4 connected to MUX8-R2 multiplexer

The MUX8-R2 multiplexer is designed for up to 128 channels (by stacking a maximum of 16 MUX8-R2's) with 2- or 3- electrode sensors or cells.

16.1 MUX8-R2 functional diagram

The figure below shows the functional diagram of the MUX multiplexer.




Figure 333 Functional diagram of the MUX8-R2 multiplexer

16.2 Specifications

The following table shows the main specifications for the MUX8-R2 multiplexer.

Table 23 MUX8-R2 specifications

number of channels	8 (up to 128 channels when daisy chained)
multiplexer	switches 8 x (WE, S, RE and CE)
on resistance for WE	1.5 Ω typical
charge injection on WE	20 pC typical
leakage current	< 20 pA (5 pA typical) at 25 °C
switching time	2 ms
compliance voltage	± 10 V
max. frequency for EIS	100 kHz when switching WE/S, RE and CE 1 MHz when switching WE/S and RE+CE combined (2-electrode configuration)

16.3 Configurations

The MUX8-R2 multiplexer is designed for use up to 128 channels with 2- or 3- electrode sensors or cells.

The multiplexer can be used with different electrode or sensor configurations:

- 1 Eight separate cells or sensors each with a working/sense, reference and counter electrode.
- 2 Eight separate cells or sensors each with a working/sense and combined reference and counter electrode.

- 3 Cell or sensor array with eight working/sense electrodes sharing one reference and one counter electrode.
- 4 Cell or sensor array with eight working/sense electrodes sharing one combined reference/counter electrode.

In all configurations the cells can be multiplexed, leaving the non-selected working electrodes either at open circuit (individually floating) or at Ground potential.

When in configurations 3 and 4, the unselected channels are switched to Ground, they will have the working electrode's potential. This is due to the fact that the active WE is always at Ground potential.

You can easily change the hardware configuration of the MUX8-R2 as part of the measurement settings in our PSTrace software or the PSTouch app for Android.




Figure 334 The MUX8-R2 multiplexer settings can be changed in the Multiplexer tab.

16.4 Connections


The table below shows which connectors are available on the MUX8-R2 for making connections to the instrument, cell or other MUX8-R2 multiplexers.

Table 24 MUX8-R2 available connectors

Connector	Function
Input	Y-cable connects to both potentiostat sensor connector and (digital) AUX
AUX	Can be used to measure auxiliary input like temperature or pH, and to switch external hardware using two digital control lines that can be set in PSTrace
Link	Connects to Input of next multiplexer, for daisy-chaining multiple multiplexers.

USB-C	For providing extra power in case more than 2 multiplexers are connected to a single instrument.
Channel 1-4	Connects to sensor cables 1-4
Channel 5-8	Connects to sensor cables 5-8


16.5 MUX8-R2 Pin-outs



The table below shows the pin-out for channel 1-4.

Table 25 Channel 1-4 connector pin-out

Pin	Function	Pin	Function	Pin	Function
1	CE4	11	RE1	21	AGND
2	RE4	12	RE_SHIELD1	22	SENSE2
3	RE_SHIELD4	13	NC	23	WE1
4	CE3	14	WE4	24	AGND
5	RE3	15	AGND	25	SENSE1
6	RE_SHIELD3	16	SENSE4		
7	CE2	17	WE3		
8	RE2	18	AGND		
9	RE_SHIELD2	19	SENSE3		
10	CE1	20	WE2		



The table below shows the pin-out for channel 5-8.

Table 26 Channel 5-8 connector pin-out

Pin	Function		Pin	Function		Pin	Function
1	CE8		11	RE5		21	AGND
2	RE8		12	RE_SHIELD5		22	SENSE6
3	RE_SHIELD8		13	NC		23	WE5
4	CE7		14	WE8		24	AGND
5	RE7		15	AGND		25	SENSE5
6	RE_SHIELD7		16	SENSE8			
7	CE6		17	WE7			
8	RE6		18	AGND			
9	RE_SHIELD6		19	SENSE7			
10	CE5		20	WE6			

16.5.1 AUX port

The diagram below shows the connections available in the auxiliary (AUX) port.




Figure 335 Auxiliary port pin out

The lead colors in the table below correspond to the standard cable PalmSens provides for use with the auxiliary port.

Table 27 Auxiliary port pin assignment and corresponding lead color

Ring	Lead color	Function
1	Red	NC
2	Black	D0 (digital I/O)
3	Yellow	Analog Input
4	White	D1 (digital I/O)
5	Green	GND
6	Blue	5V

16.6 Stacking multiple MUX8-R2

Each multiplexer has a Link connector which can be used to daisy chain to another MUX8-R2 multiplexer, expanding the number of channels. A maximum of 16 multiplexers can be connected in a daisy chain, giving a maximum of 128 channels.

The PSTrace software detects automatically how many multiplexers are daisy chained and shows the available number of channels in the user interface.

When in configurations 3 and 4, the RE/CE be connected to the CH1 of the first Multiplexer.



Figure 336 The MUX8-R2 has magnetic feet and magnets in the top for easy stacking.

17 PalmSens1 and PalmSens2



Figure 337 PalmSens1 and PalmSens2 have the same appearance.

17.1 Description

The PalmSens1 (PS1) and PalmSens2 (PS2) are both a hand-held battery powered instrument for use with electrochemical sensors or electrochemical cells. The instrument is a low noise and low-current potentiostat and galvanostat which controls the potential or current applied to the sensor and measures the current or potential response.

The suffix “1” or “2” refers to its generation and is not present in the instrument label.

The PalmSens1 and 2 are very similar. The PalmSens1 can be identified having a serialnumber between PS0201 to PS04200. The main visible difference compared to the PalmSens1 is the 5V DC-input port (barrel connector) present on the PalmSens2 for charging, whereas the PalmSens1 is charged through a mini-DIN connector.

17.2 Operating

The instrument connects to a PC using its serial RS-232 connector which is found behind the lid. See next section for more information about using a serial to USB adapter cable.



Do not use a screwdriver or tool to open the lid!

Open the lid at the left-hand side of the PalmSens by pressing the hinge upwards and then turning the lid in a single action.




Figure 338 Opening the lid of PalmSens1 or 2

Before the PalmSens is used without the adapter, the batteries must be charged. Connect the adapter to the miniDIN-connector (PS1) or dc-in connector (PS2). Switch on the PalmSens by pressing the power key until the display shows: "Selftest". During the self-test the voltage range of the instrument is tested and shown. The normal range for PalmSens is approx. -2.035 V to $+2.047\text{ V}$.

After the test the display will show:

> PalmSens	Vs. #.#
E= #.### V	
FsChg	

Figure 339 Default contents of the PalmSens1 and 2 LCD

For every reading during idle mode, the < > sign in the upper left corner flips to indicate the instrument is actively running.

If the display shows the text "FsChg", the batteries are being charged.

As soon as the display shows "PwrOK", the batteries are full and the adapter can be disconnected. The batteries have to be recharged as soon as "LowBat" is shown and the corresponding beep is heard.

The PalmSens can be used with the adapter connected and the batteries charged. The adapter however might increase the noise level.

The backlight of the LCD is switched on and off by using the key ▲. Please note that the backlight will reduce the battery-lifetime of the PalmSens from approx. 8 hours to not more than 6 hours.

The PalmSens is switched off by pressing and holding the power key for a few seconds until the display shows "Shutting down".

17.2.1 Cell connections

For more information about making a connection to the cell, see section: [Connecting a cell to the potentiostat](#) on page 25.


17.2.2 Keypad

During the conditioning, deposition, and equilibration stage two keys of the keypad are active. When the stop button is pressed the measurement is ended. The skip key is used to step to the next stage. During a measurement only the ESC or ■ key is active.

17.3 Using a USB to serial converter

When using a USB – RS232 adapter, ensure that the adapter's driver is installed in accordance with the adapter's manual. The COM number used can be identified in Windows' Device Manager.

It is crucial to note that a null-modem adapter (female-female) is required between the cable and PalmSens when using this adapter. Attempting to connect the PS1 or PS2 directly to a standard USB to Serial COM adapter (or with a direct pin-to-pin adapter) will not be successful.

To open the Device Manager quickly, simultaneously press the Windows key  and the Pause/Break key on the keyboard of your computer. Click the Hardware tab, and then click Device Manager. The installed COM port and its number are listed under Ports. Specify the corresponding COM port in PSTrace before connecting.




Figure 340 Make sure the right COM port is used in PSTrace

17.4 Auxiliary port pin-out

The following picture and table show the PalmSens1 and PalmSens2 auxiliary port pin-out and function descriptions.




Figure 341 Front view of female port (mini-DIN)

Table 28 PalmSens1 and PalmSens2 auxiliary port pin-out

Pin	Function
1	PalmSens1 (serial numbers PS0201 to PS04200): Vadapter (6 V- 1.5 A) PalmSens2 (serialnumbers PS04201 and higher): d3 digital output


2	analog ground
3	d0 digital output *
4	d1 digital output *
5	d2 digital output *
6	d0 digital input *
7	auxiliary analog input 0 V - 4.096 V
8	auxiliary analog output 0 V - 4.096 V > 10 kOhm advised resistive load on output
9	5 V digital power line (maximum current 100 mA)
shield	digital ground

* The digital input has a maximum voltage rating of -0.3 V to 5.3 V. The digital outputs have a maximum current rating of 5 mA, source or sink.

17.5 Cell connector pin-out

The following schematics show the PalmSens1 and PalmSens2 cell connector pin-out.

Front view of male plug on cable:



Front view of female connector on instrument:




Table 29 Cell connector pins

Pin	Function	Connector color
1	Working Electrode (WE)	Red
2	Analog Ground (AGND)	Green
3	Working Electrode 2 / BiPot (WE2) (if BiPot module is installed)	Yellow
4	Counter / Auxiliary Electrode (CE)	Black
5	Reference Electrode (RE)	Blue

17.6 Specifications

The following tables show the main specifications of the PalmSens1 and PalmSens2.

Table 30 PalmSens1 and PalmSens2 specifications

Potentiostat (controlled potential mode)	
dc-potential range	± 2.000 V

compliance voltage	± 8.0 V
dc-potential resolution	1 mV
dc-offset error	2 mV
accuracy	≤ 0.2 %
ac-potential amplitude	1 mV to 250 mV
current ranges	1 nA to 10 mA (8 ranges)
maximum current	± 10 mA
current resolution	0.1 % of current range 1 pA on lowest current range
accuracy	≤ 0.2 % of current range at 100 nA to 1 mA ≤ 0.5 % at 10 nA and ≤ 1 % at 1 nA all with additional 0.2 % offset error

Galvanostat (controlled current mode)	
current ranges	1 µA to 1 mA
dc-current range	-2 to + 2 times selected current range
dc-current resolution	0.1 % of selected current range
dc-offset error	≤ 0.2 %
current accuracy	≤ 0.4 %
maximum output voltage	± 8 V

General	
electrometer amplifier input	> 100 Gohm // 4 pF
rise time	approx. 50 µs
Keypad	▲▶▼◀ENTER ESC and Power (7 keys)
Display	4 lines of 16 characters with backlight
Dimensions	155 mm x 85 mm x 35 mm
Temperature range	0° C to + 40° C
Weight	0.43 kg
Power	2 AA cells NiMH 2500 mAh for > 8 hours operation. Battery charger included (6 V- 1500 mA)
Interfacing	Serial RS232 Default serial mode: 57600 baud, 8 bits data, no parity, 1 stopbit.

17.7 BiPot specifications

The table below shows the specifications of the optional Bipotentiostat module for PalmSens1 and PalmSens2.

Table 31 PalmSens1 and PalmSens2 BiPot specifications

General	
dc-potential range	± 2.000 V
dc-potential resolution	1 mV
dc-offset error	3 mV
accuracy	≤ 0.3 %
current ranges	1 nA to 10 mA (8 ranges)
maximum measured current	± 10 mA
current resolution	0.1 % of current range 1 pA on lowest current range
accuracy	≤ 0.3 % of current range at 1 μ A to 100 μ A
	≤ 0.5 % at 100 nA and ≤ 1 % at 10 and 1 nA
	all with additional 0.2 % offset error
connection	Use requires a cable with additionally a (yellow) connector for WE2

17.8 Offset calibration

When connected to a PalmSens1 or 2, a button for calibrating offsets can be found In PSTrace in the menu: 'Tools' → 'Instrument Settings' → button: 'Calibrate Offsets'.

Offset calibration can be done when the measured potential differs more than 0.3% +/- 0.003 V from the specified potential (when the cell is switched on in E control mode).

17.9 Battery replacement



This section does not apply to the PalmSens3 or PalmSens4!
A PalmSens1 or PalmSens2 has a GREEN display.



Necessary tools:
Size 14 spanner
Phillips or TORX T-10 screwdriver
Flat screwdriver
Pointy tool

Step 1



Untighten the bolt very carefully with the spanner, without scratching the surface of the housing.

Screw it off using your fingers.

Step 2



If present remove the two protective caps using the pointy tool by sticking it between the cap and housing and levering it out.

Step 3



Remove the sensor side of the housing.

Step 4



Pull out the battery connector using the tongs.

Replace the battery and connect it to the printed circuit board.

Step 5 Closing the housing

1. Place the ring back over the sensor connector and screw the sensor side of the housing back on.



Do not use too much force while screwing, this may cause the screw to slip or break!

2. Place the nut back onto the sensor connector and screw it using your fingers as fixed as you can. Then tighten the nut by turning it one quarter more using the spanner. Tightening it too much will damage the housing.
3. If available put the protective caps in place.
4. New batteries may be empty, so first connect PalmSens to the ac-adapter and charge the batteries completely.

18 PalmSens3



Figure 342 PalmSens3

18.1 Description

The PalmSens3 is a hand-held battery powered instrument for use with electrochemical sensors or electrochemical cells. The instrument is a low noise and low-current potentiostat and galvanostat which controls the potential or current applied to the sensor and measures the current or potential response.

The PalmSens3 was on the market between 2013 and 2016. The main difference compared to the PalmSens1 and 2 are much faster ARM microcontroller and FRA capabilities for Electrochemical Impedance Spectroscopy (EIS). The EIS capability was offered as a separate option.

18.2 Operating

The instrument connects to a PC using its mini-USB port which is found underneath the lid.



Do not use a screwdriver or tool to open the lid!

Open the lid at the left-hand side of PalmSens by pressing the hinge upwards and then turning the lid.




Figure 343 Instructions for opening the lid at the side of the PalmSens3

	<p>Power button</p> <p>PalmSens3 is switched on by pressing and holding the power button. If a PalmSens3 Bluetooth extension is used, make sure the Bluetooth extension is connected to the PalmSens3 auxiliary port BEFORE turning it on. PalmSens3 will detect any extensions during start-up. PalmSens3 is switched off by pressing the power button twice. If PalmSens3 is not responding, press and hold the power button for a few seconds to force it to switch off.</p>
	<p>Backlight button</p> <p>The backlight of the LCD is switched on by using the backlight button next to the power button. Please note that the backlight will reduce the battery-lifetime of PalmSens3. The backlight is turned off automatically after a while.</p>
	<p>Charging icon</p> <p>PalmSens3 is always charging when connected to a USB port. When charging the red battery icon will light up.</p>
	<p>Stop button</p> <p>Pushing the stop button will abort any running measurements immediately. If the instrument is in 'fast mode' it will return to normal mode.</p>
	<p>Start button</p> <p>No use in PSTrace.</p>
	<p>Skip button</p> <p>This button can be used to skip pretreatment stages and proceed to the next stage. See also section. Standard measurement sequence on page 32.</p>



The PalmSens3 cell connector is a push-pull connector.
Do not try to twist it!

18.2.1 Cell connections

For more information about making a connection to the cell, see section: [Connecting a cell to the potentiostat](#) on page 25.

18.3 Auxiliary port pin-out

The following picture and table show the PalmSens3 auxiliary port pin-out and function descriptions.




Figure 344 The front view of the female port (DE-15)

Table 32 PalmSens3 auxiliary port pins

Pin	Function
1	d0 digital output (5V)
2	d3 digital output (5V)
3	auxiliary analog input (>0.5 Mohm input impedance)
4	Rx (TTL default, or RS232 see pin 12)
5	Tx (TTL default, or RS232 see pin 12)
6	d1 digital output (5V)
7	d0 digital input (5V)
8	I out (V in current range)
9	5 V digital power line (max. 50mA)
10	digital ground
11	d2 digital output (5V)
12	Connect to pin 10 (DGND) for RS232 comm on pin 4 and pin 5
13	E out (-5 to +5V)
14	analog ground
15	DAC out (0-3V)
Connector housing	digital ground

18.4 Cell connector pin-out

The following schematics show the PalmSens3 cell connector pin-out.

Front view of male plug on cable:

```

  4 •   • 1
  3 •   • 2

```

Front view of female connector on instrument:

```

  1 •   • 4
  2 •   • 3

```

Table 33 Cell connector pins of the PalmSens3

Pin	Function	Connector color
1	Reference Electrode (RE)	Blue
2	Counter / Auxiliary Electrode (CE)	Black
3	Working Electrode 2 / BiPot (WE2) (if BiPot module is installed)	Yellow
4	Working Electrode (WE)	Red
Connector housing	Analog Ground (AGND)	Green

18.5 Specifications

The following tables show the main specifications of the PalmSens3.

Table 34 PalmSens3 specifications

Potentiostat (controlled potential mode)	
dc-potential range	± 5.000 V
compliance voltage	± 8.0 V
dc-potential resolution	0.15 mV
applied potential accuracy	≤ 0.2 % with max. 2 mV offset error
current ranges	100 pA to 10 mA (9 ranges)
maximum measured current	± 30 mA (typical)
current resolution	0.01 % of current range
accuracy	≤ 1 % of current range at 1 nA (≤ 5 % at 100 pA) ≤ 0.5 % at 10 nA ≤ 0.2 % at 100 nA to 1 mA ≤ 0.5 % at 10 mA all with max. 0.2 % offset error
max. acquisition rate	200 000 data points/s

Galvanostat (controlled current mode)	
current ranges	1 μ A to 10 mA
dc-current range	\pm 3.000 times selected current range
dc-current resolution	0.01 % of selected current range
max. dc-offset error	\leq 0.2 %
current accuracy (deviation)	\leq 0.4 %
maximum output voltage	\pm 8 V

Impedance measurements	
frequency range	100 μ Hz to 50 kHz
ac- amplitude range	1 mV to 0.25 V (rms)

General	
electrometer amplifier input	> 100 Gohm // 4 pF
rise time	programmable from min. 0.5 μ s

Other	
keypad	run, skip, abort, backlight and power
housing	aluminium: 155 mm x 85 mm x 35 mm
weight	430 g
temperature range	0° C to + 40° C
power supply	USB or internal Li-ion battery
battery time	>9 hours idle time with Bluetooth extension
communication	USB

18.6 BiPot specifications for PalmSens3

The table below shows the specifications of the optional Bipotentiostat module for PalmSens3.

Table 35 *PalmSens3 BiPot specifications*

General	
dc-potential range	\pm 5.000 V
dc-potential resolution	0.15 mV
dc-offset error	3 mV
accuracy	\leq 0.3 %
current ranges	1 nA to 10 mA (8 ranges)
maximum measured current	\pm 10 mA
current resolution	0.01 % of current range 0.1 pA on lowest current range

accuracy	≤ 0.3 % of current range at 1 µA to 100 µA
	≤ 0.5 % at 100 nA and ≤ 1 % at 10 and 1 nA
	all with additional 0.2 % offset error
connection	Use requires a cable with additionally a (yellow) connector for WE2

18.7 Replacement of PalmSens1 or PalmSens2 battery



This section does not apply to the PalmSens3 or PalmSens4!
A PalmSens1 or PalmSens2 has A GREEN display.



Necessary tools:

- Size 14 spanner
- Phillips or TORX T-10 screwdriver
- Flat screwdriver
- Pointy tool

Step 1



Untighten the bolt very carefully with the spanner, without scratching the surface of the housing.


Screw it off using your fingers.

Step 2



If present remove the two protective caps using the pointy tool by sticking it between the cap and housing and levering it out.

Step 3



Remove the sensor side of the housing.

Step 4



Pull out the battery connector using the tongs.

Replace the battery and connect it to the printed circuit board.

Step 5 Closing the housing

5. Place the ring back over the sensor connector and screw the sensor side of the housing back on.



Do not use too much force while screwing, this may cause the screw to slip or break!

6. Place the nut back onto the sensor connector and screw it using your fingers as fixed as you can. Then tighten the nut by turning it one quarter more using the spanner. Tightening it too much will damage the housing.

7. If available put the protective caps in place.
8. New batteries may be empty, so first connect PalmSens to the ac-adapter and charge the batteries completely.

18.8 EIS Calibration for PalmSens3



Follow the instructions of this section only if your PalmSens3 has EIS capabilities.

After the PalmSens3 has been restored it is advised to run the EIS calibration procedure.

To do this open PSTrace. Then go to menu: 'Tools' → 'Instrument settings...' and click the "Calibrate for EIS" button.

Follow the instructions on the screen.

A normal calibration procedure takes about 5 minutes.

18.9 Offset calibration

When connected to a PalmSens3, a button for calibrating offsets can be found in PSTrace in the menu: 'Tools' → 'Instrument Settings' → button: 'Calibrate Offsets'.

Offset calibration can be done when the measured potential differs more than 0.3% +/- 0.003 V from the specified potential (when the cell is switched on in E control mode).

18.10 Performing a hard reset



Please use this procedure only if your PalmSens3 shows a blue screen without characters when turning on:



18.10.1 Follow these steps to restore your PalmSens3:

1. Open the Update Firmware window from the main menu in PSTrace:
'Tools' → 'Instrument Settings...' → 'Update Firmware'
2. Make sure the USB cable is connected and PalmSens3 shows in the list next to the 'Connect' button.

- Click the 'Connect' button:




Figure 345 Connect button

- The program will say "Could not connect properly." Ignore this.
- Click the 'Browse' button.




Figure 346 Browse button

6. Select a .HEX file with the signature PalmSens3_firmware_v##.hex. Where ## is the version number.
7. Click the button 'Start update firmware'.
8. The firmware will be restored:




Figure 347 Firmware update progress window



Normally it takes only a few minutes to download new firmware to the PalmSens3.

If the update procedure does not start or something is still wrong contact PalmSens BV by email: info@palsens.com

18.11 Battery replacement

Required tools:

- Torx T10 screw driver + small flathead screw driver
- Blue Lemo tool (DCH.91.121.PA-ND) for the sensor connector nut:



Figure 348 Lemo tool



Using another tool (i.e. pliers) will probably scratch the cell cable connector.

Step 1



Remove the LEMO nut using the special LEMO tool.

Step 2




Remove both Torx screws

Step 3



Replace the battery.
The black tape around the end works as a strain relief.

Step 4

Make sure the grounding rings are placed correctly before placing the side panel back on.

Step 5

Screw the sensor side panel to the housing. Make sure to tighten the screws. Tighten the LEMO nut firmly, but do not use too much force.

19 PalmSens4



Figure 349 PalmSens4

19.1 Description

The PalmSens4 is a hand-held battery powered instrument for use with electrochemical sensors or electrochemical cells. The instrument is a low noise and low-current potentiostat and galvanostat which controls the potential or current applied to the sensor and measures the current or potential response.

The PalmSens4 is on the market since 2016. Main new features added to the PalmSens4, compared to the PalmSens3, includes a color display, internal storage and improved FRA for Electrochemical Impedance Spectroscopy (EIS).

19.2 Operating

Switch on the PalmSens4 by pressing the power button for a couple of seconds.

A full charge of the PalmSens4 takes about 8 hours.

The PalmSens4 can be used either on battery using a wireless Bluetooth connection or by connecting it via its USB-C port to a PC or Android device. The Bluetooth name of the instrument can be found on the bottom of the display when the instrument is idle. The Bluetooth identifier always starts with the letters PS, followed by the last 4 characters of its MAC address. For example: PS-6AEF.

When connecting to the USB port of a PC, Bluetooth will switch off automatically unless there is a Bluetooth connection active. This allows the instrument to be charged via the PC while a Bluetooth connection is present.

In case a pairing code is requested by the host (PC or Android device), use the pairing code 1234. If you experience issues with the Bluetooth connection, see also section: [Connecting using Bluetooth](#) on page 14.

Switch on the PalmSens4 by pressing the power button for a couple of seconds.

A full charge of the PalmSens4 takes about 8 hours.

19.2.1 Cell connections

For more information about making a connection to the cell, see section: [Connecting a cell to the potentiostat](#) on page 25.



The PalmSens4 cell connector is a push-pull connector.
Do not try to twist it!

19.3 Display and power indicator

The PalmSens4 has a color display and illuminated white ring around the power button to indicate the instrument is powered on.




Figure 350 PalmSens4 display when in idle mode.

In idle mode (not running a measurement), the display shows the following contents:

- Firmware version
- Connection status (Bluetooth or USB)
- Battery status
- Cell status (ON/OFF)
- Potential controlled mode (PSTAT) or current controlled mode (GSTAT)
- Potential readings (E)
- Current readings (I)
- Bluetooth name
- Time

For every reading during idle mode, the < > sign in the bottom left corner flips to indicate the instrument is actively running.

The time is updated whenever the instrument connects to PSTrace (Windows) or PStouch (Android). The time is used together with the date for timestamping data on the internal storage.

19.4 iR Compensation

The PalmSens4 iR Compensation module is an optional accessory available for the PalmSens4 potentiostat. It needs to be ordered concurrently with the new instrument, or alternatively, it can be installed later at the factory.

This module operates through Positive Feedback, achieved using a 16-bit MDAC in the module which scales the output of the current follower to provide a positive feedback voltage that is proportional to the current through the cell. The compensation voltage is added to the summing point before the control amplifier and thus increases (or decreases) the applied potential to counteract the iR drop.




Figure 351 Simplified schematical representation of iR compensation circuitry in the PalmSens4

Positive feedback allows for fast scan rates up to 10 V/s, depending on the characteristics of the cell. If the potential error to compensate for becomes close to the value set for E applied, the system might become unstable. Using iR compensation limits the measurement bandwidth to 10 kHz.

See also section: [Ohmic \(iR\) compensation](#) on page 46.

19.5 Auxiliary port pin-out

The following schematics and table show the PalmSens4 auxiliary port pin-out and pin functions.




Figure 352 The front view of the female port (DE-15) on the PalmSens4.

Table 36 PalmSens4 auxiliary port pin functions.

Pin	Function
1	d0 digital output (5V)
2	d3 digital output (5V)
3	auxiliary analog input -10 to +10 V, 18 bit, >0.5 M Ω input impedance
4	RESERVED
5	RESERVED
6	d1 digital output (5V)
7	d0 digital input (5V)
8	i monitor given as -V in active current range.
9	5V digital power line (max. 300mA)
10	digital ground
11	d2 digital output (5V)
12	RESERVED (NC)
13	E out (-10 to +10V)
14	analog ground
15	analog out (0 to 10 V at 12 bit) , 1 k Ω output impedance
Connector housing	digital ground

19.6 Cell connector pin-out

The following schematics show the PalmSens4 cell connector pin-out and pin functions.




Figure 353 Front view of the female Cell connector on the PalmSens4.

Table 37 Cell connector pin functions.

Pin	Function	Connector color
1	Reference Electrode (RE)	Blue
2	Reference Electrode Shield	N/A
3	Counter / Auxiliary Electrode (CE)	Black
4	Working Electrode 2 / BiPot (WE2) (if BiPot module is installed)	Yellow
5	Working Electrode (WE)	Red
Connector housing	Analog Ground (AGND)	Green



The shield of the cable must make contact with the metal case of the sensors connector.

19.7 Self Diagnostics

PalmSens4 supports Self Diagnostics from firmware v1.8 and on. Check the PalmSens4 LCD or the Instrument Settings window in PSTrace to see which firmware version the instrument is running. The Self Diagnostics tool generates a report that gives a detailed overview the instrument's functioning. This report can exported to a Word document and sent by email (info@palmsens.com) to PalmSens BV in case there's doubt about the functioning of your PalmSens4.




Figure 354 The Self Diagnostics window

To run the Self Diagnostics tool go to the menu: 'Tools' → 'Instrument Settings...' and click the button Self Diagnostics.




Figure 355 Menu 'Instrument Settings' found in 'Tools'





Figure 356 The 'Self Diagnostics' button found in 'Instrument Settings'

19.8 Performing a hard reset

The instructions in this section should be followed if the PalmSens4 is in a state where the LCD remains black after switching on and the instruments does not show up in the list with instruments in PSTrace or the Firmware Update program as “PalmSens4” or as “Programmable Device”.



Follow these instructions only if the PalmSens4 does not show as a “Programmable Device” in the device list.



If it shows as “Programmable Device”, continue with Step 5.

Required tools:

- Screwdriver: Torx 8
- Paperclip

Step 1: Remove battery lid and battery



Figure 357 PalmSens4 with the battery lid removed.

Step 2: Connect the PalmSens4 to USB

Connect the PalmSens4 to your PC with the USB cable. Place the PalmSens4 up-side-down on your desk.

Step 3: Short the two 'Rst' pads

Bend the paperclip in such a way that you can short the two 'Rst' pads in the battery compartment, as shown in the photo below.




Figure 358 Shorting the two 'Rst' pads on the PalmSens4 for recovery.

Step 4: Unplug and plug back in

Remove the USB cable from the PalmSens4, wait a few seconds and then connect the PalmSens4 back to the PC again.

Step 5: Upload new firmware

- Open PSTrace and go to menu: Tools → Instrument Settings... and click button "Update firmware".




Figure 359 PalmSens4 listed as 'Programmable Device'

- You should see "Programmable Device" selected next to the "Connect" button.
- Click "Connect"
- Click "Browse" and select file "PalmSens4_firmware_v##.hex"




Figure 360 Verifying the loaded firmware and the button to start the update procedure.

- Make sure the text shows it is firmware suitable for PalmSens4.



Make sure the text shows the selected firmware is suitable for PalmSens4!

- Click the button “Start update firmware”.




Figure 361 The firmware upload progress window.

- Wait until the Uploader gives the message “UPDATE COMPLETE” and close the window.

Step 6: Place the battery

Connect the 3000 mAh battery and place it in the battery compartment with the wires folded at the side of the battery.

Step 7: Close housing

Put the lid and the gasket back in place and close the housing. Do not use excessive force when tightening the screws.

19.9 Battery replacement

Required tools:

- Screwdriver: Torx 8

Step 1: Open the battery compartment

Unscrew the four screws and remove the battery lid and gasket. Then gently pull on the battery connector to disconnect it.



Figure 362 PalmSens4 with the battery lid removed.

Step 2: Place new battery

Connect the new 3000 mAh battery as received from PalmSens BV and place it in the battery compartment with the wires folded at the side of the battery.

Step 3: Close housing

Put the lid and the gasket back in place and close the housing. Do not use excessive force when tightening the screws.

19.10 Specifications

The following tables show the main specifications of the PalmSens4.

Table 38 PalmSens4 specifications

General			
dc-potential range	version	PS4.F#.05	PS4.F#.10
		±5 V	±10 V

compliance voltage	±10 V
maximum current	±30 mA (typical)
max. acquisition rate	150000 points/s

Potentiostat (controlled potential mode)

applied dc-potential resolution	75 µV
applied potential accuracy	≤ 0.1% ±1 mV offset
current ranges	100 pA to 10 mA (9 ranges)
measured current accuracy	≤ 0.1% at Full Scale Range
measured current resolution	0.006% of current range (5 fA on 100 pA range)

Galvanostat (controlled current mode)

current ranges	1 nA to 10 mA (8 ranges)
applied dc-current range	± 6 times applied current range
applied dc-current resolution	0.005% of applied current range
measured dc-potential resolution	75 µV at ±10 V (no gain) 7.5 µV at ±1 V (gain 10) 0.75 µV at ±0.1 V (gain 100)
measured dc-potential accuracy	≤ 0.05% or ±1 mV (for E < ±9 V) ≤ 0.2% (for E ≥ ±9 V)

FRA / EIS (impedance measurements)

frequency range	PS4.F1.##	PS4.F2.##
	10 µHz to 100 kHz	10 µHz to 1 MHz
ac-amplitude range	1 mV to 0.25 V rms, or 0.6 V p-p	

* Instruments with part number PS4.F0.## do not include EIS capability.

Electrometer	
electrometer amplifier input	> 1 TΩ // 10 pF
bandwidth	1 MHz

Other

housing	aluminium body with rubber sleeve: 15.7 x 9.7 x 3.5 cm ³
weight	+/- 500 g
temperature range	0 °C to +50 °C
power supply	USB or internal LiPo battery
communication	USB and Bluetooth

battery	11.1 Wh capacity > 16 hours idle time (> 5 hours with BiPot installed) > 4 hours with cell on at max. current Extendible by means of power bank
internal storage space	8 GB (or >100 million data points)

19.10.1 EIS contour accuracy plot




Figure 363 EIS contour accuracy plot for the PalmSens4

Note: The accuracy contour plot was determined under lab conditions and should be used for reference purposes. Please note that the true limits of an impedance measurement are influenced by all components in the system, e.g. cables, the environment, and the cell.

20 EmStat instruments series



Figure 364 The EmStat instrument series

This section describes the housed models of the EmStat-series. The EmStat-series is also available in the form of *modules* (bare PCBs) for OEM purposes.

- For more information about our products for OEM, please see the OEM section in our website.
- See also section: [EmStat Pico Based Instruments](#) on page 347.

20.1 EmStat-based instruments overview

The EmStat instruments are used with electrochemical sensors or electrochemical cells with three or four electrodes: working electrode, reference electrode and auxiliary or counter electrode and sense electrode. In the so-called two-electrode configuration the counter and reference electrodes can be combined to a single electrode. In that case the counter and reference connectors are tied together.

For more information about making a connection to the cell, see section: [Connecting a cell to the potentiostat](#) on page 25.

The EmStats are used in different hardware configurations.

20.1.1 EmStat Basic models



Figure 365 EmStat3 and EmStat4S

The basic models allow for powering and connecting via USB only. They are connected to the cell using a 1-meter cable with a high-end LEMO connector. A single multi-color LED gives an indication of the status of the instrument.

20.1.2 EmStat Blue and EmStat4X



Figure 366 EmStat3 Blue and EmStat4X LR

The EmStat Blue and the newer EmStat4X models allow to connect via USB or Bluetooth. They have an integrated battery and an additional auxiliary port which can be used with a MUX8-R2 multiplexer for example. The connection to the cell is either by means of a standard 1-meter cable with high-end LEMO connector or using the connector compatible with wide collection of screen-printed electrodes (SPEs).

The EmStat4X also allows for iR compensation.

20.1.3 EmStat4R



Figure 367 EmStat4R

The EmStat4R model has a ruggedized design, for use in the field. It allows to connect via USB or Bluetooth and has an integrated battery. The modular front-end can be exchanged for use with either a standard 1-meter cable with high-end LEMO connector or a connector compatible with a wide collection of screen-printed electrodes (SPEs). This instrument is also available for OEM with options for a tailored front-end to match a specific electrode form factor.

20.1.4 Multiplexer-integrated EmStat models



Figure 368 EmStat3-MUX8

The EmStat-modules up to 3 (and 3+) were offered in combination with the MUX8-R2 multiplexer integrated in a single housing as “EmStat MUX”. This instrument is powered and communicates via its USB port. Cell connections are made with eight 1-meter high-end cell cables divided in two sets of 4-cables which connect to two DSUB-25 connectors.

Refer to the respective Potentiostat and/or Multiplexer model(s) for more details.

20.1.5 Multi-channel EmStat models

The EmStats can also be configured as MultiEmStat where a single instrument contains 4, 8 or 12 individual EmStat modules. These instruments can be controlled simultaneously or individually in a single overview with MultiTrace for Windows. Only one channel at a time can be controlled in a PSTrace instance.

21 EmStat1 and EmStat2



Figure 369 EmStat1 (right) and EmStat2 (left)

21.1 Description

The EmStat1 and 2 are a USB-powered potentiostat-only devices. It means that these devices are specifically designed for potentiostatic applications and do not possess the capability to control currents, distinguishing them from galvanostats. The EmStat-series (Embedded potentiostat) are available in a standard shielded housing or as bare PCB. The EmStat communication protocols are well described and made available to allow users to embed the potentiostat as part in another system. The EmStat1 was enclosed in a blue/gray plastic enclosure, while the EmStat2 was enclosed in an aluminium enclosure. The EmStat1 and EmStat2 were sold from 2005 to 2012.

With their reliable communication protocols and versatile form factors, the EmStat-series devices are well-suited for a variety of laboratory and research applications.

21.2 Operating

The EmStat1 and 2 are powered and controlled directly by means of a standard mini-USB cable. They do not have a power switch or any other type of switch. The cell cable connects to a TRIAD connector.

21.2.1 Cell connections

For more information about making a connection to the cell, see section: [Connecting a cell to the potentiostat](#) on page 25.

21.3 Status LED indicators

The EmStat1 and 2 have one or two single color LEDs. The first batch of the EmStat2 was produced without a cell LED. Later batches came with a power and cell LED. The red cell LED is

switched on automatically whenever the feedback loop is closed and current can flow through the cell. The green LED is on when the device is powered.

21.4 Cell connector pin-out

The following schematics show the EmStat1 and EmStat2 cell connector pin-out and pin functions.




Figure 370 Front view of male plug on the cell cable.




Figure 371 Front view of female plug on the EmStat1 and EmStat2.

Table 39 Cell connector pins of the EmStat1 and EmStat2

Pin	Function	Connector color
1	Working Electrode (WE)	Red
2	Analog Ground (AGND)	Green
3	Reserved	N/A
4	Counter / Auxiliary Electrode (CE)	Black
5	Reference Electrode (RE)	Blue

21.5 Specifications

The following tables show the main specifications of the EmStat1 and EmStat2.

Table 40 EmStat1 and EmStat2 specifications

General	
dc-potential range	$\pm 2.000\text{ V}$
compliance voltage	$\pm 4.5\text{ V}$
maximum output current	$\pm 20\text{ mA}$
Potentiostat (controlled potential mode)	
dc-potential resolution	1 mV
applied dc-potential accuracy	$\leq 0.2\%$
max. dc-offset error	2 mV
current ranges	1 nA to 100 μA (6 ranges) 1 nA to 10 mA (8 ranges) for EmStat2

current resolution	0.1 % of current range
	1 pA on lowest current range
current accuracy	≤ 0.3 % of current range at 1 µA to 100 µA
	≤ 0.5 % at 100 nA and ≤ 1 % at 10 and 1 nA
	all with additional 0.2 % offset error

Electrometer	
electrometer amplifier input	> 100 GΩ // 4 pF
rise time	approx. 200 µs

Other	
dimensions	EmStat1: 6.2 cm x 4.6 cm x (1.7 to 2.8 cm) EmStat2: 6.2 cm x 4.8 cm x (1.7 to 2.6 cm)
power	5 V / 60 mA from USB connector (or ac-adapter ¹)
communication	USB

21.6 Offset calibration

When connected to a EmStat1 or EmStat2, a button for calibrating offsets can be found In PSTrace in the menu: 'Tools' → 'Instrument Settings' → button: 'Calibrate Offsets'.

Offset calibration can be done when the measured potential differs more than 0.3% +/- 0.003 V from the specified potential (when the cell is switched on in E control mode).

22 EmStat3 and EmStat3+



Figure 372 EmStat3(+)

22.1 Description

The EmStat3 and EmStat3+ are a USB-powered potentiostat-only devices. This means that the EmStat3 and 3+ are not capable of controlling currents. The EmStat-series (Embedded potentioStat) are available in a standard shielded housing or as bare PCB. The EmStat communication protocols are well described and made available to allow users to embed the potentiostat as part in another system. The EmStat3+ allows for applying and measuring higher voltages and currents, compared to the EmStat3.

The EmStat3 and EmStat3+ were sold from 2013 to 2022.

With their reliable communication protocols and versatile form factors, the EmStat-series devices are well-suited for a variety of laboratory and research applications.

22.2 Operating

The EmStat3 and EmStat3+ instruments are powered and controlled directly by means of a standard mini-USB cable. They do not have a power switch or any other type of switch. The cell cable connects to the high-end LEMO connector.



**The EmStat3/3+ cell connector is a push-pull connector.
Do not try to twist it!**

22.2.1 Sense lead

The Sense lead of the EmStat3+ is only active when the 100 mA range is switched in. In case low current measurements are performed, it is advised to leave the Sense connector disconnected.

22.3 Status LED indicators

The EmStat3 and 3+ have two single color LEDs. The red cell LED is switched on automatically whenever the feedback loop is closed and current can flow through the cell. The blue LED is on when the device is powered.

22.4 Cell connector pin-out

The following schematics show the EmStat3 and EmStat3+ cell connector pin-out and pin functions.



Figure 373 Front view of male plug on the cell cable




Figure 374 Front view of female plug on the EmStat3(+)

Table 41 Cell connector pins of the EmStat3(+)

Pin	Function	Connector color
1	Reference Electrode (RE)	Blue
2	Counter / Auxiliary Electrode (CE)	Black
3	Working Electrode Sense (Sense) (active in 100 mA only)	Yellow
4	Working Electrode (WE)	Red
Connector housing	Analog Ground (AGND)	Green

22.5 Specifications

The following tables show the main specifications of the EmStat3 and EmStat3+.

Table 42 EmStat3 and EmStat3+ specifications

	EmStat3	EmStat3+
dc-potential range	± 3.000 V	± 4.000 V
compliance voltage	± 5 V	± 8 V
applied dc-potential resolution	0.1 mV	0.125 mV
applied potential accuracy	≤ 0.2 % with max. 2 mV offset error	≤ 0.3 % with max. 3 mV offset error
current ranges	1 nA to 10 mA (8 ranges)	1 nA to 100 mA (9 ranges)

maximum measured current	± 20 mA typical and ± 15 mA minimum	± 100 mA typical
---------------------------------	--	----------------------

EmStat 3 and 3+ Potentiostat (controlled potential mode)	
current resolution	0.1 % of current range
	1 pA on lowest current range
current accuracy	≤ 1 % of current range at 1 nA
	≤ 0.5 % at 10 nA
	≤ 0.2 % at 100 nA to 100 μ A
	≤ 0.5 % at 1 mA, 10 mA and 100 mA
	all with additional 0.2 % offset error

Electrometer / other	
electrometer amplifier input	> 100 Gohm // 4 pF
rise time	approx. 200 μ s

EmStat 3 and 3+ regular model (not Blue) other specifications	
sensor connection	shielded cable with circular connector for WE, RE, CE and Sense (100 mA for ES3+ only)
housing	anodized aluminium: 6.7 cm x 5.0 cm x (1.9 to 2.8 cm)
weight	85 g
power supply	5 V, min. 130 mA (ES3) or 500 mA (ES3+) via USB
communication	USB

The EmStat3 Blue and EmStat3+ Blue models have an auxiliary port to control external accessories. See next section for more information.

22.6 Offset calibration

When connected to a EmStat3 or EmStat3+, a button for calibrating offsets can be found In PSTrace in the menu: 'Tools' → 'Instrument Settings' → button: 'Calibrate Offsets'.

Offset calibration can be done when the measured potential differs more than 0.3% +/- 0.003 V from the specified potential (when the cell is switched on in E control mode).

23 EmStat3 Blue and EmStat3+ Blue



Figure 375 EmStat3(+) Blue

23.1 Description

The EmStat3 Blue and 3+ Blue incorporate the EmStat3 and 3+ potentiostat modules, respectively, with additional hardware for supporting Bluetooth connections, an auxiliary port and battery charging.

The EmStat3 Blue and EmStat3+ Blue were sold from 2015 to early 2023.

23.2 Operating

The EmStat3 Blue and 3+ Blue are battery powered and therefore have a power button. Push and hold the button for two seconds to turn it on. Push and hold for 4 seconds to turn it off. Connection to the EmStat3 and 3+ Blue can either be made by means of the mini-USB port or by means of a Bluetooth connection. Press the power button briefly to turn Bluetooth off or to turn Bluetooth on again.

The name of the Bluetooth device consists of the letters PS and the last four characters of the MAC address of the device. For example: PS-CA6F. The MAC address can be found on the battery lid of the instrument.

The cell cable connects to the high-end LEMO connector.



The EmStat3 and 3+ Blue cell connector is a push-pull connector. Do not try to twist it!

23.2.1 Cell connections and Sense lead

The Sense lead of the EmStat3+ is only active when the 100 mA range is switched in. In case low current measurements are performed, it is advised to leave the Sense connector disconnected.

See also section: [Connecting a cell to the potentiostat](#) on page 25.




23.2.2 Connecting

When connecting to the USB port of a PC, Bluetooth will switch off automatically unless there is Bluetooth a connection active. This allows the instrument to be charged via a PC while a Bluetooth connection is present.

23.3 Status LED indicators

The EmStat3 Blue and 3+ Blue has a keypad with icons and a power button and a red cell LED. The follow table explains the different states for the indicators.

Table 43 Indicators on the EmStat3(+) Blue

	Power button Steady: instrument is on
	Bluetooth icon Blinking: ready to connect Steady: connected
	Battery icon Blinking red: battery low Blinking green: charging battery Steady green: battery fully charged
Cell	Steady red: cell is on

23.4 Cell connector pin-out

The following schematics show the EmStat3 Blue and EmStat3+ Blue cell connector pin-out and pin functions.



Figure 376 Front view of male plug on the cell cable




Figure 377 Front view of female plug on the EmStat3(+) Blue

Table 44 Cell connector pins of the EmStat3(+) Blue

Pin	Function	Connector color
1	Reference Electrode (RE)	Blue
2	Counter / Auxiliary Electrode (CE)	Black

3	Working Electrode Sense (Sense) (active in 100 mA only)	Yellow
4	Working Electrode (WE)	Red
Connector housing	Analog Ground (AGND)	Green

23.4.1 SPE connector

An additional connector for Screen-Printed Electrodes is present at the front of the instrument. This connector supports the most common Screen-Printed electrodes that meet the following specifications:

Table 45 Specifications of compatible SPEs

Connection pads	RE, WE, CE
Pitch	2.54 mm (0.1")
Sensor thickness	0.1 to 0.8 mm
Max. sensor width	11 mm



Ceramic substrates can have very sharp edges that might damage the SPE connector.

23.5 Auxiliary port pin-out

The EmStat3 Blue and EmStat3+ Blue auxiliary connector has the following pin-out and pin functions.




Figure 378 The front view of the female port (DE-15)

Table 46 EmStat3(+) Blue auxiliary port pins

Pin	Function	EmStat PCB pin (CON-PIN)
1	d0 digital output (5V)	1-4
2	d3 digital output (5V)	1-1
3	auxiliary analog input (0 – 4.095 V)	2-9
4	Rx (TTL comm)	2-2
5	Tx (TTL comm)	2-1
6	d1 digital output (5V)	1-3
7	d0 digital input (5V)	1-4
8	RESERVED	-

9	5 V digital power line (max. 50mA)	2-12
10	digital ground	2-11
11	d2 digital output (5V)	1-2
12	RESERVED	-
13	RESERVED	-
14	analog ground	1-9
15	DAC out (0-3V), 0.5Ω typical	2-10
Connector housing	digital ground	-

23.6 Specifications

The EmStat Blue models have similar specifications as the regular EmStat3 and 3+, but they are extended with Bluetooth connectivity, a battery and auxiliary port.

See: [EmStat3 and EmStat3+ specifications](#) on page 320

Table 47 *EmStat 3 Blue and 3+ Blue specifications*

housing	anodized aluminium: 100 mm x 60 mm x (27 to 34 mm)
weight	85 g
power supply	USB or internal Li-Po battery 5 V, min. 130 mA (ES3) or 500 mA (ES3+)
battery time	> 8 hours with cell off > 5 hours with continuous cell on at 1uA current full charge takes approx. 3 hours
communication	USB or Bluetooth

24 EmStat4S



Figure 379 EmStat4S LR (left) and EmStat4S HR (right)

24.1 Description

The EmStat4S LR and HR are a USB-powered potentiostat / galvanostat and impedance analyzers. The EmStat-series (Embedded potentioStat) are available in a standard shielded housing or as bare PCB. The EmStat communication protocols are well described and made available to allow users to embed the potentiostat as part in another system.

With their reliable communication protocols and versatile form factors, the EmStat-series devices are well-suited for a variety of laboratory and research applications.

The EmStat4S LR and HR are the market since 2021. They are the successor of the EmStat3 and EmStat3+ models.

The EmStat4S comes in two different versions; the Low Range (LR) and High Range (HR) version.

- Low Range: for lower currents and potentials. This version can be identified by its blue bottom shell.
- High Range: for higher currents and potentials. This version can be identified by its black bottom shell.

24.2 Operating

The EmStat4S is powered and controlled directly by means of a standard USB-C cable. They do not have a power switch or any other type of switch. The cell cable connects to the high-end LEMO connector.



**The EmStat4S cell connector is a push-pull connector.
Do not try to twist it!**

24.2.1 USB Power Considerations for EmStat4S HR

During measurements in the 100 mA range, the EmStat4S HR may demand more current than what is officially provided by standard USB 2.0 ports (5V at 500 mA). In cases where the available USB port cannot deliver 5V at 900 mA (or more), or if there is uncertainty about the power output of the USB port, it is essential to utilize the provided USB Y-cable with an additional USB port. This special cable configuration allows the instrument to draw additional power from a second USB port, ensuring optimal performance and accurate measurements.

24.2.2 Sense lead

The EmStat4 HR version comes with an additional Sense lead. In contrast with the EmStat3+, the Sense lead of the EmStat4S HR is active for every current range and should always be connected. This allows the Sense lead to be used in four-electrode measurements.

For more information about making a connection to the cell, see section: [Connecting a cell to the potentiostat](#) on page 25.

24.3 Status LED indicator

The EmStat4S has a multi-color status indicator. The following table explains the meaning for each color.

Table 48 *EmStat4S status led indicator*

Color	State
Steady green	Booting
Steady blue	Cell off (at open circuit)
Steady red	Cell on
Blinking green	Uploading new firmware
Blinking orange	Error state
Blinking red	Fatal error (not responding)

24.4 Cell connector pin-out

The EmStat4S LR and HR Lemo cell connector has the following pin-out and pin functions.




Figure 380 *Front view of male plug on the cell cable*




Figure 381 *Front view of female plug on the EmStat4S*

Table 49 Cell connector pin functions



Pin	Function	Connector color
1	Reference Electrode (RE)	Blue
2	Reference Electrode Shield	N/A
3	Counter / Auxiliary Electrode (CE)	Black
4	Working Electrode Sense (Sense)	Yellow
5	Working Electrode (WE)	Red
Connector housing	Analog Ground (AGND)	Green

24.5 Specifications

The Emstat4S comes in two versions: LR (Low Range) and HR (High Range). The version of your EmStat4S can easily be recognized by the color of the bottom shell. The HR as a black bottom shell, and the LR has a blue bottom shell.

The following table shows the main differences between the LR and HR versions.

Table 50 EmStat4S LR and EmStat4S HR specifications

		
	EmStat4S LR	EmStat4S HR
Potential range	±3 V	±6 V
Max. compliance voltage	±5 V	±8 V
Current ranges	1 nA to 10 mA (8 ranges)	100 nA to 100 mA (7 ranges)
Max. current	±30 mA	±200 mA
Electrode connections	WE, RE, CE, and ground, 2 mm banana plugs	WE, RE, CE, S, and ground, 2 mm banana plugs

See the EmStat4S product page on our website for more detailed specifications.

25 EmStat4X



Figure 382 The EmStat4X LR (left) and EmStat4X HR (right)

25.1 Description

The EmStat4X LR and HR are on the market since 2023. They are the successor of the EmStat3 Blue and 3+ Blue. Like the EmStat4S, the EmStat4X comes in two different versions: the Low Range (LR) and High Range (HR) version.

- Low Range: for lower currents and potentials. This version can be identified by its blue bottom shell.
- High Range: for higher currents and potentials. This version can be identified by its black bottom shell.

25.2 Operating

The EmStat4X is powered via USB-C or by its internal battery. The cell cable connects to the high-end LEMO connector.

25.2.1 Cell connections and Sense lead

The EmStat4 HR version comes with an additional Sense lead. In contrast with the EmStat3+, the Sense lead of the EmStat4S HR is active for every current range and should always be connected. This allows the Sense lead to be used in four-electrode measurements.

For more information about making a connection to the cell, see section: [Connecting a cell to the potentiostat](#) on page 25.



**The EmStat4X cell connector is a push-pull connector.
Do not try to twist it!**

25.2.2 Switching on and off

A single touch button next to the display can be used to switch the instrument on and off again. Put your finger on the button (no force needed) and hold it for a few seconds to turn in the instrument on. Do the same for turning the instrument off again.

25.3 USB power and charging

The EmStat4X can be used with any USB port. However not all USB-ports can provide enough power to prevent discharging the battery when using the EmStat4X HR at the highest current range. In order to prevent the EmStat4X HR battery to discharge when using the 100 mA range, please make sure to use the USB-C cable as provided connected to a USB-C port of the PC.



**The EmStat4X requires to be connected to a
USB-C charger or host for charging with the maximum current of
1A.**

When the EmStat4X is connected to a USB-C charger or host, charging the battery from depleted to full takes less than 3 hours. With the instrument being powered off or connected to another type of USB port, the charge current is maximum 500 mA.

25.4 Display and status LED indicator

The EmStat4X has a small display and a multi-color status indicator. The following table explains the meaning for each color.

Table 51 *EmStat4X status led indicator*

Color	State
Steady green	Booting
Steady blue	Cell off (open circuit)
Steady red	Cell on
Blinking green	Uploading new firmware
Blinking orange	Error state
Blinking red	Fatal error (not responding)

25.4.1 Display

The EmStat4X display shows the following information:

- Battery status
- Bluetooth identifier
- Connection status
- Firmware version during boot

25.5 iR Compensation

The EmStat4X iR Compensation module works using Positive Feedback. This is achieved using a 12-bit MDAC in the module which scales the output of the current follower opamp to provide a positive feedback voltage that is proportional to the current through the cell. The compensation voltage is added to the summing point before the control amplifier and thus increases the applied potential to counteract the iR drop.




Figure 383 Simplified schematical representation of iR compensation circuitry in the EmStat4X

Positive feedback allows for fast scan rates up to 10 V/s, depending on the characteristics of the cell. If the potential error to compensate for becomes close to the value set for E applied, the system might become unstable. Using iR compensation limits the measurement bandwidth to 160 kHz.

See also section: [Ohmic \(iR\) compensation](#) on page 46.

25.6 Auxiliary port pin-out

The following schematics and table show the EmStat4X auxiliary port pin-out and pin functions.




Figure 384 The front view of the female port (DE-15)

Table 52 EmStat4X auxiliary port pins

Pin	Function
1	d0 digital output (3.3 V) *
2	d3 digital output (3.3 V) *
3	auxiliary analog input -10 to +10 V, 16 bit, >0.5 MOhm input impedance
4	RESERVED
5	RESERVED
6	d1 digital output (3.3 V) *
7	d0 digital input (3.3 V) *
8	i monitor given as V in active current range. (-5 to +5 V for LR and -8 to +8 V for HR) The following ranges give an output that is divided by 10. LR: 1 mA, 10 uA, 100 nA, 1 nA HR: 1 mA, 10 uA, 100 nA
9	5V digital power line (max. 300mA)
10	digital ground
11	d2 digital output (3.3 V) *
12	d4 digital input / output (3.3 V) * configurable through MethodSCRIPT
13	E monitor, scaled to a voltage between 0 and 3 V, centered around 1.5 V, where: LR: $VRE = -2.1 * (VAUX - 1.5)$ HR: $VRE = -4.2 * (VAUX - 1.5)$
14	analog ground
15	auxiliary analog out (0 to 6 V at 12 bit)
Connector housing	digital ground

* All digital outputs have a 470 Ohm series resistance.

25.7 Cell connector pin-out

The EmStat4X LR and HR Lemo cell connector has the following pin-out and pin functions.




Figure 385 Front view of male plug on the cell cable




Figure 386 Front view of female plug on the EmStat4X

Table 53 Cell connector pin functions

Pin	Function	Connector color
1	Reference Electrode (RE)	Blue
2	Reference Electrode Shield	N/A
3	Counter / Auxiliary Electrode (CE)	Black
4	Working Electrode Sense (Sense) *	Yellow
5	Working Electrode (WE)	Red
Connector housing	Analog Ground (AGND)	Green

* Available for EmStat4X HR only

25.8 SPE connector

An additional connector for Screen-Printed Electrodes is present at the front of the instrument, for the version LR only. This connector supports the most common Screen-Printed electrodes that meet the following specifications:

Table 18 Specifications of compatible SPEs

Connection pads	RE, WE, CE
Pitch	2.54 mm (0.1")
Sensor thickness	0.1 to 0.8 mm
Max. sensor width	10 mm



Ceramic substrates can have very sharp edges that might damage the SPE connector.

25.9 Specifications

The EmStat4X LR and HR models have similar specifications as the regular EmStat4S LR and HR, but they are extended with a display, Bluetooth connectivity, a battery and auxiliary port.

See the EmStat4X product page on our website for more detailed specifications.

Table 54 *EmStat4X specifications*

EmStat4X specific		
model	LR	HR
▪ electrode connections	WE, RE, CE, and ground, with 2 mm banana plugs	WE, RE, CE, S and ground, with 2 mm banana plugs
▪ power consumption	Typical: 1W (idle) Max: 1.6W (cell on at 30 mA)	Typical: 1.5W (idle) 1.6W (cell on at 10 mA) Max: 4.6W (cell on at 200 mA)
▪ power source	USB-C or internal LiPo battery	
▪ battery	11.1 Wh capacity 80% charge in 2.5 hours, full charge in 3 hours.	
▪ communications	USB-C or Bluetooth (4.0 - Dual-Mode)	
▪ housing	aluminum body: 11.4 x 8.0 x 4.5 cm	
▪ weight	+/- 500 g	
▪ internal storage space	500 MB, equivalent to >15M datapoints	

26 Sensit series



Figure 387 The Sensit Smart (left) and Sensit BT readers (right) for use with Screen-Printed Electrodes

26.1 Description

The Sensit series comprises two models: Sensit Smart and Sensit BT. Both models are equipped with the EmStat Pico module, sharing identical main specifications and available techniques.

26.1.1 Sensit Smart

This model is a straightforward configuration featuring the EmStat Pico module housed with USB-C and SPE connectors. It serves as a compact, ready-to-use potentiostat that can be directly inserted into a smartphone or tablet. Control is achieved through the Android app PStouch. Alternatively, the Sensit Smart can be connected to a classic USB port on a PC using the USB-C Female to USB-A cable, allowing control through our PC software PSTrace.

26.1.2 Sensit BT

This model includes all features of Sensit Smart and extends functionality with additional components. It incorporates a battery, Bluetooth connection, internal memory, power switch, trigger button, and indicative lights. The Sensit BT connects to smartphones or tablets via Bluetooth, controlled through the Android app PStouch. Charging or PC connection is facilitated through the USB-C port, enabling control via PC software PSTrace.

26.2 Operating the Sensit Smart

The Sensit Smart is communicating and is powered via the USB-C connector


The LEDs indicate in which state the device is in. The follow table shows which LED behaviour is linked to which device state.

Table 55 Sensit Smart status LEDs indicators

State	Blue LED	Red LED
Powered and cell off	ON	OFF
Cell on	ON	ON
Fatal error	Blinking 5x per second	Blinking 5x per second
Warning	Blinking 1x per second	Blinking 1x per second

The Sensit Smart can be used with Screen-Printed Electrodes with the following properties:

Table 56 Specifications for compatible SPEs

Connection pad pitch	2.54 mm
Electrode connections	
Allowed sensor thickness	between 0.1 mm and 0.8 mm
Maximum sensor width	11 mm

For more information about making a connection to the cell, see section: [Connecting a cell to the potentiostat](#) on page 25.

26.3 Operating the Sensit BT



The Sensit BT comes in two different versions: SPE and SNS.

The Sensit BT.SPE can be used for running sequential measurements on two different Screen-Printed Electrodes (SPE's) each with their own Reference, Counter and Working electrodes. The second channel can also be used in Bipotentiostat mode, functioning as second Working Electrode versus the Reference and Counter electrode of channel 1. Both Working electrodes are recorded simultaneously in the Bipotentiostat mode. A connector for a SPE with 2WEs is available under request.

The Sensit BT.SNS has a lead connected to the WE of channel 2 and can be used out-of-the-box for BiPotentiostat measurements. The second Working Electrode (WE2) can either be set at a

potential offset with respect to WE1 or at a fixed potential with respect to RE1. The Bipotentiostat mode is supported in Low Speed mode for all techniques, excluding EIS and OCP.

Table 57 Sensit BT specifications

Sensit BT.SPE		Sensit BT.SNS	
			
Sensor pitch	2.54 mm	Cable length	40 cm
Electrode connections	2x CE, WE, RE	Connectors	2 mm banana
Allowed sensor thickness	between 0.1 and 0.8 mm	Electrode connections	RE, WE, WE2, CE
Maximum sensor width	11 mm		

26.4 Sensit BT trigger button

The Sensit BT trigger button is not supported natively in PSTrace, but can be used with custom MethodSCRIPTs. The button can be configured in such a way that each push starts a measurement and stores the measured data on the on-board storage.




Figure 388 Trigger button on Sensit BT

See also section: [MethodSCRIPT Sandbox](#) on page 107.

See our MethodSCRIPT trigger button example found on:

www.github.com/palmsens.

26.5 Sensit Smart and Sensit BT specifications

The following table shows the specifications of the Sensit Smart and Sensit BT for each operational mode.

Table 58 Sensit Smart and Sensit BT specifications

	Low Speed Mode	High Speed Mode	Max Range Mode
dc-potential range	-1.2 to +2 V	-1.7 to +2 V	-1.7 to +2 V
Dynamic dc-potential range *	2.2 V	1.2 V	2.6 V
Compliance voltage	-2.0 to 2.3 V		
Maximum current	± 3 mA		
Max. acquisition rate (datapoints/s)	100	1000	100
Supports FRA/EIS	NO	YES	NO

* The dynamic range is the range that can be covered during a single scan within the full potential range. For example, a linear scan can start at -1.5 V and end at 1.1 V or vice versa, covering 2.6 V dynamic range.

See the corresponding product pages on our website for detailed technical specifications.

26.5.1 Compliance voltage limitation

The Sensit series are based on the EmStat Pico modules which have a limited compliance voltage. The compliance voltage is the maximum voltage that can be applied between the working and counter electrode. Another name could be the maximum cell potential.



The Sensit series are based on the EmStat Pico modules which have a limited compliance voltage.

The higher the current at the working electrode, the higher the current required at the counter electrode. To increase the current the potentiostat increases the potential difference between the working and counter electrode. At some point, the maximum cell potential is reached and cannot be increased further.

This means the working electrode's potential might not be correct anymore, because the correct current cannot flow through the counter electrode anymore. The potentiostat ran into its compliance voltage. The potential of the working electrode will stay at the last value it could reach.

During DPV or SWV this means the differential is 0 and thus the curve shows a sudden drop in current and stays at 0. In a Cyclic Voltammogram reaching the compliance voltage is clearly visible. The constant potential will lead to a current following the Cottrell equation:




Figure 389 Typical CV showing an exceeded maximum compliance voltage

26.6 Calibration

The Sensit series has an internal automatic calibration routine which can be started in PSTrace. When connected to a Sensit-series device or EmStat Pico-based device, you can navigate via the menu: 'Tools' → 'Instrument Settings...' → button 'Calibrate'.



Your device is already calibrated. It is not recommended to run the calibration procedure, unless specifically advised by PalmSens BV.

See also section: [EmStat Pico Calibration](#) on page 347.

27 EmStat Pico based instruments



Figure 390 EmStat Pico potentiostat module

The EmStat Pico is a small potentiostat module which is available as separate product but is also at the heart of housed instruments like the Sensit Smart and Sensit BT. EmStat Pico modules with correct license can be used with PSTrace.

The EmStat Pico is built together with our partner Analog Devices Inc.

27.1 EmStat Pico specifications

See the OEM section on our website for detailed EmStat Pico specifications.

27.2 Calibration

The EmStat Pico (and Sensit-series) has an internal automatic calibration routine which can be started in PSTrace. When connected to a Sensit-series device or EmStat Pico-based device, you can navigate via the menu: 'Tools' → 'Instrument Settings...' → button 'Calibrate'.



Your device is already calibrated. It is not recommended to run the calibration procedure, unless specifically advised by PalmSens BV.

The following parts of the signal chain will be calibrated:

- ADC offset and gain;
- LPTIA gain (100 nA is calibrated externally);
- HSTIA offset and gain (100 nA and 1 uA are calibrated externally);
- HSDAC offset and gain.

The auto-calibration will not affect the 1 μ A and 100 nA calibrations values, which require an external reference. The complete procedure takes around 60 seconds.

28 Troubleshooting

28.1 Running self-diagnostics

The PalmSens4 has a built-in diagnostics tool which can be executed in PSTrace. This tool generates a detailed report showing any issues which might be present in the hardware.

See for more information section: [Self Diagnostics](#) on page 307.

28.2 Verifying your potentiostat

Your instrument can be tested by using the test sensor or dummy cell supplied with the instrument.




Figure 391 PalmSens Dummy Cell supplied with the EmStat4R SNS version.

The easiest way to verify the functioning of your instrument is to use the “WE B” circuit, which consists of a resistor with a value of 10 kΩ with a max deviation of 0.1%.

In case you are using the EmStat4R SPE version, connect the SPE dummy cell. Keep in mind that this SPE has a simplified Randles Circuit with an extra 560 Ω resistor in series with the 10 kΩ resistor.



Figure 392 The SPE dummy cell supplied with the EmStat4R SPE version.

The WE lead is connected to one side and both RE and CE to the other side of the resistor.



If your instrument has a Sense lead, make sure to connect the Sense lead to the WE lead as well.

Any of the electrochemical techniques can be applied. The current response obtained with a resistor with value R is equal to the applied potential or potential pulse divided by the value of R . So, if a potential of 0.5 V is applied on a resistor of 10 k Ω , the obtained current should be $0.5 \text{ V} / 10 \text{ k}\Omega = 50 \text{ }\mu\text{A}$.

Contact PalmSens BV if the problems are found: info@palmsens.com and report the problems as detailed as possible.

28.3 Noise

Our instruments are designed with hardware noise suppression filters to reject noise from internal and external sources. If a higher level of noise is your issue, the solving strategies are rather numerous, but the sources for noise are also numerous. Here we describe the most successful and common methods for noise reduction.

To determine the noise levels for your instrument, please refer to section “Measuring the noise level of the instrument” of the PSTrace Manual.

28.3.1 Power grid

Your power grid is usually using an alternating current. This undulating current influences the measured currents. PSTrace and PStouch have a filter for this mains frequency. In PSTrace, check in the ‘Tools’ menu under ‘General Settings’ if the mains frequency is set correctly.

28.3.2 Electrical fields

Our environment is filled with electrical fields. Some of them are created by devices around us as side effects or in case of wireless communication on purpose. Although it is a bad idea to measure directly next to an electric arc furnace, it is usually not possible to have a workspace free of electrical fields, especially not during point-of-care measurements. A Faraday cage is usually sufficient to create a field-free environment. A metal box or cage out of metal mesh is a good Faraday cage. Even a shield out of aluminum foil can help. Place your electrochemical cell inside the Faraday cage and connect the cage to the ground lead (green) of the potentiostat.

The cable delivered with your EmStat or PalmSens has an inbuilt shield and should protect your signal outside the Faraday cage. This is one of the most effective methods to reduce noise.

28.3.3 Cables

Cables should not be unnecessarily long, since they act as antennas for noise, but the cable delivered with your EmStat or PalmSens has an inbuilt shield and as long as you use the original cable, there is little reason to worry about cable induced noise.

28.3.4 Grounding

Ground your measurement equipment. The best way to connect your equipment is star-shaped, that is all parts are connected with the ground at the same point. In an electrochemical lab that point is usually one small space of the faraday cage. This way earth loops that induce noise are avoided.

28.3.5 Contacts

Check if the contacts are corroded. If so, remove the stains, for example with sandpaper.

29 Updating firmware

The firmware which runs on your instrument might need to be updated when a new version of PSTrace is installed. PSTrace will prompt a message and take care of the firmware update if this is necessary.

29.1 Firmware update

If a specific version of the firmware should be downloaded the firmware can be updated using the firmware update window found in PSTrace via the menu 'Tools' → 'Instrument settings'...

Normally this should not be necessary because PSTrace will prompt a message and take care of the firmware update if this is necessary.




Figure 393 Instrument settings menu

The 'Update firmware' button is found at the right top of the window:




Figure 394 Instrument settings window

1. If not already connected; connect the instrument to the PC and in case applicable for your instrument; please make sure to connect the ac-adapter as well.



It is essential that the instrument is not switched off or disconnected during the process of updating the firmware.

In case of a PalmSens1 or PalmSens2:

Be sure that the batteries are full or that PalmSens is connected to the ac-adapter.

1. Now click the 'Browse...' button and select the appropriate HEX file.
2. Now press the 'Start update firmware' button.
3. In case of a PalmSens1 or PalmSens2 only: The window now tells you to press the four keys 'Enter + Esc + < + >' simultaneously.

DO THIS ONLY ONCE! The updating process will start immediately after pressing the four keys.




Figure 395 Instruction showing in case of a PalmSens1 or PalmSens2 to press the four keys simultaneously.

4. If the instrument is in download mode, the status bar of the firmware update is shown. If this screen does not appear, the program can be closed and started again.
5. After completing the update, the program will show "UPDATE COMPLETE".

30 Maintenance and compliance

30.1 Lithium-Polymer batteries

The typical estimated life of a Lithium-Polymer battery is about two to three years or 300 to 500 charge cycles, whichever occurs first. One charge cycle is a period of use from fully charged, to fully discharged, and fully recharged again. Consider a two-to-three-year life expectancy for batteries that do not run through complete charge cycles.

Rechargeable Lithium-Polymer batteries have a limited life and will gradually lose their capacity to hold a charge. This loss of capacity (aging) is irreversible. As the battery loses capacity, the length of time it will power the product (run time) decreases. Lithium-Ion batteries continue to slowly discharge (self-discharge) when not in use or while in storage.

30.2 Temperature compliance

Our instruments are designed for indoor use at ambient temperatures between 0 °C and 45 °C. All the components of PalmSens products (except their batteries) are rated to the industrial temperature standard of -40 °C to +85 °C.

The battery of the potentiostat is rated -20 °C to +60 °C when discharging 0 °C to +45 °C when charging.

30.3 Humidity compliance

PalmSens instruments have not been tested in high humidity environments.

Elevated humidity however may cause measurement errors if condensation forms on the electronics. This affects measurements in the nA ranges or lower. Prolonged exposure to a condensing environment may severely decrease the life expectancy of the instrument and void its warranty.

30.4 Temperature drift

PalmSens instruments are calibrated at 21 °C. The most sensitive components of the instrument have temperature drift of 50 ppm. For instance at 1 °C or 41 °C, measurement drift of up to 0.1% may be experienced.

30.5 Atmospheric pressure

PalmSens instruments are not intended for use in safety-critical applications. Consequently, the power supplies utilized are not selected based on a specific pressure rating.

30.6 Cleaning

Make sure to disconnect your instrument from any cell or power source, if applicable, prior to cleaning. Use a cloth lightly dampened with either clean water or water containing a mild detergent to clean the outside of the instrument. Alternatively, you can use isopropyl alcohol. Avoid using a wet rag and prevent any fluids from entering the instrument. It is crucial not to immerse the instrument in any cleaning solution.

30.7 Periodic calibration and preventive maintenance

PalmSens instruments are designed in a way that eliminates the need for periodic calibration. While not mandatory, PalmSens does provide a calibration service for users with specific demands such as QC/ISO purposes. This service includes a new calibration certificate.

It's important to note that PalmSens instruments do not require preventive maintenance, further simplifying their use and reducing the overall maintenance demands on users.

30.8 Service and repair

Except for the battery in some models, your PalmSens instrument contains no user-serviceable parts internally. Any service or maintenance needs should be directed to a qualified service technician employed by PalmSens BV. Attempting to access or modify internal components without proper expertise may result in additional damage to the instrument and void warranties. It is recommended to rely on authorized service personnel for any required maintenance or repairs.

30.9 RoHS Compliance

All instruments from PalmSens have been built using lead free components and lead-free solder. They are in compliance with the European RoHS initiative.

Appendix A. CE Declarations of conformity



EU DECLARATION OF CONFORMITY

Certificate number: PSDOC-PS4-B



Manufacturer: PalmSens BV
Address: Vleugelboot 22,
3991 CL Houten,
The Netherlands

This declaration is valid for the following product:

PalmSens4 (PS4 V1.4 and higher)

Portable electrochemical analyser

- USB power and communications
- Battery power
- Bluetooth communication
- Screen display
- 1 meter cell cable

The object of the declaration described above is in conformity with the Radio Emissions Directive 2014/53/EU (RED) and applicable standards listed below:

Health & Safety

- EN 62479
- EN 61010-1

EMC

- EN 61326-1
- EN 301 489 parts 1 & 17

Efficient Usage of Radio Spectrum

- EN 300 328

This declaration is issued under the sole responsibility of PalmSens BV.

Date: 29th of November 2023

A handwritten signature in black ink, appearing to read 'C.J. van Velzen', with a long horizontal stroke extending to the right.

C.J. van Velzen, CTO

**EU DECLARATION OF CONFORMITY**

Certificate number: PSDOC-ES4S-B



Manufacturer: PalmSens BV
Address: Vleugelboot 22,
3991 CL Houten,
The Netherlands

This declaration is valid for the following product:

EmStat4S (HR/LR), Portable electrochemical analyser

- USB power and communications
- 1 meter cell cable

The object of the declaration described above is in conformity with the Electromagnetic-Compatibility Directive 2014/30/EU (EMCD) and applicable standards listed below:

EMC

- EN 61326-1

This declaration is issued under the sole responsibility of PalmSens BV.

Date: 29th of November 2023

A handwritten signature in black ink, appearing to read "C.J. van Velzen", with a long horizontal stroke extending from the end.

C.J. van Velzen, CTO



EU DECLARATION OF CONFORMITY

Certificate number: PSDOC-ES4R-B



Manufacturer: PalmSens BV
Address: Vleugelboot 22,
3991 CL Houten,
The Netherlands

This declaration is valid for the following product:

EmStat4R: (V1.0 and higher) Portable electrochemical analyser.

- USB power and communications
- Battery power
- Bluetooth communication

The object of the declaration described above is in conformity with the Radio Emissions Directive 2014/53/EU (RED) and applicable standards listed below:

Health & Safety

- EN 62479
- EN 61010-1

EMC

- EN 61326-1
- EN 301 489 parts 1 & 17

Efficient Usage of Radio Spectrum

- EN 300 328

This declaration is issued under the sole responsibility of PalmSens BV.

Date: 29th of November 2023

A handwritten signature in black ink, appearing to read "C.J. van Velzen", with a long horizontal stroke extending from the end.

C.J. van Velzen, CTO

**EU DECLARATION OF CONFORMITY**

Certificate number: PSDOC-ES4X-B



Manufacturer: PalmSens BV
Address: Vleugelboot 22,
3991 CL Houten,
The Netherlands

This declaration is valid for the following product:

EmStat4X, Portable electrochemical analyser

- USB power and communications
- Battery power
- Bluetooth communication
- Screen display
- 1 meter cell cable

The object of the declaration described above is in conformity with the Radio Emissions Directive 2014/53/EU (RED) and applicable standards listed below:

Health & Safety

- EN 62479
- EN 61010-1

EMC

- EN 61326-1
- EN 301 489 parts 1 & 17

Efficient Usage of Radio Spectrum

- EN 300 328

This declaration is issued under the sole responsibility of PalmSens BV.

Date: 29th of November 2023

A handwritten signature in black ink, appearing to read "C.J. van Velzen", with a long horizontal stroke extending to the right.

C.J. van Velzen, CTO



EU DECLARATION OF CONFORMITY

Certificate number: PSDOC-SSMT-B



Manufacturer: PalmSens BV
Address: Vleugelboot 22,
3991 CL Houten,
The Netherlands

This declaration is valid for the following product:

Sensit Smart, Portable electrochemical analyser

- USB power and communications

The object of the declaration described above is in conformity with the Electromagnetic-Compatibility Directive 2014/30/EU (EMCD) and applicable standards listed below:

EMC

- EN 61326-1

This declaration is issued under the sole responsibility of PalmSens BV.

Date: 29th of November 2023

A handwritten signature in black ink, appearing to read "C.J. van Velzen", with a long horizontal stroke extending to the right.

C.J. van Velzen, CTO

**EU DECLARATION OF CONFORMITY**

Certificate number: PSDOC-SBT-B



Manufacturer: PalmSens BV
Address: Vleugelboot 22,
3991 CL Houten,
The Netherlands

This declaration is valid for the following product:

Sensit BT, Portable electrochemical analyser

- USB power and communications
- Battery power
- Bluetooth communication
- 1 meter cell cable

The object of the declaration described above is in conformity with the Radio Emissions Directive 2014/53/EU (RED) and applicable standards listed below:

Health & Safety

- EN 62479
- EN 61010-1

EMC

- EN 61326-1
- EN 301 489 parts 1 & 17

Efficient Usage of Radio Spectrum


- EN 300 328

This declaration is issued under the sole responsibility of PalmSens BV.

Date: 29th of November 2023

A handwritten signature in black ink, appearing to read "C.J. van Velzen", with a long horizontal stroke extending to the right.

C.J. van Velzen, CTO

Appendix B. EU Waste Electrical and Electronic Equipment (WEEE) Directive

The pictogram shown above, located on the product(s) and / or accompanying documents means that used electrical and electronic equipment (WEEE) should not be mixed with general household waste. For proper treatment, recovery and recycling, please take this product(s) to designated collection points where it will be accepted free of charge.

Alternatively, in some countries, you may be able to return your products to your local retailer upon purchase of an equivalent new product. Disposing of this product correctly will help save valuable resources and prevent any potential negative effects on human health and the environment, which could otherwise arise from inappropriate waste handling. Please contact your local authority for further details of your nearest designated collection point. Penalties may be applicable for incorrect disposal of this waste, in accordance with your national legislation.

31 Index

.pssession 205, 235
 ac-voltammogram 126
 Add notes 184
 amperometric measurement 32
 Analog In 209
 Analog Out 209
 Analytical mode 214
 Auto peak search window 214
 auto ranging 34
 auto-ranging 34
 autosave 36
 Autosave settings window 36
 Auxiliary port (PalmSens) 283
 auxilay input 8
 backlight 289
 BiPot 8, 44
 Blank 45
 blank curve 35
 Bluetooth 14
 Calibration curve 213, 217
 capacitor 153
 CDC circuit editor 140
 Cell On 209
 Cell volume 214
 chi-squared 165
 Chronopotentiometric Stripping 96
 circuit components 137
 circuit editor 137

Circuit menu 145
 Constant Phase element 153
 Control E or i 209
 Corrosion 223
 Corrosion current 223
 Corrosion potential 223
 Corrosion rate 223
 CSV 205
 Current Level or Peaks 39
 current range 34, 117, 209
 Curve calculations 184, 203
 Curve calculations window 203
 curve operations 197
 Cycles 80, 93, 101
 Cyclic Polarization 221
 Data tab 176
 DC current 63
 decimal separator 10
 differintegral 192, 195
 Digital Input 209
 Digital Output 209
 double layer capacitance 153
 Dummy Cell 17
 E ac 63, 119
 E begin 54, 63, 65, 72, 119
 E cond1 222
 E cond2 222
 E dc 75, 82, 84, 101, 119
 E end 54, 63, 65, 72, 119
 E equilibration 82
 E level [n] 80
 E peak 214
 E pulse 65, 84
 E start 57, 61
 E step 54, 57, 61, 63, 65, 72, 119
 E step 1 98

E step 2	98	Levels.....	39, 80, 93
E vertex1	57, 61	Linear Polarization	221, 223, 229
E vertex2.....	57, 61	linear regression.....	203
E1(measure).....	86	Linear Sweep Potentiometry	91
E2.....	86	Linear Sweep Voltammetry	54, 223
E3.....	86	Mains Frequency	8
Edit Mode.....	136	Manual Control'.....	209
EIS Plot tab	175	Manual level.....	184
EIS Spectrum Analyzer	129	manually mark peaks	184
Emstat4S	333	Max. frequency.....	119
Enable options for E corrections.....	10	Maxwell	165
equivalent circuits	135, 229	Measure vs OCP.....	119
Excel	207	Measurement tab	117
Exporting curves	207	Min. frequency	119
Extrapolate baseline.....	184	Minimum Peak Height.....	39
fast mode	289	Minimum Peak Width	39
fit.....	135	mode	7, 84
Fit Mode	136	Multiplexer Channel	209
Fixed baseline.....	184	Multiplexer tab	39
fractional derivative order.....	195	n frequencies.....	119
Free baseline	184	noise level.....	51
Frequency	63, 119	Non-linear baseline	184
Frequency type	119	Notes.....	29, 117
Gerischer element.....	153	Number of scans.....	57
i ac.....	123	OCP	222
i applied.....	89, 101	Open Circuit Potential.....	33, 222
i begin	91, 123	Origin	131, 207
i dc.....	123	Overlapping peaks	192
i end	91, 123	overload	22
i forward/reverse	69	Override parameter'	255
i level [n].....	93	overstep	22
i step.....	91, 123	Oxidation or Reduction	39
inductor	153	peak separation window.....	192
Integration	203	Peak window.....	214
IUPAC	10, 179	Peaks	39, 190
ladder	165	Plot tab	175

Polarization resistance	223	t pulse	65, 72, 84
Polarographic convention.....	10, 179	t run	75, 82, 84, 86, 119, 123
Potential or Current	209	t1	86
potential range.....	34, 209	t2	86
printing	179	t3.....	86
Pseudo Polarography	255	Tafel slopes.....	223
real E step.....	29	techniques list	29
Recalculate.....	218	TestSensor.....	17
resistor	153	Three points	184
Sample volume.....	214	time domain	119, 123
scaling	179	underload	22
Scan rate	54, 57, 61, 63, 65, 72	Voigt.....	165
Scan type.....	119, 123	Voltage overload	22
script	243	Voltammetric Analysis.....	213
Settings	6	voltammetric measurement.....	32
simulate.....	135	Warburg element	153
smart scaling	179	warning message (method)	29
Smooth Window.....	39	WE2 Current Range	209
Solution number	214	WE2 Potential.....	209
solution resistance	153	ZView.....	10, 130
Stability.....	119		
Standard addition.....	213, 217		
status.....	22		
status bar.....	22		
stirrer	8		
Stirrer On	209		
stirrer speed	209		
Stripping chronopotentiometry	96		
Stripping Current	96		
sweep rate.....	54		
t [n].....	80, 93		
t 1.....	98		
t 2	98		
t cond1.....	222		
t cond2	222		
t interval	75, 80, 82, 84, 119		