


EmStat^{pico}™

Built with  ANALOG
DEVICES

Electrochemical interface module



Contents

EmStat Pico: Electrochemical Interface Module.....	3
Dual-channel and Bipotentiostat functionality	4
Main Specifications	4
Module pin-out	8
MethodSCRIPT™: EmStat Pico Scripting Language	9
Integrate Electrochemistry into Your Own Applications	10
EmStat Pico Development Board.....	11
PSTrace: research software for Windows	12

➤ See for more information:
www.palmsens.com/pico


EmStat Pico: Electrochemical Interface Module

The EmStat Pico is a joint development by PalmSens BV and Analog Devices Inc. PalmSens BV is known for introducing the first commercially available handheld potentiostat. Over the last decade these have evolved to become smaller and more versatile. Together with Analog Devices, PalmSens now proudly presents the world smallest potentiostat module available on the market.



With a footprint of just 18x30 mm and a height of 3 mm the EmStat Pico can be embedded into virtually any design.

Different power modes and an ultra-low power sleep mode allow the EmStat Pico to be used in wearable applications or for long term monitoring at remote sites.



The EmStat Pico supports MethodSCRIPT™, a scripting language that allows you to let the EmStat Pico do what you want without the hassle of reading into spec sheets.

Supported Techniques

The following electrochemical techniques are supported by the EmStat Pico module.

Voltammetric techniques:

- | | |
|----------------------------------|-----|
| ▪ Linear Sweep Voltammetry | LSV |
| ▪ Cyclic Voltammetry | CV |
| ▪ Square Wave Voltammetry | SWV |
| ▪ Differential Pulse Voltammetry | DPV |
| ▪ Normal Pulse Voltammetry | NPV |

The above techniques can also be used for stripping voltammetry

Techniques as a function of time:

- | | |
|---------------------------------|-----|
| ▪ Chronoamperometry | CA |
| ▪ Pulsed Amperometric Detection | PAD |
| ▪ Open Circuit Potentiometry | OCP |
| ▪ MultiStep Amperometry | MA |

Electrochemical Impedance Spectroscopy

- | | |
|------------------------------------|-----|
| ▪ Scanning or fixed frequency mode | EIS |
|------------------------------------|-----|

MethodSCRIPT™ allows for developing custom techniques. See page 9 for more information.



MethodSCRIPT™
by PalmSens

Dual-channel and Bipotentiostat functionality

The second channel of the EmStat Pico can be used for running sequential measurements on two different cells each with their own Reference, Counter and Working electrodes.

The second channel can also be used in bipotentiostat mode, functioning as second Working Electrode versus the Reference and Counter electrode of channel 1. Both channels are recorded simultaneously in the bipotentiostat mode.

The second Working Electrode (WE2) can either be set at a potential offset with respect to WE1 or at a fixed potential with respect to RE1.

The bipotentiostat mode is supported in Low-Speed mode (see table below) for all techniques, excluding EIS and OCP.

Main Specifications

The module works in three different modes;

Low-Speed mode: for scan rates up to 1 V/s or a bandwidth of 100 Hz.

High-Speed mode: for high scan rates and frequencies.

Max Range mode: a combination of the Low and High-Speed modes for optimal dynamic dc-potential range

General			
	Low-Speed mode	High-Speed mode	Max Range mode
▪ full dc-potential range	-1.2 to +2 V	-1.7 to +2 V	-1.7 to +2 V
▪ dynamic dc-potential range ¹	2.2 V	1.2 V	2.6 V
▪ compliance voltage	-2.0 to +2.3 V ²		
▪ maximum current	±3 mA		
▪ max. acquisition rate (datapoints/s)	100	1000	100
▪ supports FRA/EIS	NO	YES	NO

Potentiostat (controlled potential mode)			
	Low-Speed mode	High-Speed mode	Max Range mode
▪ applied dc-potential resolution	537 μ V	395 μ V	932 μ V
▪ applied potential accuracy	< 0.2% ±1 mV offset	< 0.5% ±1 mV offset	< 0.5% ±1 mV offset
▪ available current ranges	100 nA, 2 μ A, 4 μ A, 8 μ A, 16 μ A, 32 μ A, 63 μ A, 125 μ A, 250 μ A, 500 μ A, 1 mA, 5 mA	100 nA, 1 μ A, 6 μ A, 13 μ A, 25 μ A, 50 μ A, 100 μ A, 200 μ A, 1 mA, 5 mA	

¹ The dynamic range is the range that can be covered during a single scan within the full potential range. For example, a linear scan can start at -1.5 V and end at 1.1 V or vice versa, covering 2.6 V dynamic range.

² The compliance voltage is the maximum potential between Working and Counter electrode and depends on the selected mode.

EmStat Pico - Electrochemical interface module

	Low-Speed mode	High-Speed mode	Max Range mode
▪ current accuracy	< 0.5% of current ±0.1% of range	< 1% of current ±0.1% of range ³	
▪ measured current resolution	0.006% of selected current range (5.5 pA on 100 nA range)		
▪ measured potential resolution (for OCP)	56 µV		

FRA / EIS (impedance measurements) in High-Speed Mode only

▪ frequency range	0.016 Hz to 200 kHz
▪ ac-amplitude range	1 mV to 0.25 V rms, or 0.708 V peak-peak (max. 64 mV for current ranges ≤ 1µA and ≥ 1mA)

Bipotentiostat

▪ modes	1. WE2 at fixed potential (E offset vs RE1) 2. WE2 scanning (E offset vs WE1)
▪ max. potential WE2	$\Delta E(WE1) + \Delta E(WE2) < 1.6 \text{ V}^4$

Electrometer

▪ electrometer amplifier input	> 1 TΩ // 10 pF
▪ bandwidth	250 kHz

Communications and peripherals

▪ module communications	UART
▪ communication with external peripherals	SPI and I ² C
▪ analog I/O	3 analog input pins
▪ digital I/O	7 general-purpose I/O pins 1 wake-up pin
▪ on-board temperature sensor	±0.25 °C (optional) ⁵


³ Channel 2 has an uncompensated series resistor (typical 110 Ω) in series with the WE2 signal. This additional resistance must be taken into account.

⁴ If your main WE1 is scanning from -0.5V to +0.5V, the WE2 can only have a maximum offset of 0.6V.

⁵ The high accurate on-board temperature sensor is standard available on modules that come with the EmStat Pico Development Kit. For separate EmStat Pico modules the temperature sensor is optional.

Other	
▪ data storage	4000 datapoints on board (optional external SD card can directly be connected to Pico for mass storage)
▪ mounting	Surface mounted with castellated pads Through hole pins (2.54 mm pitch)
▪ dimensions	18 x 30 x 2.6 mm
▪ operation temperature range	-40 °C to +85 °C

EIS Accuracy Contour Plot




Note

The accuracy contour plot was determined for the first channel under lab conditions and should be used for reference purposes. Please note that the true limits of an impedance measurement are influenced by all components in the system, e.g. cables, the environment, and the cell.

Limits of some technique specific parameters for EmStat Pico		
Normal Pulse and Differential Pulse Voltammetry	Scan rate:	0.02 mV/s (0.280 mV step) to 5 V/s (10 mV step)
	Pulse time:	1ms to 300ms
Square Wave Voltammetry	Frequency:	1 Hz to 500 Hz
Linear Sweep and Cyclic Voltammetry	Scan rate:	0.02 mV/s (0.280 mV step) to 5 V/s (10 mV step)
Pulsed Amperometric Detection	Interval time:	1 ms to 10 s
	Pulse time:	1 ms to 1 s
	Maximum run time:	1000000 s
ChronoAmperometry and Open Circuit Potentiometry	Interval time:	1 ms to 300 s
	Maximum run time:	1000000 s (> 10 days at 300 s interval)
Multistep Amperometry	Interval time:	1 ms to 300 s
	Level switching overhead time:	±10 ms
	Number of levels:	1 to 255
	Number of cycles:	1 to 20000
	Maximum run time:	1000000 s per level

Module pin-out



All logic levels at 3.3V

MethodSCRIPT™: EmStat Pico Scripting Language


The EmStat Pico potentiostat module works with the new MethodSCRIPT™ scripting language. This language allows developers to program a human-readable script directly into the Pico module by means of a serial (TTL) connection. The simple script language allows for running electrochemical techniques supported by EmStat Pico and makes it easy to combine different measurements and other tasks.

Example MethodSCRIPT for EIS measurement on a test circuit

```
e
#Declare variables
var h
var r
var j
#Initialize device
set_pgstat_mode 3
#Set starting current range
set_cr 1m
#Turn cell on for measurement
cell_on
#Start EIS scan from 200kHz to 2 Hz in 41 steps
meas_loop_eis h r j 10m 200k 2 41 0
#Send results of measurement loop step
pck_start
#Send frequency
pck_add h
#Send Z real
pck_add r
#Send Z imaginary
pck_add j
pck_end
#Continue with next step of EIS scan
endloop
#Turn cell off after measurement
cell_off
```

Scripts can easily be generated and executed in PSTrace for Windows.

See page 12.



Actual measured result on dummy cell ran in Python



MethodSCRIPT is parsed on board the instrument. No DLLs or other type of code libraries are required for using MethodSCRIPT™

MethodSCRIPT features include:

- Use of variables
- (Nested) loops and conditional logic support
- User code during a measurement iteration
- Exact timing control
- Simple math operations on variables (add, sub, mul, div)
- Digital I/O, for example for waiting for an external trigger
- Logging results to internal storage or external SD card
- Reading auxiliary values like pH or temperature
- and many more...



MethodSCRIPT™

➤ See for more information, tutorials and documentation:
www.palmsens.com/methodscript

Integrate Electrochemistry into Your Own Applications

Seamless Instrument Control

- Access all PalmSens potentiostats (single- and multi-channel) through our SDKs.
- Full control of measurement techniques, data acquisition, and real-time analysis.

Cross-Platform Support

- **Python SDK**
Script and automate experiments across platforms.
- **Windows .NET SDK**
Easily integrate in C#, VB.NET, or any .NET language.
- **Android & iOS SDKs**
Build mobile apps to run PalmSens instruments in the field.
- **LabVIEW & MATLAB examples**
Quick start for engineers and researchers.



Accelerate Development

- Pre-built code sample
- Clear documentation & active support
- Sample apps to get started within minutes



PalmSens SDKs
put you in control
from the lab to the field



MATLAB




LabVIEW



➤ See for more information:
www.palmsens.com/dev


EmStat Pico Development Board

The EmStat Pico Development board allows to run your experiments conveniently in our PSTrace software for electrochemistry.



STEP 1

Connect the EmStat Pico Development Board to a PC running PSTrace



STEP 2

Fine-tune your electrochemistry for optimal use of the EmStat Pico module

```
set_pgstat_chan 0
set_pgstat_mode 3
set_max_bandwidth 200
set_pot_range 0 0
set_cr 850n
set_autoranging 850n 850n
cell_off
meas_loop_ocr o 500m 3
pck_start
pck_add o
pck_end
endloop
cell_on
core_var b 0 ab
var b o
```



STEP 3

Generate the MethodSCRIPT™ snippet for running your measurement on the EmStat Pico



STEP 4

Use the MethodSCRIPT™ snippet to run the exact same measurement on the embedded EmStat Pico in your product


Comes with code examples for:



> See for more information:
www.palmsens.com/picodevkit


PSTrace: research software for Windows

The EmStat Pico Development Board can be used directly with the PSTrace software for Windows. PSTrace automatically sets the EmStat Pico in the optimal mode based on the user specified method parameters.



Other functions in PSTrace 5

- Method validation
- Automatic peak search
- Equivalent Circuit Fitting
- Open your data in Origin and Excel with one click of a button
- Save all available curves, measurement data and methods to a single file
- Dynamic feedback on method parameters



Integration with third party software:

- Excel
- Origin
- Matlab
- ZView



Minimum System Requirements

- Windows 7, 8, 10 or 11
- 1 GHz or faster 32-bit (x86) or 64-bit (x64) processor
- 2 GB RAM (32-bit) or 4 GB RAM (64-bit)
- Screen resolution of 1280 x 800 pixels

➤ See for more information:
www.palmsens.com/pstrace

PalmSens BV has more than 50 distributors around the world.
Please contact us at info@palsens.com or go to our website to
get in touch with a distributor in your region.



Please do not hesitate to contact PalmSens for more details: info@palsens.com

PalmSens BV
The Netherlands
www.palsens.com

DISCLAIMER

Changes in specifications and typing errors reserved.
Every effort has been made to ensure the accuracy of this document. However, no rights can be claimed by the contents of this document.

