



- 1) 2*Counter Electrode
- 2) 6*O ring
- 3) 1*PTFE ring
- 4) 2*ThreadedAdapter for T24/40 Joint
- 5) 1*Reference Electrode
- 6) 1*Reference Electrode BridgeTube
- 7) 5*Reference electrode bridge core
- 8) 1*Purge Tube
- 9) 1*Electrode Mounting Rod
- 10) 1*Screw nut
- 11) 1*PTFE pad
- 12) 1*Working Electrode Hoder
- 13) 5*Replacement Teflon Gaskets
- 14) 1*Working Electrode tip
- 15) 1*1000ml Corrosion Flask
- 16) 1*Ball clip
- 17) 1*Flat Specimen Holder Kit

IS-CCS Corrosion Cell Set

17) Flat Specimen Holder Kit

