


EmStat4M™ DEVELOPMENT BOARD TOP VIEW



CD-ES4MDEV-11/21-A

EmStat4M™ DEVELOPMENT BOARD BOTTOM VIEW

Arduino Headers



Compatible with
Arduino MKR Series

AREF	+5V
DAC0/A0	VIN
A1	+3V3
A2	GND
A3	RESET
A4	TX
A5	RX
A6	SCL
~0	SDA
~1	MISO
~2	SCK
~3	MOSI
~4	7
~5	6

MKR_Trig

SW6

ON

1 2

MKR to ES

RX



TX

D4

CON10

JST PH (2.0mm)
Li-Ion Battery Connector

Battery Power



RTC Battery Compartment

CR1220 3V
Real Time Clock
backup source

Development Board Manual

A manual for the EmStat4M Development Board is available for more detailed information and instructions. It can be found in the Download section on:
www.palmsens.com/es4m-dev

