



**DETAIL A**  
SCALE 2 : 1



SPE Requirements:  
 Max. width: 11 mm  
 Thickness: between 0.1 mm and 0.8 mm  
 Connection pads pitch: 2.54 mm



UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN MILLIMETERS	DO NOT SCALE DRAWING	REVISION	17-04-2024
	PalmSens BV		
	TITLE: MUX8-R2 SPE breakout		
	DWG NO.	Sheet1	A4
	SCALE:1:1	SHEET 1 OF 1	