


# DETAIL A

SCALE 2 : 1



RE WE CE

SPE Requirements:

Max. width: 11 mm

Thickness: between 0.1 mm and 0.8 mm

Connection pads pitch: 2.54 mm

UNLESS OTHERWISE SPECIFIED:  
DIMENSIONS ARE IN MILLIMETERS

DO NOT SCALE DRAWING

REVISION 22-04-2024

PalmSens BV

TITLE:

MUX8-R2 SPE breakout

DWG NO.

Sheet1

A4

SCALE:1:1

SHEET 1 OF 1