



# OG Carbon 2.0

Screen Printed Carbon Electrodes

Droplet placed here

Plug in here



Counter Electrode  
Reference Electrode  
Working Electrode



*Disposable and high quality screen printed carbon electrodes with good reproducibility. Perfect for biosensor research and development.*

- Carbon (graphite) counter and working electrode
- Silver/silver chloride reference electrode
- 4 mm working electrode diameter
- Flexible PET support material
- 10 mm x 33 mm x 0.25 mm
- Compatible with a range of potentiostats and connectors


Easily cuttable



Available in rows of 25



Cyclic voltammometry of 3 OG Carbon 2.0 electrodes in  $\text{Fe}(\text{CN})_6^{3-/4-}$



Standard deviation of 3 OG Carbon 2.0 electrodes vs 3 competitor electrodes

More information available at [ogcarbon.co.uk](http://ogcarbon.co.uk)